Prostar XL

Percutaneous Vascular Surgical System

PPL2106464 (2015-09-18)

CE
0086

INSTRUCTIONS FOR USE

English	2	Latvian (Latviski)	74
Bulgarian (Български)	7	Lithuanian (Lietuviškai)	79
Czech (Česky)	14	Norwegian (Norsk)	85
Danish (Dansk)	20	Polish (Polski)	90
Dutch (Nederlands)	25	Portuguese (Português)	96
Estonian (Eesti keel)	32	Romanian (Română)	102
Finnish (Suomi)	37	Russian (Русский)	109
French (Français)	43	Slovakian (Slovensky)	115
German (Deutsch)	49	Slovenian (Slovenščina)	121
Greek (Ελληνικά)	55	Spanish (Español)	127
Hungarian (Magyar)	62	Swedish (Svenska)	133
Italian (Italiano)	68	Turkish (Türkçe)	139

ENGLISH

**PROSTAR XL
PERCUTANEOUS VASCULAR SURGICAL (PVS) SYSTEM
INSTRUCTIONS FOR USE**

TO ENSURE PROPER DEPLOYMENT AND USE OF THIS DEVICE AND TO PREVENT INJURY TO PATIENTS, READ ALL INFORMATION CONTAINED IN THESE INSTRUCTIONS FOR USE.

CAUTION

This device should only be used by physicians (or allied healthcare professionals, authorized by, or under the direction of such physicians) who are trained in diagnostic and therapeutic catheterization procedures and who have been trained by an authorized representative of Abbott Vascular.

Prior to use, the operator must review the Instructions for Use and be familiar with the deployment techniques associated with the use of this device.

DEVICE DESCRIPTION

The Prostar XL Percutaneous Vascular Surgical (Prostar XL PVS) device is designed to deliver polyester suture(s) to close femoral artery puncture sites following catheterization procedures. The Prostar XL device has two sutures and four needles.


The Prostar XL device consists of a sheath, which contains two pairs of sutured needles, a needle guide that precisely controls the placement of needles around the puncture site and a rotating barrel that receives deployed needles. The Prostar XL sheaths have J-shaped tips at the distal end.

A marker lumen is contained within the barrel of the device with the intra-arterial port of the lumen positioned in the needle guide.

Proximally, the marker lumen exits from the hub of the device. The marker lumen allows a pathway for back-bleeding from the femoral artery and ensures proper device positioning. The barrel rotates independently from the central core and is designed to prepare the subcutaneous track. Barrel rotation is accomplished by depressing the interlocks exiting from the hub. The Prostar XL devices track over a standard 0.038" (0.97 mm) (or smaller) guide wire.

The Prostar XL PVS System consists of a 10F Prostar XL device and a Perclose Knot Pusher. The Perclose Knot Pusher is designed to advance the tied suture knot to the arteriotomy.

The Prostar XL device and the Perclose Knot Pusher are depicted in **Figure 1**.

Figure 1**Prostar XL Percutaneous Vascular Surgical System**

- A. Prostar device
B. Perclose Knot Pusher

INDICATIONS FOR USE

The Prostar XL Percutaneous Vascular Surgical Systems are intended for percutaneous delivery of sutures for closing the common femoral artery access site, reducing the time to hemostasis and ambulation (patient walks ten feet/ three meters) for patients who have undergone catheterization procedures.

The Prostar XL PVS System reduces the time to discharge in patients having undergone diagnostic catheterization

procedures without complicating clinical conditions (refer to PRECAUTIONS and SPECIAL PATIENT POPULATIONS).

The Prostar XL PVS System is designed for use in conjunction with 8.5 to 24F sheaths.

CONTRAINDICATIONS

There are no known contraindications to the use of this device. Attention is drawn to the Warnings, Precautions, and Special Patient Populations.

WARNINGS

Do not use the Prostar XL PVS device or accessories if the packaging or sterile barrier has been previously opened or damaged, or if the components appear to be damaged or defective.

DO NOT RESTERILIZE OR REUSE. The Prostar XL PVS device and accessories are intended for single use only.

Do not use the Prostar XL PVS System if the sterile field has been broken where bacterial contamination of the sheath or surrounding tissues may have occurred, since such a broken sterile field may result in infection.

Do not use the Prostar XL PVS System if the puncture site is located above the most inferior border of the inferior epigastric artery (IEA) and / or above the inguinal ligament based upon bony landmarks, since such a puncture site may result in a retroperitoneal hematoma. Perform a femoral angiogram to verify the location of the puncture site.

Do not use the Prostar XL PVS System if the puncture is through the posterior wall or if there are multiple punctures, since such punctures may result in a retroperitoneal hematoma.

Do not use the Prostar XL PVS System if the puncture site is located in the superficial femoral artery or the profunda femoris artery, since such puncture sites may result in a pseudoaneurysm, intimal dissection, or an acute vessel closure (thrombosis of small artery lumen). Perform a femoral angiogram to verify the location of the puncture site.

PRECAUTIONS

- The Prostar XL PVS System is provided sterile and non-pyrogenic in unopened undamaged packages. Products are sterilized with ethylene oxide and intended for single use only. Do not resterilize. Store in a cool, dry place.
- Prior to use, inspect the Prostar XL PVS System to ensure that the sterile packaging has not been damaged during shipment. Examine all components prior to use to verify proper function. Exercise care during device handling to reduce the possibility of accidental device breakage.
- As with all catheter-based procedures, infection is a possibility. Observe sterile technique at all times when using the Prostar XL PVS System. Employ appropriate groin management, as per hospital protocol, post procedure and post hospital discharge to prevent infection.
- Use a single wall puncture technique. Do not puncture the posterior wall of the artery. Avoid posterior wall suture placement.
- Do not insert the Prostar XL PVS device into the femoral artery at an angle greater than 45 degrees to the longitudinal plane of the artery.
- An experienced operator should use accepted surgical technique as warranted by surgical circumstances to provide adequate knot security.
- There are no reaccess restrictions if previous arteriotomy repairs were achieved with an Abbott Vascular Suture Mediated Device.
- Do not advance or withdraw the Prostar XL device against resistance until the cause of that resistance has been determined** (see Prostar XL Device Placement section). **Excessive force used to advance or torque the Prostar XL PVS device should be avoided, as this may lead to significant vessel**

damage and / or breakage of the device, which may necessitate intervention and / or surgical removal of the device and vessel repair.

9. If excessive resistance in advancing the Prostar XL PVS device is encountered, withdraw the Prostar XL device over a 0.038" (0.97 mm) (or smaller) guide wire and reinsert the introducer sheath, or use conventional compression therapy.
10. In the event suture breakage occurs after an initial knot has been tied, care should be taken to avoid excessive force if the reintroduction of another Prostar XL device or introducer sheath is required. Any resistance to introduction should result in advancement of an introducer sheath small enough to be introduced without undue force.
11. If significant blood flow is present around the barrel of the Prostar XL device, do not deploy the needles. Remove the Prostar XL device over a 0.038" (0.97 mm) (or smaller) guide wire and insert an appropriately sized introducer sheath.
12. Remove the Prostar XL sheath before tightening the suture. Failure to remove the sheath prior to tightening the suture may result in detachment of the tip of the sheath.
13. In using this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing damage due to application of surgical instruments such as forceps or needle holders.
14. Do not attempt to re-deploy Prostar XL needles after the needles have been "backed-down" into the sheath (refer to the **TECHNIQUE FOR NEEDLE BACK-DOWN** section).
15. If using the "pre-close" technique (placement of the Prostar XL device sutures prior to dilating the access site beyond 10F), sutures should be gradually tightened as the introducer sheath is removed to maintain hemostasis.
16. Use conventional compression methods in the event bleeding from the femoral access site persists, after the use of the Prostar XL PVS device.
17. The Prostar XL device should be advanced using fluoroscopic guidance to prevent possible vessel damage.

SPECIAL PATIENT POPULATIONS

The safety and effectiveness of the Prostar XL PVS System have not been established in the following patient populations:

- Patients with introducer sheaths < 8.5F or > 24F during the catheterization procedure.
- Patients with ipsilateral arterial access sites punctured and compressed within 48 hours of closure.
- Patients having arterial access in vessels other than the common femoral artery.
- Patients having a hematoma, pseudoaneurysm or arteriovenous fistula present prior to sheath removal.
- Patients with femoral artery calcium, which is fluoroscopically visible at access site.
- Patients with small femoral arteries (< 5 mm in diameter).
- Patients with severe claudication, iliac or femoral artery diameter stenosis greater than 50% or previous bypass surgery or stent placement in the vicinity of access site.
- Patients with access sites in vascular grafts.
- Patients with prior intra-aortic balloon pump at access site at any time prior.
- Patients with ipsilateral femoral venous sheath during the catheterization procedure.
- Patients with whom there is difficulty inserting the introducer sheath or greater than one ipsilateral arterial puncture at the start of the catheterization procedure.

- Patients with intra-procedural bleeding around access site.
- Patients receiving glycoprotein IIb/IIIa inhibitors before, during, or after the catheterization procedure.
- Patients younger than 18 years of age.
- Patients who are pregnant or lactating.
- Patients with bleeding diathesis or coagulopathy.
- Patients who are morbidly obese where less than one third of the access needle is above the skin line.
- Patients with active systemic or cutaneous infection or inflammation.
- Patients with access sites above the most inferior border of the inferior epigastric artery (IEA) and / or above the inguinal ligament based upon bony landmarks.
- Patients with antegrade punctures.

CLINICAL PROCEDURE

The following instructions provide technical direction but do not obviate the necessity of formal training in the use of Prostar XL PVS System. The techniques and procedures described below are not intended as a substitute for the operator's experience and judgment in treating any specific patients.

Examination and Selection of Products

1. The outer pouch packaging of the Prostar XL PVS System provides the sterile barrier. After careful inspection of the packaging to ensure there is no damage to the sterile barrier, remove the device from the package.
2. Exercise care during device handling to reduce the possibility of accidental device breakage.
3. Exercise care when using additional surgical instruments such as forceps, hemostats or needle holders during device handling, to reduce the possibility of accidental device breakage.
4. Verify marker port patency by flushing the lumen with saline until the saline exits from the marker port. **Do not use the Prostar XL PVS device if the marker lumen is not patent.**

Arterial Puncture Considerations

1. Puncture the anterior wall of the common femoral artery at an angle of approximately 45 degrees.
2. Avoid sidewall or posterior wall femoral artery punctures.
3. Puncture locations are ideally located in the common femoral artery below the level of the inguinal ligament and above the common femoral artery bifurcation.

Prostar XL Device Placement

The following instructions detail the deployment sequence to close the access site of a catheterization procedure performed through an applicable sized sheath.

When utilizing a "pre-close" technique (placement of the Prostar XL device sutures prior to dilating the access site beyond 10F) refer to step 7 in the Suture Management section for more information.


1. Prior to Prostar XL device placement, perform a femoral angiogram through the introducer sheath to verify that the access site is in the common femoral artery.
2. Evaluate the femoral artery site for size, calcium deposits, and tortuosity, to avoid posterior wall suture placement and possible ligation of the anterior and posterior walls of the femoral artery.
3. Re-prep the access site, placing clean towels around the access site, and wearing new sterile gloves prior to handling the device and proceeding with the closure procedure.

While the introducer sheath remains in place, use a scalpel to slightly extend the incision and forceps to dilate subcutaneous tissue.

5. Place a 0.038" (0.97 mm) (or smaller) guide wire through the introducer sheath. Remove the introducer sheath while applying pressure on the groin to maintain hemostasis.
6. Carefully back-load the Prostar XL device over the guide wire until the guide wire exit port is just above the skin line. Remove the guide wire. Continue to advance the Prostar XL device until the barrel is at skin level.
7. Unlock the hub by depressing the interlocks with thumb and forefinger. Once the hub is unlocked, rotate the hub while gently advancing the barrel at a **45-degree angle**, or less.
8. A steady, continuous drip of blood from the dedicated marker lumen occurs when the Prostar XL device is properly positioned (see **Figure 2**). Marking from the lumen(s) containing suture(s) may occur but should not be used as the indicator for proper positioning and needle deployment.

Do not clamp the suture lumen with a hemostat or other instrument. Doing so will prevent suture deployment.

Figure 2


9. Lock the hub back in place.
 - 9.1. If a continuous drip of blood (luminal marking) from the dedicated marker lumen is not apparent, withdraw the device to expose the marker port. Flush the marker lumen to verify patency, and then continue to gently advance the Prostar XL device while rotating the barrel.
 - 9.2. If continuous marking is still not achieved, remove the Prostar XL device and follow conventional compression protocol, or replace the Prostar XL device with an appropriately sized introducer sheath.
 - 9.3. Do not deploy needles until a continuous drip of blood is evident from the dedicated marker lumen.

Needle Deployment


1. Confirm that the interlocks are re-engaged and correctly aligned with the locking indents in the hub (locked).
2. With the left hand, hold the hub of the device in position at an **angle of 45 degrees** (or less).
3. With the right hand, rotate the handle counterclockwise, to unlock the handle.
4. Ensure blood marking is maintained.
5. Pull the handle away from the hub to deploy the needles (see **Figure 3**).

Figure 3

6. If resistance to rotating the handle is encountered, **do not attempt to deploy the needles**. Significant resistance is an indication that the hub is not properly positioned. Refer to the **TECHNIQUE FOR NEEDLE BACK-DOWN** section to ensure proper hub positioning.
7. **Do not clamp the suture lumen with a hemostat or other instrument. Doing so will prevent suture deployment.**
8. Continue to pull the handle until the needle tips emerge at the top of the barrel.
9. While steadily holding the device in position, confirm that all four needles are visible in the hub.
10. If significant resistance is encountered prior to needle tips emerging at the top of the barrel, or if all four of the needles are not deployed, terminate deployment. Refer to the **TECHNIQUE FOR NEEDLE BACK-DOWN** section to perform needle back-down procedure.
11. Using a hemostat, remove the posterior needle(s) followed by the anterior needle(s) using the funnel-shaped hub as a fulcrum to facilitate needle removal (see **Figure 4**).


Figure 4

Suture Management

1. Once needles have been removed from the hub, remove slack in the suture(s) by pulling the suture ends to evenly matched lengths and tensioning until resistance is felt. Cut the suture ends close to the needles. Dispose of needles in accordance with hospital policy.
2. Withdraw the Prostar XL device while maintaining access to the sutures exiting the hub during this step. Create a "bow string" effect with the exposed suture(s) by bending the Prostar XL device sheath away from the operator and applying tension to the suture ends exiting the hub (see **Figure 5**). This ensures that the sutures are not tied around the sheath.

Figure 5

3. Grab the anterior suture ends (one white and one green) adjacent to the topside of the sheath and pull the suture ends through the distal end of the barrel. Place this pair of suture ends toward the patient's head.
4. Grab the posterior suture ends (one white and one green) adjacent to the underside of the sheath and pull the suture ends through the distal end of the barrel. Ensure that this pair of suture ends comes under the device and towards you. Place this pair of suture ends toward the patient's feet.
5. Identify the two ends of a single suture by color, and then tension the ends using a gentle seesawing motion.
6. Continue to withdraw the Prostar XL device until the guide wire exit port exits the skin line. Reinsert a 0.038" (0.97 mm) (or smaller) guide wire into the visible guide wire exit port so that guide wire access can be maintained until hemostasis is verified. If a guide wire straightener is used to reinsert the guide wire, care must be taken to avoid disturbing the hemostasis valve located in the device just distal to the guide wire exit port.

Access sites larger than 10F require the use of the "pre-close" technique.

7. When utilizing a "pre-close" technique, the Prostar XL device is exchanged for an appropriately sized introducer sheath or an additional Prostar XL device if additional sutures are to be placed around the same access site.
 - Sliding, self-locking surgical knots should be tied prior to introducer sheath removal according to the Knot Advancement section below.

Access sites larger than 10F may require the use of an additional Prostar XL device.

8. Use of an additional Prostar XL device:
 - 8.1 If additional sutures are to be placed, position an additional Prostar XL device rotated 45 degrees relative to the first Prostar XL device to establish alternating suture positions and repeat all the steps above beginning at **Prostar XL Device Placement**.
 - 8.2 Identify and secure the sutures from the first device separately from the sutures of the additional device.

Knot advancement

1. Identify the green suture ends. Place one green suture end in each hand.
2. Gently tension the suture to create a longer suture end. (The short end will be your rail limb of suture, for the green knot).
3. Tie a sliding, self-locking surgical knot using the green suture.
4. Lay the green suture aside.
5. Identify the white suture ends. Place one white suture end in each hand.
6. Gently tension the suture to create a longer suture end. The short end will be your rail limb of suture, for the white knot.
7. Tie a sliding, self-locking surgical knot using the white suture.
8. Securely wrap the white rail suture around your left index finger. Saturate the suture with saline. Gently pull on the white rail suture first, keeping the suture coaxial to the tissue tract. Completely remove the Prostar XL device from the artery, leaving the guide wire in the artery.

Note: If using the "pre-close" technique, sutures should be gradually tightened as the introducer sheath is removed to maintain hemostasis.

- 8.1 Do not compress the femoral access site while withdrawing the Prostar XL device from the tissue.

tract. Applying excessive force to the Prostar XL sheath, due to site compression, during removal of the device may result in a needle guide break. Please refer to **NEEDLE GUIDE BREAK PROCEDURE** section for further instruction.

- 8.2 Do not tighten the suture around the Prostar XL sheath. Tightening the suture around the Prostar XL sheath during device removal may result in a needle guide break. Please refer to **NEEDLE GUIDE BREAK PROCEDURE** section for further instructions.

9. Saturate the sutures with saline. Load the rail suture into the Knot Pusher and continue to advance the knot forward to the arteriotomy (see **Figure 6**).

Figure 6


10. With the rail suture limb securely wrapped around the left forefinger, place the Knot Pusher under the left thumb to assume a single-handed position to complete knot advancement.
11. Do not tension the non-rail suture to tighten the knot.
12. Remove the Knot Pusher from the tissue tract and from the suture without using the thumb knob.
13. Lay the white suture down, being sure to separate the rail and non-rail ends.
14. Repeat steps 8–13 with the green suture.
15. Repeat steps 1–14 for second set of sutures if additional device is used.
16. If hemostasis has been achieved, remove the guide wire from the patient. Continue with step 18.
17. If hemostasis has not been achieved:
 - 17.1 Repeat steps 9–14.
 - 17.2 Exchange the Prostar XL device with an appropriately sized introducer sheath. Care should be taken to avoid excessive force if the reintroduction of another introducer sheath is required. To avoid resistance, use an introducer sheath small enough to be introduced without undue force, but large enough to maintain hemostasis.
 - 17.3 Apply conventional compression therapy.
18. Once again, securely wrap the white rail suture around your left index finger.
19. Load the white rail suture into the Knot Pusher and advance to the level of the arteriotomy.
20. With the rail suture limb securely wrapped around the left forefinger, place the Knot Pusher under the left thumb to assume a single-handed position and to complete knot advancement.
21. With the Knot Pusher in place, tighten the knot by gently pulling on the white non-rail suture.
22. Remove the Knot Pusher from the tissue tract and from the suture without using the thumb knob.
23. Repeat these same steps 18–22 with the green suture.
24. If hemostasis is not complete, repeat the single-handed position for 20 seconds. After 20 seconds gently pull on the non-rail suture to tighten the knot. (Repeat first with the white suture and next with the green suture). Do not apply excessive pressure to the Knot Pusher or suture.
25. Once hemostasis is achieved, trim the sutures below the skin.

TECHNIQUE FOR NEEDLE BACK-DOWN

The following describes a safety feature ("needle back-down") that permits the physician to return the needles into the sheath. This feature provides the option of exchanging the Prostar XL device with another Prostar XL device or an introducer sheath so that the patient may be treated with conventional compression therapy.

1. If resistance to counterclockwise rotation of the handle is felt, **do not attempt to deploy the needles**. Check that the interlocks are correctly aligned with the locking indents in the hub. The star on the edge of the hub should be centered directly between the interlocks.
2. If the needles are not easily deployed, back the needles down into the sheath prior to removal of the device.
3. Manually remove the coiled sections of the suture lumens to uncover the suture loops.
4. Use a hemostat to grasp the pull rod close to the hub. Advance the pull rod 1 cm into the core (see **Figure 7**) and then gently pull the slack out of the uncovered suture loops.

Figure 7

5. Repeat step 4 until the handle snaps into position on the proximal end of the device.
6. Pull back on the suture loops exiting the suture lumens to ensure that all slack in the sutures has been removed.
7. Prior to device removal, use fluoroscopy to verify the needle has been returned to the needle guide. **The needle tips should be as close as possible to the proximal edge of the radiopaque sheath ring before removal of the Prostar XL device.**
8. Perform the same "needle back-down" procedure in the event that all needles are not deployed. **DO NOT REMOVE ANY DEPLOYED NEEDLES.**
9. Do not attempt to re-deploy the Prostar XL device after the needles have been "backed-down". Replace the Prostar XL device with another Prostar XL device, an introducer sheath or use conventional compression therapy.

SUTURE BREAKAGE

1. If suture breakage occurs prior to completing the initial knot, discard the suture material and remove the device over the guide wire. Use another Prostar XL PVS device to complete the procedure or replace the introducer sheath.
2. If the sutures are inadvertently tangled or removed prior to knot tying, discard the suture material and remove the Prostar XL device over the guide wire. Use another Prostar XL device to complete the procedure.
3. If suture breakage occurs after an initial knot has been tied, another Prostar XL device may be introduced or the introducer sheath may be replaced. If replacement of the introducer sheath is attempted, care should be taken to avoid excessive force during introduction. Any resistance to introduction should result in advancement of an introducer sheath small enough to be introduced without undue force, but large enough to maintain hemostasis.

NEEDLE GUIDE BREAK PROCEDURE

The following describes a safety procedure that allows the physician to remove the Prostar XL device from the patient without surgical intervention, in the event that the Needle Guide of the Prostar XL device breaks.

Procedure for Needle Guide Break pre needle deployment:

1. Do not advance or withdraw the Prostar XL device against resistance until the cause of that resistance has been determined. **Advancing or withdrawing the Prostar XL device against resistance may cause the device to break.**
2. **Note:** though the guide will appear to be broken, the sheath remains attached to the device via the holder stop. The Prostar XL device has a redundant system that prevents the sheath from becoming completely detached from the device.
3. To remove a Prostar XL device with a Needle Guide Break that occurred during device insertion, evaluate if the Guide Wire Exit Port can be visualized above the skin.
4. If the Guide Wire Exit Port cannot be visualized:
 - 4.1. Manually remove the coiled sections of the suture lumens to uncover the suture loop.
 - 4.2. Pull back on the suture loops exiting the suture lumens to ensure that all slack in the sutures has been removed.
 - 4.3. Clamp the suture lumens with hemostats.
 - 4.4. Remove the device until the Guide Wire Exit Port is visualized.
5. Insert a 0.038" (0.97 mm) (or smaller) guide wire into the Guide Wire Exit Port and replace the Prostar XL device with another Prostar XL device, an introducer sheath, or use conventional compression therapy. If a guide wire straightener is used to reinsert the guide wire, care must be taken to avoid disturbing the hemostasis valve located in the device just distal to the guide wire exit port.

Procedure for Needle Guide Break post needle deployment:

1. Do not compress the femoral access site while withdrawing the Prostar XL device from the tissue tract. **Applying excessive force to the Prostar XL sheath, due to site compression, during removal of the device may result in a needle guide break.**
2. Do not tighten the suture around the sheath. **Tightening the suture around the sheath during Prostar XL device removal may result in a needle guide break.**
3. Do not advance or withdraw the Prostar XL device against resistance until the cause of that resistance has been determined. **Advancing or withdrawing the Prostar XL device against resistance may cause the device to break.**
4. **Note:** though the guide will appear to be broken, the sheath remains attached to the device via the holder stop. The Prostar XL device has a redundant system that prevents the sheath from becoming completely detached from the device.
5. To remove the Prostar XL device when a Needle Guide Break has occurred during device removal, evaluate if the Guide Wire Exit Port can be visualized above the skin.
 - 5.1. If the Guide Wire Exit Port cannot be visualized, apply similar and simultaneous backward tension to both the circular handle and the device hub.
 - 5.2. Remove the device until the Guide Wire Exit Port is visualized.

- Insert a 0.038" (0.97 mm) (or smaller) guide wire into the Guide Wire Exit Port and replace the Prostar XL device with another Prostar XL device, an introducer sheath, or utilize conventional compression therapy. If a guide wire straightener is used to reinsert the guide wire, care must be taken to avoid disturbing the hemostasis valve located in the device just distal to the guide wire exit port.

POST PROCEDURE PATIENT MANAGEMENT

- Apply an appropriate dressing to the puncture site.
- Assess the insertion site as per hospital protocol.

COMPLICATIONS

Potential complications resulting from procedures associated with the use of the Prostar XL PVS System include, but are not limited to:

- Localized vessel wall trauma which may lead to surgical repair
- Arterial thrombus
- Deep vein thrombosis
- Pseudoaneurysm
- Hematoma
- Local infection
- Nerve injury
- Local pulse deficits or ischemia
- Blood loss that may lead to blood transfusion
- Local discomfort
- Wound dehiscence.

RECOMMENDATION FOR PATIENT AMBULATION AND DISCHARGE

- Patients may be ambulated after Prostar XL device procedures with normalized (150–180 seconds) Activated Clotting Times (ACTs).
 - Patients who have undergone Prostar XL device procedures may be eligible for earlier discharge from the hospital relative to conventional compression (i.e., manual or mechanical methods).
 - Before considering early discharge, assess for the following clinical conditions:
 - Conscious sedation
 - Anticoagulation, thrombolytic or antiplatelet therapy
 - Unstable cardiac status
 - Hematoma at the closure site
 - Hypotension
 - Pain while walking
 - Bleeding at the closure site
 - Any comorbid condition requiring observation.
- The presence of any of the above factors has generally led to the deferral of early discharge recommendations.

PRODUCT INFORMATION DISCLOSURE

Abbott Vascular Inc. has exercised reasonable care in the manufacture of this device. Abbott Vascular Inc. excludes all warranties, whether expressed or implied, by operation of law or otherwise, including but not limited to, any implied warranties of merchantability or fitness, since handling and storage of this device, as well as factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond the control of Abbott Vascular Inc. directly affect this device and the results obtained from its use. Abbott Vascular Inc. shall not be liable for any incidental or consequential loss, damage, or expense, directly or indirectly arising from the use of this device. Abbott Vascular Inc. neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device.

HOW SUPPLIED

Prostar XL Percutaneous Vascular Surgical System

Each System Includes:

- One (1) Prostar XL Percutaneous Vascular Surgical device
- One (1) Perclose Knot Pusher

The Prostar XL PVS device and accessories are provided sterile and non-pyrogenic in unopened undamaged packages. Products are sterilized with ethylene oxide and intended for single use only. This single use device cannot be reused on another patient, as it is not designed to perform as intended after the first usage. Changes in mechanical, physical, and / or chemical characteristics introduced under conditions of repeated use, cleaning, and / or resterilization may compromise the integrity of the design and / or materials, leading to contamination due to narrow gaps and / or spaces and diminished safety and / or performance of the device. Absence of original labeling may lead to misuse and eliminate traceability. Absence of original packaging may lead to device damage, loss of sterility, and risk of injury to the patient and / or user. Do not resterilize. Store in cool, dry place.

PATENTS AND TRADEMARKS

This product and / or its use may be covered by one or more of the following United States Patents: 5,417,699; 5,527,322; 5,613,974; 5,746,755; 5,779,719; 5,792,152; 5,797,929; 5,860,991; 5,902,311; 5,921,994; 6,036,699; 6,117,145; 6,206,893; 6,355,050; 6,517,553; 7,029,480; 7,445,626; EP 673228; EP 727965; EP 910288; EP 1236437. Other U.S. patents pending. Foreign patents issued and pending.

Prostar and Perclose are trademarks of the Abbott Group of Companies.

Please refer to the end of this booklet for Graphical Symbols for Medical Device Labeling.

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©2015 Abbott