
# Perclose ProGlide 6F

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RELEASED

## Perclose ProGlide 6F Suture-Mediated Closure (SMC) System

RELEASED

### INSTRUCTIONS FOR USE

#### Table of Contents

1.0	CAUTION
2.0	DEVICE DESCRIPTION
	Figure 1: Perclose ProGlide Suture-Mediated Closure System
3.0	HOW SUPPLIED
4.0	INDICATIONS
5.0	CONTRAINDICATIONS
6.0	WARNINGS
7.0	PRECAUTIONS
8.0	SPECIAL PATIENT POPULATIONS
9.0	POTENTIAL ADVERSE EVENTS
10.0	THE PERCLOSE PROGLIDE SMC SYSTEM CLINICAL PROCEDURE
10.1	Examination and Selection of Products
10.2	Arterial Site and Puncture Considerations
10.3	SMC Device Placement 5F–8F Sheath, including Optional Pre-Close and Maintaining Wire Access Techniques
10.3.1	Optional: Maintaining Wire Access During Knot Advancement (Closing over the wire)
10.3.2	Optional: Pre-Close Technique
10.4	SMC Device Placement 8.5F–21F Sheath, Utilizing Pre-Close and Maintaining Wire Access Techniques
10.5	Suture Breakage
10.6	Post-Procedure Patient Management
10.7	Recommendation for Patient Ambulation and Discharge
11.0	PRODUCT INFORMATION DISCLOSURE
12.0	PATENTS AND TRADEMARKS

### TO ENSURE PROPER DEPLOYMENT AND USE OF THIS DEVICE AND TO PREVENT INJURY TO PATIENTS, READ ALL INFORMATION CONTAINED IN THESE INSTRUCTIONS FOR USE.

#### 1.0 CAUTION

This device should only be used by physicians (or allied healthcare professionals, authorized by, or under the direction of, such physicians) who are trained in diagnostic and / or interventional catheterization procedures and who have been trained by an authorized representative of Abbott Vascular. Prior to use, the operator must review the Instructions for Use and be familiar with the deployment techniques associated with the use of this device.


#### 2.0 DEVICE DESCRIPTION

The Perclose ProGlide Suture-Mediated Closure (SMC) System is designed to deliver single monofilament polypropylene suture to close femoral artery puncture sites following diagnostic or interventional catheterization procedures. This Perclose ProGlide SMC device is composed of a plunger, handle, guide, and sheath. The Perclose ProGlide tracks over a standard 0.038" (0.97 mm) (or smaller) guide wire. A hemostasis valve restricts the blood flow through the sheath with or without the guide wire in place. The guide houses the needles, and the foot, and precisely controls the placement of these needles around the puncture site. The handle is used to stabilize the device during use. The plunger advances the needles and is used to retrieve the suture. A marker lumen is contained within the guide, with the intraluminal port of the lumen positioned at the distal end of the guide. Proximally, the marker lumen exits from the body of the device. The marker lumen allows a pathway for back-bleeding (obtaining mark) from the femoral artery to ensure proper device positioning.

Knot pusher accessory (Perclose Suture Trimmer) is included, and is designed to position the tied suture knot to the top of the arteriotomy. The Perclose Suture Trimmer is also designed to trim the trailing limbs of suture.

The Perclose ProGlide 6F SMC System is designed for use in 5F to 21F access sites. **The Perclose ProGlide SMC System is depicted in Figure 1.**

Figure 1: Perclose ProGlide Suture-Mediated Closure System

#### 3.0 HOW SUPPLIED

The Perclose ProGlide SMC device and accessories are provided sterile, non-pyrogenic in unopened, undamaged packages. Products are packaged with ethylene oxide and intended for single use only. This single use device cannot be reused on another patient, as it is not designed to perform as intended after the first usage. Changes in mechanical, physical, and / or chemical characteristics introduced under conditions of repeated use, cleaning, and / or resterilization may compromise the integrity of the design and / or materials, leading to contamination due to narrow gaps and / or

spaces and diminished safety and / or performance of the device. Absence of original labeling may lead to misuse and eliminate traceability. Absence of original packaging may lead to device damage, loss of sterility, and risk of injury to the patient and / or user. Do not resterilize. Store in a cool, dry place.

#### Perclose ProGlide 6F SMC System includes:

One (1) Perclose ProGlide 6F SMC device

One (1) Perclose Suture Trimmer

#### 4.0 INDICATIONS

The Perclose ProGlide SMC System is indicated for the percutaneous delivery of suture for closing the common femoral artery access site of patients who have undergone diagnostic or interventional catheterization procedures using 5F to 21F sheaths. For sheath sizes greater than 8F, at least two devices and the pre-close technique are required.

#### 5.0 CONTRAINDICATIONS

There are no known contraindications to the use of this device. Attention is drawn to sections 6.0 WARNINGS and 7.0 PRECAUTIONS.

#### 6.0 WARNINGS

Do not use the Perclose ProGlide SMC device or accessories if the packaging or sterile barrier has been previously opened or damaged or if the components appear to be damaged or defective.

DO NOT RESTERILIZE OR REUSE. The Perclose ProGlide SMC device and accessories are intended for single use only.

Do not use the Perclose ProGlide SMC System if the sterile field has been broken where bacterial contamination of the sheath or surrounding tissues may have occurred, since such a broken sterile field may result in infection.

Do not use the Perclose ProGlide SMC System if the puncture site is located above the most inferior border of the inferior epigastric artery (IEA) and / or above the inguinal ligament based upon bony landmarks, since such a puncture site may result in a retroperitoneal hematoma. Perform a femoral angiogram to verify the location of the puncture site. **NOTE:** This may require both a Right Anterior Oblique (RAO) and Left Anterior Oblique (LAO) angiogram to adequately visualize where the sheath enters the femoral artery. Do not use the Perclose ProGlide SMC System if the puncture is through the posterior wall or if there are multiple punctures, since such punctures may result in a hematoma or retroperitoneal bleed.

Do not use the Perclose ProGlide SMC System if the puncture site is located in the superficial femoral artery or the profunda femoris artery, or the bifurcation of these vessels, since such puncture sites may result in a pseudoaneurysm, intimal dissection, or an acute vessel closure (thrombosis of small artery lumen). Perform a femoral angiogram to verify the location of the puncture site. **NOTE:** This may require both a Right Anterior Oblique (RAO) and Left Anterior Oblique (LAO) angiogram to adequately visualize where the sheath enters the femoral artery.

#### 7.0 PRECAUTIONS

- The Perclose ProGlide SMC System is provided sterile and non-pyrogenic in unopened undamaged packages.
- Prior to use, inspect the Perclose ProGlide SMC System to ensure that the sterile packaging has not been damaged during shipment. Examine all components prior to use to verify proper function. Exercise care during device handling to reduce the possibility of accidental device breakage.
- As with all catheter-based procedures, infection is a possibility. Observe sterile technique at all times when using the Perclose ProGlide SMC System. Employ appropriate groin management, as per hospital protocol, post procedure and post hospital discharge to prevent infection.
- Use a single wall puncture technique. Do not puncture the posterior wall of the artery.
- Do not insert the Perclose ProGlide SMC device into the femoral artery at an angle greater than 45 degrees.
- There are no reaccess restrictions if previous arteriotomy repairs were achieved with Abbott Vascular SMC devices.
- If significant blood flow is present around the Perclose ProGlide SMC device, do not deploy needles. Remove the Perclose ProGlide SMC device over a 0.038" (0.97 mm) (or smaller) guide wire and insert an appropriately sized introducer sheath.
- When pushing the plunger assembly to advance the needles, stabilize the device to ensure the device does not twist or move forward during deployment. Twisting the device could lead to needle deflection resulting in a cuff miss. Do not use excessive force or repeatedly push the plunger assembly. Excessive force on the plunger during deployment could potentially cause breakage of the device, which may necessitate intervention and / or surgical removal of the device and vessel repair.
- Do not apply excessive force to the lever when returning the foot to its original position (marked #4) down to the body of the device. Do not attempt to remove the device without closing the lever. Excessive force on the lever of the device or attempting to remove the device without closing the lever could cause breakage of the device and / or lead to vessel trauma, which may necessitate intervention and / or surgical removal of the device and vessel repair.
- **Do not advance or withdraw the Perclose ProGlide SMC device against resistance until the cause of that resistance has been determined** (see section 10.3 SMC DEVICE PLACEMENT). **Excessive force used to advance or torque the Perclose ProGlide SMC device should be avoided, as this may lead to significant vessel damage and / or surgical removal of the device and vessel repair.**
- If excessive resistance in advancing the Perclose ProGlide SMC device is encountered, withdraw the device over a 0.038" (0.97 mm) (or smaller) guide wire and reinsert the introducer sheath or use manual compression.
- Remove the Perclose ProGlide sheath before tightening the suture. Failure to remove the sheath prior to tightening the suture may result in detachment of the tip of the sheath.

- In using this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing damage due to application of surgical instruments such as clamps, forceps or needle holders.
- For closure of arteriotomy sites using a 5F–8F procedural sheath, use manual compression in the event bleeding from the femoral access site persists after the use of the Perclose ProGlide SMC device.
- For closure of arteriotomy sites using an 8.5F–21F procedural sheath, depending on the type of interventional catheterization procedure, use manual compression methods, compression assisted devices, and / or other appropriate treatment methods in the event bleeding from the femoral access site persists after the use of Perclose ProGlide SMC devices.

#### 8.0 SPECIAL PATIENT POPULATIONS

The safety and effectiveness of the Perclose ProGlide SMC devices have not been established in the following special patient populations:

- Patients with introducer sheaths < 5F or > 21F during the catheterization procedure.
- Patients with small femoral arteries (< 5 mm in diameter).
- Patients with access sites above the most inferior border of the inferior epigastric artery (IEA) and / or above the inguinal ligament based upon bony landmarks.
- Patients having arterial access in vessels other than the common femoral artery.
- Patients having a hematoma, pseudoaneurysm or arteriovenous fistula present prior to sheath removal.
- Patients with femoral artery calcium which is fluoroscopically visible at access site.
- Patients with severe claudication, iliac or a femoral artery diameter stenosis greater than 50% or previous bypass surgery or stent placement in the vicinity of access site.
- Patients with access sites in vascular grafts.
- Patients with prior intra-aortic balloon pump at access site at any time prior.
- Patients with ipsilateral femoral venous sheath during the catheterization procedure.
- Patients with ipsilateral arterial access sites punctured and compressed within 48 hours of closure. **NOTE:** The previous / initial puncture site may have the potential to re-bleed due to an unstable clot and / or anticoagulants, even if the new puncture site is successfully closed with Perclose ProGlide SMC device.
- Patients where there is difficulty inserting the introducer sheath or greater than one ipsilateral arterial puncture at the start of the catheterization procedure.
- Patients with antegrade punctures.
- Patients receiving glycoprotein IIb/IIIa inhibitors before, during, or after the catheterization procedure.
- Patients who are pregnant or lactating.
- Patients with bleeding diathesis or coagulopathy.
- Patients younger than 18 years of age.
- Patients who are morbidly obese (Body Mass Index  $\geq 40$  kg/m<sup>2</sup>).
- Patients with active systemic or cutaneous infection or inflammation.

Before considering early discharge, assess the patient for the following clinical conditions:

- Conscious sedation
- Anticoagulation, thrombolytic, or antiplatelet therapy
- Unstable cardiac status
- Hematoma at the closure site
- Hypotension
- Pain while walking
- Bleeding at the closure site
- Any comorbid condition requiring observation

The presence of any of the above factors has generally led to the deferral of early discharge recommendations.

#### 9.0 POTENTIAL ADVERSE EVENTS

Potential complications associated with the use of suture-mediated closure devices may include, but are not limited to, the following:

- Allergic reaction or hypersensitivity to device components
- Anemia
- Arterial stenosis / occlusion
- Arteriovenous fistula
- Bleeding / hemorrhage
- Bruising / hematoma
- Death
- Deep vein thrombosis
- Device entrapment
- Device failure / malfunction / misplacement
- Diminished pulses distal to closure site
- Embolism
- Extended hospitalization / delayed time to ambulation
- Infection / sepsis
- Inflammation
- Intimal tear / dissection
- Ischemia distal to closure site
- Nerve injury
- Numbness
- Pain
- Perforation
- Retroperitoneal hematoma / bleeding
- Surgical exposure / closure of common femoral artery
- Thrombus formation
- Vascular injury

- Vasoconstriction / vasospasm
- Vasovagal episode
- Wound dehiscence

#### 10.0 THE PERCLOSE PROGLIDE SMC SYSTEM CLINICAL PROCEDURE

The following instructions provide technical direction but do not obviate the necessity of formal training in the use of the Perclose ProGlide SMC System. The techniques and procedures described below are not intended as a substitute for the operator's experience and judgment in treating any specific patients.

##### 10.1 Examination and Selection of Products

1. Select the Perclose ProGlide SMC device(s) for closure of 5F to 21F introducer sheath access sites.
2. After carefully inspecting the packaging of the Perclose ProGlide SMC device and accessories for damage to the sterile barrier, remove the device from the package.
3. Exercise care when using additional instruments, such as clamps, forceps or needle holders during device handling, to reduce the possibility of accidental device breakage or damage to the suture.
4. Verify marker lumen patency by flushing the marker lumen with saline until the saline exits the marker port. **Do not use the Perclose ProGlide SMC device if the marker lumen is not patent.**

##### 10.2 Arterial Site and Puncture Considerations

1. An extremely deep tissue tract can influence needle trajectory, which prevents the Perclose ProGlide SMC device needles from engaging the cuffs, or secure knot tying, as the Suture Trimmer may not be able to advance the knot to the arterial wall for complete apposition before locking the knot. An extremely deep tissue tract may require a long access needle and / or upon inserting the Perclose ProGlide device, requires compression of the subcutaneous tissue (with the handle of the device) to be able to obtain pulsatile flow.
2. Before inserting the access needle, use of ultrasound guidance to visualize the common femoral artery or fluoroscopy to visualize the femoral head is recommended. When using the femoral head as a reference point, target the middle of the femoral head as the puncture site. Performing a femoral angiogram through the introducer sheath (or procedural sheath) to verify that the access site is in the common femoral artery is recommended before anticoagulants are given.
3. Puncture the anterior wall of the common femoral artery at an angle of approximately 45 degrees. Avoid side wall or posterior wall femoral artery punctures.
4. Prior to deployment of the Perclose ProGlide SMC device, perform a femoral angiogram to evaluate the femoral artery site for vessel size, calcium deposits, tortuosity, and for disease or dissections of the arterial wall to avoid device cuff misses (device needles not engaging with the cuffs) and / or posterior wall suture placement and possible ligation of the anterior and posterior walls of the femoral artery. The puncture should be proximal to the bifurcation of the superficial femoral artery and the profunda femoris branch and distal to the inferior margin of the inferior epigastric artery.
5. There are **no** femoral artery reaccess restrictions after using Abbott Vascular closure devices.

##### 10.3 SMC Device Placement 5F–8F Sheath, including Optional

###### Pre-Close and Maintaining Wire Access Techniques

The following instructions detail the deployment sequence to close the access site of a catheterization procedure performed through a 5F to 8F sheath size.

1. Place a 0.038" (0.97 mm) (or smaller) guide wire through the procedural (or introducer) sheath. Remove the procedural sheath while applying pressure on the groin to maintain hemostasis.
2. Backload the device over the guide wire until the guide wire exit port of the device sheath is just above the skin line. Remove the guide wire before the exit port crosses the skin line.
3. Continue to advance the device just until brisk pulsatile flow of blood is evident from the marker lumen. Position the device at a 45-degree angle. Deploy the foot by lifting the lever (**marked #1**) on the top of the handle. **Do not deploy the foot unless brisk pulsatile flow of blood ("mark") is evident from the marker lumen.**
4. Gently pull the device back to position the foot against the arterial wall. If proper position of the foot has been achieved, tactile sensation will be felt **AND** blood marking will cease or be significantly reduced to a slight drip. If marking does not stop or significantly change, evaluate the angiogram for femoral artery size, calcium deposits, tortuosity, disease and for location of the puncture (ensure footplate is not in bifurcation or side branch). Reposition the device to stop blood marking or reinsert the wire, remove the device to hold manual compression or insert a new sheath.
5. While maintaining device position, stabilize the device with your free hand (the one not used to deploy the device) to maintain the gentle retraction and to ensure the device does not twist or move forward during deployment. Use your other hand to deploy needles by pushing on the plunger assembly (**in the direction marked #2**) until you visually confirm that the collar of the plunger makes contact with the proximal end of the body. Do not use excessive force or repeatedly push the plunger assembly. After visually confirming contact with the body of the device only **one time**, this step is complete. Using your thumb as a fulcrum on the handle, gently disengage the needles by pulling the plunger assembly back (**in the direction marked #3**) and completely remove the plunger and needles from the body of the device. One suture limb will be attached to the anterior needle. The posterior needle will be free of suture. Pull back on the plunger until the suture is taut, which confirms that all the sutures have been fully retracted from the body of the device.
7. Do not attempt to redeploy the needles if the suture limb is not attached to the anterior needle. Reinsert the guide wire, and remove the device over the wire. Insert a new Perclose ProGlide SMC device to complete the procedure.

RELEASED

8. Use the QuickCut suture-trimming mechanism located on the handle to cut the suture from the anterior needle distal of the link. Use of a new, sterile scalpel or scissors is optional.
9. Relax the device and then return the foot to its original position by pushing the lever (marked #4) down to the body of the device. Do not use the lever to advance the device without closing the lever.
10. Withdraw the Perclose ProGlide SMC device until the guide wire port exits the skin line.
11. Grasp the suture adjacent to the device sheath and pull the suture ends through the distal end of the proximal guide. The rail suture limb is blue and is the longer of the two suture limbs. This rail suture limb will be used to advance the knot. The shorter, non-rail limb is white tipped and will be used to lock the knot.
12. Continue with Knot Advancement (Step 13). If closing over the wire, continue with section 10.3.1.

**10.3.1 Optional: Maintaining Wire Access During Knot Advancement (Closing over the wire)**

If the operator chooses to maintain wire access, reinsert the guide wire after exposing the guide wire port at skin level AND after the sutures have been harvested from the distal guide but before removing the device for knot advancement. Steps A-H are necessary when closing over the wire.

- A. There should be sufficient guide wire exposed out of the guide wire exit port before removing the device.
- B. Wrap the rail limb (long blue limb) of suture around your left index finger, low, close to skin level.
- C. Remove the Perclose ProGlide SMC device with the right hand, while maintaining an adequate length of guide wire inside the artery. This allows placement of another Perclose ProGlide SMC device or a sheath in the event that hemostasis is not obtained.
- D. While removing the device with the right hand, simultaneously advance the knot to the arteriotomy by applying slow, consistent increasing tension to the rail suture limb, keeping the suture coaxial to the tissue tract. (Do not advance the knot with the Suture Trimmer until the wire has been completely removed from the patient.)
- E. Assess the site for adequate hemostasis. If bleeding is controlled, the operator should then remove the guide wire. Next, use the Suture Trimmer, as described in step 13, to advance and tighten the knot until complete closure is achieved. Confirm the security of the knot by having the patient cough and / or bend his / her leg. Additionally, patients may be able to move freely in bed without head of bed or leg restrictions if the close is successful.
- F. If the wire is still in place and suture breakage occurs during knot advancement, or if hemostasis is not achieved, another Perclose ProGlide SMC device can be used to complete the procedure or a sheath can be inserted. Either remove the broken suture limbs or cut the suture limbs close to the knot (using the Suture Trimmer in step 13 or a new, sterile scalpel or scissor).
- G. Care should be taken to avoid excessive force if the insertion of another Perclose ProGlide SMC device or introducer sheath is required. To avoid resistance, use an introducer sheath small enough to avoid undue force.
- H. If hemostasis cannot be achieved after the wire is removed, apply manual compression.

**10.3.2 Optional: Pre-Close Technique**

The Perclose ProGlide suture can be placed around the arteriotomy at the beginning of the procedure and knot advancement can be placed on hold until the procedure is complete. Steps A-D are necessary when using a pre-close technique.

- A. After completing steps 1-10 of Device Placement as described in section 10.3, slightly rotate the device until you can see the two suture limbs in the bend of the distal guide. Grasp the sutures adjacent to the sheath. While holding the two suture limbs together, gently pull both suture ends through the distal end of the proximal guide.
  - B. Immediately place a shodded hemostat or clamp to hold the two suture limbs together at the distal end of the non-rail suture limb (shortest blue limb with white tip). To prevent knot advancement or locking of the knot, care must be taken not to pull on the individual suture limbs until the clamp is securely holding the two limbs together.
  - C. After securing the suture limbs and before inserting the procedural sheath, gently pull on the clamp until the suture is taut to remove any suture slack from the tissue tract.
  - D. Place the clamped suture under a sterile towel during the procedure. **NOTE:** The monofilament suture can be damaged by opening and closing the clamp. If you want to attach the suture to the drape, it is recommended that you use a second clamp with the tip placed through the handle of the first clamp and attach the second clamp to the drape.
  - E. At the end of the catheterization, reinsert the guide wire into the procedural sheath.
  - F. Knot advancement to close the arteriotomy will resume starting with section 10.3.1 (In step C, the procedural sheath will be removed instead of a Perclose ProGlide device).
13. Knot Advancement Using the Suture Trimmer
- A. Securely wrap the rail (longer, blue) limb of the suture around your left forefinger, low, close to the skin.
  - B. Do not tighten the suture around the sheath. Completely remove the device or the arterial sheath (if the device was deployed at the beginning of the catheterization procedure) from the artery, while simultaneously pulling gently on the rail limb. Always keep the suture coaxial to the tissue tract.

- C. With the rail (longer, blue) suture limb securely wrapped around your left forefinger, place the rail limb into the Suture Trimmer utilizing the following steps: Retract the Thumb Knob on the Suture Trimmer with the right hand. Place the Suture Trimmer under the suture limb making an "X" or a "cross" between the suture limb and the mid-point of the Suture Trimmer. Slide the Suture Trimmer back to load the suture into the Suture Gate located at the distal end of the Suture Trimmer. Keeping the Thumb Knob retracted, turn the Suture Trimmer coaxial to the suture and then release the Thumb Knob to capture the suture in the Suture Gate. Releasing the Thumb Knob before the suture is coaxial to the Suture Trimmer can cause the suture to be caught within the sliding mechanism at the distal tip and damage the suture. Once the suture is loaded correctly, the Suture Trimmer should slide easily.
- D. With the suture limb and Suture Trimmer coaxial to the tissue tract, move the knot to the arterial surface by advancing the Suture Trimmer with the right hand while placing slow, consistent increasing tension on the rail suture with the left forefinger. Avoid quick or jerking type movements with the suture limbs. The Suture Trimmer and suture limbs should always remain coaxial to the tissue tract. The Thumb Knob should be at 12 o'clock (facing the ceiling) and the Suture Trimmer should not be rotated.
- E. With the rail (longer, blue) suture limb securely wrapped around the left forefinger, place the Suture Trimmer under the left thumb to assume a single-handed position and complete knot advancement with slow, consistent increasing tension until the suture is taut (guitar string tightness).
- F. With the Suture Trimmer in place and suture taut, tighten the knot by gently pulling the non-rail (shorter, white tipped) suture limb, keeping it coaxial to the tissue tract.
- G. Hemostasis of the access site is achieved when the knot is fully advanced to the arterial surface, the slack is gently pulled from the knot with the non-rail limb while the Suture Trimmer holds tension on the rail limb of the suture, and the tissue is in complete apposition.
- H. Remove the Suture Trimmer from the tissue tract, relax tension on the suture, and test for hemostasis by having the patient cough and / or bend his / her leg. If hemostasis has not been achieved, assume the single-handed position for 20 seconds, or until hemostasis is achieved. Secure the knot again by gently pulling on the white tipped non-rail suture limb. Do not apply excessive pressure to the suture.
- I. Confirm the security of the knot by having the patient cough and / or bend his / her leg. Additionally, patients may be able to move freely in bed without head of bed or leg restrictions if the close is successful.
- J. Once hemostasis is achieved, use the Suture Trimmer to trim the sutures below the skin. While holding both suture limbs together and pulled taut, load both suture limbs into the Suture Trimmer (as described above in step 13 C) and advance the Suture Trimmer to the arterial surface. Trim the sutures by pulling back on the red Trimming Lever. Keep the Trimming Lever pulled back while removing the Suture Trimmer and trimmed suture limbs from the tissue tract. If only one suture limb has been loaded and trimmed, repeat the same technique on the other suture limb.
- K. If hemostasis cannot be achieved, apply manual compression until hemostasis is achieved.

**10.4 SMC Device Placement 8.5F-21F Sheath, Utilizing Pre-Close and Maintaining Wire Access Techniques**

The following introductions detail the deployment sequence for closing the access site of an interventional catheterization procedure performed through 8.5F to 21F sheath size. The pre-close technique using at least two devices must be used when closing sheath sizes from 8.5F to 21F.

1. Place a 0.038" (0.97 mm) (or smaller) guide wire through the introducer sheath. Remove the introducer sheath while applying pressure on the groin to maintain hemostasis.
2. Backload the SMC device over the guide wire until the guide wire exit port of the device sheath is just above the skin line. Remove the guide wire before the exit port crosses the skin line.
3. Place the first of two Perclose ProGlide SMC devices over the guide wire. Continue to advance the device just until brisk pulsatile flow of blood is evident from the marker lumen. The device lever (marked #1) and logo should be facing the ceiling (12 o'clock).
4. Rotate the device approximately 30 degrees towards the patient's right side (approximately 10 o'clock). Position the device at a 45-degree angle. Deploy the foot by lifting the lever (marked #1) on top of the handle. **Do not deploy the foot unless brisk pulsatile flow of blood ("mark") is evident from the marker lumen.**
5. Gently pull the device back to position the foot against the arterial wall. If proper position of the foot has been achieved, tactile sensation will be felt **AND** blood marking will cease or be significantly reduced to a slight drip. If marking does not stop or significantly change, evaluate the angiogram for femoral artery size, calcium deposits, tortuosity, disease and for location of the puncture site (ensure footplate is not in bifurcation or any side branch vessel). Reposition the device to stop blood marking (maintain the 30 degree rotation) or reinsert the wire and evaluate the situation before continuing the procedure.
6. While maintaining the device position, stabilize the device with your free hand (the one not used to deploy the device) to maintain the gentle retraction and to ensure the device doesn't twist or move forward during deployment. Use your other hand to deploy needles by pushing on the plunger assembly (in the direction marked #2) until you visually confirm that the collar of the plunger makes contact with the proximal end of the body. Do not use excessive force or repeatedly push the plunger assembly. After visually confirming contact with the body of the device only **one time**, this step is complete.

7. Using your thumb as a fulcrum on the handle, gently disengage the needles by pulling the plunger assembly back **(in the direction marked #3)** and completely remove the plunger and needles from the body of the device. One suture limb will be attached to the anterior needle. The posterior needle will be free of suture. Pull back on the plunger until the suture is taut. The device is designed so that the suture has been fully retracted from the body of the device.
8. Do not attempt to redeploy the needles if the suture limb is not attached to the anterior needle. Reinsert the guide wire, and remove the device over the wire. Insert a new Perclose ProGlide SMC device to complete the procedure.
9. Use the QuickCut suture-trimming mechanism located on the handle to cut the suture from the anterior needle distal of the link. Use of a new, sterile scalpel or scissors is optional.
10. Relax the device and then return the foot to its original position by pushing the lever **(marked #4)** on top of the device, down to its original position. Do not attempt to remove the device without closing the lever.
11. Withdraw the Perclose ProGlide SMC device until the guide wire port exits the skin line.
12. Slightly rotate the device until you can see the two suture limbs in the bend of the distal guide. Grasp the sutures adjacent to the sheath. While holding the two suture limbs **together**, gently pull both suture ends through the distal end of the proximal guide.
13. **Immediately** place a shodded hemostat or clamp to hold the two suture limbs together at the distal end of the non-rail suture limb (shortest limb with white tip). To prevent knot advancement or locking of the knot, care must be taken not to pull on the individual suture until the clamp is securely holding the two suture limbs together.
14. Gently pull on the clamp until the suture is taut to **remove any suture slack from the tissue tract**. Place the clamped suture on the right side of the patient under a sterile towel. It is important to identify which suture is deployed first as this is the knot that needs to be advanced first at the end of the procedure. **NOTE:** The monofilament suture can be damaged by opening and closing the clamp. If you want to attach the suture to the drape, it is recommended that you use a second clamp with the tip placed through the handle of the first clamp and attach the second clamp to the drape. It is important to remember this suture was placed first when suture tying occurs at the end of the procedure.
15. Reinsert the guide wire. There should be sufficient guide wire inside the vessel and exposed out of the guide wire exit port for device exchange.
16. Remove the Perclose ProGlide SMC device, while holding compression above the puncture site and maintaining adequate length of guide wire inside the artery. This allows placement of another Perclose ProGlide SMC device.
17. Repeat steps 2–13 with the second Perclose ProGlide device. **NOTE:** In step 4, the second device should be rotated approximately 30 degrees towards the patient's left side (approximately 2 o'clock).
18. After removing any excess suture slack from the tissue track, place the clamped suture for the second device on the left side of the patient under a sterile towel. It is important to identify which suture was placed first and which suture was placed second. At the completion of the procedure, the first knots will be advanced in the order they were placed. The knot from the first device placed on the right side (10 o'clock) of the patient would be advanced, followed by the knot from the second device placed on the left side (2 o'clock) of the patient.
19. **Knot advancement will be placed on hold at this point and the sutures set aside under sterile towels until the operator proceeds with the procedure.**
20. When utilizing a pre-close technique, the Perclose ProGlide SMC device is exchanged for an appropriately sized introducer sheath.
21. After completing the procedure, a hydrophilic or general purpose guide wire should be advanced into the artery. Maintain adequate length of guide wire in both the vessel and exposed out of the guide wire exit port to ensure guide wire access is maintained until hemostasis is achieved.
22. Heavily irrigate the secured Perclose ProGlide sutures with heparinized saline to remove any dry blood.
23. Remove the clamp from the first suture (patient right side / 10 o'clock). The rail suture limb is blue and is the longer of the two suture limbs. This rail suture limb will be used to advance the knot. The shorter, non-rail suture limb is white tipped and will be used to lock the knot.
24. Wrap the rail suture limb securely around your left forefinger, low, close to skin level, and hold the suture coaxial to the tissue tract. **While maintaining guide wire access**, carefully remove the entire sheath system and simultaneously pull the rail limb with slow, consistent increasing tension. Avoid quick or jerking type movements with the suture limbs. Manual pressure should be applied proximal to the puncture site for hemostasis, while the sheath is removed and during initial suture advancement.
25. **Do NOT lock the knot at this point.** Due to the size of the arteriotomy, use of Suture Trimmer (section 10.3, step 13) may be needed to approximate the tissue edges. However, **Do NOT lock or excessively tighten the knot while the guide wire is still in the vessel.** Again place the suture limbs on the right side of the patient for easy identification as the first suture deployed.
26. Remove the clamp from the second suture (patient left side / 2 o'clock) and advance the knot using the same technique and **maintaining guide wire access. Again, do NOT lock the knot.** Place the suture limbs on the left side of the patient for easy identification as the second suture placed.
27. Assess for hemostasis. If brisk bleeding is observed, advance the first (patient right side / 10 o'clock) suture again and then advance the second (patient left side / 2 o'clock) suture again. Multiple knot advancements are common when closing larger sheath sizes. However, **DO NOT lock or excessively tighten the knot while the guide wire is still in the vessel. Until the wire is removed, some bleeding will be visible, but it should not be pulsatile blood flow.**
28. If acceptable hemostasis is not observed, additional Perclose ProGlide SMC devices may be deployed at this point. Repeat steps 2–11 with the next Perclose ProGlide device. **NOTE:** In step 4, the third device should not be rotated. The device will be deployed in a straight crania / caudad position (facing the ceiling / 12 o'clock). After deployment of this device, advance the knot in the same fashion. **DO NOT lock the knot or excessively tighten the knot while the guide wire remains in the vessel.**
29. Assess the site for adequate hemostasis. **If bleeding is controlled, the operator should then remove the guide wire.** With the rail suture limb (longer, blue limb) securely wrapped around the left forefinger, again advance the first suture (patient right side / 10 o'clock) and then place the Suture Trimmer under the left thumb to assume a single-handed position and complete knot advancement with slow, consistent increasing tension until the suture is taut (guitar string tightness). With the Suture Trimmer in place and the suture taut, tighten the knot by gently pulling the non-rail (shorter, white tipped) suture limb keeping it coaxial to the tissue tract.
30. **Do NOT cut the suture.** Follow the same steps to advance the second suture (patient left side / 2 o'clock), coaxial to the tissue tract and lock the knot but **do NOT cut the suture.** If applicable, advance the additional sutures and lock the knots in the order that they were placed (10 o'clock, 2 o'clock, 12 o'clock). Assess bleeding. If hemostasis is deemed adequate, cut the suture tails below the surface of the skin using the Suture Trimmer (section 10.3, step 13 B) or a new, sterile scalpel or scissors.
- 10.5 Suture Breakage**
1. If suture breakage occurs before the knot is tightened, as a wire is still in place, discard the suture material, and use another Perclose ProGlide SMC device to complete the procedure.
  2. If suture breakage occurs after a knot has been advanced and / or tightened, and a wire is still in place, another Perclose ProGlide SMC device can be used to complete the procedure or a sheath can be inserted. Care should be taken to avoid excessive force if the reintroduction of another device or introducer sheath is required. To avoid resistance, use an introducer sheath small enough to be introduced without undue force.
  3. In all cases, if another Perclose ProGlide SMC device or introducer sheath cannot be inserted, use manual compression to obtain hemostasis.
  4. To prevent suture breakage, always pull on the suture limbs with slow, consistent increasing tension. Avoid quick or jerking type movements with the suture limbs.
  5. To prevent damage to the suture and subsequent suture breaks, the Suture Trimmer and suture limbs should always remain coaxial to the tissue tract. The Thumb Knob should be at 12 o'clock (facing the ceiling) and the Suture Trimmer should not be rotated. When loading the suture into the Suture Trimmer, keep the Thumb Knob retracted until the suture and Suture Trimmer are coaxial, then release the Thumb Knob to capture the suture in the Suture Gate.
- 10.6 Post-Procedure Patient Management**
1. Apply an appropriate dressing to the access site.
  2. Assess the access site as per hospital standard of care.
- 10.7 Recommendation for Patient Ambulation and Discharge**
- Patients who have undergone a diagnostic or interventional procedure using 5F–8F sheaths may be ambulated two hours after the Perclose ProGlide SMC device procedures.
- For patients who have undergone an interventional catheterization procedure using 8.5F–21F sheaths, post-procedure ambulation and patient care is at the discretion of the physician.
- In determining whether to ambulate or discharge an individual patient, it is important to consider all clinical factors including, but not limited to, anticoagulation regimen, antiplatelet and thrombolytic agents administered, oozing or bleeding from the access site, venous access site hemostasis, the general cardiovascular condition of the patient, anesthetic levels, and the overall clinical condition of the patient.
- 11.0 PRODUCT INFORMATION DISCLOSURE**
- Abbott Vascular Inc. has exercised reasonable care in the manufacturing of this device. Abbott Vascular Inc. excludes all warranties, whether expressed or implied, by operation of law or otherwise, including but not limited to, any implied warranties of merchantability of fitness, since handling and storage of this device, as well as factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond the control of Abbott Vascular Inc. directly affect this device and the results obtained from its use. Abbott Vascular Inc. shall not be liable for any incidental or consequential loss, damage, or expense, directly or indirectly arising from the use of this device. Abbott Vascular Inc. neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device.
- 11.0 PATENTS AND TRADEMARKS**
- This product and / or its use may be covered by one or more of the following United States Patents: 6,136,010; 6,190,396; 6,358,258; 6,558,399; 6,730,102; 6,746,457; 6,964,668; 7,001,400; 7,029,487; 7,048,747; 7,094,246; 7,147,646 ; 7,201,762; 7,235,087; EP 721313; EP 1093346. Other US patents pending. Foreign patents issued and pending.
- QuickCut, Perclose, and ProGlide are trademarks of the Abbott Group of companies.

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©2004, 2013, 2015 Abbott