RELEASED

English

Deutsch / German

Français / French

Español / Spanish

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Svenska / Swedish

ROTATING HEMOSTATIC VALVE
0.115" (2.92 mm or 8.8F)
(Tooly-Burst type adapter)

Reading this manual carefully prior to use. Failure to observe all warnings and precautions may result in complications.

HOW SUPPLIED

Sterile. Sterilized with ethylene oxide gas. Non-pyrogenic.

This single use device cannot be reused on another patient, it is not designed to perform as intended after the first use. Changes in mechanical, physical, and/or chemical characteristics introduced under conditions of repeated use, cleaning, and/or resterilization may compromise the integrity of the design and/or materials, leading to contamination due to narrow gaps and/or a space and dimensionality safety and/or performance of the device. Absence of original labeling may lead to misuse and ultimately hazardous. Absence of original packaging may lead to device damage, loss of sterility, and risk of injury to the patient and/or user.

Contents. One (1) Rotating Hemostatic Valve 0.115" (2.92 mm or 8.8F).

Storage. Store in a dry, dark, cool place.

INDICATION

The Hemostatic Valve is recommended for maintaining a fluid-tight seal around vascular catheters or related equipment during invasive vascular procedures. The inside diameter of the smallest portion of the valve is 0.115" (2.92 mm or 8.8F).

WARNINGS

This device is intended for single use only.

Do not inject any fluid if air bubbles are visible within the valve. First, aspirate the valve to remove the air, then flush the valve with normal saline as described in the DIRECTIONS FOR USE section.

If using a balloon dilatation catheter, during insertion and removal through the adapter, flush and fill the assembly with normal saline. To fill the valve sections, open one valve, place one finger over the luer fitting, and continue to fill the assembly.

PRECAUTION

Rotating Hemostatic Valve should not undergo pressures > 400 psi (2758 kPa).

PREPARATION FOR USE

Inspect the device prior to use to verify that no damage has occurred during handling. Flush and fill the assembly with normal saline. To fill the valve sections, open one valve, place one finger over the luer fitting, and continue to fill the assembly.

DIRECTIONS FOR USE

- Connect the hemostatic valve to the catheter or sheath. Aspirate the system to remove any trapped air and flush the assembly with normal saline. Remove blood from the valve section by opening the valve while continuing to flush the assembly and then closing the valve when the blood has been purged.
- Open the side valve, insert the catheter or equipment through the valve, and advance an appropriate distance into the vasculature. Close the valve around the shaft of the catheter or equipment so that the valve forms a liquid-tight seal around the shaft, just distal to initial movement of equipment.

WARNING: It is important that the valve be closed tightly enough to prevent blood leakage around the dilatation catheter shaft, yet not so tight as to restrict the flow of fluid (i.e., contrast medium) through the catheter.

- Perform the intended operation by following the recommended procedure of the equipment.
- Withdraw the equipment to the point that only the distal end of the equipment remains 10 to 20 cm inside the vasculature. Open the valve of the hemostatic valve and withdraw the equipment.
- Disconnect the hemostatic valve from the guiding catheter or sheath.

Follow the recommended procedures for performing post-dilatation angiography of the treated vessel.

ROTIERENDES HÄMOSTASEVENTIL
2,92 mm (0,115 in. oder 8,8 F)
(Tooly-Burst-Adapter)

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Speicherung: Store in a dry, dark, cool place.

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