TORQUE DEVICE

CAUTION

CAREFULLY READ ALL INSTRUCTIONS PRIOR TO USE. FAILURE TO OBSERVE ALL WARNINGS AND PRECAUTIONS MAY RESULT IN COMPLICATIONS.

DESCRIPTION

The torque device is designed to be secured to the proximal portion of a guide wire with a diameter from 0.009" to 0.018" (0.23 mm to 0.46 mm).

This device should be used only by physicians trained in invasive vascular procedures that utilize interventional and / or diagnostic devices (e.g., balloon dilatation catheters, atherectomy devices, stent delivery systems, intravascular ultrasound devices).

HOW SUPPLIED

Sterile. Sterilized with ethylene oxide gas. Non-pyrogenic.

This single use device cannot be reused on another patient, as it is not designed to perform as intended after the first usage. Changes in mechanical, physical, and / or chemical characteristics introduced under conditions of repeated use, cleaning, and / or resterilization may compromise the integrity of the design and / or materials, leading to contamination due to narrow gaps and / or spaces and diminished safety and / or performance of the device. Absence of original labeling may lead to misuse and eliminate traceability. Absence of original packaging may lead to device damage, loss of sterility, and risk of injury to the patient and / or user.

Contents. Ten (10) Torque Devices per box

Storage. Store in a dry, dark, cool place.

PRECAUTIONS

This device is intended for single-use only; do not reuse. Do not resterilize.

Refer to the appropriate device labeling for intended use, contraindications, and potential complications associated with the use of vascular interventional and / or diagnostic devices.

INSTRUCTIONS FOR USE

Inspect all product prior to use. Do not use if the package is open or damaged, or if product is damaged.

To utilize the Torque Device:

- Loosen the cap of the Torque Device.
- Insert the proximal end of the guide wire into the funnel-shaped hole on the distal end of the Torque Device. Once positioned at the desired location, tighten the cap to secure the Torque Device to the guide wire.
- Rotate the Torque Device to steer the guide wire to the desired position.
- To move the Torque Device to a new position, loosen the cap, slide the device along the guide wire to the desired position and tighten the cap.

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©2010 Abbott