20/30 PRIORITY PACK

Accessory Kit

CAUTION

Carefully read all instructions prior to use. Failure to observe all warnings may result in complications.

DESCRIPTION
The 20/30 PRIORITY PACK Accessory Kit is comprised of four disposable devices. (1) 20/30 INDEFLATOR Inflation Device – The 20/30 INDEFLATOR Inflation Device is an inflatable device with an analog pressure gauge that registers pressures from 140 to 30 atmospheres (VAC to 441 PSi/3041 kPa). The 20/30 INDEFLATOR Inflation Device employs a threaded plunger and locking watch that allow for the generation and sustaining of pressure. A three-way stopcock is included for use during preparation of the device and in conjunction with interventional devices. (2) Rotating hemostatic valve (RHV) – The rotating hemostatic valve is a Tuohy-Borst type adapter that has been designed to maintain a fluid-tight seal around vascular interventional and / or diagnostic devices. (3) Guide wire introducer – The guide wire introducer is a device designed to assist in the introduction of 0.009" (0.23 mm) to 0.018" (0.46 mm) guide wires used in conjunction with vascular interventional and / or diagnostic devices. (4) Torque device – The torque device is designed to be used with 0.009" (0.23 mm) to 0.018" (0.46 mm) guide wires to facilitate steering of the guide wires within the vascular anatomy.

HOW SUPPLIED

Sterile. Sterilized with ethylene oxide gas. Non-proprietary.

This single use device cannot be reused on another patient, as if so designed to prevent an intended after the first usage. Changes in mechanical, physical, and / or chemical characteristics introduced under conditions of repeated use, cleaning, and / or restoration may compromise the integrity of the design and / or materials, leading to contamination due to narrow gaps and / or spaces and diminished safety and / or performance of the device. Absence of original labeling may lead to misuse and eliminate traceability. Absence of original packaging may lead to device damage, loss of sterility, and risk of injury to the patient and / or user.

CONTENTS: (1) 20/30 INDEFLATOR Inflation Device, (1) Rotating Hemostatic Valve, (1) Guide Wire Introducer, (1) Torque Device, (1) Stopcock

Storage: Store in a dry, dark, cool place.

INDICATIONS:
The 20/30 PRIORITY PACK Accessory Kit is recommended for use during vascular procedures in conjunction with interventional and / or diagnostic devices (e.g., balloon dilatation catheters, atherectomy devices, stent delivery systems, intravascular ultrasound devices). The 20/30 INDEFLATOR Inflation Device is recommended for use in conjunction with balloon catheters.

WARNINGS

Use only with radiopaque inflation medium recommended by the interventional device manufacturer.

The in vivo interventional device inflation pressure should not exceed the maximum inflation pressure recommended by the manufacturer of the device.

This device is intended for single use only. Do not reuse. Do not sterilize, as this can compromise device performance and increase the risk of cross contamination and reprocessing.

It is important that the knurled knob on the rotating hemostatic valve (RHV) be closed tightly to prevent blood leakage around the shaft of the device, but do not so tight to restrict the flow of fluid (i.e., contrast media) through the device.

Do not inject any fluid if air bubbles are visible within the rotating hemostatic valve.

During insertion and removal of inflatable devices through the valve, the balloon must be completely deflated to enable it to safely and easily traverse the rotating hemostatic valve.

PRECAUTIONS

Note the "Use By" date specified on the package.

Use this product prior to use. Do not use if the package is open or damaged, or if product is damaged.

Refer to the appropriate device labeling for indications, contraindications, and potential complications associated with the use of interventional and / or diagnostic devices. This pack should be used by physicians trained in invasive vascular procedures that utilize interventional and / or diagnostic devices such as balloon dilatation catheters (e.g., balloon dilatation catheters, atherectomy devices, stent delivery systems, intravascular ultrasound devices).

Rotating Hemostatic Valve should not under pressures > 400 psi (2758 kPa).

CLINICIAN USE INFORMATION

I. To utilize the 20/30 INDEFLATOR Inflation Device:

- Connect the 20/30 INDEFLATOR Inflation Device to a 10 cc syringe filled with contrast medium.
- To fill the 20/30 INDEFLATOR Inflation Device with contrast medium:
 - Open the stopcock and rotate the locking switch counterclockwise to its unlocked position. Ensure that the ↕ icon located on the locking switch is centered within the syringe portion of the 20/30 INDEFLATOR Inflation Device.
 - Aspirate contrast medium by slowly withdrawing the plunger.
 - Evacuate air in the 20/30 INDEFLATOR Inflation Device by slowly injecting the necessary volume of 20/30 INDEFLATOR Inflation Device back into the syringe.
 - Repeat steps 2b through 2c until all air is removed from the 20/30 INDEFLATOR Inflation Device and 6 cc of contrast remains.
- To create a closed system, rotate the plunger locking switch clockwise into its locked position. Ensure that the ↶ icon on the locking switch is centered within the syringe portion of the 20/30 INDEFLATOR Inflation Device. Close the stopcock and disconnect the syringe.
- Connect the 20/30 INDEFLATOR Inflation Device stopcock to the inflation port of the interventional device.
- When balloon inflation is desired, verify that the plunger locking switch is in its locked position ↶ on the locking switch). Rotate the plunger handle clockwise to increase pressure or counterclockwise to decrease pressure.
- Pressure can be quickly released by rotating the lock switch one-half turn counterclockwise, to the unlocked position (↕ position on the locking switch) and then retracting the plunger handle. To maintain this position, turn the locking switch one-half turn clockwise to the locked position (↶ icon on the locking switch).
- To quickly release from the vacuum position, turn the locking switch one-half turn counterclockwise, to the unlocked position (↕ icon on the locking switch).

II. To utilize the rotating hemostatic valve:

- Flush and fill the rotating hemostatic valve with normal saline.
 - Connect the rotating hemostatic valve to the guiding catheter or sheath. Aspirate the system to remove any trapped air and flush the assembly with normal saline. Remove blood from the valve section by opening the knurled knob while continuing to flush the rotating hemostatic valve. Close the knob when the blood has been purged from the rotating hemostatic valve.
 - Open the valve to insert the interventional and / or diagnostic device through the valve. Close the knurled knob of the rotating hemostatic valve so that the valve forms a liquid-tight seal around the shaft and does not inhibit movement of the equipment or restrict the flow of fluid (e.g., contrast medium) through the device.
- III. To utilize the guide wire introducer:
- Insert the guide wire introducer, shaft end first, into the proximal opening of the rotating hemostatic valve that is attached to the end of the guiding catheter or sheath. Advance the guide wire introducer until its tip is distal to the seal in the rotating hemostatic valve.
 - Adjust the knurled knob of the rotating hemostatic valve to seal around the guide wire introducer in order to minimize blood loss.
 - Carefully insert the distal tip of the guide wire through the guide wire introducer and into the guiding catheter or sheath.
 - After the guide wire has been positioned at the desired location, the guide wire introducer may be removed.
- IV. To utilize the torque device:
- Unscrew the cap of the torque device.
 - Insert the proximal end of the guide wire into the funnel-shaped hole on the distal end of the torque device cap. Once it is positioned at the desired location, tighten the cap to secure the torque device to the guide wire.
 - Rotate the torque device to steer the guide wire to the desired position.
 - To remove the torque device to the new position, loosen the cap, slide the device along the guide wire to the desired position, and tighten the cap.

TRADEMARKS

INDEFLATOR is a trademark of the Abbott Group of Companies.

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