# ***Inflation Device***

<b>D</b> <b>Inflations Vorrichtung</b> (Einmalgebrauch)	<b>E</b> <b>Inflation Device</b> (Single Use)
<b>F</b> <b>Dispositif</b> <b>d'inflation/déflation</b> (Usage unique)	<b>I</b> <b>Dispositivo per</b> <b>Insufflazione</b> (Monouso)

- D** Gebrauchsanweisung
- E** Instruction Manual
- F** Manuel d'utilisation
- I** Manuale d'Uso


CE 0124



STERILE EO

C  
Contrast

	<b>DEUTSCH</b>	<b>ENGLISH</b>	<b>FRANÇAIS</b>	<b>ITALIANO</b>
	Autoklavierbar	Autoclavable	Autoclavable	Autoclavable
	Kontrastmittel injizierbar	Contrast Injectable	Liquide de contraste	Liquido di contrasto
	Zerlegbar	Detachable	Détachable	Disassemblabile
	Bespülbar	Flushable	Rinçable	Risciacquabile
	Hochfrequenz-Anschluss	High Frequency	Haute fréquence	Alta Frequenza
	Luer-Lock-Anschluss	Luer-Lock-Connector	Connecteur Luer-Lock	Connettore Luer-Lock
	Rotierbar	Rotatable	Rotatif	Rotante
	Empfohlener Führungsdraht	Recommended Guide Wire	Fil guide recommandé	Con guida
	Empfohlener Führungsdraht	Recommended Guide Wire	Fil guide recommandé	Con guida
	Empfohlener Führungsdraht	Recommended Guide Wire	Fil guide recommandé	Con guida
	Empfohlener Führungsdraht	Recommended Guide Wire	Fil guide recommandé	Con guida

**D****Warnung !**

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**E****Warning !**

The **Inflation Device** is for **SINGLE USE ONLY** and has to be disposed after use.  
**Reprocessing, resterilization and/or reuse of this product is prohibited!**

**Intended Purpose**

The **Inflation Device** instrument is intended for use during balloon angioplasty procedures to inflate, monitor the pressure and deflate the balloon.

**Warning !**

The user of this equipment should be thoroughly trained in the applicable procedure. Furthermore, failure to read and thoroughly understand the contents of this instruction manual may result in serious injury to the patient and/or user. It is essential to follow the instructions enclosed, as well as instructions found in other instruction manuals, infection control requirements and current hospital protocols for usage, cleaning and sterilization. Possible injuries related to endoscopic procedures may include perforation, electrical burns and shock, hemorrhage, infection, explosion, etc.

**! Failure to follow these instructions may result in malfunction and/or damage of the instrument or endoscope!**

**Storage**

Please store the instrument at a temperature of 10°C to 30°C and a humidity of 30 to 60%. Keep dry and avoid exposure to sunlight and chemical fume.

**General Safety Instruction**

- Wearing protective garments (i.e., gloves, protective eyewear, gowns, etc.) will assist in the protection against risk of cross-contamination and injury.
- Do not use an instrument if it does not function properly. If the device has any physical defect(s), replace the instrument with a new and/or functioning instrument.
- Check the endoscope's working channel diameter prior to introducing the instrument into the endoscope. To ensure the appropriate fit of the instrument the diameter of the endoscope's working channel needs to be at least "0.2 mm" larger than the instrument's outer diameter.  
**! Never use products which do not comply with the technical specifications !**
- Do not force an instrument down the working channel of the endoscope. This could result in endoscope and/or instrument damage.
- Do not under any circumstances make any modifications to the instrument or components described in this manual.
- If applicable check the active cord for cracks and/or damages on its insulation.

## **Inspection and Testing**

- Note any damage/open seal of the packaging.
  - ! *If damage or an open seal is detected, sterilization may be compromised and the instrument should be discarded and replaced.*
- Check date of sterilization.
  - ! *Never use products with non-valid sterilization date.*

## **Description**

Manufacturer's instructions for use should be read and understood for the ***Inflation Device*** and/or components connected to the ***Inflation Device***.

The disposable ***Inflation Device*** contains:

- One-piece 30 cc syringe with lock design
- Relating palm-grip
- Manometer
- High pressure connecting tube with a male rotating adapter

The manometer measure pressures ranging from vacuum to 30 ATM with an accuracy of 3%. The gauge is marked 1.0 ATM increments. There is also a scale of comparable PSI measurements.

## **Preparation**

- Prepare a reservoir of contrast media per balloon catheter or media manufacturer's recommendation.
- Hold the device so the manometer is upside down and insert the end of the connecting tube into the contrast solution.
- Push the lock release button and aspirate enough solution to fill the syringe.
- Hold the device upright to purge the air from the syringe and the connecting tube. Inspect the syringe and the connecting tube to ensure that the device has been completely purged and release the button.

### **Warning !**

***Check manometer needle position before each inflation.  
! It should be on zero +/- 1 ATM.***

## **Application**

- The piston is unlocked by pushing the lock release button forward, towards to the manometer.
- In this position you can freely pull the piston back for aspiration or push it for injection.
- To lock the piston in position, release the button. Turn the palm grip clockwise to infuse the contrast mixture into the balloon.
- The pressure in the balloon is registered on the manometer.

## **Balloon Inflation and Deflation**

- To inflate the balloon, engage the lock lever.
- Turn the palm grip clockwise slowly until the desired inflation pressure is reached.
- The lock device maintains the increasing pressure.
- To deflate the balloon, push the lock button forward, releasing the piston and pull back to aspirate. Release the button to lock.
- After use, dispose the ***Inflation Device*** according to legal requirements and your facility's infection control protocols.

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Medi-Globe GmbH