

D Polypektomieschlinge (Einmalgebrauch)	GB Polypectomy Snare (Single Use)
F Anse à polypectomie (Usage unique)	I Ansa da polipectomia (monouso)

D Gebrauchsanweisung
GB Instruction for Use
F Manuel d'utilisation
I Istruzioni per l'uso

No.	Deutsch	English	Français	Italiano
1	HF - Anschluss	HF-Connector	Connecteur HF	Connettore HF
2	Fingerschlitten	Fingerslide	Poignée Coulissante	Manico scorrevole
3	Teflontubus	Teflon Tube	Tubulure en Téflon	Tubatura in teflon

D

Achtung !

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GB

Warning!

The polypectomy snare is for SINGLE USE ONLY and has to be disposed after use. Reprocessing, resterilization and/or reuse of this product is prohibited!

Intended Purpose

The polypectomy snare is intended to remove polyps and small tumors from the gastrointestinal tract.

Contraindications

The following contraindications can be found in the medical literature with regard to the devices polypectomy snares, papillotomes, biopsy forceps and coagulation probes:

- Missing patient consent and lack of patient cooperation
- Non-fasting patient
- Inadequate preparation
- Intestinal wall fragility: e.g. highly florid inflammation of the colon (e.g. ulcerative colitis, diverticulitis, radiation colitis, toxic megacolon)
- Peritonitis, acute abdomen e.g. intestinal perforation, ileus
- Sepsis
- Comorbidity, e.g. severe cardiopulmonary diseases and decompensation
- Non-controllable haemorrhagic diatheses (prothrombin time below 50 %, partial thromboplastin time PTT prolonged by more than two-fold, thrombocytes below 50 000)
- Older pacemakers or pacemakers that cannot be set to a fixed frequency due to parasystolic rhythm
- Hip endoprosthesis or other metal implants
- Pregnancy
- Recently created gastrointestinal anastomosis

With regard to individual endoscopic indications for instrument application, the qualified user is ultimately responsible to assess contraindications and immanent risks of instrument application vs. its benefits for each specific medical situation.

Warning !

The user of this equipment should be thoroughly trained in the applicable procedure. Furthermore, failure to read and thoroughly understand the contents of this instruction manual may result in serious injury to the patient and/or user. It is essential to follow the instructions enclosed, as well as instructions found in other instruction manuals, infection control requirements and current hospital protocols for usage, cleaning and sterilization. Possible injuries related to endoscopic procedures may include perforation, electrical burns and shock, hemorrhage, infection, explosion, etc.

! Failure to follow these instructions may result in malfunction and/or damage of the instrument or endoscope.

Storage

Please store the instrument at a temperature of 5°C to 30°C and a humidity of 30 to 60%. Keep dry and avoid exposure to sunlight and chemical fume.

General Safety Instructions

- Wearing protective garments (i.e., gloves, protective eyewear, gowns, etc.) will assist in the protection against risk of cross-contamination and injury.
- Do not use an instrument if it does not function properly. If the device has any physical defect(s), replace the instrument with a new and/or functioning instrument.
- Check the endoscope's working channel diameter prior to introducing the instrument into the endoscope. To ensure the appropriate fit of the instrument the diameter of the

endoscope's working channel needs to be at least "0.2 mm" larger than the instrument's outer diameter.

- ! Never use products which do not comply with the technical specifications !
- Do not force an instrument down the working channel of the endoscope. This could result in endoscope and/or instrument damage.
- Do not under any circumstances make any modifications to the instrument or components described in this manual.
- If applicable check the active cord for cracks and/or damages on its insulation

Inspection and Testing

- Note any damage/open seal of the packaging.
- ! If damage or an open seal is detected, sterilization may be compromised and the instrument should be discarded and replaced.
- Check date of sterilization.
- ! Never use products with non-valid sterilization date.
- Check polypectomy snare prior to usage. Make sure the instrument is operating smoothly by advancing the Fingerslide **(2)** to open the snare.
- ! Should any irregularities occur, dispose the instrument and replace it with a new one.
- Connect your HF-Generator with the active cord to the connecting pin **(1)** on the fingerslide **(2)**.
- ! The polypectomy snare is compatible with generators of the following manufacturers: **ERBE, MARTIN, OLYMPUS und VALLEYLAB.**

Application

- Insure that the snare loop is fully retracted into the plastic sheath before slowly advancing the snare into the endoscopes working channel.
- ! Should any resistance be encountered when advancing the polypectomy snare through the curved sections of the scope, straighten the endoscope as much as possible to ease passage
- Open the snare loop by advancing the finger slide forward. Direct the snare loop over the polyp to be excised from the gastrointestinal wall and slowly retract the Fingerslide **(2)** to close the snare loop around the polyp at the desired location.
- ! Don't use too much force when closing the snare, it could get damaged!
- Insure that the proper connection is made from the polypectomy snare to the high frequency generator. According to the application refer to the table below for specific H.F. generator settings.

Application	Maximum permissible generator power	Maximum permissible voltage
Cutting	200 W	1300 Vpp – 1500 Vpp
Coagulating	120 W	3100 Vpp – 3200 Vpp

- ! Do not touch the metal working channel outlet of the endoscope with the metallic components of a (HF) device while the device is being energized via the electrosurgical generator! Metal components of instruments such as polypectomy snare loops can result in damage to the device such as short circuit, breakage and/or melting of the instrument.
- Apply short impulses of high frequency current until the polyp has been completely excised from the gastrointestinal wall. Repeat polypectomy procedures as necessary for removal of additional polyp(s).
- Retract the snare loop into the sheath and gently remove the snare from the endoscope. Proceed with polyp retrieval as applicable.
- ! Make sure that the snare loop remains closed while removing the snare.
- After use, dispose the entire instrument according to legal requirements and your facility's infection control protocols.

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Medi-Globe GmbH, Medi-Globe-Str. 1-5, 83101 Achenmühle - Germany
Medi-Globe GmbH