<p>D Biliärer und Pankreas Stent (Einmalgebrauch)</p>	<p>E Biliary and Pancreatic Stent (Single Use)</p>	<p>F Stents biliaires et pancréatiques (Usage unique)</p>
---	--	---

- D** Gebrauchsanweisung
- E** Instruction Manual
- F** Manuel d'utilisation

CE 0123

STERILE EO

.035"
Guide Wire

	DEUTSCH	ENGLISH	FRANÇAIS
	Autoclavierbar	Autoclavable	Autoclavable
	Kontrastmittel injizierbar	Contrast Injectable	Liquide de contraste
	Zerlegbar	Detachable	Détachable
	Bespülbar	Flushable	Rinçable
	Hochfrequenz-Anschluss	High Frequency	Haute fréquence
	Luer-Lock-Anschluss	Luer-Lock-Connector	Connecteur Luer-Lock
	Rotierbar	Rotatable	Rotatif
	Empfohlener Führungsdraht	Recommended Guide Wire	Fil guide recommandé
	Empfohlener Führungsdraht	Recommended Guide Wire	Fil guide recommandé
	Empfohlener Führungsdraht	Recommended Guide Wire	Fil guide recommandé
	Empfohlener Führungsdraht	Recommended Guide Wire	Fil guide recommandé

ART.NO	DEUTSCH	ENGLISH	FRANÇAIS
GBS-01-...S – GBS-09-...S	Stent Sets	Stent Sets	Sets de stent
GBS-01-... – GBS-09-...	Stents	Stents	Stents
GBS-31-...	Stent-Einführungs-Sets	Introducer Sets	Set d'introduction de stent
GBS-41-...	Führungskatheter	Guiding Catheter	Cathéter guide
GBS-42-...	Positionierer	Pusher	Poussoir

D**Warnung !**

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E**Warning !**

All Stents are for **SINGLE USE ONLY** and have to be disposed after use.
Reprocessing, resterilization and/or reuse of this product is prohibited!

Intended Purpose

The  **Medi-Globe** Stents and related components are intended to aid the drainage of obstructed ducts within the pancreatic and biliary trees. The stents are designed to be used during endoscopic procedures. It is important to follow your facility's current infection protocols, antimicrobial agents/processes and accessory manufacturer's instructions and recommendations. Failure to do so may result in patient/self injury, instrument malfunction and or endoscope damage.

Warning !

The user of this equipment should be thoroughly trained in the applicable procedure. Furthermore, failure to read and thoroughly understand the contents of this instruction manual may result in serious injury of the patient and/or user. It is essential to follow the enclosed instructions, as well as instructions found in other instruction manuals, infection control requirements and current hospital protocols for usage, cleaning and sterilization. Possible injuries related to endoscopic procedures may include perforation, electrical burns and shock, hemorrhage, infection, explosion, etc.

! Failure to follow these instructions may result in malfunction and/or damage of the instrument or endoscope!

Storage

Please store the instrument at a temperature of 5°C to 30°C and a humidity of 30 to 60%. Keep dry and avoid exposure to sunlight and chemical fume.

General Safety Instruction

- Wearing protective garments (i.e., gloves, protective eyewear, gowns, etc.) will assist in the protection against risk of cross-contamination and injury.
- Do not use an instrument if it does not function properly. If the device has any physical defect(s), replace the instrument with a new and/or functioning instrument.
- Check the endoscope's working channel diameter prior to introducing the instrument into the endoscope. To ensure the appropriate fit of the instrument the diameter of the endoscope's working channel needs to be at least "0.2 mm" larger than the instrument's outer diameter.
! Never use products which do not comply with the technical specifications !
- Do not force an instrument down the working channel of the endoscope. This could result in endoscope and/or instrument damage.
- Do not under any circumstances make any modifications to the instrument or components described in this manual.
- If applicable check the active cord for cracks and/or damages on its insulation.

Inspection and Testing

- Note any damage/open seal of the packaging.
 - ! *If damage or an open seal is detected, sterilization may be compromised and the instrument should be discarded and replaced.*
- Check date of sterilization.
 - ! *Never use products with non-valid sterilization date.*
- Check the Stent prior to each use.
 - ! *Should any irregularities occur, dispose the instrument and replace it by a new one.*
 - ! *Using a damaged instrument with physical or functional abnormalities may result in endoscope damage, patient and/or user injury.*

Application

Stent Placement (Set): (all Stents except Pigtail Stent)

- Cannulate the duodenal papilla and position a guide wire in the biliary duct/ pancreatic duct.
- Load the Guiding Catheter over the Guide Wire and into the biliary duct at the location where the Stent is to be positioned. For instillation of contrast medium (as required) remove the Guide Wire and instill via the Luer-Lock Connecting Piece at the proximal end of the Guiding Catheter. Remove the Luer-Lock Connecting Piece from the Guiding Catheter after use or if not required.
 - ! *A Guiding Catheter is not used with Stents of 7.0 or 5.0 French size. When using Stents of 7.0 or 5.0 French the Guide Wire is left in position for Stent Placement.*
- Introduce the tapered end of the Stent over the Guiding Catheter and place the Stent Positioner (short clear tube) over the Stent and proximal retention flap. Advance the Stent over the Guiding Catheter and into the working channel of the endoscope. After the proximal retention flap has entered the endoscope working channel, remove the Stent Positioner from the Stent and Guiding Catheter.
- Introduce the Pusher over the Guiding Catheter and into the endoscope working channel. Gently advance the Stent with the Pusher through the working channel and into the desired Stent placement position within the biliary duct/ pancreatic duct.
- After Stent placement, carefully remove the Pusher, Guiding Catheter and/or the Guide Wire from the endoscope.

Pigtail Stent Placement (Set):

- Follow procedures in step (see above) for Guide Wire positioning.
- Advance the Pigtail Stent over the Guide Wire and into the endoscope working channel. Stent Positioner is required for Stents with Flaps. (see above)
- Introduce the Pusher over the Guide Wire and into the endoscope working channel. Gently advance the Pusher through the working channel and into the desired Stent placement position within the biliary duct/ pancreatic duct.
- After Stent placement, carefully remove the Pusher and Guide Wire from the endoscope.

Stents: (purchased individually)

Stents purchased individually from  **Medi-Globe** should follow all user instructions contained in this manual. We further suggest that the following components are available and used as recommended in this instruction manual. A copy of this manual should be kept on file for reference when the Stents are purchased individually.

- **Guide Wire:** Used for cannulation of the biliary and/or pancreatic ducts allowing for instrument exchange and Stent placement.
- **Guiding Catheter:** Used to assist and guide the Stent into position and to allow for instillation of contrast medium as necessary. The Platinum band on the distal tip of the Guiding catheter is radiopaque and improves visualization under fluoroscopic examination.
- **Stent Positioner:** Placed over the proximal end of the Stent to minimize resistance of the Stent retention flap when inserting the Stent into the endoscope working channel.
- **Pusher:** Used to advance the Stent through the endoscope and into position within the obstructed duct.

Additional Information

- Verification of Stent placement is obtained by multiple magnified fluoroscopic examination. Cholangiograms may be necessary to detect and/or rule out other possible strictures within the affected duct. Multiple Stents may be required to assist in drainage of additional strictures within the duct.
- Stents may remain in the body for 29 days. If a Stent(s) become occluded over time it should be replaced with a new Stent.
- Although Stent migration is a rare occurrence, it remains a possibility. Full distal migration occurs when a Stent becomes dislodged from its originally position within the duct and moves into the duodenum. If a Stent has fully migrated, it will normally pass via the stool. It is possible that the Stent could cause internal injury to the duodenal wall and may require removal using various instruments used for foreign body retrieval prior to replacement with a new Stent.
- After use, dispose the stent and the components according to legal requirements and your facility's infection control protocols.

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Medi-Globe GmbH