
# p64 Flow Modulation Device

## INSTRUCTIONS FOR USE

2	DEUTSCH GEBRAUCHSANWEISUNG	29	DANSK BRUGSANVISNING	56	БЪЛГАРСКИ ИНСТРУКЦИИ ЗА УПОТРЕБА
5	ENGLISH INSTRUCTIONS FOR USE	32	NORSK BRUKSANVISNING	59	EESTI KASUTUSJUHEND
8	FRANÇAIS MODE D'EMPLOI	35	POLSKI INSTRUKCJA UŻYTKOWANIA	62	HRVATSKI UPUTE ZA UPOTREBU
11	ITALIANO ISTRUZIONI PER L'USO	38	SLOVENSKY NÁVOD NA POUŽITIE	65	SLOVENŠČINA NAVODILA ZA UPORABO
14	ESPAÑOL INSTRUCCIONES DE USO	41	ČESKY NÁVOD K POUŽITÍ	68	LIETUVIŲ K. NAUDOJIMO INSTRUKCIJA
17	PORTUGUÊS INSTRUÇÕES DE UTILIZAÇÃO	44	ROMÂNĂ INSTRUCȚIUNI DE UTILIZARE	71	TÜRKÇE KULLANIM TALİMATLARI
20	NEDERLANDS GEBRUIKSAANWIJZING	47	MAGYAR HASZNÁLATI ÚTMUTATÓ	74	РУССКИЙ ЯЗЫК ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
23	SVENSKA BRUKSANVISNING	50	ΕΛΛΗΝΙΚΑ ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	77	اللغة العربية إرشادات الاستخدام
26	SUOMI KÄYTTÖOHJEET	53	LATVISKI LIETOŠANAS PAMĀCĪBA		

## Hersteller:

phenox GmbH  
Lise-Meitner-Allee 31  
44801 Bochum, Deutschland  
Telefon: +49 234 36 919 0  
Telefax: +49 234 36 919 19

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1 x p64 Flow Modulation Device

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1 x p64 Flow Modulation Device

**Product description**


Fig. 1: p64 and delivery system (with distal wire tip) in peel-away sheath

Fig. 2: Detached delivery system and deployed p64 implant

The p64 Flow Modulation Device is a tubular vascular implant and consists of 64 interwoven Nitinol wires ⑨. Two wires, which are located opposite one another, are wrapped by platinum spirals ⑦ and ensure visibility under x-ray fluoroscopy. In addition, a platinum marker is located on each of the eight ends ⑤ on the proximal end of the implant.

To connect the implant to the delivery system, the eight implant ends ⑤ are flush-mounted into a slotted ring on the delivery wire ①. This slotted ring is covered by a polymer tube (detachment tube ③), which extends to the proximal end of the delivery system. The polymer tube is fixed to the proximal end of the delivery system using a detachable torquer ②.

The delivery system of the p64 has an x-ray marker on the slotted ring and another one at the distal wire tip ⑧ to allow the operator to determine its position.

The product is stored in a peel-away sheath ④ and is transferred into a microcatheter with an inside diameter of 0.027 inches. This sheath is peeled away during the insertion of the p64 to enable complete passage through the microcatheter.

The implant self-expands as it is pushed out of the microcatheter. Even after it is fully deployed in the target vessel, the implant can be completely recovered back into the microcatheter to allow it to be repositioned or removed.

The p64 is always deployed by means of a coordinated movement, whereby the microcatheter is withdrawn and the delivery system is advanced to avoid any movement of the distal implant end by shortening effect. Due to the shortening effect, the distal delivery wire tip ⑧ of the p64 moves distally during implant deployment.

The p64 is always recovered by means of a coordinated movement, whereby the microcatheter is advanced and the delivery system is withdrawn.

All manipulations are carried out under x-ray visualization.

After a final check of the expansion and position, the implant is manually detached, whereby the white torquer ② at the proximal end of the delivery system is unlocked and replaced by any standard torquer (being compatible with 0.014 inch or 0.016 inch microguidewire), which is then locked approximately 1 cm proximally to the end of the detachment tubing. The polymer tube ③ is then withdrawn towards the torquer while countering that force by pushing slightly on the torquer and delivery wire. The polymer tube ③ has an additional handle ⑩ at its proximal end for easier detachment.

**Intended use**

The p64 Flow Modulation Device is a self-expanding, tubular vascular implant and allows the controlled and selective modulation of blood flow in extra- and intracranial arteries.

In addition, the physical properties of the p64 straighten the target vessel slightly and strengthen it.

These properties aid the endovascular reconstruction of diseased arteries along their cervical and intracranial course.

**Indications**

The p64 Flow Modulation Device is a self-expanding, tubular implant and is used in the endovascular treatment of vascular diseases such as

- saccular and fusiform aneurysms and pseudoaneurysms
- vascular dissections in the acute and chronic phases
- vascular perforations and AV fistulae

**Contraindications**

The treatment of patients is contraindicated:  
- In whom antiplatelet and/or anticoagulation therapy is contraindicated or antiplatelet therapy

did not start in a timely manner prior to treatment,

- In whom angiography demonstrates the anatomy is not appropriate for endovascular treatment, due to conditions such as severe vessel tortuosity or stenosis.

**Compatibility**

All p64 models are compatible with microcatheters which have an internal diameter of 0.027 inches (0.69 mm). We recommend the

- ReVerse 27 (Reverse Medical/Covidien/Medtronic),
- EXCELSIOR XT-27 (150/6 cm, not Flex; Stryker Neurovascular),
- HEADWAY 27 (Microvention),
- MARKSMAN (ev3/Covidien/Medtronic).

In its relaxed state, the diameter of the p64 is approx. 0.3 or 0.5 mm (depending on the model) larger than the nominal diameter. The length specifications on the packaging describe the clinically usable length without the proximal ends.

The p64 must be used according to the specifications regarding the minimum and maximum target vessel diameters which are stated on the packaging.

**Recommended procedure**

1. Gather and document as complete a case history as possible, especially regarding the history of the present illness, comorbidities, previous interventions and current medication.
2. As far as possible, inform the patient and their relatives and document the patient's consent to go ahead with the intervention, pointing out possible complications and their potential consequences (disability, care dependency, death). In the case of patients who cannot give consent themselves, their relatives should be asked, wherever possible, what they think the will of the patient would be. Otherwise, the rules of emergency care for incapable patients apply.
3. Promptly take all measures necessary for a suitable pretreatment with drugs to ensure the **inhibition of platelet aggregation**. Based on current knowledge, dual antiplatelet therapy when implanting the p64 and similar products is suitable for the prevention of thrombus formation caused by the implant. For this purpose, 100 mg of aspirin and 75 mg of clopidogrel can be administered orally every day for at least 3 days before the planned intervention. Alternatively, 500 mg of aspirin and 600 mg of clopidogrel can be administered orally as one-off doses on the day before treatment. The safety of the treatment is presumably increased if the effective inhibition of platelet function is verified by means of an appropriate test (e.g. Multiplate, VerifyNow) before the intervention. Regarding substitutes in case of resistance to clopidogrel and the use of Gp IIb/IIIa antagonists, we refer you to the respective current specialist literature.
4. A CT or MRI examination of the cranium and, where necessary, the throat is advisable beforehand, in order to ensure a comprehensive preliminary diagnosis.
5. The diagnostic angiography and endovascular treatment should be carried out under general anesthesia with neuromuscular relaxation and invasive hemodynamic monitoring. During anesthesia, aim to maintain suitable systolic blood pressure values.
6. After preparing both groins, a 6F or 8F catheter is inserted, preferably into the right femoral artery.
7. Then moderate **heparinization** should begin, which should also last for the duration of the intervention. An intravenous dose of between 3000 and 5000 units of Heparin has proven suitable in practice. Where available, determining the ACT ("activated clotting time") is advisable.
8. Angiographic visualization of the internal and external carotid arteries on both sides and of the vertebral artery on at least one side is recommended, along with the respective dependent vessels. Enlarged images and, where necessary, oblique images of the affected vessel(s) are recommended.
9. The target vessel(s) for the endovascular treatment must be defined.
10. A 6F guide catheter or the combination of an 8F guide catheter and a suitable extension catheter is inserted into the afferent cervical vessel, avoiding vasospasm.
11. Measure the diameter of the target vessel, where the distal and proximal ends of the p64 are to be released, as accurately as possible.  
It is important that the implant is only inserted into target vessels of suitable size. In this connection, observe the specifications regarding the minimum and maximum vessel diameters on the packaging, as well as the instructions regarding selection of a model of the correct size (see Information on size selection). The length of the p64 must be selected so that the implant overlaps the lesion at the distal and proximal ends by at least a few millimeters.
12. Insert a suitable microcatheter with a corresponding microguidewire into the target vessel using a hemostatic valve and pressurized irrigation. Here, the use of so-called "road map" technology is advisable. Never probe against resistance! Aim to position the tip of the microcatheter 10-15 mm distal to the treatment target.  
Once the target vessel treatment segment has been reached, carefully pull on the microcatheter in order to remove any surplus catheter length and straighten the catheter.
13. Remove the micro-wire from the microcatheter under x-ray fluoroscopy.
14. Device preparation: Flush the device inside the peel-away sheath by deploying the implant out of the sheath in saline solution. Keep the tip of the sheath inside the saline solution and recapture the p64 again into the sheath.
15. With the aid of a tight-closing hemostatic valve and under continuous pressurized irrigation with heparinized physiological saline solution, the p64 is transferred from its peel-away sheath into the microcatheter. For this purpose, the hemostatic valve is opened. The peel-away sheath of the p64 is inserted through the open valve. The hemostatic valve is closed carefully and the peel-away sheath of the p64 is flushed by retrograde entry of the irrigation fluid.

16. Once the peel-away sheath of the p64 is completely flushed in this manner, it is advanced until it reaches the distal end of the hub adapter of the microcatheter. The peel-away sheath is held fixed in this position. The p64 is then advanced from the peel-away sheath into the microcatheter using the delivery wire to which the implant is fixed. This process is continued until approx. 60 cm of the delivery wire is inside the microcatheter.

17. The peel-away sheath is then pulled back proximally as far as the torquer. **Since the torquer must not be detached until the final detachment of the p64**, the peel-away sheath is removed fully by carefully peeling it away. For this purpose, a short loop is located on the proximal end of the peel-away sheath, which is pulled back proximally over the handle at the detachment tube and thereby tears open the peel-away sheath longitudinally. In the process, it is imperative to avoid buckling the delivery wire.

The process of inserting the p64 generally corresponds to that of inserting other similar implants. Should you encounter particular resistance which can only be overcome with effort, the implant and possibly also the microcatheter must be removed and the vessel accessed once more.

18. The p64 is slowly advanced to the tip of the microcatheter under continuous fluoroscopy. The distal tip of the implant should reach the tip of the microcatheter. **Never push the p64 beyond the distal tip of the microcatheter.** This can lead to a dissection or perforation of the target vessel.

19. Release the implant fully by carefully and very slowly withdrawing the microcatheter. **Please note that the distal delivery wire tip moves distally during the implant deployment.** Once the distal end of the implant is fully expanded and anchored in the distal vessel, continue to deploy the implant by continuously pushing on the delivery wire in order to facilitate the expansion of the p64. In order to ensure optimal wall apposition, the deployment must be a coordinated effort of continuous pushing of the delivery wire and adjustments (advancing or withdrawing) of the microcatheter so that the microcatheter is centralized longitudinally along the vessel. The release of the p64 should take place under fluoroscopy in order to ensure that the implant is properly deployed and the distal end has not moved.

20. The p64 is self-expanding and, when properly deployed, apposes itself against the vascular wall. The implant may over-expand slightly at the neck of the aneurysm due to the increased diameter at that point. Correct deployment can be verified using the helical, longitudinal platinum markers.

21. Injecting approx. 6-10 ml of x-ray contrast medium through the guide catheter allows one to check whether the aneurysm/dissection/target vessel has been satisfactorily covered by the deployment and release of the p64.

When using DSA systems with a digital detector and CT technology ("flat panel detector CT", e.g. DynaCT [Siemens], XperCT [Philips]), the implant can be visualized on the sectional image. This has proven particularly effective in the evaluation of the deployment and apposition to the wall.

22. If the radial deployment of the p64 is insufficient or the position or the model size selected is unsuitable, the p64 can be recovered into the microcatheter in order to allow the implant to be repositioned, redeployed or removed altogether. To do so, the microcatheter is advanced while the delivery system is slowly withdrawn.

23. If the position and deployment of the p64 are satisfactory, the implant is detached immediately.

**Detachment: Loosen the torquer at the proximal end of the delivery wire by unscrewing the torquer cap while the torquer handle is not rotated. Remove the torquer. Place another torquer about 1 cm proximally to the end of the detachment tubing and tighten the torquer to the delivery wire. While holding the delivery wire in a fixed position, gently pull the proximal end of the polymer tube towards the torquer.**

The polymer tube provides an additional handle at its proximal end for easier handling of the tube.

The slot of the implant ends of the delivery system is thus exposed and the ends can deploy. Due to the radial expansion of the proximal end, a contraction of the implant takes place.

24. Remove the delivery system by gently withdrawing the delivery wire.

25. Insufficient deployment of the p64 can be improved by means of a subsequent balloon dilatation. As far as possible, the p64 should appose against the vascular wall.

26. After the entire p64 is detached, if a subsequent telescoping device is required, **gently** advance the microcatheter through the p64. When the microcatheter tip is distal to the p64, gently retract the wire tip into the microcatheter and remove the delivery system completely out of the microcatheter. The microcatheter is now in position for a subsequent p64 to be advanced and deployed.

27. Injecting approx. 6-10 ml of x-ray contrast medium through the guide catheter allows one to check once more, if necessary, whether the target vessel has been sufficiently covered by the application of the p64. This check should be repeated 10 to 15 minutes later where necessary.

28. Take steps to ensure adequate inhibition of platelet aggregation. Proven medications following implantation include a 100 mg oral dose of aspirin every day on an ongoing basis and a 75 mg oral dose of clopidogrel every day for at least 12 months, but longer when necessary or on an ongoing basis.

#### Information on the selection of patients and lesions

If compliance to the antiplatelet medication described above cannot be guaranteed following the implantation of a p64, thrombotic closure of the implant and the vessel around it may occur within just a few days. Patients who cannot comply with the prescribed medication may not be suitable for treatment with a p64.

From the time of implantation of a p64, several weeks or months may pass before an aneurysm is no longer a risk. In this period, no complete protection from a (fresh) rupture/bleed can be guaranteed. Therefore, patients who are in the acute phase post aneurysm rupture should be treated by options that offer greater protection from re-rupture/bleed.

#### Information on size selection

- Select the implant diameter so that the deployed diameter comes as close as possible to the target vessel diameter, in order to achieve proper vessel wall apposition of the p64.
- Do not use the implant in target vessels whose diameter is not within the range of application specified on the packaging. Caution: Slight oversizing (selection of a p64 with a range of application above the diameter of the target vessel) is possible in order to reduce the mesh density of the implant.
- Caution: Substantial oversizing (selection of a p64 with a range of application considerably above the diameter of the target vessel) poses the risk of mal-deployment (incomplete expansion).
- Caution: Undersizing (selection of a p64 with a range of application below the diameter of the target vessel) leads to insufficient fixation of the p64 within the vessel and allows blood to flow around the outside of the implant (a so-called "endoleak"). The implant is then unstable, is subject to migration and hemodynamically ineffective.
- Ensure that the implant overlaps the lesion distally and proximally. If the selected product is too short or too long, it can be removed and replaced with a suitable one.
- Ensure that the implant does not end proximally in a narrow vessel curve because this may constrain a full proximal expansion. Choose an implant length which results in a complete coverage of the proximal vessel curve by the p64.

#### Precautions

- The p64 may be deployed up to three (3) times in the target vessel.
- For flushing deploy the p64 out of the peel-away sheath in saline solution and recapture the implant again into the sheath. Thorough flushing of the peel-away sheath is essential in order to remove any trapped air bubbles.
- All manipulations must be carried out under fluoroscopic visualization.
- If the implant is advanced beyond the distal end of the microcatheter, the vessel may be dissected or perforated.
- Do not pull the deployed implant back through the vessel into the microcatheter. Instead, push the microcatheter over the p64 while simultaneously fixating the wire to reposition and redeploy the implant if necessary.
- The p64 is a delicate implant and requires careful handling. Never push the microcatheter onto the p64 against resistance. Never twist the delivery system. If necessary, remove the p64 along with the microcatheter.
- Forceful pulling or pushing or twisting on the delivery wire or premature loosening of the torquer may inadvertently detach the p64 from the delivery system. In such a case, recovery with a foreign body retrieval device (e.g. Micro-Snare, ev3/Covidien/Medtronic) is recommended.
- Extreme stress to the delivery system tip could lead to separation of some of its parts. In such cases, recovery with a foreign body retrieval device (e.g. Alligator Device, ev3/Covidien/Medtronic) is recommended.
- If there is the risk that the patient is unresponsive or not fully responding to dual antiplatelet therapy, the time between implant deployment and detachment is supposed to be as short as possible to avoid any insufficient deployment of proximal implant end by bonding of braiding wires caused by blood ingredients (e.g. fibrin). Any additional treatment (e.g. coiling of the aneurysm while the associated microcatheter is "jailed" by the deployed p64) is supposed to be done after p64 detachment.
- Do not detach the implant in case of poor wall apposition or if the implant is in a stretched condition. Improve expansion by re-deploying or removing/replacing the device.
- If all eight ends of the p64 do not emerge out of the slotted ring of the delivery wire, careful movement of the delivery system helps to loosen them.
- Insufficient deployment of the proximal implant end can be improved by manipulation with the delivery system, the microcatheter, or means of a subsequent balloon dilatation. If the proximal end of the p64 does not open, recovery with a foreign body retrieval device (e.g. Micro-Snare, ev3/Covidien/Medtronic) is recommended.
- If vasospasm is suspected in the affected vascular region, all necessary measures, e.g. medication, should be used to aid regression prior to implantation.
- For implant detachment, use a separate torque device locked to the delivery wire. A suitable torquer is compatible with microguidewires that are 0.014 to 0.016 inches in diameter.
- If the torquer cannot be removed from the detachment tubing (by unscrewing and pulling it to the proximal direction), unscrew the torquer device completely. A higher resistance may be felt prior to complete separation of torquer handle and cap. Remove the torquer handle while the torquer cap and inlay stay in place. Then continue using the separate torque device for the detachment process.
- If the implant cannot be detached (or is not supposed to be detached) after torquer release and needs to be removed, tighten the torquer again to the detachment tubing while the detachment tube is compressed to create an initial load. This load avoids unintended p64 detachment during microcatheter advancement and p64 withdrawal in the microcatheter.
- Patients with known hypersensitivity to nickel-titanium materials may suffer an allergic reaction to the implant.
- Certification: The implant may only be used by specialized and appropriately trained physicians. Completion of a product training course run by phenox GmbH is a prerequisite

for the use of the p64. At least three (3) interventions using the p64 must be supervised by a physician or other qualified person delegated by phenox GmbH, and their course and outcome must be documented.

**General information**

- Keep away from heat. Store in a cool, dry place.
- Use is only permissible prior to the expiration date, as sterility is not guaranteed otherwise.
- Do not use damaged packaging, as sterility cannot be assumed otherwise.
- Intended for single use only. The implant must not be withdrawn into its peel-away sheath, resterilized or reprocessed for use in other patients since it cannot be cleaned reliably.
- The delivery system and, where necessary, packaging components must be disposed of appropriately in marked containers.

**Complications**

The following complications, amongst others, may arise during or after implant of the p64:

- Air embolism, embolism in distal vessels, thrombosis and cerebral ischemia
- Perforation, rupture, dissection and other arterial lesions
- Occlusion of side branch/perforator
- (Transient) stenosis of target vessel
- Vasospasm, occurrence of a pseudoaneurysm, intracranial hemorrhage
- Bleeding at the puncture site, allergic reaction, infection
- Space-occupying infarction, neurological deficit including the consequences of a stroke
- Pain and hemiplegia, disorders of movement and/or sensation, death

**Magnetic resonance imaging**

Non-clinical tests have shown that the p64 is suitable for MRI at a magnetic flux density of 3T. Under clinical conditions, 1.5T has proven to be unproblematic for the implant.

**Symbols and their meaning**

	Caution
	Consult instructions for use
	Do not reuse
	Do not resterilize
<b>LOT</b>	Batch code
<b>STERILE EO</b>	Sterilized using ethylene oxide
	Do not use if package is damaged
	Non-pyrogenic
	Use-by date
<b>REF</b>	Catalog number
<b>CONT</b>	Content
	Keep dry. Keep away from sunlight.
	This product was brought onto the market in accordance with Directive 93/42/EEC concerning medical devices.
	Manufacturer

**Limitation of liability**

phenox GmbH shall not be liable for damage caused by use other than that intended of the product.

phenox and p64 are registered trademarks of phenox GmbH in the Federal Republic of Germany and other states.

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Le dispositif de modulation du débit p64 est un implant vasculaire tubulaire constitué de 64 fils de Nitinol entrelacés ©. Deux fils situés à l'opposé l'un de l'autre et entourés de spires en platine ⑦ garantissent la visibilité de l'implant sous contrôle fluoroscopique. En outre, un marqueur au platine est situé sur chacune des huit extrémités ⑤ de l'extrémité proximale de l'implant.
phenox GmbH décline toute responsabilité concernant les dommages causés par une mauvaise utilisation du produit.
phenox et p64 sont des marques déposées de la société phenox GmbH en Allemagne et dans d'autres pays.