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ENGLISH

INSTRUCTIONS FOR USE

pITA Rapid Exchange Neurovascular PTA Balloon Catheter

Manufacturer

phenox GmbH
Lise-Meitner-Allee 31
44801 Bochum, Germany
Phone: +49 234 36 919 0
Fax: +49 234 36 919 19

Intended Purpose:

pITA is intended for angioplasty of intra- and extracranial brain supplying arteries, harbouring previously untreated stenoses or recurrent stenoses after angioplasty and/or stenting or flow diverter deployment. pITA is also to be used for dilatation of insufficiently adapted stents and flow diverters.

Indications

pITA is indicated for

1. Angioplasty of intracranial arteries
 - with symptomatic, previously untreated stenoses, with a degree of stenosis of at least 50%, and refractory to best medical therapy OR
 - with symptomatic stenoses, causing hemodynamic impairment and an instable neurological condition OR
 - with recurrent stenoses after angioplasty and/or stent implantation or stenoses after flow diverter deployment or stent assisted treatment of aneurysm, all of which causing recurrent clinical symptoms or having insufficient collaterals OR
 - for the treatment of acute stroke with occlusions or high grade stenoses caused by atherosclerosis or dissection OR
 - in order to enable another endovascular treatment.
2. Angioplasty of extracranial, brain supplying arteries
 - with symptomatic, previously untreated stenoses, with a degree of stenosis of at least 50% (NASCET). The consideration of surgery as an alternative treatment option has to be made on an individual basis. OR
 - with asymptomatic, previously untreated stenoses, with a degree of stenosis of at least 70% (NASCET). The consideration of surgery as an alternative treatment option has to be made on an individual basis. OR
 - with recurrent stenoses after angioplasty and/or stent implantation or flow diverter treatment or surgery, with a degree of stenosis of at least 50% (NASCET) OR
 - for the treatment of acute stroke with occlusions or high grade stenoses caused by atherosclerosis or dissection OR
 - in order to enable another endovascular treatment.
3. Dilatation of stents and flow-diverters with poor wall apposition in intra- and extracranial brain supplying arteries.

Contraindications

The balloon catheter is contraindicated for use in patients

- in whom antiplatelet therapy and/or anticoagulation therapy is contraindicated,
- with severe allergy or other contraindications for contrast media,
- with contraindication for whichever accompanying medication is necessary,
- who are judged to have lesions that prevent effective angioplasty.

Package contents

- 1 Rapid Exchange PTA Balloon Catheter
- 1 Flush cannula
- 2 PTA storing clips
- 1 Compliance data sheet

The catheter is only sterile and non-pyrogenic if the package is not opened, damaged or broken. Sterilized with ethylene oxide. Read instructions prior to using this product. For single use only.

WARNING: Re-use of single-use products carries a potential risk for the patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the product may lead to injury, illness or death of the patient.

Product description

pITA PTA-Catheter is a balloon-tipped, rapid exchange catheter. The distal shaft (40 cm) consists of two lumina. One lumen enables inflation and deflation of the balloon, the other permits the insertion of a guidewire.

Two radio-opaque markers indicate the length of the cylindrical portion of the balloon. The diameter 1.25 mm/1.5 mm balloon catheters have one centrally located x-ray marker only. The balloon is protected with a removable sheath, which protects the initial profile. The balloon is designed to provide an inflatable segment of known diameter and length at recommended pressures.

The two markers on the shaft of the catheter (115 cm and 125 cm from distal tip) help estimate the introduced catheter length while positioning the catheter. The proximal section of the catheter is a single-lumen, stainless steel hypotube. The distal portion has a double lumen for separate insertion of the guidewire. A cannula is provided for flushing of the distal guidewire lumen.

The product has a hydrophilic-coated surface at the distal catheter portion, which provides optimal guiding properties during the intervention.

WARNING: Do not expose the catheter to organic solvents, e.g. alcohol and DMSO.

WARNING: The enclosed clips are suitable to store the catheter safely and space-saving in the sterile field. The clips should only be used at the proximal shaft (hypotube). The clips should not be used at the distal shaft.

Material required to perform PTA with pITA Rapid Exchange PTA catheter

- Rapid Exchange pITA PTA catheter
- Guiding catheter (min. inner diameter 0.055"/1.40 mm)
- Manometer-controlled dilatation syringe
- 20 ml syringe for de-airing of the balloon catheter
- Introducer sheath with dilator
- 0.014" (0.36 mm) x 200 cm guidewire
- Contrast medium
- Sterile physiological saline solution
- Y-connector with hemostatic adapter

Recommended procedure

Only physicians familiar with and trained in angioplasty of cervical and intracranial brain supplying arteries may use the pITA balloon catheter. It is at the responsibility of the treating physician to keep informed with the current stage of development, indications and limitations of PTA of supraaortic vessels and to act accordingly.

PTA of supraaortic vessels may only be performed after sufficient premedication with antiplatelet agents and under systemic heparinization. Over- or underdose has to be avoided since it may increase the risk of hemorrhage or stroke.

Prior to use, it is essential that all the instruments, catheters and individual components be carefully checked for integrity and prepared according to the manufacturer's recommendations.

Selection of the balloon-tipped catheter

The nominal diameter must not exceed the diameter of the artery proximal and distal to the stenosis. For the treatment of intracranial vessels, undersizing is strongly recommended to increase safety.

If the stenosis cannot be crossed with the desired balloon catheter, predilatation with a smaller diameter balloon may facilitate passage of a more appropriately sized balloon catheter.

Preparation of the balloon-tipped catheter

1. Remove the balloon catheter from the sterile packaging. Do not use or re-use if the sterile packaging is damaged or any part of the device is bent or kinked.
2. The balloon filling fluid is a mixture of 50:50 contrast medium and sterile physiological saline solution. Higher concentrations of contrast medium may lead to increased deflation time or improper deflation. Never use air or gaseous inflation media.
3. Prior to use, the balloon must be completely de-aired. Draw up approximately 3 ml of the contrast medium mixture into a 20 ml syringe. Connect the syringe to the luer adapter with the syringe hub pointing downwards. Aspirate for 15-20 seconds and then release the plunger. Remove the syringe and evacuate all air from the barrel. Then reconnect it and repeat the process. Several repetitions may be necessary to completely de-air the balloon catheter.
4. Remove catheter from dispenser. Remove protective sheath. Reassure integrity of the balloon by visual inspection.
5. Irrigate the guidewire lumen of the monorail catheter with physiological saline solution using the accessory cannula.
6. Hydrate the hydrophilic coating at the distal catheter portion for at least 30 seconds.

Deployment technique

7. The procedure starts by inserting a guiding catheter fitted with a hemostatic adapter using standard technique.
8. Carefully introduce the guidewire and place tip distal to the stenosis. Difficult anatomic situations may require microcatheter-supported catheterization.
9. Carefully introduce the balloon-tipped catheter over the proximal end of the guidewire and into the hemostatic adapter of the guiding catheter. The balloon should always be deflated and folded around the shaft like an umbrella before usage, to maintain its narrow profile.
10. Advance the PTA catheter forward under x-ray and road-map control. Care must be taken to avoid uncontrolled guidewire movements. Never torque the device during use.
11. Position the balloon across the target lesion.
12. The manometer-controlled syringe must now be prepared and connected. First, the manometer-controlled syringe must be de-aired.
13. Verify that a meniscus of contrast medium appears on both the catheter balloon port and the manometer-controlled dilatation syringe.
14. Connect the manometer-controlled dilatation syringe and the catheter balloon port.
15. Slowly inflate the balloon under x-ray conditions to the appropriate pressure. Refer to compliance data table at the end of this manual. Balloon pressure must not exceed the rated burst pressure. The use of a manometer-controlled inflation syringe is mandatory in supraaortic PTA procedures. Consider premedication with parasympatholytic drugs prior to PTA of the extracranial carotid artery to prevent bradycardia. If difficulty is experienced during balloon inflation, do not continue the inflation. Remove the product and do not attempt to reuse it.

Removal procedure

16. After dilating under x-ray, completely deflate the balloon. Before withdrawing the product, visually confirm complete balloon deflation. If difficulty is experienced while deflating, connect a large-bore syringe and manually deflate the balloon.
17. Pull the balloon-tipped catheter back into the guiding catheter and check the results of dilatation angiographically.

18. If necessary, the procedure can be repeated with the same or different size **pITA** balloon catheter.

Complications

- Arterial perforation, rupture
- Arterio-venous fistulas
- Cardiac arrhythmia/cardiac arrest
- Death
- Dissection
- Distal embolism (caused by air, tissue, or thrombotic emboli)
- Hematoma
- Hemorrhage (intracranial and elsewhere)
- Hypotension/Hypertension
- Infection
- Ischemia (intracranial and elsewhere)
- Pain
- Pseudoaneurysm
- Restenosis
- Side effects of antiplatelet agents and contrast medium or other medication given during the procedure
- Stroke/TIA (including severe strokes with permanent disability and dependence from nursing care)
- Thrombembolism
- Thrombus formation
- Total occlusion of the artery
- Vascular complications leading to a prolonged procedure or additional endovascular or surgical procedures
- Vasospasm

Precautions

- Monitoring of vital signs is necessary throughout the procedure. For intracranial procedures, general anaesthesia is recommended.
- The risk of ischemic or hemorrhagic events may be increased in patients with hematologic diseases or hemorrhagic diathesis. The decision about the use of pITA shall be made weighing the treatment risk on an individual basis. Medical measures to decrease the risk shall be taken prior to treatment.
- PTA of supraaortic vessels may only be performed after sufficient premedication with platelet-inhibiting drugs and under heparinization.
- The pITA balloon was designed and intended for single use only. Do not resterilize. Do not reinsert. Do not use after expiration of the indicated expiry date or if damaged. Only the inside of the pouch and the content of the pouch can be considered sterile. Therefore, only remove the content prior to use.
- Prior to beginning of the angioplasty procedure, all equipment to be used, including the balloon catheter should be carefully examined to verify their functionality and integrity.
- Use only appropriate liquid balloon inflation medium. Failing to have the appropriate fluid (e.g. having a higher concentration of contrast dye) may lead to prolonged in- and deflation times. Never use air or gaseous inflation media. Do not expose the device to alcohol or DMSO.
- To achieve optimal performance and to maintain lubricity of the hydrophilic coating, continuous flow of appropriate flush solution between **pITA** and guiding catheter, and **pITA** and any intraluminal device should be maintained.
- Select the balloon size according to the target vessel size diameter. Avoid oversizing. For intracranial vessels, undersizing is recommended.
- To decrease the risk of air embolism, careful attention must be paid to ensure that tight catheter connections are maintained, that all devices are thoroughly flushed and prepared according to the manufacturer's recommendations.
- If thrombus apposition at the target lesion is detected during diagnostic angiography, the risk of thromboembolism is increased. The intervention should be postponed unless the medical condition of the patient requires an urgent intervention. Protective devices (e.g. filters) should be used if possible.
- During procedures in supraaortic, vessels never advance the catheter or other devices if resistance is felt. Determine the cause of the resistance and take appropriate action. If the cause cannot be determined, withdraw the catheter or device. Otherwise, damage of the device or vessel injury may occur.
- The product should only be manipulated in the body while under fluoroscopic observation with high quality radiographic imaging equipment.
- Balloon pressures during the intervention must not exceed the rated burst pressure. The product is tested that with 95% confidence 99.9% of all catheters, respectively catheters connection sites, are resistant to indicated maximum pressure (rated burst pressure) without failure or leakages. Higher pressures may cause balloon rupture or damage the catheter or the balloon joints respectively. Should this occur, deflate the balloon and remove it.
- Manometer-controlled inflation of the balloon is mandatory. Inflation and deflation **must be** performed slowly. The balloon must be deflated completely prior to removal.
- Do not reuse the balloon once deflated.
- If there is suspicion of vasospasm in the affected vascular region, all measures required for reversion should be taken, e.g. with appropriate vasodilating drug therapy. Prolonged severe vasospasm may cause a stroke due to thrombus formation or hypoperfusion.

- Blood pressure should be controlled during the procedure and up to seven days in patients at high risk for reperfusion hemorrhage (high-grade stenosis, poor collaterals, hemodynamic impairment, hypertension, old age), unless a concomitant medical condition requires blood pressure elevation.
- Patients should be instructed to take their antiplatelet medication regularly after treatment and be informed about possible consequences of non-compliance (stent thrombosis).

Storage

The packaged product should be protected from direct sunlight. If stored correctly, the product can be used up to the expiration date stated on the package label. The packaged products should not be exposed to temperatures below +10° C and above +40° C.

Disposal instructions

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

Subject to technical alterations

Symbols Glossary

	Manufacturer
	Consult instructions for use
	Caution
	This product is intended for single use only. Do not reuse!
	Batch number (Lot number)
	This product was sterilized using ethylene oxide.
	This product is non-pyrogenic.
	Do not resterilize
	Expiration date
	Catalog Number
	Contents
	Keep dry! Keep away from sunlight.
	Do not use if the package is damaged!
	This product was brought onto the market in accordance with Directive 93/42/EEC concerning medical devices.
	Temperature limitation: +10°C to +40°C

Limits of Liability

phenox GmbH is not liable for damage caused by other than intended use or reuse of the product.

phenox and **pITA** are registered trademarks of phenox GmbH in Germany and other countries.

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phenox et **pITA** sont des marques déposées de phenox GmbH en Allemagne et dans d'autres pays.