pORTAL®

INSTRUCTIONS FOR USE

2	DEUTSCH GEBRAUCHSANWEISUNG	29	DANSK BRUGSANVISNING	56	БЪЛГАРСКИ ИНСТРУКЦИИ ЗА УПОТРЕБА
5	ENGLISH INSTRUCTIONS FOR USE	32	NORSK BRUKSANVISNING	59	EESTI KASUTUSJUHEND
8	FRANÇAIS MODE D'EMPLOI	35	POLSKI INSTRUKCJA UŻYTKOWANIA	62	HRVATSKI UPUTE ZA UPORABU
11	ITALIANO ISTRUZIONI PER L'USO	38	SLOVENŠČINA NAVODILA ZA UPORABO	65	SLOVENČINA NÁVOD NA POUŽITIE
14	ESPAÑOL INSTRUCCIONES DE USO	41	ČEŠTINA NÁVOD K POUŽITÍ	68	LIETUVIŲ K. NAUDOJIMO INSTRUKCIJA
17	PORTUGUÊS INSTRUÇÕES DE UTILIZAÇÃO	44	ROMÂNĂ INSTRUCȚIUNI DE UTILIZARE	71	TÜRKÇE KULLANIM TALİMATLARI
20	NEDERLANDS GEBRUIKSAANWIJZING	47	MAGYAR HASZNÁLATI ÚTMUTATÓ	74	РУССКИЙ ЯЗЫК ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
23	SVENSKA BRUKSANVISNING	50	ΕΛΛΗΝΙΚΑ ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	77	اللغة العربية إرشادات الاستخدام
26	SUOMI KÄYTTÖOHJEET	53	LATVISKI LIETOŠANAS PAMĀCĪBA		

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Package contents

The table below shows the package contents for the different REF numbers:

	pPORTAL Steerable Hydrophilic Guidewire	pPORTAL EXT Extension Wire
Single Pack	REF: PORT-14-200-1 Contents: 1 Instructions for Use 1 sterile pouch with: 1 Guidewire 1 Torquer 1 Insertion Tool	REF: PORT-14-115-1-EX Contents: 1 Instructions for Use 1 sterile pouch with: 1 Extension Wire
5-Pack	REF: PORT-14-200-5 Contents: 1 Instructions for Use 5 sterile pouches with: 1 Guidewire 1 Torquer 1 Insertion Tool	REF: PORT-14-115-5-EX Contents: 1 Instructions for Use 5 sterile pouches with: 1 Extension Wire

Product Description

pPORTAL is a guidewire with a diameter of 0.36 mm (0.014"). pPORTAL may be used in standard length (205 cm), comprising

and a core wire (G). The distal segment has a hydrophilic coating (C, length 36 cm). The distal tip (E) of the guidewire is provided straight, but shapeable. The proximal portion is coated with Polytetrafluoroethylene (PTFE). pPORTAL Steerable Hydrophilic Guidewire's proximal end is equipped with an extension dock (F), which can be used for the connection to pPORTAL EXT Extension Wire.

Description of the component pPORTAL EXT Extension Wire:

pPORTAL EXT Extension Wire may extend pPORTAL Steerable Hydrophilic Guidewire to a total length of 318 cm, in order to simplify the exchange of over-the-wire catheters or other interventional devices. pPORTAL EXT Extension Wire has a diameter (a) of 0.36 mm (0.014"). The total length (b) is 115 cm. It consists of a stainless steel core wire (d) with a coating made of Polytetrafluoroethylene (PTFE). A Nitinol hypotube (c, 45 mm long) is attached to pPORTAL EXT Extension Wire's distal end to allow for connection with pPORTAL Steerable Hydrophilic Guidewire's proximal extension dock (F).

Intended Use

The intended use of pPORTAL is to facilitate the selective introduction and positioning of catheters and other interventional devices in extra- and intracranial brain-supplying arteries and in extra- and intracranial brain-draining veins. The intervention using pPORTAL is usually done through femoral artery access using a sheath and a guide catheter. For transvenous access a sheath is introduced in a suitable venous vessel (e.g. femoral vein). The duration of use is transient, thus, the device is under normal circumstances intended for continuous use for less than 60 minutes.

Indications

pPORTAL, consisting of pPORTAL Steerable Hydrophilic Guidewire and optional pPORTAL EXT Extension Wire, is intended for use in extra- and intracranial brain-supplying arteries and in extra-


Fig. 1

of the component pPORTAL Steerable Hydrophilic Guidewire only (fig.1, left). pPORTAL may in addition be extended to 318 cm, by connecting pPORTAL EXT Extension Wire (fig. 1, right) to pPORTAL Steerable Hydrophilic Guidewire.

Description of the component pPORTAL Steerable Hydrophilic Guidewire:

pPORTAL Steerable Hydrophilic Guidewire is a guidewire with a diameter (A) of 0.36 mm (0.014") and a length (B) of 205 cm. It consists of a radiopaque, distal platinum coil (D, length 5 cm)

and intracranial brain-draining veins. It can be used to facilitate the selective introduction and positioning of catheters and other interventional devices within target vessels.

Contraindication

The device is not intended for use in coronary arteries. The device is contraindicated for use in patients with contraindications for whichever accompanying medication is necessary.

Compatibility

pPORTAL Steerable Hydrophilic Guidewire is compatible with **pPORTAL EXT** Extension Wire.

pPORTAL Steerable Hydrophilic Guidewire and **pPORTAL EXT** Extension Wire may only be combined with each other and may not be combined with any other guidewire or extension wire.

When using other ancillary devices commonly used in intravascular procedures, verify compatibility to **pPORTAL** Steerable Hydrophilic Guidewire and **pPORTAL EXT** Extension Wire.

pPORTAL is compatible with the following products:

- microcatheters with an inner diameter of at least 0.0165" (0.42 mm) and at most 0.027" (0.69 mm), and
- other catheters, balloon catheters or balloon-expandable stents, whose IFU or label states compatibility to guidewires with a diameter of 0.014" (0.36 mm),

and which additionally fulfill the following properties:

- the maximum usable length is 170 cm, in case it is not a rapid-exchange catheter.
- their indication includes the use in extra- and intracranial brain-supplying arteries and in extra- and intracranial brain-draining veins.

Recommended Procedure

pPORTAL used in standard length of 205 cm (**pPORTAL** Steerable Hydrophilic Guidewire) may facilitate the delivery of catheters and interventional devices (described below in steps 1 to 12).

If **pPORTAL** Steerable Hydrophilic Guidewire is placed inside the target vessel, it may be extended using the optional **pPORTAL EXT** Extension Wire for the exchange of catheters/over the wire devices (described below in steps 13 to 21).

Preparation for use of **pPORTAL** Steerable Hydrophilic Guidewire:

1. Remove the dispenser (containing the guidewire) from the sterile packaging.
2. Before removing the guidewire from the dispenser, fill it with sterile saline and allow the guidewire to hydrate for at least 30 seconds.
3. Carefully withdraw the guidewire from the dispenser. If resistance is felt, repeat flushing until the guidewire can be easily removed from the dispenser.
4. Inspect the guidewire thoroughly to ensure that it is not damaged.
5. If tip shaping is desired, gently shape the guidewire's flexible distal tip according to standard practice using the insertion tool.

Directions for use of **pPORTAL** Steerable Hydrophilic Guidewire:

6. Before inserting the guidewire in a catheter/interventional device, flush the lumen.
7. Insert the hydrated guidewire through the guidewire port of the catheter/interventional device. The guidewire insertion tool may be used to facilitate insertion of the guidewire distal tip into catheter hubs or hemostatic valves.
8. Ensure that at least 5 cm of the guidewire extends from the most proximal hub. This will prevent the guidewire from slipping down inside the catheter/the interventional device.
9. If necessary, attach a torquer. To use the torquer, slip the device over the proximal end of the guidewire. Secure the torquer in place by tightening the rotating knob. Loosening and retightening the rotating knob will allow for the repositioning of the torquer.
10. Select the target vessel by carefully advancing or rotating the guidewire under fluoroscopy using standard techniques.
11. Maintain the position of the wire while catheters or interventional devices are guided to the target vessel.

12. The guidewire can be extended to a total length of 318 cm with the extension wire, in order to facilitate the exchange of catheters/interventional devices, if necessary.

Preparation for use of **pPORTAL EXT** Extension Wire:

13. Carefully remove the extension wire from the dispenser.
14. Inspect the extension wire thoroughly to ensure that it is not damaged.

Directions for use of **pPORTAL** Steerable Hydrophilic Guidewire in connection with **pPORTAL EXT** Extension Wire:

15. Remove any devices, such as a torquer, from the proximal end of **pPORTAL** Steerable Hydrophilic Guidewire.
16. While attaching and removing the extension wire from the guidewire, monitor the position of the guidewire's distal tip and the catheter/interventional device under fluoroscopy.
17. Insert the proximal end of **pPORTAL** Steerable Hydrophilic Guidewire into the distal tip of **pPORTAL EXT** Extension Wire. Hold the extension wire and the guidewire as closely as possible to the ends to avoid damage. Maintain the position of the guidewire during attachment. Push the extension wire and the guidewire firmly together to ensure a good connection.
18. In order to verify secure attachment, carefully pull on the guidewire as well as the extension wire.
19. Open the hemostatic Y-connector on the proximal end of the inserted catheter. Withdraw the catheter/interventional device via the guidewire connected to the extension wire.
20. Insert the replacement catheter/interventional device via the extension wire connected to the guidewire. Maintain the position of the guidewire tip.
21. After exchanging the catheter/interventional device, carefully remove the extension wire from the guidewire. Firmly hold the extension wire and the guidewire and pull the extension wire away from the guidewire. During removal, the position of the guidewire distal tip should not be changed.

Precautions

- Select the most appropriate guidewire to facilitate the treatment of the lesion.
- Read the instructions for use carefully before using **pPORTAL**.
- Do not use **pPORTAL**, if label is not legible.
- **pPORTAL** should only be used by physicians who are thoroughly trained in percutaneous intravascular techniques and procedures.
- If another device is used with **pPORTAL**, then refer to that product's labeling for intended use, contraindications and potential complications associated with the use of that device.
- Before every use, inspect guidewire and extension wire for surface irregularities, bends, and kinks or other damage.
- Do not use any damaged or irregular devices.
- Do not use the guidewire if the guidewire's tip cannot be shaped as desired.
- Handle the guidewire and extension wire carefully.
- **pPORTAL** may only be manipulated inside the body under fluoroscopy. Do not move the wire without monitoring the resulting movement of the distal tip.
- Other devices can be damaged or vessels injured by unintentional movement of the guidewire tip.
- Never advance or withdraw **pPORTAL** when resistance is felt. When resistance is felt, stop the procedure. Determine the cause of the resistance and take appropriate action. Excessive force against resistance may result in separation of the guidewire tip, damage to the catheters/interventional devices or vessel perforation.
- Appropriate premedication with antiplatelet agents and systemic heparinization are to be considered according to the overall intervention.

pORTAL®

- If there is suspicion of vasospasm, all measures for reversion should be taken, e.g., appropriate therapy with vasodilating medications. Prolonged severe vasospasm may cause a stroke due to thrombus formation or hypoperfusion.
- Do not expose **pORTAL** Steerable Hydrophilic Guidewire or **pORTAL** EXT Extension Wire to organic solvents (such as disinfection agents or DMSO).
- The guidewire must always remain wet for at least 30 seconds prior to use in order to properly hydrate the hydrophilic coating.
- Neither the guidewire insertion tool nor the torquer should be inserted into the body.
- A tight catheter may cause abrasion of the hydrophilic coating due to the variation in catheter size. It is advisable to stop using such catheters.
- Securely fasten the torquer onto the wire to prevent slippage of the torquer and to avoid damaging the device.
- **pORTAL** EXT Extension Wire must not be used alone. It must be securely connected with **pORTAL** Steerable Hydrophilic Guidewire in order to be used.
- The extension wire must not be inserted in the body.
- Do not bend the extension wire repeatedly at the same place. This can lead to deformation, breakage, or separation of the extension wire.
- Do not manipulate or rotate the extension wire while it is attached to the guidewire. This can cause the extension wire to detach from the guidewire.
- If a **pORTAL** Steerable Hydrophilic Guidewire or **pORTAL** EXT Extension Wire are removed for use from a package with five devices, ensure that the remaining sterile devices remain in the original packaging so that the integrity of the unused sterile bag is protected.
- Before use, check the packaging and devices carefully for transport damage.
- Do not use if the sterile packaging is opened or damaged.

General information

- Keep away from heat. Store in a cool, dry place.
- The device must be used before the expiration date, as sterility cannot otherwise be guaranteed.
- Do not use damaged packaging, as sterility cannot be assumed otherwise.
- Intended for single use only. **pORTAL** Steerable Hydrophilic Guidewire and **pORTAL** EXT Extension Wire must not be resterilized or reprocessed for use in other patients, since sterility cannot be ensured and the coating may degrade otherwise.
- The device, and if required, the packaging components must be disposed of following your institution's procedure for disposal of hazardous waste in appropriately marked containers.

Complications

The following malfunctions of the device and complications, among others, may arise during or after use of the **pORTAL**:

- allergic reaction
- death
- embolism
- failure to cross a lesion
- hematoma at the puncture side
- hemorrhage
- infection
- ischemia
- neurological/intracranial sequelae
- pseudo aneurysm
- seizure
- stroke

- thrombus formation
- transient ischemic attack
- vascular occlusion
- vasospasm
- vessel or aneurysm damage, including vessel perforation and dissection

Symbols and their Meaning

	Caution
	Consult instructions for use
	Do not reuse
	Do not resterilize
LOT	Batch code
STERILE EO	Sterilized using ethylene oxide
	Do not use if package is damaged
	Non-pyrogenic
	Use-by date
REF	Catalogue number
Contents	Content
	Keep dry. Keep away from sunlight.
	This product was brought onto the market in accordance with Directive 93/42/EEC concerning medical devices.
	Manufacturer

Limitation of Liability

phenox GmbH shall not be liable for damages caused by use other than that intended or reuse of the device.

phenox and **pORTAL** are registered trademarks of phenox GmbH in the Federal Republic of Germany and other countries.

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phenox et **pPORTAL** sont des marques commerciales déposées de phenox GmbH en Allemagne et d'autres États.