Set NeoPiston

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Révision : 2

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Sterile Single Use **Air/Water Valve, Suction Valve, Biopsy Valve, Water Jet Adapter and Cleaning Valve.**

- Eliminate manual cleaning and reprocessing of reusable valves
- Help create consistent practices
- Reduce the potential for errors

【Intended Use】

Air/Water Valve is used to control the air / water function of an endoscope during GI endoscopic procedures.

Suction Valve is used to control the suction function of an endoscope during GI endoscopic procedures.

Biopsy Valve is used to cover the opening to the biopsy/suction channel of Olympus®/Fujinon®/Pentax® endoscopes. The **Biopsy Valve** provides access for endoscopic device passage and exchange, helps maintain insufflation, minimizes leakage of biomaterial from the biopsy port through the endoscopic procedure and provides access for irrigation.

Water Jet Adapter is used in conjunction with the irrigation tubing (tubing and accessories to accommodate various endoscopes and irrigation pumps) and are intended to provide irrigation via irrigation fluids such as sterile water supplied during gastrointestinal endoscopic procedures when used in conjunction with an irrigation pump (or electrosurgical unit).

Cleaning Valve is intended as part of the precleaning process, to help clear the air/water channel of Olympus®/Fujinon®/Pentax® GI endoscopes post procedure.

【Warnings and Precautions】

- Guidelines for Preventing Cross Contamination:
 - These devices are intended to be used with Olympus®/Fujinon®/Pentax® Endoscopes, and follow the Olympus®/Fujinon®/Pentax® guidelines for air/water/suction use in the Olympus®/Fujinon®/Pentax® endoscopes.
 - Yangzhou Fartley did not design these devices to be reprocessed or reused, and therefore cannot verify that reprocessing can clean and/or sterilize or maintain the structural integrity of these devices to ensure patient and/or user safety.
 - These disposable medical devices are not intended for reuse. Any institution, practitioner, or third party, who reprocesses, refurbishes, remanufactures, resterilizes, and or reuses these disposable medical devices must bear full responsibility for their safety and effectiveness.
 - Endoscopic procedures should only be performed by persons having adequate training and familiarity with endoscopic techniques.

- Consult the medical literature relative to techniques, contraindications, complications and hazards prior to any endoscopic procedure.
- When handling these devices use appropriate aseptic technique.

Air/Water Valve and Suction Valve:

Make sure the holes in the air/water and suction valves are not clogged or blocked. (This could cause continuous air/CO₂ flow, or inadequate suction.)

Biopsy Valve:

- Do not use a sharp or pointed object to prime the valve prior to use.
- Do not leave a device hanging from the valve. Doing so can cause the creation of a larger valve slit/hole that may cause leakage.
- If the lid of the valve is opened while attached to the endoscope during a procedure, scope suction will be compromised and leakage may occur.
- If leakage occurs, a sterile gauze should be used to cover the valve.

Water Jet Adapter:

- Never attach or reattach the Water Jet Adapter to an endoscope that is/was inside the patient and has not been reprocessed.
- Exposure to bodily fluids may occur during the connection or disconnection of this device. The User is responsible for adhering to Body Substance Isolation protocols.
- Do not insert this device into the patient. The device should be attached to the auxiliary water port of a flexible gastrointestinal endoscope and an irrigation source. Use proper techniques to avoid risk of aspiration of a patient.
- Always prime the auxiliary water channel or biopsy channel prior to inserting the endoscope into the patient.

Cleaning Valve:

- Do Not Use the Disposable Cleaning Valve during patient procedures. Continuous air flow could cause patient injury.
- Endoscopic cleaning procedures should only be performed by persons having adequate training with endoscope cleaning techniques.
- Do not attempt to reuse, reprocess, refurbish, remanufacture or resterilize this device. Fartley did not design this device, nor is it intended to be reused, reprocessed, refurbished, remanufactured, or resterilized. Performing such activities on this disposable medical device presents a safety risk to patients

(i.e. compromised device integrity, cross-contamination, infection).

Contraindications

Contraindications include those specific to any endoscopic procedure.

【Instructions for Use】

- 1) Prior to clinical use, familiarize yourself with the device and read all the Instructions for Use.
- 2) Inspect the package for shipping or handling damage. If damage is evident, DO NOT USE THESE DEVICES, SAVE THEM FOR RETURN, AND CONTACT YOUR LOCAL PRODUCT SPECIALIST.

Suction Valve

- 3) Remove the suction valve from the tray.
- 4) Insert the suction valve into suction cylinder of an Olympus®/Fujinon®/Pentax® endoscope (see Figure 1).


Figure 1

- 5) Twist valve slightly to seat in proper orientation within the cylinder.
Note: The valve should not be able to be rotated once it is fully inserted into the cylinder.
- 6) Actuate the valve by pressing and releasing a few times to confirm smooth function.
Note: Dipping in water may facilitate ease of use.

Air/Water Valve

- 7) Remove the air/water valve from the tray.
- 8) Insert the air/water valve into the air/water cylinder of an Olympus®/Fujinon®/Pentax® endoscope (see Figure 1).
- 9) Confirm that the valve fits securely and that there is no gap between the valve end cap and the endoscope.

- 10) Actuate the valve by pressing and releasing a few times to confirm smooth function.
Note: Dipping in water may facilitate ease of use.

Biopsy Valve

- 11) If using the 3 piece kit (See Figure 2), attach the biopsy valve.


Figure 2

- 12) Securely place the biopsy valve onto the biopsy/suction channel opening of an Olympus®/Fujinon®/Pentax® endoscope.

Note: This biopsy valve does not accommodate a luer-lock syringe without an irrigating adaptor.

Water Jet Adapter

- 13) If using the 4 piece kit (See Figure 3), attach the Torrent® irrigation scope connector.


Figure 3

- 14) Attach the Water Jet Adapter to the irrigation tubing.
15) Attach the Water Jet Adapter to the auxiliary water port by screwing it onto the port. When attaching the scope connector to the port, take care to ensure that the scope connector is aligned to the port and avoid excessive tightening. Misalignment or excessive tightening of the scope connector could lead to the possibility of leakage.
16) Once the Water Jet Adapter is attached to the auxiliary water port of the endoscope, attach irrigation tubing to the Water Jet Adapter. The Water Jet Adapter is

manufactured with a one-way valve to minimize the risk of cross-contamination of the irrigation system.

- 17) Place distal tip of endoscope into water filled receptacle.
- 18) Press and hold the suction valve for a few seconds and then release.
- 19) Confirm that endoscope is suctioning water.
- 20) Confirm that air is exiting the tip of scope, which is evident by the production of air bubbles.
- 21) Remove distal tip of endoscope from water filled receptacle.
- 22) Press the air/water valve, and hold for 5-10 seconds.
- 23) Confirm that water exits the air/water nozzle on the end of the endoscope.
- 24) Remove and discard the air/water valve, suction valve, biopsy valve and the Water Jet Adapter after each patient use.

Cleaning Valve

- 25) If using the Cleaning Valve, see figure 4

- 26) Remove the device from the package.
- 27) Turn off the airflow regulator on the light source after initial rinsing of scope (per facility guidelines).
- 28) Remove the air/water valve from the scope.
- 29) Insert the Disposable Cleaning Valve into the air/water cylinder on scope handle.
- 30) Place the distal end of the scope into a basin of water.
- 31) On the light source, turn on the airflow regulator and set it to a maximum output of “high” .
- 32) Depress the Disposable Cleaning Valve for the recommended time per the scope manufacturer’ s
- 33) reprocessing instructions. Note: This flushes the air/water channel with water.
- 34) Release the Disposable Cleaning Valve.
- 35) Keep the distal end of the scope in the basin of water.

- 36) Flush air for the recommended time per the scope manufacturer’ s reprocessing instructions. Note: Air bubbles will be seen in the basin of water.
- 37) Turn off the air/flow regulator on the light source.
- 38) Remove distal end of scope from basin of water.
- 39) Remove the Disposable Cleaning Valve from scope after precleaning step is completed.

40) Discard the Disposable Cleaning Valve in accordance with institutional guidelines.

Warning :

After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

【Caution】

- ✓ Sterility is guaranteed unless the sterile pouch is damaged or seal is broken.
- ✓ Disposable device for single use.
- ✓ Sterilized by ETO Sterilization.
- ✓ Do NOT re-use or re-sterile.
- ✓ Check the production and expiry date printed on the sterile pouch.
- ✓ Do not store in direct sunlight, at extreme temperature, or in high humidity.
- ✓ Re-use may cause cross-infection.

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Symbol	Explanation	Symbol	Explanation
	Manufacturer		Keep away from sunlight
	Authorized representative in the European Community		Keep dry
	Date of manufacture		Temperature limit
	Use-by date		Do not re-use
	Batch code		Consult instructions for use
	Sterilized using ethylene oxide		Medical device
	Do not use if package is damaged		Notified Body, and product meet the basic requirements of MDD93/42/EEC

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