DONJOY®

ENGLISH

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USER PROFILE

The intended user should be a licensed medical professional, the patient, the patient's caretaker, or a family member providing assistance. The user should be able to:

- Read, understand and be physically capable to perform all the directions, warnings and cautions provided in the information for use.

INTENDED USE/INDICATIONS

- Post-traumatic, pre-op or post-operative use, functional rehabilitation use, return to activities.
- Light to moderate ACL/PCL injuries, moderate to severe lateral ligament injuries (MCL, LCL), meniscal injuries, knee surgery.

CONTRAINDICATIONS

None

WARNING & PRECAUTIONS

- We recommend that the first application be made under supervision of a health care professional.
- If pain, swelling, changes in sensation or other unusual reactions occur while using this product, you should contact your doctor immediately.

APPLICATION INSTRUCTIONS

The knee brace has 2 lateral hinges blocked in extension. It also has a removable posterior stay along the popliteal area. This stay is pre-shaped, but can be moulded according to patient's needs or it can be removed.

With the patient lying on their back, open the brace.

1 - Position the patient's leg on the brace, with the knee centred between the 2 bilateral hinges. Close the anti-slip strap on the calf using the hook and loop closure.

2 - Secure the brace by closing the thigh cuff and calf cuff using the hook and loop closures.

3 - Adjust the 4 straps to the required length and close using the buckle/clip closure. Tighten the brace according to level of compression required.

Hinges

- To adjust flexion/extension, open the cap of each hinge, position the flexion/extension stops at the desired angle, and close the hinge cap.
- To avoid damaging the hinges, the flexion/extension stops must be positioned at the same angle on both hinges.

CARE

- Hand wash in lukewarm water (30°C) with a mild detergent
- Air dry (do not place on radiators/heaters or other heat sources)
- Do not use bleach or chlorine, do not iron.

COMPOSITION

Knee brace: Aluminium, polyester, polyamide, monofilament, polyurethane, cotton, elastomer, resin, acetal (POM), polypropylene, polyethylene

NOT MADE WITH NATURAL RUBBER LATEX

INTENDED FOR USE ON A SINGLE PATIENT

WARRANTY DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale.

NOTICE: WHILE EVERY EFFORT HAS BEEN MADE IN STATE-OF-THE-ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.

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