INTENDED PURPOSE

The Penumbra Delivery Microcatheters are intended to assist in the delivery of diagnostic agents and therapeutic devices to the vasculature.

DEVICE DESCRIPTION

The Penumbra Delivery Microcatheters are single lumen medical devices designed to aid a physician in accessing distal vasculature when used in conjunction with a guide catheter and a micro guidewire. Penumbra Delivery Microcatheters are supplied with a peel away introducer sheath, annealed stainless steel mandrels that can be used to shape the distal tip as desired, and a rotating hemostasis valve (RHV).

The Penumbra Delivery Microcatheters includes the following:

- Penumbra PX Slim™ Microcatheter
- Penumbra Velocity® Microcatheter
- Penumbra LANTERN™ Microcatheter
- ACCESS25™ Microcatheter

The Microcatheters may be available in the following shapes:

- Straight
- 45 degree
- 90 degree
- 130 degree (J)

The Penumbra Delivery Microcatheters are supplied sterile for use in a sterile clinical field.

INDICATION FOR USE

The Penumbra Delivery Microcatheters are intended to assist in the delivery of diagnostic agents, such as contrast media, and therapeutic devices, such as occlusion coils to the peripheral and neuro vasculature.

INTENDED PATIENT POPULATION

The intended patient population of the Penumbra Delivery Microcatheters is patients who require diagnostic or therapeutic interventional procedures.

CONTRAINDICATIONS

There are no known contraindications.

WARNINGS

- Do not advance or withdraw the Penumbra Delivery Microcatheters against resistance without careful assessment of the cause using fluoroscopy. If the cause cannot be determined, withdraw the device. Moving or torquing the device against resistance may result in damage to the vessel or device.
- The devices are intended for single use only. Do not resterilize or reuse. Resterilization and/or reuse may compromise the structural integrity of the device and/or lead to device failure which in turn may result in patient injury, illness, or death.

PRECAUTIONS

- The Penumbra Delivery Microcatheters should only be used by physicians who have received appropriate training in interventional techniques.
- Do not use kinked or damaged devices. Do not use open or damaged packages. Return all damaged devices and packaging to the manufacturer/distributor.
- Use prior to the “Use By” date.
- Use the Penumbra Delivery Microcatheters in conjunction with fluoroscopic visualization.
- As in all fluoroscopy procedures, consider all necessary precautions to limit patient radiation exposure by using sufficient shielding, reducing fluoroscopy times, and modifying radiation technical factors whenever possible.
- Maintain a constant infusion of an appropriate flush solution.
- If flow through the device becomes restricted, do not attempt to clear the lumen by infusion. Remove and replace the device.
- Avoid using alcohol, antiseptic solutions, or other solvents to pre-treat the device because it may cause unpredictable changes in the coating which could affect the device safety and performance.

POTENTIAL ADVERSE EVENTS

Possible complications include, but are not limited to, the following:

- |                                                                             |                                              |                                                                                                                    |
|-----------------------------------------------------------------------------|----------------------------------------------|--------------------------------------------------------------------------------------------------------------------|
| • Acute vessel occlusion                                                    | • Emboli including foreign body embolization | • Neurological deficits including stroke                                                                           |
| • Allergic reaction and anaphylaxis from contrast media or device materials | • Hematoma or hemorrhage at access site      | • Pseudoaneurysm                                                                                                   |
| • Air embolism                                                              | • Hemorrhage                                 | • Radiation exposure that may lead to cataracts, skin reddening, burns, alopecia, or neoplasia from x-ray exposure |
| • Aneurysm rupture                                                          | • Infection                                  |                                                                                                                    |
| • Arteriovenous Fistula                                                     | • Inflammation                               |                                                                                                                    |
| • Death                                                                     | • Intracranial hemorrhage                    |                                                                                                                    |
|                                                                             | • Ischemia                                   |                                                                                                                    |

- Renal impairment or acute renal failure from contrast media
- Vessel dissection
- Vessel perforation
- Vessel spasm
- Vessel thrombosis

\*If a serious incident related to the device occurs, contact your Penumbra, Inc. representative and the competent authority in your respective country/region.

COMPATIBILITY INFORMATION

- Guide Catheter: 0.041in [1.04mm] ID minimum (or 0.038in [0.97 mm] guidewire-compatible diagnostic catheter)
- Guide Wire: 0.024in [0.61mm] maximum

FLOW RATE INFORMATION

| PX SLIM, Velocity, and LANTERN Delivery Microcatheters |                 |                       |                    |                    |                    |                    |
|--------------------------------------------------------|-----------------|-----------------------|--------------------|--------------------|--------------------|--------------------|
| Total Catheter Length (cm)                             | Dead Space (ml) | Flow Rate (ml/sec) at |                    |                    |                    |                    |
|                                                        |                 | 150 PSI (1034 kPa)    | 300 PSI (2068 kPa) | 450 PSI (3102 kPa) | 600 PSI (4136 kPa) | 750 PSI (5171 kPa) |
| 115 cm                                                 | 0.53 ml         | 1.9                   | 2.8                | 3.5                | 4.4                | 5.1                |
| 135 cm                                                 | 0.59 ml         |                       |                    |                    |                    |                    |
| 150 cm                                                 | 0.64 ml         |                       |                    |                    |                    |                    |
| 160 cm                                                 | 0.67 ml         |                       |                    |                    |                    |                    |

| ACCESS25 Delivery Microcatheter |                 |                                          |                       |          |
|---------------------------------|-----------------|------------------------------------------|-----------------------|----------|
| Total Catheter Length (cm)      | Dead Space (ml) | Flow Rate (ml/sec) at 300 PSI (2068 kPa) |                       |          |
|                                 |                 | Saline                                   | 50/50 Saline/Contrast | Contrast |
| 150 cm                          | 0.67 ml         | 2.9                                      | 2.4                   | 0.7      |
| 160 cm                          | 0.70 ml         |                                          |                       |          |
| 170 cm                          | 0.73 ml         |                                          |                       |          |

DEVICE PREPARATION AND USE

1. Refer to Warnings, Precautions, and Potential Adverse Events prior to use.
2. Remove microcatheter from packaging and inspect for damage or kinks.
3. Hydrate the outer surface of the microcatheter.
4. If desired, the microcatheter distal tip may be shaped using the provided shaping mandrel. Bend to the desired shape and hold mandrel/microcatheter assembly above a sterile steam source. Cool tip in sterile saline. Remove mandrel from microcatheter and retain for additional shaping use, or discard.
5. Connect syringe with heparinized saline to microcatheter hub and flush to facilitate guidewire insertion.
6. Prepare guidewire according to manufacturer's instructions.
7. Insert guidewire into microcatheter hub according to manufacturer's instructions.
8. Attach a rotating hemostasis valve (RHV) to microcatheter hub. Attach stopcock to RHV sidearm.
9. Connect RHV to flush solution pressure bag. Maintain pressure above arterial pressure to prevent retrograde blood flow into microcatheter.

DEVICE USE

1. Place appropriate guide catheter in the desired vascular location according to manufacturer's instructions. Attach a rotating hemostasis valve (RHV) to guide catheter hub. Attach stopcock to RHV sidearm. Connect RHV to flush solution pressure bag. Maintain pressure above arterial pressure to prevent retrograde blood flow into guide catheter.
2. Use peel away introducer to insert microcatheter/guidewire assembly into the RHV connected to the proximal hub of the guide catheter.
3. Advance microcatheter/guidewire assembly to the distal end of the guide catheter. Remove the peel away introducer by pulling the tab away from the microcatheter.
4. Use conventional catheterization techniques under fluoroscopic guidance to advance the microcatheter to the desired anatomy. Remove the guidewire.

**WARNING: Microcatheter advancement beyond the guidewire may result in vessel trauma.**

5. After completing the procedure, withdraw the microcatheter.

**CLINICAL BENEFITS AND PERFORMANCE CHARACTERISTICS**  
The Penumbra Delivery Microcatheters assist in the delivery of diagnostic agents and therapeutic devices to the vasculature.

**SAFETY AND CLINICAL PERFORMANCE**  
A copy of the Summary of Safety and Clinical Performance can be viewed by searching the device name on the EUDAMED website:  
<https://ec.europa.eu/tools/eudamed/eudamed>  
Basic UDI DI: 081454801PDMET

**MATERIALS**  
The Penumbra Delivery Microcatheters are not made with natural rubber latex.

**STORAGE AND DISPOSAL**

- Store in a cool, dry place.
- Dispose of device in accordance with hospital, administrative, and/or local government policy.

|                                                                                     |                                   |                                                                                                   |                                    |
|-------------------------------------------------------------------------------------|-----------------------------------|---------------------------------------------------------------------------------------------------|------------------------------------|
|    |                                   | Attention, see Instructions For Use                                                               |                                    |
| <div>Rx Only</div>                                                                  |                                   | Prescription only – US Federal Law restricts this device to use by or on the order of a physician |                                    |
|    |                                   | Nonpyrogenic                                                                                      |                                    |
|    |                                   | Sterilized using ethylene oxide. Single sterile barrier system with protective packaging outside. |                                    |
| <div>MD</div>                                                                       | Medical Device                    |                  | Date of Manufacture                |
|    | Manufacturer                      | <div>EC REP</div>                                                                                 | EU Representative                  |
| <div>LOT</div>                                                                      | Lot number                        |                 | Do not reuse                       |
|  | Use by                            | <div>REF</div>                                                                                    | Catalogue number                   |
|  | Do not resterilize                |                | Not made with natural rubber latex |
|  | Do not use if package is damaged. |                                                                                                   |                                    |

**WARRANTY**  
Penumbra, Inc. (Penumbra) warrants that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular use. Handling, storage, cleaning, and sterilization of this device as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond Penumbra’s control directly affect the device and the results obtained from its use. Penumbra’s obligation under this warranty is limited to the repair or replacement of this device and Penumbra shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this device. Penumbra neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. Penumbra assumes no liability with respect to devices reused, reprocessed, or resterilized and makes no warranties, expressed or implied, including but not limited to merchantability or fitness for intended use, with respect to such device.

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