
INSTRUCTIONS FOR USE / BRUGSANVISNING /
GEBRAUCHSANWEISUNG / KASUTUSJUHEND /
KÄYTTÖOHJEET / INSTRUCTIONS D’UTILISATION /
ISTRUZIONI PER L’USO / 사용 설명서 /
NAUDOJIMO INSTRUKCIJOS / LIETOŠANAS
INSTRUKCIJA / GEBRUIKSAANWIJZING /
INSTRUKSJONER FOR BRUK / BRUKSANVISNING

LEAD CAP KIT / LEDNINGSHÆTTESÆT /
ELEKTRODENKAPPEN-SET / ELEKTROODIKORGI
KOMPLEKT / JOHDON SUOJATULPPASARJA /
KIT DE CAPUCHON DE LA SOND E / KIT
CAPPUCCIO PER ELETTROCATETERE /
전극선 캡 키트 / LAIDO GAUBTO
RINKINYS / ELEKTRODA VĀCIŅU KOMPLEKTS /
LEADDOPSET / SETT MED LEDNINGSHETTER /
ELEKTRODFÖRSLUTNINGSSATS

REF 6623

Contains: / Indeholder: / Inhalt: / Komplektis: / Sisältö: / Contient: / Contiene: / 내용물: / Komplektą sudaro: / Satur: / Inhoud: / Inneholder: / Innehåller:	
	(2) 4,75-mm Lead cap / 4,75 mm ledningshætte / 4,75-mm-Elektrodenkappe / 4,75 mm elektroodikork / 4,75 mm johdon suojatulppa / Capuchon de sonde de 4,75 mm / Cappuccio per elettrocatetere da 4,75 mm / 4.75mm 전극선 캡 / 4,75 mm laido gaubtas / 4,75 -mm elektroda vāciņš / 4,75-mm leaddop / 4,75 mm ledningsshette / 4,75 mm elektroförslutning
	(2) IS-1* Lead cap / IS-1* ledningshætte / IS-1*-Elektrodenkappe / IS-1* elektroodikork / IS-1* johdon suojatulppa / Capuchon de sonde IS-1* / Cappuccio per elettrocatetere IS-1* / IS – 1 * 전극선 캡 / IS-1* laido gaubtas / IS-1* elektroda vāciņš / IS-1* leaddop / IS-1*-ledningsshette / IS-1* elektroförslutning
	(2) DF-1** Lead cap / DF-1** ledningshætte / DF-1**-Elektrodenkappe / DF-1** elektroodikork / DF-1** johdon suojatulppa / Capuchon de sonde DF-1** / Cappuccio per elettrocatetere DF-1** / DF – 1 * * 전극선 캡 / DF-1** laido gaubtas / DF-1** elektroda vāciņš / DF-1** leaddop / DF-1**-ledningsshette / DF-1** elektroförslutning

* IS-1 ISO 5841-3:2013

** DF-1 ISO 11318:2002

ENGLISH

This literature is intended for use by professionals trained or experienced in device implant and/or follow-up procedures.

WARNING: For single patient/single procedure use only. Do not reuse, reprocess, or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness, or death. Reuse, reprocessing, or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

CAUTION: The accessory is sterilized with ethylene oxide (EO) gas before final packaging. When the accessory is received, it is sterile provided the container is intact. If the container is wet, punctured, opened, or otherwise damaged, return the accessory to Boston Scientific.

CAUTION: Do Not use with other Boston Scientific lead terminals.

STERILE

Sterilized using ethylene oxide.

INTENDED USE

Use to protect and isolate terminals of an implanted Boston Scientific lead when the lead is not connected to a pulse generator.

The overall clinical benefit for this accessory is derived from the device system with which it is intended to be used.

DESCRIPTION

The lead caps are made of silicone rubber. 100% of the patient-contacting material is silicone rubber.

Surface Area

The patient contacting surface area of each lead cap in this kit is approximately 5 cm².

Expected Lifetime

The expected lifetime of each lead cap in this kit is a minimum of 10 years.

INSTRUCTIONS FOR USE

Select the lead cap that corresponds to the type of lead terminal. Place the lead cap over the terminal end of a Boston Scientific lead. Tie a suture around the lead cap groove to secure the lead cap to the terminal.

For additional information on how and when to use this accessory (including any warnings, precautions, and other safety information related to its use for the specific procedure being performed), see the literature accompanying the lead to be capped.

Any serious incident that occurs in relation to this device should be reported to Boston Scientific and the relevant local regulatory authority.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

For customers in the European Union, use the accessory name found in the labeling to search for the device’s Summary of Safety and Clinical Performance, which is available on the European database on medical devices (Eudamed) website:

https://ec.europa.eu/tools/eudamed

EXPLANT AND DISPOSAL

Explant

See the physician’s lead manual regarding the explantation and return of lead caps after lead explantation.

Disposal

CAUTION: Dispose of all items used during a procedure using standard biohazard handling techniques since they are considered biohazardous.

All items used during a procedure, such as accessories, consumables, and packaging, may be contaminated with infectious substances. Consider the following to minimize the risk of infection, microbial hazards, or physical harm:

- Biohazardous waste should be disposed in a biohazard container that is labeled with the biological hazard symbol and taken to a designated facility for biohazardous waste for proper treatment in accordance with hospital, administrative, and/or local government policy.

- Biohazardous waste should be treated with an appropriate thermal or chemical process.

NOTE: *Untreated biohazardous substances should not be disposed of in the municipal waste system.*

WARRANTY DISCLAIMER

Boston Scientific disclaims all expressed and implied warranties for this product, including without limitation any implied warranties of merchantability or fitness for a particular purpose. Buyer assumes all risk of loss or damages arising from use of this product.

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All rights reserved.
Boston Scientific Corporation
Boston Scientific (Australia) Pty Ltd
© 2021 Boston Scientific Corporation or its affiliates.