Flexiva™ Pulse ID TracTip™

MAX POWER: 50 W

MAX FREQUENCY: 80 HZ

Flexiva™ Pulse ID

MAX POWER: 50 W / 100 W

MAX FREQUENCY: 80 HZ

High Power Single-Use Laser Fiber

en	Instructions for Use	2
es	Instrucciones de uso	6
fr	Mode d'emploi	11
de	Gebrauchsanweisung	16
it	Istruzioni per l'uso	21
nl	Instructies voor gebruik	26
pt-EU	Instruções de Utilização	31

TABLE OF CONTENTS

REUSE WARNING 3

DEVICE DESCRIPTION 3

 User Information..... 3

 Table 1. Product Information 3

 Figure 1. Flexiva Pulse ID laser fiber..... 3

 Figure 2. Flexiva Pulse ID TracTip laser fiber 3

 Contents 3

 Operating Principle..... 3

 Materials..... 3

INTENDED USE/INDICATIONS FOR USE 3

 Clinical Benefits Statement..... 3

CONTRAINDICATIONS..... 3

WARNINGS..... 3

PRECAUTIONS 4

ADVERSE EVENTS 4

HOW SUPPLIED 4

 Handling and Storage..... 4

OPERATIONAL INSTRUCTIONS 5

 Additional Required Items..... 5

 Preparation..... 5

 Procedure 5

 Intra-operative Fiber Removal and Inspection 5

 Disposal 5

 Post-Procedure..... 5

PATIENT COUNSELING INFORMATION 5

WARRANTY..... 5

SYMBOL DEFINITIONS..... 5

 Table 2. Symbol Definitions 5

Flexiva™ Pulse ID TracTip™

MAX POWER: 50 W

MAX FREQUENCY: 80 HZ

Flexiva™ Pulse ID

MAX POWER: 50 W/100 W

MAX FREQUENCY: 80 HZ

High Power Single-Use Laser Fiber

Rx ONLY

REUSE WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient serious bodily injury or illness. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to serious bodily injury or illness of the patient. After use, see *Disposal* to properly dispose of the device.

DEVICE DESCRIPTION

Flexiva Pulse ID and Flexiva Pulse ID TracTip laser fibers are high power single-use fiber optic laser energy delivery devices consisting of a SMA-905 connector, strain relief, and a silica core fiber jacketed with ethylene tetrafluoroethylene (ETFE). The Flexiva Pulse ID laser fibers are equipped with a polished, flat output tip. The Flexiva Pulse ID TracTip laser fibers are equipped with a polished and reinforced ball-shaped output tip designed to allow passage through a fully deflected scope. These fibers may be used in a variety of laser-based surgical cases as an integral part of laser systems. Fiber is recognized by open and closed systems. A RFID (Radio-frequency identification) tag enables read/write data storage for compatible RFID-equipped laser systems (closed systems), to recognize the connected laser fiber, control the maximum energy delivered based on fiber capability, and limit the fiber to a single use.

User Information

This device is intended to be used by, or under the supervision of a skilled physician with experience in urological laser surgery. Physicians are assisted by surgical technicians and nursing teams. For use with Ho:YAG laser systems with a standard SMA-905 connector that have been cleared for surgical use.

Recommended for use with compatible open and compatible Lumenis manufactured Ho:YAG laser systems including*:

- Lumenis Pulse™ 120H
- Lumenis Pulse™ 50H
- Lumenis Pulse™ 20
- Lumenis Pulse™ 100H
- Lumenis Pulse™ 30H
- VersaPulse™ PowerSuite™

*This is not an exhaustive list.

Please refer to the laser system's User Manual for complete information regarding applications, contraindications, precautions, and warnings. Before using the Flexiva Pulse ID and Flexiva Pulse ID TracTip laser fibers, the physician should fully understand the use of the holmium (Ho:YAG) laser, safety considerations, tissue interactions, and proper technique specific to the treatment for which the physician intends to use this product.

Table 1. Product Information

Product Description	Core Size (µm)	Cladding Outside Diameter (µm)	Jacket Outside Diameter (µm)	Blade Diameter (µm)	Maximum Power	
					Hertz (Hz)	Watts (W)
Flexiva Pulse ID TracTip	242 750 (TracTip Outside Diameter)	288	480 +/- 20	406 (0.016 in)	80	50
Flexiva Pulse ID 242	242	288	480 +/- 20	406 (0.016 in)	80	50
Flexiva Pulse ID 365	365	409	600 +/- 50	508 (0.020 in)	80	100
Flexiva Pulse ID 550	550	614	1000 +/- 50	711 (0.028 in)	80	100
Flexiva Pulse ID 910	910	1014	1400 +/- 50	1117 (0.044 in)	80	100


Figure 1. Flexiva Pulse ID laser fiber


Figure 2. Flexiva Pulse ID TracTip laser fiber

1. Laser Input
2. Fiber Connector Hub
3. Fiber with Jacket
4. Distal Output Tip (Flat Tip (Fig. 1) or TracTip (Fig. 2))

Contents

- (1) High power single-use laser fiber with SMA-905 connector; or,
- (5) High power single-use laser fibers with SMA-905 connectors.

Operating Principle

Flexiva Pulse ID and Flexiva Pulse ID TracTip laser fibers deliver holmium (Ho:YAG) laser energy endoscopically from a laser console to soft tissues or urinary calculi. The energy travels within the fiber's glass core and exits through the distal end of the fiber. On impact the holmium energy fragments urinary calculi and/or ablates soft tissue.

Materials

Flexiva Pulse ID and Flexiva Pulse ID TracTip laser fibers are constructed of medical grade materials, a silica core fiber jacketed with ethylene tetrafluoroethylene (ETFE), tested to be biocompatible for their intended use.

INTENDED USE/INDICATIONS FOR USE

Flexiva Pulse ID and Flexiva Pulse ID TracTip laser fibers are intended to be used as a device that transmits Ho:YAG laser energy from cleared laser consoles to urological anatomy. Flexiva Pulse ID and Flexiva Pulse ID TracTip laser fibers are indicated for urologic applications for which the laser systems are cleared, limited to endoscopic procedures involving vaporization, ablation, hemostasis, coagulation, excision, resection, incision of soft tissue, and lithotripsy of urinary calculi. The fiber is designed for use with a standard SMA-905 connector and has been cleared for surgical use.

Clinical Benefits Statement

Clinical benefits include stone clearance, improved urodynamics, and improved quality of life.

CONTRAINDICATIONS

Use of the Flexiva Pulse ID and Flexiva Pulse ID TracTip laser fibers is contraindicated for treatment of patients for whom endoscopic procedures are contraindicated including the following:

- Uncontrolled bleeding disorders and coagulopathy
- Untreated urinary tract infection (UTI)
- Carcinoma of the prostate

WARNINGS

- Carefully read all instructions before use. Observe all warnings and precautions noted throughout these and other instructions relevant to the procedure. Failure to do so may result in complications.
- Use of opened, damaged, or expired product could expose the patient, user, procedure room personnel, or a combination of them, to a non-sterile laser fiber, or damaged product which could cause infection or injury.
- Improper use of the device or use of a damaged device may result in severe eye or tissue damage, fire in the procedure room and accidental laser exposure to the procedure room personnel or patient.
- Do not bend fiber at sharp angles. Visible light seen leaking from the fiber indicates fiber damage.
- Immediately discontinue use if breaks or fractures appear in the laser fiber. These breaks or fractures may allow undirected emission of laser energy, rendering the distal tip useless and causing harm to surrounding tissues.
- Do not use this laser fiber in the presence of flammable anesthetics, combustible materials, or in an oxygen enriched environment due to risk of fire or explosion.
- Gently wrap, hold, or coil the length of the laser fiber to avoid tangling. The fiber should not be clamped with forceps or other securing instruments as it could damage the fiber.

- The accumulation of blood and tissue residue can overheat the fiber creating a patient burn risk, and can result in fiber breaks, sharp edges, low or ineffective laser radiation and fiber tip degradation. In this case make sure the laser system is in "STANDBY" mode then carefully remove the laser fiber from the scope. Use a soft, sterile cloth to clean the glass tip, avoiding side pressure that may damage the tip. Do not scrape the laser fiber with sharp instruments as damage may result. The laser fiber tip can also be shortened with a compatible tool.
- Perforation may occur while executing laser treatment near sensitive areas (blood vessels, intestine). The laser energy should only be applied for the time to achieve the desired effect.
- Never leave large detached fiber pieces inside a patient as they may cause stone formation, obstruct the urinary system, require stent placement, or a combination of these issues.
- Use caution when using the laser fiber in pediatric patients and pregnant patients as the device has not been studied in these populations.
- The patient and all personnel within the procedure room must wear protective laser glasses appropriate for Ho:YAG lasers throughout the procedure. Failure to do so may result in ocular injury. Refer to the laser system's User Manual for requirements concerning eye protection.
- Use irrigation throughout the procedure to absorb any heat produced and to clear debris and particulates from the field of view. Avoid excessive irrigation that could cause the patient to absorb too much fluid.
- For tip cleaning, use care when removing the fiber from the scope to avoid contact with unsterile areas or tools, or damaging the fiber.
- After use, this device and any detached fiber pieces can pose a biological and sharps hazard. Take precautions to handle and properly dispose of these directly into a sharps disposal container labeled with a biological hazard symbol. Sharps and biological hazard waste should be safely disposed by using available sharps and biological waste channels in accordance with hospital, administrative, and/or local government policy.

If intra-operative fiber removal and inspection is performed:

- Whenever the fiber is removed from the scope, use a sterile technique to prevent compromising fiber sterility and compatible tools to avoid damaging the fiber.

PRECAUTIONS

- During all steps of the procedure, including disposal, handle the fiber with care.
- In soft tissues, increased laser exposure time will result in a deeper and wider zone of necrosis.
- If tissue or blood appears to be stuck to the lasing tip, damage to the laser fiber can result. Put the laser console in "STANDBY" mode then carefully remove the tip from the target tissue. Remove the laser fiber from the scope. Use a soft, sterile cloth to clean the glass tip, avoiding side pressure that may damage the tip. Do not scrape the laser fiber with sharp instruments as damage may result. The laser fiber tip can also be shortened with a compatible tool.
- Care should be exercised with the glass tip to avoid severe impact or side stresses that may fracture the tip.
- To prevent the tip from sticking to soft tissue, move the distal tip back and forth or side to side during each lasing interval.
- When handling the fiber, accessories, tools, disposables, or materials that were in contact with the patient, take protective measures to prevent cross-contamination.
- Observe all safety information provided in the laser system's User Manual. Failure to thoroughly understand and follow all instructions, precautions, and warnings may result in injury to the patient and/or user, and/or may result in damage to or malfunction of devices.
- The general regulations and information for the safe handling of laser energy must be applied to avoid injury to the patient and/or user, and/or may result in damage to or malfunction of devices. Use the laser fiber only under operating room conditions.
- When using the fiber in contact surgery procedures, avoid forceful movements on the tissue surface and lateral pressure on the fiber tip.
- Do not exceed the recommended maximum power levels and frequency limits associated with the selected fiber size.
- Ensure the laser fiber tip and target stone or soft tissue are properly visualized when applying laser energy.
- Ensure that the fiber tip and aiming beam are in clear view and at a safe distance from the tip of the endoscope when laser energy is being applied to avoid scope damage.
- Avoid advancing or removing the Flexiva Pulse ID TracTip laser fiber through a deflected scope when:
 - scope deflection is greater than 270°
 - scope bend diameter is less than 1.9 cm (0.748 in)
 - scope working channel is less than 3.6F
- Do not pass a laser fiber with a flat tip through deflected scopes to avoid excessive pressure on the scope liner and potential scope and/or fiber damage.

- Avoid direct laser energy contact with the endoscope and endoscopic accessories. Failure to do so may damage the endoscope and endoscopic accessories.
- If resistance is encountered during advancement or withdrawal, stop. Do not continue without first determining the cause of resistance and taking remedial action.
- Avoid repassing the Flexiva Pulse ID TracTip laser fiber through a deflected scope after energy has been applied to the treatment area as the TracTip may not maintain its original shape and may cause scope damage.
- Take caution. In the event of no output energy, the laser fiber connector hub may be hot to touch.

ADVERSE EVENTS

Potential adverse events associated with Ho:YAG laser and fibers may include but are not limited to:

- Adhesions
- Allergic Reaction
- Bladder Spasms
- Burns
- Discomfort
- Dysuria
- Edema
- Embolism
- Erectile Dysfunction
- Extravasation
- Fever
- Fistula
- Fluid Overload (potential sequelae include congestive heart failure, edema, electrolyte imbalance)
- Healing Impaired
- Hematoma
- Hematuria
- Hemorrhage
- Hydronephrosis
- Hypertension
- Incontinence
- Infection
- Inflammation
- Laceration
- Nerve Injury
- Obstruction
- Ocular Injury
- Pain
- Perforation
- Reflux, Genitourinary (GU)
- Restenosis
- Retrograde Ejaculation
- Stenosis, Non-vascular
- Sterility/Infertility
- Therapeutic Response, Altered
- Tissue Damage
- Unretrieved Device Fragments
- Urgency
- Urinary Retention
- Vasovagal Response

HOW SUPPLIED

The Flexiva Pulse ID and Flexiva Pulse ID TracTip laser fibers are supplied sterile as single-use products not intended or validated for reprocessing or reuse.

- Do not use if package is damaged or unintentionally opened before use.
- Do not use if labeling is incomplete or illegible.
- Do not use past the use by date displayed on product label. See *Disposal* to properly dispose of this device.

Handling and Storage

This product has no special handling or storage requirements.

OPERATIONAL INSTRUCTIONS

Before activating the laser, review the laser system’s User Manual for instructions for proper set-up and operation of the laser system.

Additional Required Items

Additional devices, accessories, equipment, or materials for the safe use of this device are:

- Ho:YAG laser systems with a standard SMA-905 connector that have been cleared for surgical use.
- Recommended for use with compatible open and compatible Lumenis manufactured Ho:YAG laser systems.
- Protective laser glasses.

WARNING

Do not use this laser fiber in the presence of flammable anesthetics, combustible materials, or in an oxygen enriched environment due to risk of fire or explosion.

WARNING

The patient and all personnel within the procedure room must wear protective laser glasses appropriate for Ho:YAG lasers throughout the procedure. Failure to do so may result in ocular injury. Refer to the laser system’s User Manual for requirements concerning eye protection.

Preparation

This single-use device is intended for use on patients under anesthesia in sterile conditions using common urology scopes.

1. Identify and retrieve the fiber model to be used in the procedure and confirm the label information is correct on the device retrieved. Ensure that the sterile packaging is not open, torn, or punctured, and the product is not damaged. If necessary, return the device to Boston Scientific for replacement.
2. Using a sterile technique, carefully open the fiber pouch to prevent compromising fiber sterility. Remove the hoop from the pouch and present to a sterile field. Use care removing the fiber connector hub from the hoop clip using an outward motion and avoid dropping or tangling it as damage may occur if the fiber is struck or bent sharply. Inspect and treat the output tip with particular care as it represents the most-delicate, easily damaged portion of the assembly.
3. Hold the fiber connector hub and remove the protective cap with care to avoid touching the exposed connector end. Do not hold the fiber from the rubber strain relief or the fiber body.
4. Make sure the laser console is “OFF” or on “STANDBY” mode before inserting the fiber into the SMA port.
5. Hold the fiber connector hub and screw it onto the fiber port on the laser console. Finger-tighten to avoid applying excessive force or over-tightening and breaking the fiber connector hub. The console’s microswitch should trigger when the optical coupling is secure.
6. If the laser console is in “OFF” mode, turn the laser console “ON” then leave in “STANDBY” mode.
7. Activate the laser aiming beam, adjust aiming beam brightness, and verify it is not dim or misshapen to confirm the laser fiber functions as intended. Hold the fiber tip to a non-reflective surface and ensure that a circular colored spot appears. If the spot is dim, not visible, or misshapen, confirm console settings are correct, and if so, do not use and return the device to Boston Scientific for replacement.
8. Using a sterile technique, manage, hold, coil, or secure the fiber length during fiber visual inspection using the aiming beam.
9. Carefully inspect the fiber for kinks, punctures, fractures, a broken tip (flat or ball), jacket/fiber or other particulates/pieces, or other visual damage. If the fiber appears damaged, do not use the device. If it is an unused device, return it to Boston Scientific for replacement.
10. Set the laser treatment parameters, as instructed in the laser system’s User Manual and within the appropriate power level for the fiber size. Note, not all console settings may be available to use.

Procedure

1. Verify the scope adaptor is in the open position for Flexiva Pulse ID fibers. Insert the fiber through the working channel of the endoscopic delivery system until the fiber tip is in clear view. Flexiva Pulse ID laser fibers with core sizes 242 (flat tip) and 365 can be advanced through a flexible scope provided the scope is straightened as much as possible during fiber insertion.
2. The Flexiva Pulse ID TracTip laser fiber can be advanced through a fully deflected flexible scope.

Precaution: Avoid advancing or removing the Flexiva Pulse ID TracTip laser fiber through a deflected scope when:

- scope deflection is greater than 270°
- scope bend diameter is less than 1.9 cm (0.748 in)
- scope working channel is less than 3.6F

3. Position the fiber tip and/or aiming beam on the target treatment area.
4. Place the laser console in “READY” mode.
5. Depress the footswitch to deliver the treatment beam.
6. During the procedure, check fiber tip quality, such as fiber tip length, sharp edges, adhesion of debris to the fiber tip, and jacket/fiber particles. If damage cannot be corrected, stop using the device.

Intra-operative Fiber Removal and Inspection

WARNING

Whenever the fiber is removed from the scope, use a sterile technique to prevent compromising fiber sterility and compatible tools to avoid damaging the fiber.

See *Table 1. Product Information*, to select tools compatible with your fiber. See *Preparation* for inspection steps. See *Disposal* to properly dispose of the device.

1. Before removing the fiber from the scope, place the console in “STANDBY” or “OFF” mode. It is recommended to straighten the scope to avoid potential scope and/or fiber damage.
2. Using a sterile technique, remove the fiber from the scope carefully to avoid compromising fiber sterility or damaging the fiber. It is recommended to straighten the scope if deflected to avoid potential scope and/or fiber damage.
3. Whenever the fiber is removed from the scope, inspect fiber tip quality. It is also recommended to complete the aiming beam check steps before re-inserting the fiber into the scope. See *Preparation* for inspection steps. If the fiber is damaged, do not use the device. See *Disposal* to properly dispose of the device.

Disposal

To minimize the risk of infection or microbial hazards after use, dispose device and packaging as follows:

After use, device may contain biohazardous substances. The device and packaging should be treated and disposed of as biohazardous waste or have them treated and disposed of in accordance with any applicable hospital, administrative, and/or local government regulations. Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system.

This device contains a sharps hazard. Take precautions to ensure that sharps are handled properly. Dispose of all sharps directly into a sharps disposal container labeled with a biological hazard symbol. Sharps waste should be safely disposed of using available sharps waste channels in accordance with hospital, administrative, and/or local government policy.

Post-Procedure

Verify the console is in “STANDBY” or “OFF” mode. If the scope is deflected, it is recommended to straighten the scope to avoid potential scope damage then remove the fiber from the scope. Disconnect the fiber connector hub from the console port. The Flexiva Pulse ID family of laser fibers are single-use devices. After use, dispose of the fiber and fiber fragments (sharps). See *Disposal* to properly dispose of the device. Any serious incident that occurs in relation to this device should be reported to Boston Scientific and the relevant local regulatory authority.

PATIENT COUNSELING INFORMATION

Brief the patient on any relevant post-procedure instructions, contraindications, warnings, precautions, and/or adverse events found within this fiber IFU, pertaining to the patient. Also, refer to the applicable laser system’s User Manual.

WARRANTY

For device warranty information, visit www.bostonscientific.com/warranty.

The following are trademarks of Boston Scientific Corporation or its affiliates: Flexiva, TracTip. All other trademarks are the property of their respective owners.

SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at www.bostonscientific.com/SymbolsGlossary. Additional symbols are defined in Table 2.

Table 2. Symbol Definitions

	Contents
---	----------

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All rights reserved.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd
© 2024 Boston Scientific Corporation or its affiliates.