Medtronic

1X8 SC	3776
	3876
1X8	3777
	3877
1X8 COMPACT	3778
	3878

Low impedance lead kit
Kit d'électrode faible impédance
Elektrodenset mit niedriger Impedanz
Equipo de electrodo de baja impedancia
Geleidingsdraadset met lage impedantie
Kit dell'elettrocattetera a bassa impedenza
Elektrodset med låg impedans
Σετ ηλεκτροδίου χαμηλής σύνθετης αντίστασης
Elektrodesæt med lav impedans
Conjunto de electrodos de baixa impedância



Implant manual • Manuel d'implantation
Implantationsanleitung • Manual de implantación
Implantatiehandleiding • Manuale per l'impianto
Implantationsmanual • Εγχειρίδιο εμφύτευσης
Implantationshåndbog • Manual de implante

USA Rx only

A13106001.fm 11/23/05 2:26 pm
UC200xxxxxx EN
4 x 8 inches (101 mm x 203 mm)

Medtronic Confidential
LeadExtTemplate.fm


MA13106A002 Rev A

Explanation of symbols on product or package labeling


Refer to the appropriate product for symbols that apply.

 Open here

  Lead length


 Do not reuse

 Sterilization: ethylene-oxide gas


 Caution, consult accompanying documents.

 Manufacturing date

 Use by

 Storage temperature

 Lot number

 Conformité Européenne (European Conformity).
This symbol means that the device fully complies
with European Directive AIMD 90/385/EEC.

 For USA audiences only

Medtronic® is a registered trademark of Medtronic, Inc.

Table of contents

Device description 5

Package contents 5

Device specifications 5

Instructions for use 8

Preparing for surgery 8

Placing a percutaneous lead 8

Testing stimulation intraoperatively 10

Anchoring the lead 11

Tunneling the percutaneous extension or the lead for test stimulation 13

Connecting the percutaneous extension (if used) to the lead 14

Tunneling completion and incision closure 15

Tunneling the lead for direct connection to a neurostimulator 16

Refer to the indications sheet for indications and related information.

Refer to the appropriate information for prescribers booklet for contraindications, warnings, precautions, adverse events summary, individualization of treatment, patient selection, use in specific populations, resterilization, and component disposal.

Refer to System Eligibility, Battery Longevity, Specifications reference manual packaged with the software application card for neurostimulator selection, battery longevity calculations and specific neurostimulator specifications.

USA Refer to the clinical summary booklet packaged with the neurostimulator for information on the clinical study results of the neurostimulation system and individualization of treatment.

Device description

The Medtronic Models 3776 / 3876 1 X 8 SC, 3777 / 3877 1 x 8, and 3778 / 3878 1 x 8 Compact Leads are part of a neurostimulation system for pain therapy.

Package contents

- Lead
- Percutaneous extension (3876, 3877, 3878 only)
- Anchors:
 - two-wing
 - three-wing
 - EZ
 - twist-lock
- Connector boots, 1 x 8 low-profile (2)
- Guide wire
- Needle, 14 gauge, 8.9 cm (epidural, modified Tuohy)
- Stylets:
 - curved at a 28° angle
 - straight, 0.25 mm diameter (inserted in the lead)
- Tunneling tools: tunneling rod, tunneling tip (narrow), tube (passing straw)
- Wrench, torque
- Product literature
- Warranty card (USA only)
- Registration form (USA only)

Device specifications

The lead has electrodes on the distal end; the proximal end fits into an eight-conductor connector. A stylet has been inserted into the proximal end of the lead to aid in positioning the lead.


Figure 1. Electrode, lead contact, and stylet.

Table 1. Device specifications^a for lead: Models 3776 / 3876, 3777 / 3877, and 3778 / 3878

Description	Model 3776 / 3876	Model 3777 / 3877	Model 3778 / 3878
Connector	Octapolar, in-line	Octapolar, in-line	Octapolar, in-line
Conductor resistance ^{b,c}	Maximum 16 Ω for all lengths	Maximum 16 Ω for all lengths	Maximum 16 Ω for all lengths
Length	20 – 100 cm	20 – 100 cm	20 – 100 cm
Diameter	1.3 mm	1.3 mm	1.3 mm
Distal end			
Number of electrodes	8	8	8
Electrode shape	Cylindrical	Cylindrical	Cylindrical
Electrode length	3.0 mm	3.0 mm	3.0 mm
Electrode spacing	1.5 mm	6.0 mm	4.0 mm
Electrode distance	34.5 mm	66.0 mm	52.0 mm
Distal tip distance	Not applicable	Not applicable	Not applicable
Radiopaque marker	No	No	No
Proximal end			
Lead contact length	1.5 mm	1.5 mm	1.5 mm
Lead contact spacing	2.8 mm	2.8 mm	2.8 mm
Lead contact distance	22.5 mm	22.5 mm	22.5 mm
Stylet handle length	29.0 mm	29.0 mm	29.0 mm
Stylet wire diameter	0.25, 0.30 mm	0.25, 0.30 mm	0.25, 0.30 mm

^a All measurements are approximate.

^b Electrical resistance of this device only.

^c Resistance is proportional to length: long lengths have higher resistance that may limit the amplitude.

Table 2. Material of components in the Models 3776 / 3876, 3777 / 3877, and 3778 / 3878 package

Component	Material	Material contacts human tissue
Lead		
Conductor wire	Silver core MP35N	No
Conductor wire insulation	Fluoropolymer	No
Electrodes	Platinum-iridium	Yes
Insulation	Polyurethane	Yes
Proximal connector	MP35N	Yes
Lead tip	Epoxy	Yes
Adhesive	Epoxy	Yes
Percutaneous extension^a		
Conductor wire	Silver core MP35N	No
Conductor wire insulation	Fluoropolymer	No
Insulation	Polyurethane, Silicone rubber	Yes
Setscrew connector block	Stainless steel	Yes
Setscrew	Titanium	Yes
Spring contacts	MP35N, Platinum-iridium	Yes
Anchor		
Two-wing, three-wing, EZ	Silicone rubber	Yes
Twist-lock	Polysulfone	Yes
Connector boot	Silicone rubber	Yes
Guide wire	Stainless steel	Yes
Needle	Stainless steel	Yes
Stylets (curved, straight)		
Handle	Polypropylene	Yes
Wire	Stainless steel with parylene coating	Yes
Tunneling tools		
Tunneling rod	Stainless steel	Yes
Tunneling tip	Stainless steel	Yes
Tube (passing straw)	Fluoropolymer	Yes
Wrench, torque		
Handle	Polymer	Yes
Shaft	Stainless steel	Yes

^a Applies to Models 3876, 3877, and 3878 only.

Instructions for use

Implanting physicians should be experienced in epidural-access procedures and should be thoroughly familiar with all product labeling.

△ Cautions:

- Do not bend, kink, or stretch the lead or extension body, which may damage the component.
- Do not use any instrument to handle the lead. The force may compress the lead, resulting in the inability to insert or withdraw the stylet.
- Do not bend or kink the stylet, which may damage the component.
- When reinserting a stylet, use care because too much pressure on the stylet could damage the lead, resulting in intermittent or loss of stimulation.
- Do not use saline or other ionic fluids at connections, which could result in a short circuit.

Preparing for surgery

Before opening the lead package, verify the model number, use-by date, lead-length, and connector type.

Placing a percutaneous lead

△ **Caution:** To reduce the risk of lead damage that may result in intermittent stimulation or loss of stimulation, requiring additional surgery to replace the lead:

- use only the modified Tuohy needle supplied in the kit.
- use a shallow needle-insertion angle (45° or less) when inserting or withdrawing the needle into or out of the epidural space.

△ **Warning:** As with any spinal procedure, the risk of serious injury to the patient (eg, hemorrhage, hematoma, or paralysis) increases as location of the selected needle insertion site progresses up the vertebral column - from a lower risk at a lumbar location to a higher risk at a cervical location. Select a vertebral location that provides the widest and easiest access to the epidural space during needle insertion to reduce the risk of serious patient injury resulting from direct trauma to the spinal cord.

1. Make an incision at the needle-entry site to the depth of the subcutaneous fascia.
2. Using a paramedial approach and under fluoroscopy, insert the needle (included in the kit) into the epidural space at the appropriate angle until you encounter resistance from the ligamentum flavum.
Note: Midline placement may cause passive damage to the lead over time because of ligament or spinous process movement.
3. Confirm needle location under fluoroscopy.
4. After rotating the needle so that the beveled edge faces cephalad, remove the needle stylet.
5. Advance the needle and confirm entry into the epidural space (eg, using the loss-of-resistance technique with air or sterile [United States Pharmacopeia—USP] water).

△ **Caution:** Do not use contrast media or a saline flush. Contrast media may obscure the field of view and a saline flush may increase the difficulty of lead placement.

6. For a second lead, repeat steps 1 – 5 noting these recommendations:
 - Implant the second lead parallel to the first lead and approximately 1 – 3 mm lateral of the physiological midline.
 - Introduce the second lead one vertebral space below the first lead to help prevent nicking or cutting the first lead and to allow sufficient space for suturing both lead anchors.
 - Stagger the lead tips or place them several vertebral spaces apart, depending on the position that produces the most effective paresthesia.
7. After inserting the guide wire through the needle, advance the guide wire no farther than 1 – 3 cm past the needle tip. Next, remove the guide wire from the needle.

Note: If the guide wire track deviates from the intended pathway, steering and manipulating the lead will be more difficult.
8. Using fluoroscopy, slowly insert the lead through the needle and advance the lead to the initial target placement site. A stylet may need to be reinserted.

Notes:

 - The curved stylet is keyed so that the flat notch of the stylet handle faces the same direction as the distal curved tip.
 - When using a curved stylet and resistance is encountered during lead advancement, exchange the curved stylet for the straight stylet and use short, firm movements to advance the lead.
9. After verifying the lead position under fluoroscopy (anterior-posterior and lateral views), compare that location with the location that has the highest probability of paresthesia coverage.

Notes:

 - To increase lead stability, insert enough lead length to extend at least three vertebral bodies into the epidural space.
 - Position the lead so that the center bipole electrodes (eg, electrodes 3 and 4 on an octapolar lead) will be active. If lead migration occurs, effective stimulation may be regained through programming adjustments instead of surgical repositioning.

Testing stimulation intraoperatively

Note: This procedure requires a snap-lid connector cable, which is available as an accessory kit.

△ **Caution:** To prevent possible uncomfortable or unexpected stimulation (jolting or shocking sensation):

- program parameter changes in small increments above the perception threshold (the amplitude value(s) at which the patient first perceives paresthesia).
- decrease the amplitude(s) to 0.0 V before:
 - changing electrode polarities.
 - connecting or disconnecting the cable to the external neurostimulator.
 - replacing the external neurostimulator batteries.
 - turning ON the neurostimulator.

1. Use the lid latch to open the snap-lid connector cable's snap-lid connector (Figure 2a).


Figure 2. Snap-lid connector cable.

2. While holding the lead in place, disconnect the stylet handle from the lead (proximal end), and partially withdraw the stylet.
3. Align the proximal end of the lead against the inside end of an 8-contact lead groove of the connector (Figure 2a).
4. Check that the lead contacts line up with the electrical contacts inside the lead groove.

Note: The lead fits only one way into the snap-lid connector.

5. Press the lead gently into the lead groove (Figure 2b). If a second lead is used, repeat steps 2 – 5.
6. Close the lid until the latch snaps shut in order to seat the lead(s) in the lead groove(s).

Note: Do not force the lid closed; it should close easily. If it does not, disassemble the components and repeat steps 1-6 in this section.

7. After connecting the cable to the patient's implanted lead, pass the plug end out of the sterile field.
8. After verifying that the external neurostimulator is turned OFF, connect the cable plug to the receptacle of the external neurostimulator.

△ **Warning:** Maintain adequate slack in the cable. If there is not enough slack and the cable is pulled, the percutaneous components may be dislodged.

9. Identify optimal stimulation parameters, beginning at a pulse width of 210 – 240 μ s and a rate of 30 Hz.

Note: Ensure that the patient can provide immediate feedback.

10. Increase the amplitude beginning from 0.0 V while asking the patient close-ended questions to identify the perception threshold (the amplitude at which the patient first perceives paresthesia), the discomfort threshold (the amplitude at which paresthesia is beyond the patient's tolerance), and the paresthesia coverage.

Note: If good paresthesia coverage is not attained, change electrode settings before repositioning the lead to confirm the direction of lead movement.

11. If two leads were placed, repeat steps 9 – 10 for the second lead; optimize paresthesia coverage using both leads.
12. In the patient's chart, document the lead position that provided appropriate stimulation coverage (ie, record the settings and patient responses and include a fluoroscopic image of the final lead position).
13. Disconnect the cable from the lead by opening the snap-lid connector and gently lifting the lead(s) from the electrical contacts in the lead groove. Leave the cable in the sterile field for additional parameter testing before closing.

△ **Caution:** Do not pull on the cable. Pulling on the cable may break a wire or dislodge the lead. A broken wire or dislodged lead may result in loss of stimulation and may require surgery to replace the lead.

Anchoring the lead

△ **Caution:** Do not tie ligatures around the lead or extension body, which may damage the insulation.

1. Prepare the anchor site by making a 5 – 7 cm (2.0 – 2.8 in) longitudinal incision around the needle shaft, dissecting down to the supraspinous ligament, and establishing hemostasis.
2. Disconnect the stylet from the lead. Partially expose, but do not completely withdraw the stylet.
3. While maintaining lead position, use minimal force to remove the needle and stylet.

△ **Caution:** Use minimal traction to remove the needle because quick or sudden removal may dislodge the lead.

4. Slide the anchor onto the proximal end of the lead and continue sliding the anchor down as close as possible to where the lead emerges from the vertebral column. Use care to maintain the lead position.

Notes:

- The anchor can be lubricated with sterile water.
- For a twist-lock anchor, ensure that the anchor is unlocked before sliding the anchor onto the lead. (An unlocked anchor is twisted so the two parts are at a 90° angle to each other.)
- An EZ anchor can be trimmed to a shorter length.

5. Use 2-0 nonabsorbable suture (such as silk or some types of braided polyester mesh) to secure the anchor.

△ **Cautions:**

- Do not use polypropylene suture material on silicone components. Polypropylene may damage the component, resulting in component failure.
- Do not overtighten ligatures on the anchor or connector boot. Ligatures that are too tight may damage the component.

For the twist-lock anchor: Suture the anchor to the supraspinous ligament or deep fascia. Next, lock the anchor to secure the lead.

Note: If minor repositioning is necessary, unlock the lead, then relock the anchor.

For the EZ anchor: Tie a ligature around the distal end of the anchor to secure the anchor to the lead. If needed, make minor adjustments to the lead position. Then, use four or more ligatures to secure the anchor to the supraspinous ligament or deep fascia. Last, suture the proximal end of the anchor.

For the winged anchor: Tie ligatures around the grooves in the middle of the anchor to secure the anchor to the lead. Next, suture the anchor to the supraspinous ligament or deep fascia.

Note: Laboratory testing has shown that injecting silicone medical adhesive between the anchor and lead after tying the ligatures increases the anchor-to-lead holding force.

6. Verify test stimulation parameters to ensure that the lead has not moved: connect the snap-lid connector cable to the lead, then connect the cable plug to the external neurostimulator. If the lead has moved, reposition it.

Tunneling the percutaneous extension or the lead for test stimulation

1. After simulating the tunneling route, mark the patient's skin at the lead-extension connection site and at the percutaneous extension or lead exit site. The exit site should be on the side opposite that intended for the neurostimulator and at least 10 cm (4 in) lateral. (The neurostimulator should be placed on the opposite side of the body from another active implanted device and should be placed preferably on the right side of the body to allow for future placement of cardiac devices on the patient's left side.)

△ **Cautions:**

- When routing the extension or lead, avoid sharp bends or kinks, which may break the wires. Broken wires may create an open circuit, resulting in loss of stimulation or component failure and requiring surgical replacement.
- When two leads are implanted, route the lead-extensions so the area between them is minimized (Figure 3). If the lead-extensions are routed in a loop and the patient is exposed to some sources of electromagnetic interference (eg, theft detectors), the patient may perceive a momentary increase in stimulation, which some patients have described as uncomfortable stimulation (jolting or shocking sensation).


Figure 3. Routing for two lead-extensions (or two leads).

2. After assembling the tunneling tools, make an incision at the exit site.
3. Begin at the exit site and tunnel subcutaneously to the lead incision.

Notes:

- Deep tunneling is not desirable.
 - Avoid the lower thoracic ribs.
 - If the tunneling tool does not extend to the lead incision site, make an intermediate incision.
4. Unscrew the tunneling tip; then, remove the tunneling tool, leaving the passing straw in place in the tunnel.
 5. Gently insert the proximal end of the extension or lead through the passing straw to the exit site.

Note: If tunneling the percutaneous extension, use blunt dissection to form a tunnel from the lead incision to the lead-extension connection site.

Connecting the percutaneous extension (if used) to the lead

△ **Caution:** Before connecting components, wipe off any body fluids and dry all connections. Fluids in the connections may result in stimulation at the connection site, intermittent stimulation, or loss of stimulation.

1. Slip the connector boot onto the lead body with the open end of the connector boot facing outward, exposing the lead contacts.

Notes:

- Irrigate the connector boot with a nonionic antibiotic solution. Dry the interior of the boot.
- A clear and white radiopaque boot are provided for use in dual-lead systems.

2. Wipe the lead body and extension setscrew connector junction with sterile gauze. If necessary, use sterile (USP) water or a nonionic antibiotic solution. Dry all connections.

3. Insert the lead connector fully into the extension setscrew connector junction.

Notes:

- Each lead contact must be aligned under each extension connector contact.
- During insertion, some resistance is typical.

4. Use the torque wrench to tighten the setscrew, which completes the electrical circuit with the lead contacts.

△ **Cautions:**

- To avoid overtightening, do not use a hex wrench to tighten an extension setscrew. Overtightening an extension setscrew may damage the lead contacts and cause an open or short circuit, resulting in intermittent or loss of stimulation.
- Discard the torque wrench after making all connections. Reusing a torque wrench may result in undertightening or overtightening and subsequently, intermittent or loss of stimulation.

5. Confirm the integrity of the connection by checking the impedance of all electrode pairs with the clinician programmer.

6. After sliding the connector boot over the lead-extension connection, tie a ligature around the wide end of the connection.

Note: Do not suture the narrow end of the connector boot.

Tunneling completion and incision closure

1. Carefully pull the lead-extension connection or lead through the tunnel.

△ **Cautions:**

- Do not pull the lead or extension taut. Allow enough slack in the lead or extension to accommodate patient movement. Pulling the lead or extension taut may result in a short or open circuit or migration of implanted components.
- Coil the excess lead body into a circular loop greater than 2 cm (0.8 in) in diameter. Do not use a U-shaped loop or bend (Figure 4). A circular loop decreases the possibility of electromagnetic interference and its effects and prevents kinking or damaging the lead body.


Figure 4. Coiling excess lead.

2. After carefully removing the passing straw from the tunnel, verify that the location of the lead and the pattern of stimulation have not changed.
3. Close the lead incision and the exit site, leaving the percutaneous extension or lead wires looped under a dressing with the connector contacts exiting the bandage.
4. Connect the snap-lid connector cable to the percutaneous extension or lead so that the cable is available for test stimulation.

Tunneling the lead for direct connection to a neurostimulator

Creating a subcutaneous pocket for the neurostimulator

1. Make an incision equal to the length of the neurostimulator.
2. Use blunt dissection to create a subcutaneous pocket.

Note: Refer to the neurostimulator implant manual for the proper subcutaneous pocket depth.

Tunneling the lead

1. After simulating the tunneling route, mark the patient's skin at the neurostimulator pocket. (The neurostimulator should be placed on the opposite side of the body from another active implanted device and should be placed preferably on the right side of the body to allow for future placement of cardiac devices on the patient's left side.)

△ Cautions:

- When routing the lead, avoid sharp bends or kinks, which may break the wires. Broken wires may create an open circuit, resulting in loss of stimulation or component failure and requiring surgical replacement.
- When two leads are implanted, route the leads so the area between them is minimized (Figure 5). If the leads are routed in a loop and the patient is exposed to some sources of electromagnetic interference (eg, theft detectors), the patient may perceive a momentary increase in stimulation, which some patients have described as uncomfortable stimulation (jolting or shocking sensation).


Figure 5. Routing for two leads.

2. After assembling the tunneling tools, begin at the neurostimulator pocket and tunnel subcutaneously to the lead incision site.

△ **Caution:** Proceed slowly when the tunneling tool approaches the lead site. If excess force is used, the patient could experience additional trauma when resistance to tunneling suddenly ceases.



Notes:

- Deep tunneling is not desirable.
 - Avoid the lower thoracic ribs.
 - If the tunneling tool does not extend to the lead incision site, make an intermediate incision.
3. Unscrew the tunneling tip; then, remove the tunneling tool, leaving the passing straw in place in the tunnel.

4. Gently insert the proximal end of the lead through the passing straw to the neurostimulator pocket.
5. After carefully removing the passing straw, verify that the location of the lead has not changed.
6. Refer to the neurostimulator implant manual for instructions on connecting the lead to the neurostimulator.

A13106001.fm 11/23/05 2:26 pm
UC200xxxxxx EN
4 x 8 inches (101 mm x 203 mm)



Medtronic Confidential
LeadExtTemplate.fm

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All Rights Reserved
Medtronic B.V.
Medtronic International Ltd.
© Medtronic, Inc. 2005