Medtronic

PISCES QUAD®	3487A
PISCES QUAD COMPACT™	3887
PISCES QUAD PLUS®	3888

Lead kit

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Elektrodenset

Equipo de electrodos

Geleidingsdraadset

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Elektrodset

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

Medtronic Confidential
NeuroLdExt_R01










198855002 Rev A

Explanation of symbols on product or package labeling

Refer to the appropriate product for symbols that apply.

Medtronic®, PISCES Quad®, PISCES Quad Compact™, and PISCES Quad Plus® are trademarks of Medtronic, Inc.

Table of contents

Device description 5

Package contents 5

Device specifications 5

Instructions for use 7

Preparing for surgery 8

Placing a percutaneous lead 8

Testing stimulation intraoperatively 9

Anchoring the lead 10

Tunneling the percutaneous extension 10

Connecting the percutaneous extension to the lead 11

Refer to the indications sheet for indications and related information.

Refer to the appropriate information for prescribers booklet for contraindications, warnings, precautions, adverse events summary, individualization of treatment, patient selection, use in specific populations, resterilization, and component disposal.

Refer to System Eligibility, Battery Longevity, Specifications reference manual packaged with the software application card for neurostimulator selection, battery longevity calculations and specific neurostimulator specifications.

USA Refer to the clinical summary booklet packaged with the neurostimulator for information on the clinical study results of the neurostimulation system and individualization of treatment.

Device description

The Medtronic Models 3487A PISCES Quad, 3887 PISCES Quad Compact, or 3888 PISCES Quad Plus Leads are part of a neurostimulation system for pain therapy.

Package contents

- Lead
- Percutaneous extension
- Anchors (3)
- Connector boot
- Guide wire
- Needle, 15 gauge (epidural, curved or modified Tuohy)
- Screener cable
- Stylets:
 - bent at a 20° angle (inserted in Models 3487A and 3887; separate for Model 3888)
 - straight (inserted in Model 3888; separate for Models 3487A and 3887)
 - short
- Tunneling tools: tunneling rod, tunneling tips (1 narrow, 1 wide), tube (passing straw)
- Wrench, hex
- Product literature
- Warranty card (USA only)
- Registration form (USA only)

Device specifications

The lead has electrodes on the distal end; the proximal end fits into a four-conductor connector. A stylet has been inserted into the proximal end of the lead to aid in positioning the lead and performing intraoperative test stimulation.

Radiopaque marker (3888 only)


Figure 1. Electrode, lead contact, and stylet. Model 3888 shown.

Table 1. Device specifications^a for lead Models: 3487A, 3887, and 3888.

Description	Model 3487A	Model 3887	Model 3888
Connector	Quadripolar, in-line	Quadripolar, in-line	Quadripolar, in-line
Shape	Straight	Straight	Straight
Conductor resistance ^{b,c}	<125 Ω (28 cm) <250 Ω (56 cm)	<125 Ω (28 cm) <250 Ω (56 cm)	<110 Ω (28 cm) <220 Ω (56 cm)
Length	10 – 110 cm	10 – 110 cm	10 – 110 cm
Diameter	1.3 mm	1.3 mm	1.3 mm
Distal end			
Number of electrodes	4	4	4
Electrode shape	Cylindrical	Cylindrical	Cylindrical
Electrode length	3.0 mm	3.0 mm	6.0 mm
Electrode spacing	6.0 mm	4.0 mm	12.0 mm
Electrode distance	30.0 mm	24.0 mm	60.0 mm
Distal tip distance	Not applicable	Not applicable	7.9 mm
Radiopaque marker	No	No	Yes
Proximal end			
Lead contact length	2.3 mm	2.3 mm	2.3 mm
Lead contact distance	15.2 mm	15.2 mm	15.2 mm
Stylet handle length	40.0 mm	40.0 mm	40.0 mm

^a All measurements are approximate.

^b Electrical resistance of this device only.

^c Resistance is shown for 2 standard models. Resistance is proportional to length: long lengths have higher resistance that may limit the amplitude.

Table 2. Material of components in the Models 3487A, 3887, or 3888 package.

Component	Material	Material contacts human tissue
Lead		
Conductor wire	MP35N	No
Conductor wire insulation	Fluoropolymer	No
Electrodes	Platinum-iridium	Yes
Insulation	Polyurethane	Yes
Proximal connector	MP35N	Yes
Percutaneous extension		
Conductor wire	Stainless steel	No
Conductor wire insulation	Fluoropolymer	No
Handle	Polycarbonate	Yes
Insulation	Silicone rubber	Yes
Setscrew connector blocks	Stainless steel	Yes
Setscrew	Titanium	Yes
Anchor		
Twist-lock	Polymer	Yes
Two-wing, three-wing	Silicone rubber	Yes
Connector boot	Silicone rubber	Yes
Guide wire	Stainless steel	Yes
Needle	Stainless steel	Yes
Screener cable		
Connector	Polycarbonate	Yes
Cable	Rubber	Yes
Stylets (bent, straight, short)		
Handle	Acetal resin	Yes
Wire	Stainless steel with fluoropolymer coating	Yes
Tunneling tools		
Tunneling rod	Stainless steel	Yes
Tunneling tips	Stainless steel	Yes
Tube (passing straw)	Fluoropolymer	Yes
Wrench, hex		
Handle	Polypropylene	Yes
Shaft	Stainless steel	Yes

Instructions for use

Implanting physicians should be experienced in epidural-access procedures and should be thoroughly familiar with all product labeling.

Cautions:

- Do not bend, kink, or stretch the lead or extension body, which may damage the component.
- Do not bend or kink the stylet, which may damage the component.
- When reinserting a stylet, use care because too much pressure on the stylet could damage the lead, resulting in intermittent or loss of stimulation.
- Ensure that the lead contacts are aligned with the stylet handle. If the lead contacts are not aligned properly, intermittent or loss of stimulation may occur.
- Use only rubber-tipped forceps on the lead body. Do not use sharp-edged instruments (eg, hemostat), which may nick or cut the insulation.
- Do not use saline or other ionic fluids at connections, which could result in a short circuit.

Preparing for surgery

Before opening the lead package, verify the model number, use-by date, lead-length, and connector type.

Placing a percutaneous lead

Cautions:

- Use only the curved or modified Tuohy needle supplied in the kit. Using other needles may damage the lead.
- Use the recommended needle-insertion angle for the vertebral location listed below. An angle other than that recommended increases the risk of damage to the lead during insertion or manipulation.

Vertebral location	Recommended needle-insertion angle
Lumbar	Less than 45° to the skin
Thoracic	Less than 30° to the skin
Cervical	Greater than 45° to the skin

- Make an incision at the needle-entry site to the depth of the subcutaneous fascia.
- Using a paramedial approach and under fluoroscopy, insert the needle (ie, the curved or modified Tuohy included in the kit) into the epidural space at the appropriate angle until you encounter resistance from the ligamentum flavum.

Note: Midline placement may cause passive damage to the lead over time because of ligament or spinous process movement.

- Confirm needle location under fluoroscopy.
- After rotating the needle so that the beveled edge faces cephalad, remove the needle stylet.
- Advance the needle and confirm entry into the epidural space (eg, using the loss-of-resistance technique with air or sterile [United States Pharmacopeia—USP] water).

Caution: Do not use contrast media or a saline flush. Contrast media may obscure the field of view and a saline flush may increase the difficulty of lead placement.

- For a second lead, repeat steps 1 – 5 noting these recommendations:
 - Implant the second lead parallel to the first lead and approximately 1 – 3 mm lateral of the physiological midline.
 - Introduce the second lead one vertebral space below the first lead to help prevent nicking or cutting the first lead and to allow sufficient space for suturing both lead anchors.
 - Stagger the lead tips or place them several vertebral spaces apart, depending on the position that produces the most effective paresthesia.
- After inserting the guide wire through the needle, advance the guide wire no farther than 1 – 3 cm past the needle tip. Next, remove the guide wire from the needle.

Note: If the guide wire track deviates from the intended pathway, steering and manipulating the lead will be more difficult.

- Using fluoroscopy, slowly insert the lead through the needle and advance the lead to the initial target placement site. A stylet may need to be reinserted.

Notes:

- The bent stylet is keyed so that the flat notch of the stylet handle faces the opposite direction of the distal bent tip.

- When using a bent stylet and resistance is encountered during lead advancement, exchange the bent stylet for the straight stylet and use short, firm movements to advance the lead.
9. After verifying the lead position under fluoroscopy (anterior-posterior and lateral views), compare that location with the location that has the highest recommended probability of paresthesia coverage.

Notes:

- To increase lead stability, insert enough lead length to extend at least three vertebral bodies into the epidural space.
- Position the lead so that the center bipole electrodes (eg, electrodes 1 and 2 on a quadripolar lead) will be active. If lead migration occurs, effective stimulation may be regained through programming adjustments instead of surgical repositioning.

Testing stimulation intraoperatively


 **Caution:** To prevent an abrupt change in stimulation, which some patients have described as uncomfortable stimulation (jolting or shocking sensation):

- Program parameter changes in small increments above the perception threshold (the amplitude value at which the patient first perceives paresthesia).
- Decrease the amplitude to 0.0 V before:
 - changing electrode polarities.
 - connecting the screener cable to the screener.
 - turning ON the neurostimulator or screener.

1. After connecting the screener cable to the patient's implanted lead, pass the plug end out of the sterile field.

Note: When testing two leads, mark the Lead 1 and Lead 2 screener cables to ensure placement in the correct receptacle of the screener.

2. After verifying that the screener is turned OFF, connect the screener cable plug to the appropriate receptacle of the screener.

 **Warning:** Maintain adequate slack in the screener cable. If there is not enough slack and the cable is pulled, the percutaneous components may be dislodged.

3. Identify optimal stimulation parameters, beginning at a pulse width of 210 – 240 μ s and a rate of 30 Hz.

Note: Ensure that the patient can provide immediate feedback.

4. Increase the amplitude beginning from 0.0 V while asking the patient close-ended questions to identify the perception threshold (the amplitude at which the patient first perceives paresthesia), the discomfort threshold (the amplitude at which paresthesia is beyond the patient's tolerance), and the paresthesia coverage.

Note: If good paresthesia coverage is not attained, change electrode settings before repositioning the lead to confirm the direction of lead movement.

5. If two leads were placed, repeat steps 1–4 for the second lead; optimize paresthesia coverage using both leads.
6. In the patient's chart, document the lead position that provided appropriate stimulation coverage (ie, record the settings and patient responses and include a fluoroscopic image of the final lead position).
7. Disconnect the screener cable from the lead. Leave the screener cable in the sterile field for additional parameter testing before closing.

 **Caution:** Do not pull directly on the cable to disconnect the connector because this may lead to wire breakage and inadequate or discontinuous stimulation.

Anchoring the lead

△ **Caution:** Do not tie ligatures around the lead or extension body, which may damage the insulation.

1. Prepare the anchor site by making a 5 – 7 cm (2.0 – 2.8 in) longitudinal incision around the needle shaft, dissecting down to the supraspinous ligament, and establishing hemostasis.
2. Disconnect the stylet from the lead. Partially expose, but do not completely withdraw the stylet.
3. While maintaining lead position, use minimal force to remove the needle and stylet.

△ **Caution:** Use minimal traction to remove the needle because quick or sudden removal may dislodge the lead.

4. Slide the anchor onto the proximal end of the lead and continue sliding the anchor down as close as possible to where the lead emerges from the vertebral column. Use care to maintain the lead position.

Notes:

- The anchor can be lubricated with sterile water.
 - For a twist-lock anchor, ensure that the anchor is unlocked before sliding the anchor onto the lead. (An unlocked anchor is twisted so the two parts are at a 90° angle to each other.)
5. Use 2-0 nonabsorbable suture (such as silk or some types of braided polyester mesh) to secure the anchor.

△ Cautions:

- Do not use polypropylene suture material on silicone components. Polypropylene may damage the component, resulting in component damage or failure.
- Do not overtighten ligatures on the anchor or connector boot. Ligatures that are too tight may damage the component.

For the twist-lock anchor: Suture the anchor to the supraspinous ligament or deep fascia. Next, lock the anchor to secure the lead.

Note: If minor repositioning is necessary, unlock the lead, then relock the anchor.

For the winged anchor: Tie ligatures around the grooves in the middle of the anchor to secure the anchor to the lead. Next, suture the anchor to the supraspinous ligament or deep fascia.

Note: Laboratory testing has shown that injecting silicone medical adhesive between the anchor and lead after tying the ligatures increases the anchor-to-lead holding force.

6. Verify test stimulation parameters to ensure that the lead has not moved: connect the short stylet handle to the lead, connect the screener cable to the lead, then, connect the screener cable plug to the screener. If the lead has moved, reposition it.

Tunneling the percutaneous extension

1. After simulating the tunneling route, mark the patient's skin at the lead-extension connection site and at the percutaneous extension exit site (exit site). The exit site should be on the side opposite that intended for the neurostimulator and at least 10 cm (4 in) lateral. (The neurostimulator should be placed on the opposite side of the body from another active implanted device and should be placed preferably on the right side of the body to allow for future placement of cardiac devices on the patient's left side.)

Cautions:

- When routing the extension, avoid sharp bends or kinks, which may break the wires. Broken wires may create an open circuit, resulting in loss of stimulation or component failure and requiring surgical replacement.
- When two leads are implanted, route the lead-extensions so the area between them is minimized (Figure 2). If the lead-extensions are routed in a loop and the patient is exposed to some sources of electromagnetic interference (eg, theft detectors), the patient may perceive a momentary increase in stimulation, which some patients have described as uncomfortable stimulation (jolting or shocking sensation).


Figure 2. Routing for two lead-extensions.

2. After assembling the tunneling tools, make an incision at the exit site.
3. Begin at the exit site and tunnel subcutaneously to the lead incision.

Notes:

- Deep tunneling is not desirable.
 - Avoid the lower thoracic ribs.
 - If the tunneling tool does not extend to the lead-extension connection site, make an intermediate incision.
4. Use blunt dissection to form a tunnel from the lead incision to the lead-extension connection site.
 5. After tying a ligature to the handle of the percutaneous extension, gently pull the extension through the passing straw to the exit site.

Connecting the percutaneous extension to the lead

Caution: Before connecting components, wipe off any body fluids and dry all connections. Fluids in the connections may result in stimulation at the connection site, intermittent stimulation, or loss of stimulation.

1. After disconnecting the short stylet from the lead, slip the connector boot onto the lead body with the open end of the connector boot facing outward, exposing the lead contacts.

Note: Irrigate the connector boot with a nonionic antibiotic solution. Dry the interior of the boot.

2. Wipe the lead body and extension setscrew connector junction with sterile gauze. If necessary, use sterile (USP) water or a nonionic antibiotic solution. Dry all connections.
3. Insert the lead connector fully into the extension setscrew connector junction.

Note: Each lead contact must be aligned under each setscrew.

4. Tighten the setscrews to complete the electrical circuit with the lead contacts.

△ **Caution:** Do not overtighten the extension setscrews. Overtightening the extension setscrews may damage the lead contacts and cause an open or short circuit, resulting in intermittent or loss of stimulation.

5. Confirm the integrity of the connection by rechecking the stimulation pattern.
6. After sliding the connector boot over the lead-extension connection, tie a ligature around the wide end of the connection.

Note: Do not suture the narrow end of the connector boot.

7. Carefully pull the lead-extension connection through the tunnel to the lead-extension connection site.

- △ **Cautions:**
- Do not pull the lead or extension taut. Allow enough slack in the lead and extension to accommodate patient movement. Pulling the lead or extension taut may result in a short or open circuit or migration of implanted components.
 - Coil the excess lead body into a circular loop greater than 2 cm (0.8 in) in diameter. Do not use a U-shaped loop or bend (Figure 3). A circular loop decreases the possibility of electromagnetic interference and its effects and prevents kinking or damaging the lead body.


Figure 3. Coiling excess lead.

8. After carefully removing the passing straw from the tunnel, verify that the location of the lead and the pattern of stimulation have not changed.
9. Close the lead incision and the percutaneous extension exit site, leaving the percutaneous extension wires looped under a dressing with the plastic handle exiting the bandage.
10. Connect the screener cable to the percutaneous extension so that the screener cable is available for trial evaluation.

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All Rights Reserved
Medtronic B.V.
Medtronic International Ltd.
© Medtronic, Inc. 2004