NANOSTIM™ Synthetic Bone Paste

Classification

NANOSTIM™ is a medical device class 3 and not a drug.

Material and composition

NANOSTIM™ contains needle shaped phasepure hydroxyapatite nanosize crystals as a suspension in water.

NANOSTIM™ is a fully synthetic, highly compatible bone substitute material. There are no substances of biological origin used in the synthesis or processing of the product. No additional preservatives or medical substances are present.

NANOSTIM™ is a highly viscous, fully degradable bone substitute.

NANOSTIM™ is used for the filling of bone defects. The chemical composition and crystalline structure are essentially identical with the mineral phase of natural bone.

NANOSTIM™ is packaged in cartridges ready to use for direct application in the operating room.

NANOSTIM™ does not harden in situ.

Responsible manufacturer

aap Biomaterials GmbH & Co. KG
Lagerstr. 11 - 15
64807 Dieburg Germany

Responsible distributor

Medtronic Sofamor Danek USA, Inc.
1800 Pyramid Place
Memphis TN 38132
Tel. 1.901.396.3133

Pharmaceutical form

NANOSTIM™ is a ready-to-use sterile paste available in a range of pack sizes.

Indications

NANOSTIM™ is used for the filling and reconstruction of bone defects in orthopaedics, traumatology and neurosurgery.

It is defined for:

- Filling of bone defects due to fractures
- Filling of bone defects due to resection of benign tumors
- Filling of bone cysts
- Filling of donor sites after removal of autografts
- Filling of metaphyseal bone defects due to osteoporosis

NANOSTIM™ can be used to stop bleeding in cancellous bone.

NANOSTIM™ is suitable as a bone graft extender of autograft as well as allograft.

A well vascularized implantation site and a correct reposition of bone in case of fractures in combination with a stable fixation is a requirement for successful treatment with NANOSTIM™ and is an important prerequisite for bone healing.

Contraindications

NANOSTIM™ should not be applied or administered in the following conditions:

- Acute or chronic infections at the implantation site
- Defects due to malignant tumors
- Defects in the area of open epiphyseal plates

Particular caution is indicated in patients with disorders of bone metabolism. If the application is considered, prior to administration a benefit/risk analysis should be performed.

Use in pregnancy or nursing mothers

No information is available in pregnant or nursing patients. Therefore prior to application of NANOSTIM™ an individual benefit and risks analysis should be performed.

Use in children or elderly patients

To date, no special precautions are known in relation to the age of patients which would make particular precautionary measures necessary.

Side effects

There are no side effects known in relation to NANOSTIM™.

Interactions

No interactions of NANOSTIM™ with other medical devices or drugs have been found in case of appropriate application.

Precautionary measures and warnings

In unstable bone defects it is absolutely imperative to use NANOSTIM™ only in combination with adequate fixation or stabilisation measures.

In order to prevent the occurrence of emboli the highly pressurised application of NANOSTIM™ into a tightly confined space with venous or arterial access, should be avoided. In patients with increased wound secretion, e. g. in patients with impaired blood coagulation a drainage system should be used.

In rare cases an increased wound secretion with a delayed healing was

observed. Therefore the application of NANOSTIM™ with halliux valgus and metatarsal arthroses is not recommended.

NANOSTIM™ is not to be re-sterilized.

Dosage and application

NANOSTIM™ is intended to remain locally in bony tissues.

The quantity of NANOSTIM™ required to fill a bone defect will depend upon the size of the defect.

Defects should be completely filled and contoured to the natural surfaces of the bone without pressing the material into the defect. Excess material is to be removed. NANOSTIM™ is supplied as a sterile paste in an application cartridge ready for use.

The content of one cartridge is for single use only.

Immediately before application of NANOSTIM™ the locking cap is removed from the cartridge and the paste is injected by continuous pressure directly into the defect.

NANOSTIM™ should only be implanted into vital bone in order to ensure rapid revascularization and new bone formation within the implant.

Special application devices for NANOSTIM™ are available in order to facilitate the implantation (see accessories).

Sterilisation procedure

NANOSTIM™ is supplied sterile in a peel off package. Sterilisation is performed by gamma radiation.

Storage

NANOSTIM™ should be protected against freezing. NANOSTIM™ should be stored at temperatures between 4 and 25 °C. Direct contact with a heating system or storage in direct sunlight should be avoided. The content of the peel off package is sterile, therefore the product must only be used if the package is undamaged.

Expiration

The expiration date is printed on the sterile peel off package and on the cardboard box. The product must not be used after this date.

Accessories for application

Screw on applicator for 5 cc. cartridges.

Application needle for 1 and 2 cc. cartridges.

Description of accessories

The application needle enables the percutaneous application of NANOSTIM™ e.g. in radius fractures, the application through screw wholes of fixation plates and provides an easy access to deep localisations of bone defects.

The flexible screw on plastic applicator for the 5 cc. cartridges is useful for an easy access into deep defect locations.



These products are only intended for the use in combination with NANOSTIM™ cartridges. They are sterilized by Ethylene oxide. The products must only be used if the package is undamaged.

Usage of the accessories

The application needle is to be firmly fixed into the conus of the 1 or 2 cc. cartridges. The screw on applicator is firmly screwed on to the 5 cc. cartridge. Prior to application the firm anchorage of the accessories should be examined.

After single use the application devices should be disposed. Prior to disposal of the needle the locking cap should be put on again.

Used symbols printed in the packages

	Sterile
	Sterilisation by Gamma-radiation
	Sterilisation by Ethylene oxide
	Order number/article number
	Pay attention to the directions for use
	Do not re-use
	Temperature limits 4 to 25 °C
	Batch number
	Expiration date
	Manufacturer
	Sign of conformity
	registered trade mark

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