DRUG ELUTING BALLOON CATHETER

AcoArt Orchid

Drug Eluting Balloon 0.035"

AcoArt Tulip

Drug Eluting Balloon 0.018"

AcoArt Litos

Drug Eluting Balloon 0.014"

EN	Drug Eluting Balloon Catheter Instructions for use	4	NO	Medikamenteluerende ballongkateter Bruksanvisning	28
DE	Medikamenten-freisetzender Ballonkatheter Gebrauchsanweisung	6	SV	Läkemedelsavgivande ballongkateter Bruksanvisning	30
FR	Cathéter à ballonnet à élution médicamenteuse Mode d'emploi	9	FI	Lääkeainetta vapauttava pallonsaajennuskatetri Käyttöohje	33
IT	Catetere a palloncino a rilascio di farmaco Istruzioni d'uso	11	CS	Balónkový katér uvolňující léčivo Návod k použití	35
ES	Catéter de balón liberador de fármaco Instrucciones de uso	14	HU	Gyógyszerelvező ballongkatéter Használati utasítások	38
NL	Geneesmiddel-eluerende ballonkatheter Gebruiksaanwijzing	17	TR	İlaç Kaplı Balon Kateter Kullanım Kılavuzu	40
DA	Lægemiddelreluerende ballonkateter Brugsanvisning	19	KO	약물 방출 풍선 카테터 사용설명서	43
PL	Cewnik balonowy uwalniający lek Sposób użycia	22	TW	“先瑞達”藥物塗層氣球擴張導管 使用說明書	45
PT	Cateter de balão com eluição de medicamento Instrução de utilização	25			

Medical Device // Medizínprodukt // Dispositif médical // Dispositivo medico // Producto sanitario // Medisch hulpmiddel //
Medicinsk udstyr // Wyrób medyczny // Dispositivo médico // Medisinsk udstyr // Medicinteknisk produkt // Lääkinnällinen laite //
Zdravotnický prostředek // Orvostechikai eszköz // Tibbi cihaz // 의료기기 // 醫療器械

Contains a medicinal substance // Enthält ein Arzneimittel // Contient une substance médicinale // Contiene una sostanza medicinale // Contiene una sustancia medicinal // Bevat een geneesmiddel // Indeholder et lægemiddel // Zawiera substancję leczniczą // Contém uma substância medicinal // Innehåller ett medicinskt stoff // Innehåller ett läkemedel // Sisältää lääkeainetta // Obsahuje léčivou látku // Gyógyszert tartalmaz // Tibbi bir maddə içərir // 약용 물질 함유 // 包含藥物

Translation/Übersetzung/Traduction/Traduzione/Traducción/Vertaling/Oversættelse/Thumaczenie/Tradução // Oversættelse // Översättning // Käänös // Preklad // Fordítás // Tercüme // 번역 // 翻譯

Balloon Nominal Diameter // Ballon Nominaldurchmesser // Diamètre Nominal du Ballon // Diametro nominale del palloncino // Globo Diámetro Nominal // Nominale diameter ballon // Ballon nominel diameter // Nominálna šírka balónika // Diamêtro nominal do balão // Ballong nominell diameter // Ballong nominell diameter // Pallon läpimitta // Nominální průměr balónku // Ballon névleges átmérője // Balon Nominal Çapı // 풍선 명목 직경 // 球囊標稱直徑

Balloon Nominal Length // Ballon Nominallänge // Longueur nominale du ballon // Lunghezza nominale del palloncino // Globo Largo Nominal // Nominale lengte ballon // Ballon Nominel længde // Nominálna dĺžka balónika // Comprimento nominal do balão // Ballong nominell lengde // Ballong nominell diameter // Pallon pituus // Nominální délka balónku // Ballon névleges hossza // Balon Nominal Uzunluğu // 풍선 명목 길이 // 球囊標稱長度

Usable Catheter Length // Nutzbare Katheterlänge // Longueur utilisable du cathéter // Lunghezza utile del catetere // Utilizable Longitud Catéter // Bruikbare katheterlengte // Brugbar kateter længde // Długość użytkowa cewnika // Comprimento Efectivo do Cateter // Den effektive lengden av kateteret // Användbar kateterlängd // Katerin käyttöpituus // Použitelná délka katetru // Felhasználható katéter hossza // Kullanabilir Kateter Uzunluğu // 카테터의 유효길이 // 導管有效長度

Recommended Introducer // Empfohlenes Einführbesteck // Tuyau de la gaine recommandé // Canale della guaina consigliato // Recomendada vaina // Aanbevolen inbrengorgaan // Anbefalede Introducer // Zalecan prowadnik // Bainha Introduitora Recomendada // Anbefalt kateter skjede // Rekommenderad Introduktör // Suositeltu sisäänviejäholkki // Doporučený zavaděč // Ajánlott bevezető // Tavsye edilen introdüser // 레코넨트 삽관기 // 推薦套管

Maximum recommended guide wire // Empfohlener Führungsdraht // Fil de guidage maximum recommandée // Filo guida massimo raccomandato // Alambre de guía recomendada máxima // Aanbevolen maximale geleidedraad // Maksimal anbefalet guidewren // Zalecana maksymalna prowadnica // Fio de Guia Máximo Recomendado // Maksimal anbefalt guidekabel // Rekommenderad maximal styrtråd // Suurin suositeltu ohjainvaijeri // Maximální doporučený vodič drát // Maximális ajánlott vezetődórt // Tavsye edilen maksimum kilavuz tel // 일정한 가이드 와이어 // 匹配導線

Nominal Pressure // Nominaldruck // Pression nominale // Pressione nominale // Presión nominal // Nominale druk // Nominelt tryk // Nominálne ciśnienie // Pressão normal // Nominelt trykk // Nominellt tryck // Nimellispaine // Nominální tlak // Névleges nyomás // Nominal Basınc // 명목 압력 // 公稱壓力

Rated Burst Pressure // Maximal zulässiger Arbeitsdruck // Pression de Rupture Nominale // Pressione di rottura nominale // Ráfaga de presión nominal // Nominale barstdruk // Nominel Sprængningstryk // Nominálne cisnienie wiazki // Pressão de Ruptura Nominal // Rangert sprengetrykk // Märksprängtryck // Nimellinen puhkeamisaine // Jmenovitý tlak protřzeni // Névleges felrepedési nyomás // Nominal Patlama Basıncı // 액정 폭파 압력 // 額定爆破壓

Date of manufacture // Herstellungsdatum // Date de Fabrication // Data di produzione // Fecha de Fabricación // Fabricage datum // Produktionsdato // Data produkcyj // Data de Fabricação // Produktionsdato // Tilverkningsdatum // Valmistuspäivä // Datum výroby // Gyártási dátum // Üretim Tarihi // 제조 일자 // 生產日期

Manufacturer // Hersteller // fabricant // Produttore // Fabricante // Fabrikant // Producent // Producent // Fabricante // Produsenten // Tillverkare // Valmistaja // Vyrobce // Gyártó // Üretici // 제조자 // 生產者

CE Marking of Conformity. Notified Body number// CE-Kennzeichnung. Konformitäts-Kennnummer // Numéro d'identification d'un organisme notifié de Marquage CE // Marchio CE di conformità Numero di certificazione // El mercado de conformidad CE Certificación No // CE-markering van conformiteit. Aangemelde instantie nummer // CE-merkning. Nummer fra det notificerede organ // Oznakowanie zgodności CE. Numer jednostki notyfikowanej // Qualificação confirmada.Código de Qualidade // CE-merking av samsvars Sertifisering nummer // CE Märkning om överensstämelse anmält organ nummer // CE-merkintä, tunnustenumero // Označení shody CE Číslo oznámeného subjektu // EK-megfelelősgéti jelölés. Tanúsított testület száma // CE Uygunluk İşareti Onaylanmış Kuruluş numarası // CE 합격 표시 인증번호 // CE 合格標誌 認證號

Sterile, sterilized using ethylene Oxide // Steril, Ethylenoxid-sterilisiert // Stérile, stérilisé avec de l'oxyde d'éthylène // Sterile, sterilizzato con ossido di etilene // Esterilizado con óxido de etileno // Steriel, gesteriliseerd met ethylenoxide // Steril, steriliseres med ethylenoxid // Sterilny, wysterylizowany za pomocą tlenku etylenu // Esteril, esterilizados con óxido de etileno // Etylenoksid sterilisering // Steril, steriliserades med etylenoxid // Steriili, steriloitu etyleenioksidilla // Sterilni, sterilizováno etylenoxidem // Steril, etilenoxydall sterilizált// Steril, Etilen oksit ile sterilize edilmiştir // 에틸렌옥사이드 멸균 // 經環氧乙烷滅菌

Do not re-use // Nicht wiederverwenden // Ne pas réutiliser // Non riutilizzare // No reutilizar // Niet hergebruiken // Du må ikke genbruge // Nie wykorzystywać ponownie // Não re-uso // Ikke gjenbruk // inte återanvändas // Ei saa käyttää uudelleen // Nepouzivat opakovaně // Ne használt újra // Tekrar kullanmayınız // 반복 사용 불가 // 不得二次使用

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ACOART ORCHID DRUG ELUTING BALLOON CATHETER 0.035"							
Usable Catheter Length: 80 cm							
ORC03020AG	ORC35020AG	ORC04020AG	ORC45020AG	ORC05020AG	ORC55020AG	ORC06020AG	ORC07020AG
ORC03030AG	ORC35030AG	ORC04030AG	ORC45030AG	ORC05030AG	ORC55030AG	ORC06030AG	ORC07030AG
ORC03040AG	ORC35040AG	ORC04040AG	ORC45040AG	ORC05040AG	ORC55040AG	ORC06040AG	ORC07040AG
ORC03060AG	ORC35060AG	ORC04060AG	ORC45060AG	ORC05060AG	ORC55060AG	ORC06060AG	ORC07060AG
ORC03080AG	ORC35080AG	ORC04080AG	ORC45080AG	ORC05080AG	ORC55080AG	ORC06080AG	ORC07080AG
ORC03100AG	ORC35100AG	ORC04100AG	ORC45100AG	ORC05100AG	ORC55100AG	ORC06100AG	ORC07100AG
ORC03120AG	ORC35120AG	ORC04120AG	ORC45120AG	ORC05120AG	ORC55120AG	ORC06120AG	ORC07120AG
ORC03150AG	ORC35150AG	ORC04150AG	ORC45150AG	ORC05150AG	ORC55150AG	ORC06150AG	ORC07150AG
ORC03200AG	ORC35200AG	ORC04200AG	ORC45200AG	ORC05200AG	ORC55200AG	ORC06200AG	ORC07200AG
ORC03250AG	ORC35250AG	ORC04250AG	ORC45250AG	ORC05250AG	ORC55250AG	ORC06250AG	ORC07250AG
ORC03300AG	ORC35300AG	ORC04300AG	ORC45300AG	ORC05300AG	ORC55300AG	ORC06300AG	ORC07300AG

ORC08020AG	ORC09020AG	ORC10020AG	ORC12020AG
ORC08030AG	ORC09030AG	ORC10030AG	ORC12030AG
ORC08040AG	ORC09040AG	ORC10040AG	ORC12040AG
ORC08060AG	ORC09060AG	ORC10060AG	ORC12060AG
ORC08080AG	ORC09080AG	ORC10080AG	ORC12080AG
ORC08100AG	ORC09100AG	ORC10100AG	ORC12100AG
ORC08120AG	ORC09120AG	ORC10120AG	ORC12120AG

ACOART ORCHID DRUG ELUTING BALLOON CATHETER 0.035"							
Usable Catheter Length: 130 cm							
ORC03020BG	ORC35020BG	ORC04020BG	ORC45020BG	ORC05020BG	ORC55020BG	ORC06020BG	ORC07020BG
ORC03030BG	ORC35030BG	ORC04030BG	ORC45030BG	ORC05030BG	ORC55030BG	ORC06030BG	ORC07030BG
ORC03040BG	ORC35040BG	ORC04040BG	ORC45040BG	ORC05040BG	ORC55040BG	ORC06040BG	ORC07040BG
ORC03060BG	ORC35060BG	ORC04060BG	ORC45060BG	ORC05060BG	ORC55060BG	ORC06060BG	ORC07060BG
ORC03080BG	ORC35080BG	ORC04080BG	ORC45080BG	ORC05080BG	ORC55080BG	ORC06080BG	ORC07080BG
ORC03100BG	ORC35100BG	ORC04100BG	ORC45100BG	ORC05100BG	ORC55100BG	ORC06100BG	ORC07100BG
ORC03120BG	ORC35120BG	ORC04120BG	ORC45120BG	ORC05120BG	ORC55120BG	ORC06120BG	ORC07120BG
ORC03150BG	ORC35150BG	ORC04150BG	ORC45150BG	ORC05150BG	ORC55150BG	ORC06150BG	ORC07150BG
ORC03200BG	ORC35200BG	ORC04200BG	ORC45200BG	ORC05200BG	ORC55200BG	ORC06200BG	ORC07200BG
ORC03250BG	ORC35250BG	ORC04250BG	ORC45250BG	ORC05250BG	ORC55250BG	ORC06250BG	ORC07250BG
ORC03300BG	ORC35300BG	ORC04300BG	ORC45300BG	ORC05300BG	ORC55300BG	ORC06300BG	ORC07300BG

ORC08020BG	ORC09020BG	ORC10020BG	ORC12020BG
ORC08030BG	ORC09030BG	ORC10030BG	ORC12030BG
ORC08040BG	ORC09040BG	ORC10040BG	ORC12040BG
ORC08060BG	ORC09060BG	ORC10060BG	ORC12060BG
ORC08080BG	ORC09080BG	ORC10080BG	ORC12080BG
ORC08100BG	ORC09100BG	ORC10100BG	ORC12100BG
ORC08120BG	ORC09120BG	ORC10120BG	ORC12120BG

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Acotec Scientific Co., Ltd
4th and 5th Floor, Building
No. 8, 16 Hongdabei Road,
BDA, 100176 Beijing, China
Tel : +86 10 67872107
Fax: +86 10 67866678

European Authorized
Representative:
MedDevConsult GmbH
Airport Center (Building C)
Flughafenstrasse 52a
22335 Hamburg
Germany

Publication Date : 2020 08 03