**INTENDED PURPOSE**

The Penumbra Access Catheter System is designed to provide a conduit for introduction of interventional devices to the peripheral and neurovascular anatomy.

**DEVICE DESCRIPTION**

The Penumbra Access Catheter System is comprised of multiple components. The BENCHMARK BMX81 Access System, Neuron MAX Access System, and BENCHMARK BMX96 Access System include a Delivery Catheter, Select™ Catheter, and Dilator. The Neuron Intracranial Access System and the Benchmark Intracranial Access System include a Delivery Catheter and a Select™ Catheter. The Delivery Catheters can be used individually with a 0.038 in [0.97 mm] guidewire or together with the Select Catheter to access the desired anatomy. The Delivery Catheter may be provided with a hemostasis valve adapter, rotating hemostasis valve, and peelable introducer.

Delivery Catheter:

The Delivery Catheter is a single lumen, variable stiffness catheter with a radiopaque markerband or radiopaque section on the distal end and a Luer hub on the proximal end. The Delivery Catheter dimensions are included on the individual device label. The Delivery Catheter is compatible with introducer sheaths appropriately sized for the outer diameter of Delivery Catheter. The Delivery Catheter has a hydrophilic coating on the distal segment of the catheter shaft. The hydrophilic coating length per catheter is noted in Table 1.

Select Catheter:

The Select Catheter is a single lumen, braid-reinforced, variable stiffness catheter with a radiopaque distal end and a Luer hub on the proximal end. The Select Catheter is available in multiple tip shapes (see Table 1). The Select Catheter tip shape and dimensions are included on the individual device label. The Select Catheter is compatible with the Delivery Catheters as noted in Table 2.

Dilator (Accessory for Sheath Configuration only):

The Dilator is a single lumen, radiopaque catheter with a tapered distal end and a Luer hub on the proximal end. The Dilator provided is compatible with the Delivery Catheter. The Dilator facilitates the percutaneous entry of the Delivery Catheter by forming an atraumatic transition from the skin through the subcutaneous tissue to the vessel.

Table 1 – Tip Shapes and Hydrophilic Coating		
Catheter	Tip Shapes	Hydrophilic Coating Length (cm)
<b>Delivery Catheter</b>		
Neuron 053	Straight, Multipurpose (MP)	11
Neuron 070		
Benchmark		9
BMX81		18
Neuron MAX		7
BMX96		9
<b>Select Catheter</b>		
5F	SIM, H1, or BER	N/A
6F	SIM, H1, BER, or SIM-V	N/A

**INDICATIONS FOR USE**

The Penumbra Access Catheter System is indicated for the introduction of interventional devices into the peripheral and neuro vasculature.

**INTENDED PATIENT POPULATION**

The intended patient population is patients who require neuro and/or peripheral intervention.

**CONTRAINDICATIONS**

There are no known contraindications.

**WARNINGS**

- The Penumbra Access Catheter System should only be used by physicians who have received appropriate training in interventional techniques.
- The risks and benefits for radial access against a transfemoral approach should be carefully weighed and considered for each patient.
- The device is intended for single use only. Do not resterilize or reuse. Resterilization and/or reuse may compromise the structural integrity of the device and/or lead to device failure which in turn, may result in patient injury, illness, or death.
- Do not advance or withdraw the Penumbra Access Catheter System against resistance without careful assessment of the cause using fluoroscopy. If the cause cannot be determined, withdraw the device. Unrestrained moving or torquing the device against resistance may result in damage to the vessel or device.
- Do not use automated high-pressure contrast injection equipment with the Penumbra Access Catheter System because it may damage the device.
- Failure to abide by the warnings in this labeling might result in damage to the device coating, which may necessitate intervention or result in serious adverse events.

**PRECAUTIONS**

- Do not use kinked or damaged devices. Do not use open or damaged packages. Return all damaged devices and packaging to the manufacturer / distributor.
- Use prior to the “Use By” date.
- Use the Penumbra Access Catheter System in conjunction with fluoroscopic visualization.
- Maintain a constant infusion of an appropriate flush solution.
- If flow through the device becomes restricted, do not attempt to clear the lumen by infusion. Remove and replace the device.
- Prior to beginning radial artery access, conduct screening, such as an Allen test, to ensure that radial access is appropriate for the patient.
- As in all fluoroscopy procedures, consider all necessary precautions to limit patient radiation exposure by using sufficient shielding, reducing fluoroscopy times and modifying radiation technical factors whenever possible.
- Avoid using alcohol, antiseptic solutions, or other solvents to pre-treat the device because this may cause unpredictable changes in the coating which could affect the device safety and performance.

**POTENTIAL ADVERSE EVENTS\***

Possible complications include, but are not limited to, the following:

- Acute vessel occlusion
- Access site complications such as hematoma, hemorrhage, spasm, occlusion, or compartment syndrome
- Air embolism
- Allergic reaction and anaphylaxis from contrast media or device material
- Arteriovenous fistula
- Death
- Emboli including foreign body embolization
- Inflammation
- Infection
- Intracranial hemorrhage
- Ischemia
- Nerve injury
- Neurological deficits including stroke
- Pseudoaneurysm
- Renal impairment or acute renal failure from contrast media
- Radiation exposure that may lead to cataracts, skin reddening, burns, alopecia, or neoplasia from x-ray exposure
- Vessel dissection
- Vessel perforation
- Vessel spasm
- Vessel thrombosis

\*If a serious incident related to the device occurs, contact your Penumbra, Inc. representative and the competent authority in your respective country/region.

Table 2 – Device Compatibility and Selection Information			
Delivery Catheter	Introducer Sheath	Guidewire	Select Catheter
Neuron 053	6F	0.038" (0.97 mm)	5F
Neuron 070	6F		5F
Benchmark	6F		5F
BMX81	7F		5F
Neuron MAX	8F		6F
BMX96	8F		6F

**DEVICE PREPARATION AND USE (GUIDE CATHETER CONFIGURATION)**

1. Refer to **Warnings, Precautions, and Potential Adverse Events** prior to use.
2. Refer to **Table 2** for device compatibility information.
3. Select an appropriately sized Delivery Catheter based on the anatomy and length.
4. If using, select appropriate Select Catheter shape based on the target vessel and its surrounding anatomy.
5. Gently remove the Delivery Catheter and packaging card from the pouch by grasping the Delivery Catheter hub and packaging card and slowly pulling them out of the pouch.
6. Remove the Delivery Catheter from the packaging card by removing the hub from the card tabs before gently removing the Delivery Catheter shaft.
7. Remove the Rotating Hemostasis Valve (RHV) from the card tabs, if provided.
8. Inspect the Delivery Catheter for kinks or other damage.
9. Connect the RHV to the hub of the Delivery Catheter. Flush the lumen and hydrate the outer surface of the Delivery Catheter with heparinized saline.
10. If using, insert the Select Catheter into the Delivery Catheter, and advance the Select Catheter until the distal tip of the Select Catheter is at the distal tip of the Delivery Catheter.
11. If a 0.038 in [0.97 mm] guidewire is used instead of the Select Catheter, insert the 0.038 in [0.97 mm] guidewire into the Delivery Catheter, and advance the guidewire until the distal tip of the wire is at the distal tip of the Delivery Catheter.
12. Place an introducer sheath appropriately sized for the Delivery Catheter's outer diameter.
13. Advance the Delivery Catheter and Select Catheter or 0.038 in [0.97 mm] guidewire simultaneously into the introducer sheath and extend the Select Catheter or 0.038 in [0.97 mm] guidewire 5 cm to 6 cm distal to the tip of the Delivery Catheter.
14. If desired, a 0.038 in [0.97 mm] guidewire may be introduced through the Select Catheter.
15. Advance the Delivery Catheter and Select Catheter or guidewire system, and then select the appropriate vessel with the Select Catheter or 0.038 in [0.97 mm] guidewire.
16. After gaining access to the desired vessel with the Select Catheter or 0.038 in [0.97 mm] guidewire, advance the Delivery Catheter over the Select Catheter or 0.038 in [0.97 mm] guidewire into the desired vessel.
17. Under fluoroscopy, advance the Delivery Catheter until the desired position is achieved.
18. If used, remove the Select Catheter, and replace with a 0.038 in [0.97 mm] or smaller guidewire.
19. Advance the Delivery Catheter and guidewire to the vascular site and remove the guidewire.

**DEVICE PREPARATION AND USE (SHEATH CONFIGURATION - BENCHMARK BMX81 ACCESS SYSTEM, NEURON MAX ACCESS SYSTEM, AND BENCHMARK BMX96 ACCESS SYSTEM ONLY)**

1. Refer to **Warnings, Precautions, and Potential Adverse Events** prior to use.
2. Refer to **Table 2** for device compatibility information.
3. Select an appropriately sized Delivery Catheter based on the anatomy and length.
4. Gently remove the Delivery Catheter and packaging card from the pouch by grasping the Delivery Catheter hub and packaging card and slowly pulling them out of the pouch.
5. Remove the Delivery Catheter from the packaging card by removing the hub from the card tabs before gently removing the Delivery Catheter shaft.
6. Gently remove the Dilator from the packaging card by removing the hub from the card tabs before gently removing the Dilator shaft.
7. Inspect the Delivery Catheter for kinks or other damage.
8. Flush the lumen and hydrate the outer surface of the Delivery Catheter with heparinized saline.
9. Insert the Dilator through the hemostasis valve of the Delivery Catheter, snapping it into place at the hub.
10. Flush the lumen of the Dilator with heparinized saline.
11. Introduce the cannula of an angiographic needle into the vessel using aseptic technique. Holding the needle in place insert the flexible end of a 0.038 in [0.97 mm] guidewire through the needle and into the vessel. Gently advance the guidewire to the desired depth.
12. Holding the guidewire in place, withdraw the needle and apply pressure to the puncture site until the Delivery Catheter is inserted into the vasculature.
13. Thread the Delivery Catheter/Dilator assembly over the guidewire, grasping the Delivery Catheter close to the skin to prevent buckling. Using a rotating motion, advance the assembly through the tissue into the vessel.
14. Under fluoroscopy, advance the Delivery Catheter until the desired position is achieved.
15. Detach the Dilator from the Delivery Catheter by releasing the snap-fit ring at the hub. Withdraw the guidewire and Dilator. To avoid damage to the Delivery Catheter, do not withdraw the Dilator until the Delivery Catheter is in the vessel.
16. Aspirate from the Hemostasis Valve Adapter (HVA) to remove any potential air. After aspiration and flushing, consider establishing a heparinized solution or suitable isotonic solution via the HVA.
17. If using, select appropriate Select Catheter shape based on the target vessel and its surrounding anatomy.
18. If using, insert the Select Catheter into the Delivery Catheter, and advance the Select Catheter until the distal tip of the Select Catheter is at the distal tip of the Delivery Catheter.
19. If a 0.038 in [0.97 mm] guidewire is used instead of the Select Catheter, insert the 0.038 in [0.97 mm] guidewire into the Delivery Catheter, and advance the guidewire until the distal tip of the Delivery Catheter.
20. Extend the Select Catheter or 0.038 in [0.97 mm] guidewire 5 cm to 6 cm distal to the tip of the Delivery Catheter.
21. If desired, a 0.038 in [0.97 mm] guidewire may be introduced through the Select Catheter.
22. Advance the Delivery Catheter and Select Catheter or guidewire system, and then select the appropriate vessel with the Select Catheter or 0.038 in [0.97 mm] guidewire.
23. If using, remove the Select Catheter, and replace with a 0.038 in [0.97 mm] or smaller guidewire.
24. Advance the Delivery Catheter and guidewire to the vascular site and remove the guidewire.
25. After the procedure, remove the Delivery Catheter when clinically indicated by placing compression on the vessel above the puncture site and slowly withdraw the Delivery Catheter. Discard the Delivery Catheter appropriately.
26. Important: Upon removing the Delivery Catheter, aspirate via the HVA to collect any fibrin that may have been deposited within or at the tip of the catheter.

**CLINICAL BENEFITS AND PERFORMANCE CHARACTERISTICS**

The Penumbra Access Catheter System aids in the use of interventional devices to the peripheral and neurovascular anatomy. The detailed listing of all known benefits is included in the Summary of Safety and Clinical Performance (SSCP).

**SAFETY AND CLINICAL PERFORMANCE**


A copy of the Summary of Safety and Clinical Performance can be viewed by searching the device name on the EUDAMED website: <https://ec.europa.eu/tools/eudamed>

Basic UDI DI: 081454801ACCESSY8

**STORAGE AND DISPOSAL**

- Store in a cool, dry place.
- Dispose of device in accordance with hospital, administrative, and /or local government policy.

**SYMBOLS GLOSSARY**

	Attention, see Instructions for Use		
	Prescription only – US Federal Law restricts this device to use by or on the order of a physician		
	Nonpyrogenic		
	Sterilized using ethylene oxide. Single sterile barrier system with protective packaging outside.		
	Do not use if package is damaged.		
	Medical Device		Date of Manufacture
	Manufacturer		EU Representative
	Lot number		Do not reuse
	Use by		Catalogue number
	Do not resterilize		Not made with natural rubber latex

**WARRANTY**

Penumbra Inc. (Penumbra) warrants that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular use. Handling, storage, cleaning, and sterilization of this device as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond Penumbra's control directly affect the device and the results obtained from its use. Penumbra's obligation under this warranty is limited to the repair or replacement of this device and Penumbra shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this device. Penumbra neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. Penumbra assumes no liability with respect to devices reused, reprocessed, or resterilized and makes no warranties, expressed or implied, including but not limited to merchantability or fitness for intended use, with respect to such devices.

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