pNOVUS 21 Microcatheter

DIRECTIONS FOR USE

DIRECTIONS FOR USE

pNOVUS 21 Microcatheter

 phenox Ltd
Kamrick Court,
Ballybrit Business Park,
Galway, Ireland H91 XY38
Phone: +353 91 740100
Email: info@phenox.ie
Website: www.phenox.ie

Contents

- 1 x pNOVUS 21 Microcatheter
- 1 x shaping mandrel
- 1 x introducer sheath

The Product Family Includes the Following:

- pNOVUS 21 150cm Microcatheter
- pNOVUS 21 160cm Microcatheter

Device Description

The pNOVUS 21 Microcatheter (pNOVUS 21) is a variable stiffness single lumen catheter. The pNOVUS 21 has a radiopaque marker on the distal end to facilitate fluoroscopic visualization, and a luer hub on the proximal end to facilitate the infusion of agents and the smooth transfer of other devices (e.g. guidewires) into the inner lumen of the microcatheter. The distal 100cm of the microcatheters external surface has a hydrophilic coating applied for increased lubricity during use. The pNOVUS 21 has a straight tip. Steam shaping of the distal tip allows for one-time customizing of the tip shape.

The microcatheter is supplied sterile in a packing hoop, within a sealed pouch and shelf carton configuration. An introducer sheath and shaping mandrel are supplied on a backing card inside the pouch. The packaging is designed to facilitate ease of handling and aseptic presentation of the device.

Table 1 along with Figure 1 presents the critical dimensions of the pNOVUS 21.


Figure 1: Schematic Drawing of pNOVUS 21

Table 1: Dimensions of pNOVUS 21

	Description	pNOVUS 21 150	pNOVUS 21 160
A	Effective length	150 cm	160 cm
B	Total length	157 cm	167 cm
C	Distal section length	6.5 cm	6.5 cm
D	Hydrophilic coating length	100 cm	100 cm
E	Inner diameter	0.021 inch / 1.6F / 0.54 mm	0.021 inch / 1.6F / 0.54 mm
F	Outer diameter distal	0.033 inch / 2.5F / 0.83 mm	0.033 inch / 2.5F / 0.83 mm
G	Outer diameter proximal	0.035 inch / 2.7F / 0.90 mm	0.035 inch / 2.7F / 0.90 mm
Distance of distal tip to radiopaque marker		1.1 mm	

Intended Purpose

The pNOVUS 21 Microcatheter is intended to support the stable access of suitable therapeutic or diagnostic devices, i.e. devices compatible with the inner diameter of the microcatheter, and the infusion of therapeutic and diagnostic agents subject to use by their IFU.

Indications for Use

The pNOVUS 21 Microcatheter is used in the neuro vasculature during diagnostic and/or therapeutic procedures.

Contraindications

The pNOVUS 21 Microcatheter is not intended for use in peripheral or coronary vasculature.

Patient Target Group

The pNOVUS 21 Microcatheter is intended for patients undergoing diagnostic and/or therapeutic procedures in the neurovasculature.

Warnings

- The pNOVUS 21 may only be used by physicians thoroughly trained in percutaneous, intravascular techniques and procedures.
- The microcatheter must only be manipulated under fluoroscopy inside the body. Limit the exposure to X-ray radiation doses to patients and physicians by using sufficient shielding, reducing fluoroscopy times, and modifying X-ray technical factors when possible.
- Do not advance or retract the microcatheter against resistance without careful assessment of the cause using fluoroscopy. If the cause cannot be determined, withdraw the microcatheter. Excessive force against resistance may result in damage of the microcatheter, damage or dislodgment of other interventional devices, or vessel injury.
- The safety and effectiveness of the device has not been established or is unknown in vascular regions other than those specifically indicated.
- This device is coated with a hydrophilic coating at the distal end. Hydrate the microcatheter for at least 30 seconds before removing from the hoop.
- Do not use if packaging is open or damaged, as sterility cannot be assumed otherwise.
- Do not exceed 150 psi injection pressure. Excessive pressure may result in catheter damage, vessel injury or thromboembolism.
- This device is intended for single use only. Discard the microcatheter after one procedure. Structural integrity and/or function may be impaired through re-use or cleaning. Microcatheters are extremely difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused.
- Endovascular procedures have a higher risk of complication for individuals with known generalized vessel wall weakness (e.g., Ehlers-Danlos syndrome). The risks and benefits must be thoroughly assessed for these individuals before proceeding.

Precautions

- Only use the device if DFU and label can be read clearly and can be fully understood and if the label properly adheres to the packaging.
- Inspect the microcatheter prior to use to verify its condition is suitable for the specific procedure. Inspect for any damage, surface irregularities, bends or kinks. Do not use kinked, damaged, or irregular microcatheters.
- Only insert the microcatheter into an appropriately sized vessel.
- The safety and effectiveness of the device has not been established for neonatal use.
- The pNOVUS 21 is supplied sterile and intended for single use only. Do not re-sterilize or re-use. Re-sterilization or Re-use may result in ineffective microcatheter coating lubrication, which may result in high friction and the inability to access the target neurovasculature.
- Sufficient irrigation of the microcatheter is necessary to remove enclosed air bubbles.
- Use prior to the "Use By" date.
- Only select compatible guide wires, guide catheters and instruments or agents to be used with the microcatheter.
- Incorrect steam shaping of the microcatheters may result in deformation or damage of the hydrophilic coating. Moreover, improper device delivery and deployment may be caused, depending on the degree of shaping and catheter deflection during device delivery.
- Do not expose the microcatheter to organic solvents, such as disinfection agents, alcohol or DMSO.
- Do not wipe the device.
- Do not excessively torque the device. The device may fail if 12 or more turns are applied.
- The shaping mandrel and introducer sheath are not intended for use in the human body.
- If there is suspicion of vasospasm in the affected vascular region, all measures required for reversion should be taken, e.g. with appropriate drug therapy.
- Appropriate premedication with antiplatelet agents and systemic heparinization are to be considered according to the type and duration of the intervention.
- In case of any unusual microcatheter behavior, stop the procedure and replace the microcatheter.
- Avoid blockage of blood flow by the microcatheter.
- If flow through the catheter is hampered, do not try to reestablish flow by forceful flushing. The catheter must be removed and replaced.
- Return all damaged packaging to the manufacturer or distributor.

Storage

Keep away from heat. Store in a cool, dry place.

Disposal

The device and if applicable, its packaging, has to be disposed as biohazard waste.

Potential Adverse Events

The following and further adverse events may occur during or despite the use of the pNOVUS 21:

- Air-, catheter-, plaque-, thrombus- emboli
- Aneurysm rupture / perforation
- Allergic or toxic reaction
- Adverse response by patient (other than above)
- Emboli (thromboembolic/device components)
- Extracranial Hemorrhage / Hematoma
- Extracranial / intracranial complications during procedure (other than above)
- Harm of 3rd person via contaminated products / infection
- Hemorrhagic stroke
- Infection
- Inflammation
- Ischemic stroke
- Neurovascular / Neurological complications
- Patient Injury
- Severe harm (permanent disability / Death)
- Transient ischemic Attack (TIA)/seizure
- Vasospasm
- Vessel occlusion/thrombosis
- Vessel perforation
- Vessel dissection
- Vessel damage (other than above)

Compatibility

Devices used with the pNOVUS 21 must be compatible with the microcatheter's dimensions according to the manufacturer's description.

The pNOVUS 21 is not compatible with glue, glue mixture or non-adhesive liquid embolic agents.

pNOVUS 21 is compatible with the following products (not supplied – selected based on compatibility and physician's experience and preference):

- Guidewires with a maximum outer diameter of 0.018 inch/0.46mm
- (Guide) Catheters and Distal Access Catheters with a minimum inner diameter of 0.053 inch/1.35mm / 4.1F
- Non-liquid interventional devices which are declared compatible to a microcatheter with an inner diameter of 0.021 inch / 0.54mm according to the respective IFU/DFU
- In particular, pNOVUS 21 is compatible with the phenox pRESET Thrombectomy Device (size 4-20, 5-40, 6-30 and 6-50) and phenox p48 MW (HPC) and p64 MW (HPC) Flow Modulation Devices
- Therapeutic and diagnostic agents:
 - Saline
 - Contrast Media
 - Agents subject to use by their IFU/DFU
- Rotating hemostatic valves (RHV), with a standard luer adaptor
- Syringes, with a standard luer adaptor

Recommended Procedure

1. Refer to Warnings, Precautions, and Potential Adverse Event's prior to use.
2. Introduce a suitable catheter sheath into an appropriate vessel according to standard practice and the respective IFU/DFU.
3. Prepare a guide catheter and/or distal access catheter for use per its DFU. Connect the RHV and continuous pressure infusion to the guide catheter.
4. Insert the guide catheter into an appropriate vessel by the help of a suitable guidewire according to standard practice and their IFU/DFU's.
5. Alternatively, in case a long sheath is used, the procedure may be conducted without an additional catheter.
6. Confirm the appropriate length pNOVUS 21 Microcatheter is selected.
7. Carefully open microcatheter packaging. Peel open the pouch from the chevron end and remove the hoop from pouch.

8. Flush the hoop with sterile saline through the luer fitting connected to the end of the packaging hoop before removing the microcatheter. Hydrate with sterile saline for at least 30 seconds.
 9. Carefully withdraw the microcatheter from the hoop.
 10. Inspect the microcatheter thoroughly to ensure that there is no damage or irregularities. If there are any irregularities or damage observed, replace the device.
 11. Flush the inner lumen of the microcatheter with sterile saline.
 12. Do not allow the microcatheter to dry once hydrated or after tip shaping.
 13. If desired, shape the microcatheter tip using the following steps (a-g):
The following instructions should strictly be followed to maintain the integrity of the inner diameter of the microcatheter
 - a. Remove the shaping mandrel from its packaging
 - b. Insert the shaping mandrel's entire working length into the microcatheter's distal lumen. Flush or dip the microcatheter tip into sterile saline to aid the mandrel's movement.
 - c. Carefully bend the shaping mandrel and distal microcatheter portion into the desired shape. A slight exaggeration may be required to account for microcatheter relaxation.
 - d. Shape the microcatheter by holding the shaped portion in front of a steam source for approximately 30 seconds. Hold the microcatheter no closer than 2.54 cm (1 inch) from the steam source.
 - e. Before removing the mandrel, allow the microcatheter tip to cool in air or saline. Remove the mandrel and discard. Flush the microcatheter lumen after removal of the mandrel.
 - f. Multiple shaping is not recommended.
 - g. Inspect the microcatheter tip for damage that may have resulted from steam shaping. If damage is found, do not use the microcatheter.
 14. Connect the RHV and continuous pressure infusion to the microcatheter.
 15. Prepare the guidewire according to the respective IFU/DFU.
 16. Carefully insert the guidewire into the microcatheter. If desired, the guidewire- microcatheter assembly may be inserted into a compatible larger catheter (min. inner diameter of 0.053 inch / 1.35mm / 4.1F).
 17. Introduce the guidewire- microcatheter assembly through the guiding catheter's RHV using the peel away introducer sheath.
 18. Remove the introducer sheath from the guidewire- microcatheter assembly by retracting the introducer sheath from the RHV and peeling it off.
- Warning:** Removing the introducer sheath from the microcatheter without a guidewire inserted may cause damage to the microcatheter.
19. Tighten the O-ring valve around the microcatheter just enough to prevent back flow, but not so tight as to inhibit microcatheter movement.
- Warning:** Excessive tightening of the RHV around the microcatheter may lead to damage of the device.
20. Maintain continuous flush in the microcatheter and support catheter.
 21. Advance the guidewire and microcatheter past the distal tip of the guide catheter.
 22. Carefully advance the guidewire and microcatheter to the intended vessel under fluoroscopy. Advance the guidewire and follow with the microcatheter in an alternate fashion. If a larger catheter is used in conjunction with the microcatheter, follow the microcatheter with the larger catheter.
 23. Remove the guide wire from the microcatheter after the target vessel has been reached.
 24. Infuse the therapeutic or diagnostic agent or advance the therapeutic or diagnostic device according to the manufacturer's description.
 25. When the procedure is complete, remove the microcatheter carefully under fluoroscopy.

NOTE: Flow rates for the given pressure and infusates are provided below as a reference and are approximate values.

Flow Rate Table			
pNOVUS 21 Microcatheter			
Length (cm)	Pressure (psi)	100% Saline (mL/sec)	50% Ultravist-370* / 50% Saline (mL/sec)
150	150	1.17	0.96
160	150	1.15	0.81

*The dynamic viscosity of Ultravist® 370 injectate at 37°C is 10 mPa·s.

Performance Characteristics

The pNOVUS 21 Microcatheter used in combination with standard guide catheters and accessories as well as compatible therapeutic and diagnostic devices and /or agents enables access/diagnosis/treatment of conditions in the neurovasculature.

Symbols Glossary as per ISO 15223-1

Reference Number or Symbol	Title of Symbol	Description of Symbol
 5.1.1	Manufacturer	Indicates the medical device manufacturer, as defined in the EU directives 90/385/EEC, 98/42/EEC and 98/79/EC.
 5.1.3	Date of manufacture	Indicates the date when the medical device was manufactured
 5.1.4	Use-by date	Indicates the date after which the medical device is not to be used.
 5.1.5	Batch code	Indicates the manufacturers batch code so that the batch or lot can be identified.
 5.1.6	Catalogue number	Indicates the manufacturers catalogue number so that the medical device can be identified.
 5.2.3	Sterilized using ethylene oxide	Indicates medical device has been sterilized using ethylene oxide.
 5.2.6	Do not resterilize	Indicates a medical device that is not to be resterilized.
 5.2.8	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened.
 5.3.2	Keep away from sunlight	Indicates a medical device that needs protection from light sources.
 5.3.4	Keep dry	Indicates medical device that needs to be protected from moisture.
 5.4.2	Do not reuse	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
 5.4.3	Consult instructions for use	Indicates the need for the user to consult the instructions for use.
 5.4.4	Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
 5.6.3	Non-pyrogenic	Indicates a medical device that is non-pyrogenic.
 0344	CE Mark with Notified Body Number/ ID	Indicates that a device has received CE approval and the number of the Notified Body.
 5.7.7	Medical Devices Symbol	Indicates the item is a medical device.
 5.2.14	Single Sterile Barrier System with Protective Packaging Outside	Indicates a single sterile barrier system with protective packaging outside. Used in conjunction with symbol 5.2.3.
 5.7.10	Unique Device Identifier	Indicates a carrier that contains Unique Device Identifier information.
Other Symbols included in this document		
 CONT	Content	Indicates the contents included in the packaging
 hydrophilic	Hydrophilic coating length	Indicates the location of the hydrophilic coating length on the device drawing
	Inner diameter of microcatheter	Indicates the inner diameter of the microcatheter

Summary of Safety and Clinical Performance

The summary of Safety and Clinical Performance (SSCP) for this device is available to view at

<https://phenox.net/international/mdr/mc/>.

Once Eudamed is available, the SSCP will also be available to view on

<https://webgate.ec.europa.eu/eudamed/landingpage> with UDI-DI 539153091 040 8 (PNOV-21-150) and UDI-DI 539153091 038 5 (PNOV-21-160).

Specification of the Clinical Benefit to be Expected

This device aids in reaching the target vessel in the neurovasculature and delivery of therapeutic and diagnostic devices and/or agents.

For Use by Professional Clinicians Trained
in the Use of these device types only

Any serious incident that has occurred in
relation to the device should be reported to
the manufacturer and the authority per
complaints guidance in your jurisdiction.

Limits of Liability

phenox Ltd. shall not be liable for damage caused by use other than that intended or reuse of the product.

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phenox Ltd, Kamrick Court, Ballybrit Business Park, Galway, Ireland H91 XY38