Endo-SPONGE®

CS	Návod k použití	-	5	LV	Lietošanas instrukcija	-	43
DA	Brugsvejledning	-	8	NL	Gebruiksaanwijzing	-	46
DE	Gebrauchsanweisung	-	11	NO	Bruksinstruksjoner	-	49
EL	Οδηγίες χρήσης	-	14	PL	Instrukcja stosowania	-	52
EN	Instructions for use	-	18	PT	Instruções de uso	-	56
ES	Instrucciones de uso	-	21	RO	Instrucțiuni de utilizare	-	59
ET	Kasutusjuhend	-	24	SK	Návod na použitie	-	62
FI	Käyttöohjeet	-	27	SL	Navodila za uporabo	-	65
FR	Mode d'emploi	-	30	SV	Användarinstruktioner	-	68
HU	Használati utasítás	-	33	TR	Kullanım talimatı	-	72
IT	Istruzioni per l'uso	-	37				
LT	Naudojimo instrukcijos	-	40				

CE 0123

02-09/23 506570

B | BRAUN

SHARING EXPERTISE


B. Braun Surgical, S.A.
Carretera de Terrassa, 121
08191 Rubí, Spain

Phone +34 93 586 62 00
Fax +34 93 588 10 96


www.bbraun.com


1


7


13


19


2


8


14


20


3


9


15


21


4


10


16


22


5


11


17


23


6

12

18

24

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(EN) Medical device / Catalogue number / Batch code / Use by date / Date of manufacture / Sterile unless package is opened or damaged. Method of sterilization: Ethylene Oxide / Single Sterile Barrier System / Consult Instructions for Use / Caution / Do not resterilize / Do not reuse / Do not use if package is damaged and consult instructions for use / Unique Device Identifier / power button

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Endo-SPONGE®

Description of the device

Endo-SPONGE®, an endoscopic vacuum therapy (EVT) device, is a minimally invasive method intended for drainage of the lower pelvic area, indicated for the treatment of anastomotic or Hartmann's stump leakages.

The Endo-SPONGE® system consists of a drainage tube with an attached open-pore polyurethane sponge, a connection system (Y-piece with an extensor tube), an application system (overtube, pusher) and a rinsing set (syringe, tip).

The drainage consists of Redon Drainage (12 CH) with the open-pore sponge attached to the end. The drainage tube in the sponge has perforated holes on its side. The sponge on the drainage tube can be cut to size for the respective application.

The Endo-SPONGE® application system consists of one coaxially arranged tube, the insertion tube made of silicone (overtube) and the more rigid inner tube with a handle (pusher). The lumen of the overtube is slightly larger than the outside diameter of the applied endoscope and is used as a guide for insertion of the sponge system. The pusher is used to push the sponge forward and position it.

Two different vacuum sources are compatible with Endo-SPONGE®: Redyrob® Trans Plus bottle and Endo-SPONGE® Pump. Depending on the vacuum source, the connection system may vary.

The Endo-SPONGE® product provides an extensor tube which has a Y-piece (with a double-drain connector) and a Luer lock connector which allows connection to the Redyrob® Trans Plus bottle.

If an Endo-SPONGE® Pump is used as a vacuum source, the connection between the Redon Drain Sponge and the pump is made by means of a Y-piece attached to the suction tube included in the Endo-SPONGE® Pump Disposable canister.

Contents

1. Endo-SPONGE®, open-pore PUR sponge (ø 3.3 x 7.5 cm) with Redon drain CH12, med. PVC, 40 cm long.

2. Pusher, ABS + PVC, CH 30, 30 cm long.

3. Overtubes in 2 sizes, depending on device and sponge size:

- Silicon tube, each 29 cm long.

- Tapered rounded tip.

- Size 1: inner diameter 13 mm, outer diameter 17 mm.

- Size 2: inner diameter 15 mm, outer diameter 19 mm.

4. Irrigation set comprised of:

- 20 ml syringe

- Dilator: component made of PP used to facilitate the rinsing of the drainage.

- Slide clamp: component made of PP used to clamp off the Redon drain while removing the Endo-SPONGE®

5. Y connecting tube with Luer lock fitting for Redyrob® Trans Plus bottle (purchased separately – Ref. 5526604).

NOT INCLUDED

1. Vacuum source:

a) Redyrob® Trans Plus bottle (Ref. 5526604) – controllable wound drainage system.

b) Endo-SPONGE® Pump (Ref. 5526650) and Endo-SPONGE® Pump Disposable canister. (Ref. 5526653) – variable speed medical vacuum pump and canister with a suction tube with a Y connector. This pump is intended to be used for internal application supplying a constant negative vacuum pressure of 150-200 mmHg.

2. Sterile hydrogel based on glycerol.

3. Standard or therapeutic flexible endoscope.

4. Additional endoscopic accessories to be used at the physician's discretion.

Materials used

Acrylonitrile Butadiene Styrene (ABS), Polyvinyl chloride (PVC), Polyurethane (PUR), Polyethylene terephthalate (PET), Hydrogel, Silicone, Polyethylene (PE), Polypropylene (PP), Isoprene rubber (IR).

Indications for use

Treatment of anastomotic or Hartmann's stump leakages following colorectal surgery in the lower pelvic area (extraperitoneal position) by means of negative pressure.

The leak must have created a drainable cavity with or without a local infection. Endo-SPONGE® is intended for use in adults. Since limited clinical evidence is available in the paediatric population, Endo-SPONGE® is not indicated for the paediatric population.

Moreover, since there is no clinical evidence for the use of Endo-SPONGE® in pregnant and breastfeeding women either, Endo-SPONGE® is not indicated for use during pregnancy or breastfeeding.

This treatment should only be performed by experienced doctors with many years of practice both in interventional treatment of the lower gastrointestinal tract using flexible endoscopy and in negative pressure wound therapy in general.

Mode of action

Endo-SPONGE® consists of an open-pore sponge connected to a drainage tube. After endoscopic insertion of the sponge into the leakage cavity, the drainage tube is then routed out through the anus and connected to a vacuum system. The application of the vacuum leads to continuous drainage of the fluid and the sponge in the cavity promotes cleaning of the surface. To achieve an effective treatment, the size of the sponge is cut to fit the cavity. Depending on the size of the leakage cavity, it may be necessary to place more than one sponge into the cavity. The sponge system is changed every 24-72 hours. To change the sponge, the vacuum is disconnected. Removal of the sponge must be done with prior irrigation with 0.9 % saline solution to remove the granulated tissue from the surface of the sponge. The sponge is removed through the anus and the size of the new sponge is adapted to the size of the leakage cavity. The Endo-SPONGE® treatment is stopped once the cavity reaches a size of less than 2 x 1 cm, because no further reduction of the sponge size is technically possible.

Contraindications

- Contraindications derived from the device's intended purpose

- Malignant tumour wound

- Necrotic tissue/gangrene

- Untreated osteomyelitis

- Sponge position directly adjacent to vessels, urinary bladder or small intestinal loops

- Non-drainable septic focus

- Clotting disorders

- Treatment with a therapeutic dose of anticoagulant or platelets aggregation inhibitor drugs

- Generalised peritonitis or sepsis

- Contraindications derived from the residual risks

- Patients with known sensitivities or allergies to its components (refer to Materials used for details about the components).
- The dimensions of the device should be taken into consideration for its use in specific patient groups (e.g., small-framed persons).

Mode of application

1. Preparation of the Endo-SPONGE®

1.1. Preparing the patient

- Prepare the patient according to the usual standards for an endoscopy. Follow all facility standards and guidelines.

1.2. Cleaning the leak cavity

- Clean the leak cavity using appropriate methods that follow facility-approved procedures, such as irrigating with Ringer's solution.

1.3. Measuring the leak cavity

- Insert the endoscope in order to provide visualisation for measuring the size of the leak cavity.

- Measure the dimensions (length and diameter) of the leak cavity.

Note: it is recommended to use radiographic imaging to help with measuring the dimensions and surroundings of the cavity (Fig. 1).

- Remove the endoscope after measuring the leak cavity.

1.4. Sizing the sponge

The sponge should fill the leak cavity without touching the intact mucosa. Follow the instructions below to re-size the sponge if needed:

- *If shortening the sponge length by 1 cm or less:* cut from the distal end of the sponge.

- *If shortening the sponge by more than 1 cm:* shorten the drainage tube as well so the sponge protrudes at least 3 mm from the end of the tube.

- To shorten the drainage tube, peel the sponge back to access the end of the drainage tube, then cut the drainage tube (Fig. 2). Ensure that the drainage tube is cut below the top of the sponge surface.

Note: particles and sponge residues are produced when the sponge is cut to size. The sponge must be cut at an appropriate distance from the patient and in a suitable environment in which particles are permitted.

- Document the length of material that was removed from the sponge. This length will affect an insertion step later in the procedure.

- Inspect the sponge for edges and points. Round off any edges or points found on the sponge. Edges and points could result in sponge particles being left in the leak cavity during removal.

- Remove all sponge residues and particles from the sponge surface by tapping the sponge.

1.5. Selecting the large or small overtube

- *If the sponge has not been resized:* use the large overtube.
- *If the sponge has been resized in length and/or width:* use the small or large overtube according to clinical choice.

- Push the endoscope into the overtube to check that the endoscope slides easily into the overtube.

Note: if the endoscope does not slide easily into the overtube, apply a sterile hydrogel to the endoscope or to the inside of the overtube.

1.6. Preparing for use with rigid rectoscopy or proctoscopy (if necessary)

Note: if using a rigid rectoscopy or proctoscope, it will also serve as the overtube.

- *If using a rectoscopy or proctoscope:* trim the pusher so it is no more than 1 cm longer than the rectoscopy or proctoscope (Fig. 3).

1.7. Inserting the endoscope into the overtube

- Push the non-tapered end of the overtube over the endoscope until the tapered end is about 3 mm from the end of the endoscope (Fig. 4).

2. Inserting the overtube

2.1. Inserting into the anal canal

- Lubricate the outside of the overtube (or rectoscopy or proctoscope) with sterile hydrogel.

- Insert the endoscope and overtube (or rectoscopy or proctoscope) into the anal canal. Gently push the endoscope and overtube until they reach the base of the leak cavity (Fig. 5).

2.2. Pushing into the leak cavity

- Slowly push the endoscope with the overtube (or rectoscopy or proctoscope) forward towards the end of the leak cavity (Fig. 6).

Note: if the entrance to the leak cavity is too small for the device to go through, the defect can be carefully enlarged using appropriate endoscopic techniques.

2.3. Pushing into the leak cavity

- Firmly stabilise the endoscope and keep it still to avoid further penetration.
- Push the overtube forward over the stationary endoscope until the tapered end is close to the distal end of the leak cavity, without touching the cavity wall (Fig. 7). Closely watch the position of the overtube while pushing it forward.

2.4. Removing the endoscope

- Firmly stabilise the overtube and keep it still to avoid further penetration.
- Gently remove the endoscope from the leak cavity and anal canal while

keeping the overtube still (Fig. 8).

3. Inserting the Endo-SPONGE®

3.1. Lubricating the sponge

Note: after performing step 3.1, continue with step 3.2 without delay. The gel may dry out if step 3.2 is not performed quickly, and the lubrication will not be effective.

- Lubricate the lumen of the overtube (or of the rectoscopy or proctoscope) and the re-sized sponge sufficiently with a sterile hydrogel based on glycerol.
- Perform step 3.2 without delay.

3.2. Inserting the sponge

- Firmly stabilise the overtube and keep it still to avoid further penetration.
- While stabilising the overtube, compress the sponge and insert it fully into the overtube (or rectoscopy or proctoscope). (Fig. 9)

3.3. Positioning the sponge

- Place the drain of the sponge inside the pusher.
- Continue stabilising the overtube and keep it still to avoid further penetration.

- Position the sponge in the overtube according to options A or B below:

a) If the length of the sponge was not previously cut: - push the sponge forward slowly and using steady force until the non-tapered end of the overtube reaches the black mark on the pusher (Fig. 10).

b) If the length of the sponge was previously shortened: push the sponge forward slowly and using steady force until the black mark on the pusher exceeds the non-tapered end of the overtube by the same length the sponge was cut (Fig. 11).

- The tip of the sponge should now be positioned at the deepest point in the leak cavity.

3.4. Expanding the sponge

- Keep the pusher stabilised and still to avoid further penetration.
- Slowly retract the overtube (or rectoscopy or proctoscope) until it touches the handle of the pusher.

- The sponge will now expand in the leak cavity.

- Remove the pusher and overtube (or rectoscopy or proctoscope) carefully without dislocating the sponge (Fig. 12).

3.5. Checking the position of the sponge

- Using the endoscope, check the position of the sponge carefully (Fig. 13). The sponge should fill the leak cavity without touching the intact mucosa. If the sponge is not properly positioned or sized in the leak cavity, the treatment will not be effective.

Note: make sure the sponge is not disturbed by the endoscope.

- If the sponge does not fill the leak cavity, a second sponge should be added to the leak cavity.

4.a. Connecting the Redyrob® Trans Plus bottle

4.a.1. Connecting the drainage tube

- Connect the drainage tube of the sponge to the open port on the Y connecting tube.

Note: if necessary, the end of the drainage tube can be cut at an angle for easier insertion.

- Insert the drainage tube to the minimum depth marked by the dashed line in Fig. 14.

- If only one sponge has been inserted into the patient, continue to step 4.a.2. If using two or three sponges, continue with the following sub-steps.

- *If connecting two or three sponges:* cut off the top section of the opposite port of the Y connecting tube as described in Fig. 15.

- Connect the second drainage tube to the other side of the Y connecting tube.

- *If using three sponges:* connect the third drainage tube to the open port of a separate (second) Y connecting tube. This Y connecting tube will eventually be connected to a second Redyrob® Trans Plus bottle.

4.a.2. Removing the red stopper

- Unscrew and remove the red stopper from the Redyrob® Trans Plus bottle (Fig. 16).

4.a.3. Connecting the Luer lock fitting

- Connect the white Luer lock fitting of the Y connecting tube to the line on the Redyrob® Trans Plus bottle (Fig. 17).

- *If using three sponges in the leak cavity:* repeat this step by connecting the second Y connecting tube to the second Redyrob® Trans Plus bottle.

4.a.4. Setting the control knob

- Set the control knob on the Redyrob® Trans Plus bottle to suction power 1 (Fig. 18).

- Observe the flow of secretion and check the quantity in the Redyrob® Trans Plus bottle.

Notes:

- *Secretion of liquid normally begins immediately. If it does not, all the connections of the system should be checked, as well as the setting of the control knob.*

- *Additional sealing is not necessary, since the patient's sphincter apparatus*

will create a tight closure.

- Additional fixing of the drainage, which is sticking out per anum, is likewise unnecessary, since the sponge is kept in the leak cavity by the application of suction.

Warnings:

Do not use suction power 2 or 3 of the Redyrob® Trans Plus bottle.

Do not use conventional high vacuum or medium vacuum wound drainage bottles, as they produce too much suction.

4.a.4. Inspection of the system

- Inspect the system regularly. If the vacuum gauge of the Redyrob® Trans Plus bottle (blue column indicator, Fig. 19) shows "Low", the vacuum bottle must be changed to a new one.

- Maximum use of the Endo-SPONGE® is 72 hours. At that point, it must be replaced, since the growth of granulation tissue into the sponge pores can make it difficult to remove the sponge.

4.b. Connecting the Endo-SPONGE® pump.

4.b.1. Setting the vacuum system

- Connect the Endo-SPONGE® Pump, Ref. 5526650 (not supplied in the set) to the Endo-SPONGE® Pump Disposable canister, Ref. 5526653 (not supplied in the set). (Fig. 20)

Notes:

- The Endo-SPONGE® Pump Disposable canister includes a suction tube with a Y connector.

- Follow the Instructions for Use of the Endo-SPONGE® Pump.

4.b.2. Connecting the drainage tube

- Connect the drainage tube of the sponge to the open port on the Y connecting tube provided with the Endo-SPONGE® Pump Disposable canister.

Note: if necessary, the end of the drainage tube can be cut at an angle for easier insertion.

- Insert the drainage tube to the minimum depth marked by the dashed line in Fig. 21.

- If only one sponge has been inserted into the patient, continue to step 4.b.3. If using two or three sponges, continue with the following sub-steps.

- If connecting two or three sponges: cut off the top section of the opposite part of the Y connecting tube as described in Fig. 22.

- Connect the second drainage tube to the other side of the Y connecting tube.

- If using three sponges: connect the third drainage tube to the open port of a separate (second) Y connecting tube. This Y connecting tube will eventually be connected to a second Endo-SPONGE® Pump.

4.b.3. Setting the negative pressure

Set a negative pressure of 125-200 mmHg following the Instructions for Use of the Endo-SPONGE® Pump.

Notes:

- Secretion of liquid normally begins immediately. If it does not, all the connections of the system should be checked.

- Additional sealing is not necessary, since the patient's sphincter apparatus will create a tight closure.

- Additional fixing of the drainage, which is sticking out per anum, is likewise unnecessary, since the sponge is kept in the leak cavity by the application of suction.

Warnings:

For more information on using the Endo-SPONGE® Pump, please refer to its IFU.

4.b.4. Inspection of the system

- Inspect the system regularly. If an alarm sounds, follow the Instructions for Use of the Endo-SPONGE® Pump to find out how to resolve it.

- Maximum use of the Endo-SPONGE® is 72 hours. At that point, it must be replaced, since the growth of granulation tissue into the sponge pores can make it difficult to remove the sponge.

5. Removal of the Endo-SPONGE®

5.1.a. Disconnecting the Redyrob® Trans Plus

- Change the position of the control knob to 0.

- Close the Y connecting tube with the slide clamp.

- Disconnect the white Luer lock fitting of the Y connecting tube from the Redyrob® Trans Plus line.

5.1.b. Disconnecting the Endo-SPONGE® Pump

Switch off the Endo-SPONGE® Pump by pressing the power button for approximately 3 seconds.

5.2. Irrigating the Endo-SPONGE® drain line

- Disconnect the Y connector from the Endo-SPONGE® drain line.

- Fill the provided syringe with 20 ml of Ringer's solution or 0.9% saline solution.

- Attach the green syringe cap to the syringe (Fig. 23).

- Flush 20 ml of Ringer's or 0.9% saline solution through the Endo-SPONGE® drain line (Fig. 24).

- Attach the slide clamp to the Endo-SPONGE® drain line and clamp the tubing closed.

- The irrigation step should be repeated 2-3 times for best removal.

5.3. Removing the Endo-SPONGE® from the patient

- Pull the drainage tube(s) out carefully, steadily and with increasing force until the sponge(s) detaches from the leak cavity.

- After the sponge(s) has detached, continue to pull with a small amount of force in order to remove the tube(s) from the anal canal.

- If removing more than one Endo-SPONGE®: remove the sponges under endoscopic vision to ensure the easiest order of removal.

a) Ensure that the sponge that is intended to be removed is identified and removed.

b) Ensure that the drainage tubes are not tangled. The most proximal sponge tube must be removed first.

5.4. Checking the integrity of the sponge

- Thoroughly inspect the integrity of the sponge. Look for breakage and missing parts.

- Using the endoscope, inspect the treated leak cavity for sponge residue. Remove any sponge residue using endoscopy grasping forceps.

- If it is necessary to insert another Endo-SPONGE®, then perform all steps in the instructions, starting in Section 1 - Preparation of the Endo-SPONGE®.

- Document the outcome of the treatment.

Warning: sponge residue left in patients could result in foreign body reactions, surgical intervention and increased operating times.

6. Subsequent treatment

6.1. Assessing the leak cavity

- Additional Endo-SPONGE® treatments are recommended until both of the following criteria have been met:

a) The length of the leak cavity is less than 2 cm.

b) The diameter is less than 1 cm.

If only one or neither of the above criteria have been met, continue performing subsequent Endo-SPONGE® treatments according to the instructions. Do not use the Endo-SPONGE® treatment for longer than 90 days.

6.2. Resizing subsequent sponges

- If additional Endo-SPONGE® treatment is needed: follow all instructions related to resizing the sponge starting in Section 1 - Preparation of the Endo-SPONGE®.

Warnings

- The Endo-SPONGE® drain sponge is made of plastic materials and is MRI compatible. See Precautions.

- The sponge may generate residual sponge particles when reshaping and/or removing.

- The residual sponge particles could cause fistula formation or foreign body reactions. Possible need for endoscopic removal.

- Endo-SPONGE® must not be used in body openings other than that indicated.

- Due to the underlying disease, most patients have a localised infection that can lead to sepsis (i.e. peritonitis, necrosis, etc.).

- Due to the underlying disease, some patients may develop chronification of sinus.

- Endo-SPONGE® can only be used in combination with the vacuum sources mentioned below:

- Redyrob® Trans Plus bottle, Ref. 5526604 (not included).

- Endo-SPONGE® Pump, Ref. 5526650 (not included) in combination with the Endo-SPONGE® Pump Disposable canister, Ref. 5526653 (not included).

Precautions

- It must be ensured that the drain sponge is not connected to the pump during the MRI exploration, and also that there is no metal clip or metal piece used to stick or fix the drain tube to the patient.

- Particles and sponge residues are produced when the sponge is cut to size. The sponge must be cut (e.g. with scissors or a scalpel) at an appropriate distance from the patient and in a suitable environment in which particles are permitted.

- After cutting the sponge, all residues and particles must be removed from the surface by tapping the sponge, which must be collected and disposed of in the usual manner.

- When reshaping the sponge, ensure that there are no sharp edges or points and round it off, otherwise these can easily break off during removal of the sponge. After cutting, no cut in the sponge or the drainage tube should exist.

- When shortening the sponge length, also shorten the drainage tube. The sponge must protrude at least 3 mm from the end of the drainage tube.

- If using the Endo-SPONGE® Pump, the prescribed negative pressure (125 to 200 mmHg) must be ensured.

- If using the Redyrob® Trans Plus bottle, set the control knob to position 1.

- The dwell time depends on the local clinical situation. If applicable, it is recommended to adjust the sponge change time according to the amount of debris and growth of granulation tissue.

- A dwell time of 48 hours is recommended; more than 72 hours must be ruled out due to the risk of granulation tissue overgrowing into the sponge,

and, as a result, the sponge possibly breaking during removal, leaving part of it in the application area and embedded in granulation tissue. If this occurs, an endoscopic loop must be used to detach the sponge from the surrounding tissue for its removal. The duration of the treatment must not exceed 90 days.

- Prior to application, an appropriate imaging procedure must be performed to rule out an abscess in the wound area, which can only be treated by a surgical or interventional procedure and not with Endo-SPONGE®.
- The flow of secretion must be monitored and the quantity checked. The drainage of secretion usually begins immediately. If this is not the case, the connection between the Endo-SPONGE® and the vacuum source must be checked, as well as its settings.
- Do not use items if the packaging is not intact.
- Be careful not to damage to the endoscope sheath.
- Damage to endoscope cables can occur in the case of excess bending of the flexible distal end within the overtube.
- Do not re-use the Endo-SPONGE®. All items are for single use except the variable speed medical pump (Endo-SPONGE® Pump).
- After each replacement, the sponge must be disposed of due to the risk of clogging or obstruction of the sponge, because as a result of this, the soaking and draining capacity of the sponge may be affected.
- The Endo-SPONGE® device must not be used after the expiry date indicated on the label.

Side effects

- Erosion of structures adjacent to sponge (vessels, urinary bladder, small intestine, colon, etc.)
- Injury to intestinal wall and perforation has been reported in some cases
- Fistula formation
- Abscess has been reported in patients undergoing Endo-SPONGE treatment
- Scarring
- Bleeding, which depending on patient condition could lead to severe bleeding
- Sponge dislocation
- Post-interventional stricture / stenosis
- Pain

Sterilisation

The Endo-SPONGE® is supplied sterile. The Endo-SPONGE® is sterilised by Ethylene Oxide and must not be re-sterilised. Only use the Endo-SPONGE® if the package is undamaged. Check for the absence of channels along the sealing of the sterile barrier system, as well as for the absence of perforations. If such defects are found, the device should be discarded in the recommended manner.

Storage

Store the Endo-SPONGE® at room temperature. Do not expose it to extreme temperatures for extended periods of time. Use the product only if it has been stored correctly and before the expiry date indicated on the packaging.

Device disposal

Once the treatment is completed, the different components of the Endo-SPONGE® must be discarded in specific containers. It is the user's responsibility to determine if disposal material is hazardous according to federal, state and local regulations. Dispose of the contents and the container in compliance with applicable local, state, national and international regulations.

Information for the User/Patient

If a serious incident occurs in relation to the Endo-SPONGE®, it should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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