INTENDED PURPOSE

The INDIGO® Aspiration System is designed to remove thrombus from the vasculature using mechanical aspiration.

Lightning Bolt™ Aspiration Tubing is a sterile aspiration tubing component of the INDIGO® Aspiration System and is designed to serve as a conduit to assist in thrombus removal and restoration of blood flow in the peripheral vasculature and for the treatment of pulmonary embolism.

DEVICE DESCRIPTION

The INDIGO® Aspiration System is comprised of several devices:

- INDIGO Aspiration Catheter
- Penumbra Aspiration Pump
- INDIGO Aspiration Pump Canister
- INDIGO Aspiration Tubing
- INDIGO Separator™

The INDIGO Aspiration System is designed to remove thrombus from the vasculature using mechanical aspiration. The INDIGO Aspiration Catheter targets aspiration from the pump directly to the thrombus. The INDIGO Separator may be used to clear the lumen of the INDIGO Aspiration Catheter should it become blocked with thrombus. The INDIGO Aspiration Catheter is introduced through a guide catheter or vascular sheath into the peripheral vasculature and guided over a guidewire to the site of the primary occlusion. The INDIGO Aspiration Catheter is used with the Penumbra Aspiration Pump to aspirate thrombus from an occluded vessel. As needed, an INDIGO Separator may be deployed from the INDIGO Aspiration Catheter to assist with thrombus removal. The INDIGO Separator is advanced and retracted through the INDIGO Aspiration Catheter at the proximal margin of the primary occlusion to facilitate clearing of the thrombus from the INDIGO Aspiration Catheter tip.

The INDIGO Aspiration Catheter may be provided with a steam shaping mandrel, rotating hemostasis valve and introducer. The INDIGO Separator may be provided with an introducer and torque device. The devices are visible under fluoroscopy.

For the aspiration source, the INDIGO Aspiration Catheter is used in conjunction with the Penumbra Aspiration Pump, which is connected using the INDIGO Aspiration Tubing and the INDIGO Aspiration Pump Canister. The Lightning Bolt Aspiration Tubing (INDIGO Aspiration Tubing) is designed to serve as a conduit to assist in thrombus removal, facilitating transfer of vacuum between the Penumbra Aspiration Pump and the INDIGO Aspiration Catheter while providing intermittent, continuous, or modulated aspiration. Modulated aspiration is provided when the Lightning Bolt Aspiration Tubing alternates between connecting the INDIGO Aspiration Catheter to the Penumbra Aspiration Pump and a sterile saline intravenous (IV) bag at ambient pressure.

INDIGO System Aspiration Catheters have a hydrophilic coating on the distal segment of the catheter shaft (see **Table 1**).

Table 1 – Hydrophilic Coating Length	
INDIGO Aspiration Catheter	Length (cm)
CAT3	75
CAT5	30
CAT6	11
CAT6X	30
CAT7/CAT7D	14
CAT8/CATD	7
CAT12	22

INDICATIONS FOR USE

INDIGO Aspiration Catheters and Separators

As part of the INDIGO Aspiration System, the INDIGO Aspiration Catheters and Separators are indicated for the removal of fresh, soft emboli and thrombi from vessels of the peripheral arterial and venous systems, and for the treatment of pulmonary embolism.

INDIGO Aspiration Tubing

As part of the INDIGO Aspiration System, the INDIGO Sterile Aspiration Tubing is indicated to connect the INDIGO Aspiration Catheters to the Penumbra Aspiration Pump.

Penumbra Aspiration Pump

The Penumbra Aspiration Pump is indicated as a vacuum source for Penumbra Aspiration Systems.

INTENDED PATIENT POPULATION

The intended patient population for the INDIGO Aspiration System is patients with soft emboli or thrombi in peripheral arteries and veins, and/or for patients with pulmonary embolism.

CONTRAINDICATIONS

Not for use in the coronaries or the neurovasculature.

WARNINGS

- The INDIGO Aspiration System should only be used by physicians who have received appropriate training in interventional techniques.
- Do not advance, retract, or use any component of the INDIGO System against resistance without careful assessment of the cause using fluoroscopy. If the cause cannot be determined, withdraw the device or system as a unit. Unrestrained torquing or forced insertion of the catheter or separator against resistance may result in damage to the device or vessel.
- Placing the guidewire too distal in the pulmonary vasculature or excessive manipulation of the aspiration/guiding catheter in the smaller, peripheral, and segmental pulmonary artery branches can result in vessel perforation.
- Do not use the INDIGO Aspiration System with a pump other than a Penumbra Aspiration Pump.
- Do not pressurize the saline IV bag connected to the Lightning Bolt Aspiration Tubing in set up.
- Use of Lightning Bolt Aspiration Tubing adjacent to other equipment should be avoided because it could result in improper operation. If such use is necessary, Lightning Bolt Aspiration Tubing and the other equipment should be observed to verify that they are functioning properly.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 12 inches (30 cm) to any part of Lightning Bolt Aspiration Tubing. Otherwise, this could result in degradation of the performance of this equipment.

PRECAUTIONS

- The device is intended for single use only. Do not resterilize or reuse. Resterilization or reuse of the device may compromise the structural integrity of the device and/or lead to device failure which in turn may result in patient injury, illness, or death.
- Do not use kinked or damaged devices. Do not use open or damaged packages. Return all damaged devices and packaging to the manufacturer/ distributor.
- Use prior to the “Use By” date.
- Use the INDIGO Aspiration System in conjunction with fluoroscopic visualization.
- Maintain a constant infusion of appropriate flush solution.
- When performing aspiration, ensure that the INDIGO Aspiration Tubing is open for only the minimum time needed to remove the thrombus. Excessive aspiration or failure to stop the INDIGO Aspiration Tubing when aspiration is complete is not recommended.
- Hemoglobin and hematocrit levels should be monitored in patients with >700 mL blood loss from the clot aspiration procedure.
- The INDIGO Separator is not intended for use as a guidewire. If repositioning of the INDIGO Aspiration Catheter is necessary during the revascularization procedure, such repositioning should be performed over an appropriate guidewire using standard catheter and guidewire techniques.
- Do not use automated high-pressure contrast injection equipment with the INDIGO Aspiration Catheter because it may damage the device.
- Do not use in the presence of a flammable anesthetic mixture with air or nitrous oxide.
- Do not use in oxygen rich environment.

POTENTIAL ADVERSE EVENTS*

Possible complications include, but are not limited to, the following:

- acute vessel occlusion
- air embolism
- allergic reaction and anaphylaxis from contrast media or device material
- anemia
- arrhythmia
- arteriovenous fistula
- cardiac injury, cardiac perforation, cardiac tamponade
- cardio-respiratory arrest
- compartment syndrome
- death
- emboli
- emergent surgery
- foreign body embolization
- hematoma or hemorrhage at access site
- hemoptysis
- hemorrhage
- hypotension / hypertension
- infarction leading to organ damage
- infection
- ischemia
- myocardial infarction
- neurological deficits including stroke
- pneumothorax
- pseudoaneurysm
- renal impairment or acute renal failure from contrast media
- residual thrombus due to inability to completely remove thrombus or control blood flow
- respiratory failure
- valvular damage
- vessel spasm, thrombosis, dissection (intimal disruption), or perforation

*If a serious incident related to the device occurs, contact your Penumbra, Inc. representative and the competent authority in your respective country/region.

PROCEDURE

- 1. Refer to **Warnings, Precautions, and Potential Adverse Events** prior to use.
- 2. Refer to **Table 2** for device compatibility information.
- 3. As each device of the INDIGO Aspiration System is used, remove the device from the packaging and inspect for damage or for kinks.
- 4. Prepare INDIGO Aspiration System devices for use by flushing the packaging and device with heparinized saline (includes INDIGO Aspiration Catheter and INDIGO Separator only).
- 5. Prepare a guide catheter or vascular sheath according to the manufacturer's Instructions for Use.
- 6. Using conventional catheterization techniques under fluoroscopic guidance, place the guide catheter or vascular sheath into the appropriate artery or vein that is proximal to the thrombus occlusion site over an appropriate guidewire.

Table 2 – Device Compatibility and Selection Information			
INDIGO Aspiration Catheter	Vascular Sheath / Guide Catheter	Guidewire	Indigo Separator
CAT3	5F / 7F	0.014"	SEP3
CAT5	6F / 8F	0.014"	SEP5
CAT6	6F / 8F	0.038"	SEP6
CAT6X	6F / 8F	0.014"	SEP5
CAT7	7F / 9F	0.038"	SEP7
CAT7D	7F / 9F	0.038"	SEP7D
CAT8	8F / 10F	0.038"	SEP8
CATD	8F / 10F	0.038"	SEPD
CAT12	12F / 14F	0.038"	SEP12

The INDIGO Separator is intended to be used only as part of the INDIGO Aspiration System.


INDIGO ASPIRATION CATHETER PREPARATION AND USE

- 1. Confirm vessel diameter and select an appropriately sized INDIGO Aspiration Catheter.
- 2. Attach the provided rotating hemostasis valve to the INDIGO Aspiration Catheter.
- 3. If using a guide catheter, insert the INDIGO Aspiration Catheter into a rotating hemostasis valve connected to the proximal hub of the guide catheter. If a vascular sheath is being used, insert the INDIGO Aspiration Catheter through the valve of the vascular sheath using the introducer (if provided). After inserting the INDIGO Aspiration Catheter through the valve, remove the introducer.
- 4. Advance the INDIGO Aspiration Catheter into the target vessel over the guidewire. Position the INDIGO Aspiration Catheter proximal to the thrombus. Remove the guidewire from the INDIGO Aspiration Catheter.
NOTE: Access to the site of occlusion may be aided by the coaxial use of a compatible inner catheter.

INDIGO LIGHTNING BOLT ASPIRATION TUBING PREPARATION AND USE

- 1. Refer to **Figure 1** for the Lightning Bolt Aspiration Tubing components.


Figure 1: Lightning Bolt Aspiration Tubing

- 2. Attach the Lightning Bolt Aspiration Tubing to the INDIGO Aspiration Canister on the Penumbra Aspiration Pump via the suction connector. Plug the Lightning Bolt Aspiration Tubing power cord into Pump power port.
- 3. Turn on the Penumbra Aspiration Pump (refer to the Penumbra Aspiration Pump *Operation Manual*). Confirm that the Pump reaches maximum vacuum level.

- 4. Connect the Lightning Bolt Aspiration Tubing saline spike to a bag of saline and suspend saline above the patient on an IV pole.
- 5. Press the PRIME button. The prime step is complete when the flow of saline from the distal tubing stops.


Figure 2: Lightning Bolt Aspiration Tubing Control Panel

Indicator	Status	
Power Indicator	Off = No Power	Green = Ready
Vacuum Indicator	Yellow = Insufficient Vacuum	
Saline Indicator	Yellow = Insufficient Saline	
PRIME Button	Blue = Prime	

- 6. Confirm that the Aspiration Catheter is positioned proximal to the thrombus per INDIGO Aspiration Catheter Preparation and Use.
- 7. Remove all devices from the lumen of the Aspiration Catheter.
- 8. Remove the rotating hemostasis valve from the proximal hub of the Aspiration Catheter.
- 9. Ensure the proximal hub of the Aspiration Catheter is filled with fluid.
- 10. Create a wet-to-wet connection by pressing the PRIME button and connecting the Lightning Bolt Aspiration Tubing to the proximal hub of the Aspiration Catheter as shown in **Figure 3**.
NOTE: Connect Lightning Bolt Aspiration Tubing directly to the proximal hub of the Aspiration Catheter.

Figure 3: Aspiration Tubing Connection to Aspiration Catheter


- 11. Press and release the GO button to begin aspiration.
- 12. The condition indicator lights on the Lightning Bolt control will illuminate to indicate the following conditions:
 - a. GO button is Green = Lightning Bolt is powered ON
 - b. Low/No Flow indicator is illuminated = Low/No Flow is detected
 - c. Sustained No Flow indicator is illuminated = Sustained No flow is detected. When Sustained No Flow is detected, Lightning Bolt will provide continuous aspiration.
- 13. To stop aspiration, press the STOP button on the control panel of the Lightning Bolt Aspiration Tubing and turn off the Penumbra Aspiration Pump.
- 14. **Intermittent and Continuous Aspiration:** Lightning Bolt Aspiration Tubing may be set to deliver intermittent and continuous aspiration by turning off modulated aspiration. When Lightning Bolt is ready for use (GO button flashing), press and hold the GO button until the Modulated Aspiration OFF indicator illuminates. To stop aspiration, press the STOP button on the Control Panel of the Lightning Bolt Aspiration Tubing.
- 15. If the red SYSTEM light is observed, stop using the device and replace as needed.

INDIGO SEPARATOR PREPARATION AND USE


- 1. Select an appropriate size INDIGO Separator corresponding to the Aspiration Catheter.
- 2. Advance the INDIGO Aspiration Catheter to gently embed the distal tip into the thrombus.
- 3. Insert the INDIGO Separator through the rotating hemostasis valve and the proximal hub of the Aspiration Catheter. If provided, an introducer may be used.
- 4. Advance the INDIGO Separator through the Aspiration Catheter until the Separator cone is located at the distal tip of the Aspiration Catheter.
- 5. Connect the Aspiration Tubing to the side port of the rotating hemostasis valve as shown in **Figure 4**.

Figure 4: Assembled Aspiration System with INDIGO Separator

- 6. Begin aspiration and advance/retract the INDIGO Separator to assist with aspiration and removal of thrombus as shown in **Figure 5**.

Figure 5: Thrombus Aspiration using the INDIGO Separator

- 7. To stop aspiration, press the STOP button on the control panel of the Lightning Bolt Aspiration Tubing and turn off the Penumbra Aspiration Pump.
- 8. Remove the INDIGO Separator.
- 9. If necessary, additional INDIGO Aspiration Catheters and Separators may be used to further remove thrombus at the discretion of the physician.
- 10. Using a 5cc or larger syringe, aspirate approximately 5cc of blood from the INDIGO Aspiration Catheter to remove any thrombus that may remain in the catheter.
- 11. Obtain a post-treatment angiogram by injecting contrast media through the guide catheter or vascular sheath.

CLINICAL EXPERIENCE

Penumbra EXTRACT-PE Trial Clinical Data Summary:

The Penumbra EXTRACT-PE trial was a prospective, multicenter, single arm trial to determine the safety and efficacy of the Indigo® Aspiration System for mechanical thrombectomy in subjects with acute pulmonary embolism (PE). The trial enrolled 119 subjects at 22 centers in the U.S. An independent imaging core lab and a clinical events committee (CEC) reviewed safety endpoints data.

INCLUSION CRITERIA

- 1. Clinical signs and symptoms consistent with acute PE with duration of 14 days or less. Evidence of PE must be from CTA.
- 2. Systolic BP ≥ 90 mmHg with evidence of dilated RV with an RV/LV ratio > 0.9
- 3. Patient is 18 years of age or older

EXCLUSION CRITERIA

- 1. tPA use within 14 days prior to baseline CTA
- 2. Systolic BP < 90 mmHg for 15 mins or the requirement of inotropic support to maintain systolic BP ≥ 90 mmHg
- 3. Pulmonary hypertension with peak PA > 70 mmHg by right heart catheterization
- 4. History of severe or chronic pulmonary hypertension
- 5. FiO2 requirement $> 40\%$ or > 6 LPM to keep oxygen saturation $> 90\%$
- 6. Hematocrit $< 28\%$
- 7. Platelets $< 100,000/\mu\text{L}$
- 8. Serum creatinine > 1.8 mg/dL

- 9. INR > 3
- 10. aPTT (or PTT) > 50 seconds on no anticoagulation
- 11. History of heparin-induced thrombocytopenia (HIT)
- 12. Contraindication to systemic or therapeutic doses of anticoagulants
- 13. Major trauma < 14 days
- 14. Presence of intracardiac lead
- 15. Cardiovascular or pulmonary surgery within last 7 days
- 16. Cancer requiring active chemotherapy
- 17. Known serious, uncontrolled sensitivity to radiographic agents
- 18. Life expectancy < 90 days
- 19. Female who is pregnant
- 20. Intracardiac Thrombus
- 21. Patients on ECMO
- 22. Current participation in another investigational study

RESULTS:

There were 119 subjects that were enrolled and included in the Intent-To-Treat (ITT) analysis population. The Modified ITT (mITT) population was a subset of the ITT population and was the primary population for all efficacy parameters. The mITT population excluded subjects who received adjunctive treatments or thrombolytics intra-procedure through 48 hours post-procedure for the purpose of reducing clot burden in the pulmonary artery. There were 110 subjects in the mITT population. Safety endpoints are evaluated based on ITT population. There were 110 patients who completed the study 30 day follow-up, 3 patients died, 1 withdrew from the study, and 5 were lost to follow-up.

DEMOGRAPHICS AND MEDICAL HISTORY:

Of the 119 subjects, the median age was 59.8 years old (range 21, 92), 44.5% (53/119) were female, a previous deep vein thrombosis (DVT) was reported in 60.5% (72/119), any cancer in 22.7% (27/119), and a previous pulmonary embolism (PE) was reported in 17.6% (21/119) of subjects.

PROCEDURAL CHARACTERISTICS:

The study procedure was initiated in 119 subjects and completed in 118 subjects; one subject had the procedure aborted prior to aspiration. Conscious sedation was used in 97.5% of subjects and 78.2% had a right femoral access site. The median procedure time (venous puncture to Indigo device removal) was 66 min; IQR [46.0, 94.0]. The median time from the first Indigo device insertion to the last Indigo device removal was 37 minutes; IQR [23.5, 60.0]. The median ICU stay was 1.0 day [1.0, 2.0]. The treatment included the following pulmonary embolism locations:

Treatment Location	All Subjects (N=119)
Main PA with Unilateral	0.8% (1/119)
Main PA with Bilateral	23.5% (28/119)
Unilateral Only	5.9% (7/119)
Bilateral Only	69.7% (83/119)

DEVICE USAGE:

There were 129 Penumbra aspiration catheters and 99 Penumbra separators used during the study procedure.

Penumbra Devices Used – By Subject	All Subjects (N=119)
Penumbra Catheters Used per Subject	
1	91.6% (109/119)
2	8.4% (10/119)
Catheters	
Aspiration Catheter 8, TORQ, 85 cm	13.4% (16/119)
Aspiration Catheter 8, XTORQ, 115 cm	83.2% (99/119)
Aspiration Catheter 8, XTORQ, 115 cm Kit	6.7% (8/119)
Penumbra Separators Used per Subject	
0	19.3% (23/119)
1	78.2% (93/119)
2	2.5% (3/119)
Separators (Separator 8, 150 cm)	80.7% (96/119)

PRIMARY EFFICACY ENDPOINT:

The primary efficacy endpoint was the reduction in RV/LV ratio from baseline to 48 hours assessed by CTA and evaluated by independent core laboratory. The lower limit of the 95% confidence interval (CI) of the change in RV/LV ratio at 48 hours was >0.20. The mean baseline RV/LV ratio was 1.46 (0.29) and it was reduced to 1.04 (0.17) at 48 hours. The endpoint was met with absolute reduction in RV/LV ratio of 0.42±0.25 (95%CI 0.37, 0.46), representing a 26.9% reduction (p<0.0001).

Primary Effectiveness Endpoint	mITT (N=110)	95% CI ^[3]
Absolute Reduction in RV/LV Ratio from Baseline to 48 Hours ^[1]		
N ^[2]	106	
Mean (SD)	0.42 (0.25)	(0.37, 0.46)
Median [IQR]	0.41 [0.25, 0.54]	
Range (Min, Max)	-0.19, 1.23	

^[1] p <0.0001 based on a one-sample, two-sided asymptotic test with a null proportion of 0.20 and the standard error computed from the sample proportion.
^[2] Data were analyzed based on mITT with evaluable 48 hours CTA
^[3] Confidence intervals have not been adjusted for multiplicity

SAFETY:

The primary safety endpoint was major adverse events, a composite of device- related death, major bleeding, and device-related SAEs (defined as clinical deterioration, pulmonary vascular injury, and cardiac Injury) within 48 hours. The safety hypothesis stated that the 48-hour rate of major adverse events would not equal 40%. The primary safety endpoint was met and the rate was 1.7% (95% CI 0.0%, 4.0%), p<0.0001). Two subjects (1.7%) experienced 3 events: 1 subject experienced Sustained Ventricular Tachycardia (Device-related Death/ Major Bleeding/Clinical Deterioration/Pulmonary Vascular Injury) and Hemoptysis (Major Bleeding/Pulmonary Vascular Injury); 1 subject experienced a groin access site bleeding (Major Bleeding).

PRIMARY SAFETY ENDPOINT (AS-TREATED POPULATION):

Primary Safety Analysis	All Subjects (N=119)	95% CI ^[2]
Composite of major adverse events within 48 Hours ^[1]	1.7% (2/119)	(0.0%, 4.0%)

^[1] p <0.0001 based on a one-sample, two-sided asymptotic test with a null proportion of 0.40 and the standard error computed from the sample proportion.
^[2] Confidence intervals have not been adjusted for multiplicity

SECONDARY ENDPOINTS (AS-TREATED POPULATION):

Secondary safety endpoints included device-related death within 48 hours, major bleeding within 48 hours, clinical deterioration within 48 hours, pulmonary vascular injury within 48 hours, cardiac injury within 48 hours, any-cause mortality within 30 days, device-related SAEs within 30 days, and symptomatic PE recurrence within 30 days. Serious Adverse Events (SAEs) and Adverse events (AEs) were collected through the 30-day follow-up. The key events are included below.

A total of three deaths occurred during the trial. One death was adjudicated as procedure and device related. Two additional non-procedure/device related deaths occurred due to progression of pre-existing diseases (cancer progression and a recent brain resection due to a previous ischemic stroke). A total of 7 subjects experienced 8 secondary endpoint events.

SECONDARY SAFETY ENDPOINTS:

Secondary Safety Analysis*	All Subjects (N=119)
Device-related Death within 48 Hours	0.8% (1/119)
Major Bleeding within 48 Hours	1.7% (2/119)
Clinical Deterioration within 48 Hours	1.7% (2/119)
Pulmonary Vascular Injury within 48 Hours	1.7% (2/119)
Cardiac Injury within 48 Hours	0.0%
Any-cause Mortality within 30 Days	2.5% (3/119)
Device-related SAEs within 30 Days	1.7% (2/119)
Symptomatic PE Recurrence within 30 Days	0.0%
Device and procedure related Adverse Events within 30 days from procedure*	
Procedure-related SAEs	2.5% (3/119)
Device-related SAEs	1.7% (2/119)
Procedure-related Non-Serious AEs	9.2% (11/119)
Device-related Non-Serious AEs	4.2% (5/119)

*One subject can meet several endpoint events and can appear in multiple categories

CLINICAL BENEFITS AND PERFORMANCE CHARACTERISTICS

Clinical benefits of the INDIGO Aspiration System are the removal or reduction of thrombi and emboli from the arterial or venous peripheral vasculature or for patients with pulmonary embolism. The detailed listing of all known benefits is included in the Summary of Safety and Clinical Performance (SSCP).

SAFETY AND CLINICAL PERFORMANCE

A copy of the Summary of Safety and Clinical Performance can be viewed by searching the device name on the EUDAMED website: <https://ec.europa.eu/tools/eudamed>
Basic UDI DI: 081454801INDIGOSYSTEMYF

LIGHTNING BOLT ASPIRATION TUBING OPERATING CONDITIONS

- Temperature: 65°F - 75°F (18°C - 24°C)
- Humidity: <75% RH
- Pressure: Sea Level – 6000 ft (0 m – 1828 m)

LIGHTNING BOLT ASPIRATION TUBING TRANSPORTATION CONDITIONS

- Temperature: 14°F - 131°F (-10°C - 55°C)
- Humidity: <85% RH

STORAGE

- Store in a cool, dry place.

DISPOSAL PROCEDURES

Dispose of devices using standard hospital biohazard procedures.

APPROVALS

- Conforms to
- AAMI Std. ES60601-1
 - IEC Std. 60601-1-6
 - ISO Std. 10079-1
 - Certified to CSA Std. C22.2# 60601-1

Lightning Bolt Aspiration Tubing complies with all current international safety and EMC requirements for medical electrical equipment as follows:

- IEC 60601-1-2 including:
 - CISPR 11 Class A
 - EN/IEC 61000-3-2
 - EN/IEC 61000-3-3
 - EN/IEC 61000-4-2
 - EN/IEC 61000-4-3
 - EN/IEC 61000-4-4
 - EN/IEC 61000-4-5

- EN/IEC 61000-4-6
- EN/IEC 61000-4-8
- EN/IEC 61000-4-11
- EN/IEC 61000-4-39
- AIM 7351731

Lightning Bolt Aspiration Tubing Electromagnetic Compatibility

Guidance and manufacturer's declaration – electromagnetic emissions			
Lightning Bolt Aspiration Tubing is intended for use in the electromagnetic environment specified below. The customer or the user of Lightning Bolt Aspiration Tubing should assure that it is used in such an environment.			
Emissions Test		Compliance	Electromagnetic Environment - Guidance
RF emissions CISPR 11		Group 1	Lightning Bolt Aspiration Tubing creates RF energy as a byproduct of its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11		Class A	Lightning Bolt Aspiration Tubing is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2		Class A	
Voltage Fluctuations/ Flicker emissions EN 61000-3-3		Complies	
Guidance and manufacturer's declaration – electromagnetic immunity			
Lightning Bolt Aspiration Tubing is intended for use in the electromagnetic environment specified below. The customer or the user of Lightning Bolt Aspiration Tubing should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/ output lines	±2 kV for power supply lines ±1 kV for input/ output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Voltage Dips 30% reduction, 25/30 periods at 0°	Voltage Dips 30% reduction, 25/30 periods at 0°	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Lightning Bolt Aspiration Tubing requires continued operation during power mains interruptions, it is recommended that Lightning Bolt Aspiration Tubing be powered from an uninterruptible power supply or a battery.
	Voltage Dips > 95% reduction, 0.5 period at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	Voltage Dips > 95% reduction, 0.5 period at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	
	Voltage Dips > 95% reduction, 1 period at 0°	Voltage Dips > 95% reduction, 1 period at 0°	
	Voltage Interruptions > 95% reduction, 250/300 periods	Voltage Interruptions > 95% reduction, 250/300 periods	

Guidance and manufacturer's declaration – electromagnetic immunity			
Lightning Bolt Aspiration Tubing is intended for use in the electromagnetic environment specified below. The customer or the user of Lightning Bolt Aspiration Tubing should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
(50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz (6 Vrms in ISM radio Bands within 150 kHz – 80 MHz 80% AM at 1 kHz)	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of Lightning Bolt Aspiration Tubing, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	3 V/m	Recommended separation distance $d = 1.2\sqrt{P}$ $d = 0.35\sqrt{P}$ 80 MHz to 800 MHz $d = 0.70\sqrt{P}$ 800 MHz to 2.7 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which Lightning Bolt Aspiration Tubing is used exceeds the applicable RF compliance level above, Lightning Bolt Aspiration Tubing should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating Lightning Bolt Aspiration Tubing.			
^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Recommended separation distances between portable and mobile RF communications equipment and Lightning Bolt Aspiration Tubing				
Lightning Bolt Aspiration Tubing is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of Lightning Bolt Aspiration Tubing can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and Lightning Bolt Aspiration Tubing as recommended below, according to the maximum output power of the communications equipment.				
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter (m)			
	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.7 GHz $d = 2.3\sqrt{P}$	710, 745, 780, 5240, 5500, 5785 MHz $d = 0.7\sqrt{P}$	385, 450, 810, 870, 930, 1720, 1845, 1970, 2450 MHz $d = 0.2\sqrt{P}$
0.01	0.12	0.23	0.07	0.02
0.1	0.38	0.74	0.21	0.07
1	1.2	2.3	0.67	0.21
10	3.8	7.4	2.1	0.70
100	12	23	6.7	2.1
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.				
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.				
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.				

Immunity to RF Wireless Communications Equipment						
Test Frequency (MHz)	Band ^{a)} (MHz)	Service ^{a)}	Modulation ^{b)}	Maximum Power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
385	380 – 390	TETRA 400	Pulse modulation ^{b)} 18 Hz	1.8	0.3	27
450	430 – 470	GMRS 460, FRS 460	FM ^{c)} ±5 kHz deviation 1 kHz sine	2	0.3	28
710	704 – 787	LTE Band 13, 17	Pulse modulation ^{b)} 217 Hz	0.2	0.3	9
745						
780						
810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation ^{b)} 18 Hz	2	0.3	28
870						
930						
1720	1700 – 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation ^{b)} 217 Hz	2	0.3	28
1845						
1970						

Immunity to RF Wireless Communications Equipment						
Test Frequency (MHz)	Band ^{a)} (MHz)	Service ^{a)}	Modulation ^{b)}	Maximum Power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
2450	2400 – 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^{b)} 217 Hz	2	0.3	28
5240	5100 – 5800	WLAN 802.11 a/n	Pulse modulation ^{b)} 217 Hz	0.2	0.3	9
5500						
5785						
Immunity to Proximity Magnetic Fields						
Test Frequency		Modulation		Immunity Test Level (A/m)		
30 kHz		CW		8		
134.2 kHz		Pulse modulation ^{b)} 2.1 kHz		65		
13.56 MHz		Pulse modulation ^{b)} 50 kHz		7.5		
^{a)} For some services, only the uplink frequencies are included.						
^{b)} The carrier shall be modulated using a 50% duty cycle square wave signal.						
^{c)} As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.						

EQUIPMENT CLASSIFICATIONS

Protection from electric shock Class II and externally powered
Degree of protection against electric shock..... Type BF Applied Part
Degree of protection against ingress of liquids or particles..... IP22
Mode of Operation..... Continuous Operation
Flow Rate 0-18 LPM

SYMBOL DEFINITIONS

	Manufacturer		European Union Authorized Representative
	Sterilized using ethylene oxide. Single sterile barrier system with protective packaging outside.		Catalog Number
	Lot number		Medical Device
	Type BF Applied Part		Use By
	Do not reuse		Refer to Instruction Manual / Booklet
	Do not resterilize		Prescription Only – US Federal Law restricts this device to use by or on the order of a physician
	Not made with natural rubber latex		Nonpyrogenic
	Do not use if package is damaged		Attention, see Instructions for Use
	High Vacuum/ Low Flow (Intermittent suction)		US and Canada Certification
	Date of Manufacture		

WARRANTY

Penumbra, Inc. (Penumbra) warrants that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular use. Handling, storage, cleaning, and sterilization of this device as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond Penumbra's control directly affect the device and the results obtained from its use. Penumbra's obligation under this warranty is limited to the repair or replacement of this device and Penumbra shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this device. Penumbra neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. Penumbra assumes no liability with respect to devices reused, reprocessed, or re-sterilized and makes no warranties, expressed or implied, including but not limited to merchantability or fitness for intended use, with respect to such devices.

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