INTENDED PURPOSE

The MIDWAY Delivery Catheters are designed to provide a conduit for introduction of interventional devices to the peripheral and neurovascular anatomy.

DEVICE DESCRIPTION

The MIDWAY Delivery Catheters include the MIDWAY 43 Delivery Catheter and MIDWAY 62 Delivery Catheter. The MIDWAY Delivery Catheters can be used with a guidewire to access the desired anatomy. The MIDWAY Delivery Catheters are provided with a rotating hemostasis valve and peelable introducer.

MIDWAY Delivery Catheter:

The MIDWAY Delivery Catheter is a single lumen, variable stiffness catheter with a radiopaque markerband on the distal end and a Luer hub on the proximal end. The MIDWAY Delivery Catheter dimensions are included on the individual device label. The MIDWAY Delivery Catheter is compatible with introducer sheaths and guide catheters appropriately sized for the outer diameter of the MIDWAY Delivery Catheter. The MIDWAY Delivery Catheter has a hydrophilic coating on the distal segment of the catheter shaft. The hydrophilic coating length per catheter is noted in Table 1.

Table 1 – Hydrophilic Coating	
Catheter	Hydrophilic Coating Length (cm)
MIDWAY 43	60
MIDWAY 62	30

INDICATIONS FOR USE

The MIDWAY Delivery Catheter is indicated for the introduction of interventional devices into the peripheral and neuro vasculature.

INTENDED PATIENT POPULATION

The intended patient population is patients who require neuro and/or peripheral intervention.

CONTRAINDICATIONS

There are no known contraindications.

WARNINGS

- The MIDWAY Delivery Catheter should only be used by physicians who have received appropriate training in interventional techniques.
- The risks and benefits for radial access against a transfemoral approach should be carefully weighed and considered for each patient.
- The device is intended for single use only. Do not resterilize or reuse. Resterilization and/or reuse may compromise the structural integrity of the device and/or lead to device failure which in turn, may result in patient injury, illness, or death.
- Do not advance or withdraw the MIDWAY Delivery Catheter against resistance without careful assessment of the cause using fluoroscopy. If the cause cannot be determined, withdraw the device. Unrestrained moving or torquing the device against resistance may result in damage to the vessel or device.
- Do not use automated high-pressure contrast injection equipment with the MIDWAY Delivery Catheter because it may damage the device.
- Failure to abide by the warnings in this labeling might result in damage to the device coating, which may necessitate intervention or result in serious adverse events.

PRECAUTIONS

- Do not use kinked or damaged devices. Do not use open or damaged packages. Return all damaged devices and packaging to the manufacturer / distributor.
- Use prior to the “Use By” date.
- Use the MIDWAY Delivery Catheter in conjunction with fluoroscopic visualization.
- Maintain a constant infusion of an appropriate flush solution.
- If flow through the device becomes restricted, do not attempt to clear the lumen by infusion. Remove and replace the device.
- Prior to beginning radial artery access, conduct screening, such as an Allen test, to ensure that radial access is appropriate for the patient.
- As in all fluoroscopy procedures, consider all necessary precautions to limit patient radiation exposure by using sufficient shielding, reducing fluoroscopy times and modifying radiation technical factors whenever possible.
- Avoid using alcohol, antiseptic solutions, or other solvents to pre-treat the device because this may cause unpredictable changes in the coating which could affect the device safety and performance.

POTENTIAL ADVERSE EVENTS*

Possible complications include, but are not limited to, the following:

- Acute vessel occlusion
- Access site complications such as hematoma, hemorrhage, spasm, occlusion, or compartment syndrome
- Air embolism
- Allergic reaction and anaphylaxis from contrast media or device material
- Arteriovenous fistula
- Death
- Emboli including foreign body embolization
- Inflammation
- Infection
- Intracranial hemorrhage
- Ischemia
- Nerve injury
- Neurological deficits including stroke
- Pseudoaneurysm
- Renal impairment or acute renal failure from contrast media
- Radiation exposure that may lead to cataracts, skin reddening, burns, alopecia, or neoplasia from x-ray exposure
- Vessel dissection
- Vessel perforation
- Vessel spasm
- Vessel thrombosis

*If a serious incident related to the device occurs, contact your Penumbra, Inc. representative and the competent authority in your respective country/region.

Table 2 – Device Compatibility and Selection Information		
Delivery Catheter	Introducer Sheath / Guide Catheter	Guidewire
MIDWAY 43	6F / 0.071" (1.80 mm)	0.038" (0.97 mm)
MIDWAY 62	6F / 0.081" (2.06 mm)	

DEVICE PREPARATION AND USE

1. Refer to **Warnings, Precautions, and Potential Adverse Events** prior to use.
2. Refer to **Table 2** for device compatibility information.
3. Select an appropriately sized MIDWAY Delivery Catheter based on the anatomy and length.
4. Gently remove the MIDWAY Delivery Catheter and packaging hoop from the pouch by grasping the MIDWAY Delivery Catheter hub and packaging hoop and slowly pulling them out of the pouch.
5. Prepare for use by flushing the packaging hoop and device with heparinized saline.
6. Remove the MIDWAY Delivery Catheter from the packaging hoop and inspect the MIDWAY Delivery Catheter for kinks or other damage.
7. Remove the Rotating Hemostasis Valve (RHV) from the card tabs.
8. Connect the RHV to the hub of the MIDWAY Delivery Catheter. Flush the lumen of the MIDWAY Delivery Catheter with heparinized saline.
9. Insert a guidewire into the MIDWAY Delivery Catheter and advance the guidewire until the distal tip of the wire is at the distal tip of the MIDWAY Delivery Catheter.
10. Place an introducer sheath or guide catheter appropriately sized for the MIDWAY Delivery Catheter's outer diameter.
11. Advance the MIDWAY Delivery Catheter and guidewire simultaneously into the introducer sheath or guide catheter and extend the guidewire 5 cm to 6 cm distal to the tip of the MIDWAY Delivery Catheter.
12. Advance the MIDWAY Delivery Catheter and guidewire system, and then select the appropriate vessel with the guidewire.
13. After gaining access to the desired vessel with the guidewire, advance the MIDWAY Delivery Catheter over the guidewire into the desired vessel.
14. Under fluoroscopy, advance the MIDWAY Delivery Catheter until the desired position is achieved.
15. Advance the MIDWAY Delivery Catheter and guidewire to the vascular site and remove the guidewire.

CLINICAL BENEFITS AND PERFORMANCE CHARACTERISTICS

The MIDWAY Delivery Catheter aids in the use of interventional devices to the peripheral and neurovascular anatomy. The detailed listing of all known benefits is included in the Summary of Safety and Clinical Performance (SSCP).


SAFETY AND CLINICAL PERFORMANCE

A copy of the Summary of Safety and Clinical Performance can be viewed by searching the device name on the EUDAMED website: <https://ec.europa.eu/tools/eudamed>
Basic UDI DI: 081454801ACCESSY8

STORAGE AND DISPOSAL

- Store in a cool, dry place.
- Dispose of device in accordance with hospital, administrative, and /or local government policy.

SYMBOLS GLOSSARY

	Attention, see Instructions for Use		
	Prescription only – US Federal Law restricts this device to use by or on the order of a physician		
	Nonpyrogenic		
	Sterilized using ethylene oxide. Single sterile barrier system with protective packaging outside.		
	Do not use if package is damaged.		
	Medical Device		Date of Manufacture
	Manufacturer		EU Representative
	Lot number		Do not reuse
	Use by		Catalogue number
	Do not resterilize		Not made with natural rubber latex

WARRANTY

Penumbra Inc. (Penumbra) warrants that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular use. Handling, storage, cleaning, and sterilization of this device as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond Penumbra’s control directly affect the device and the results obtained from its use. Penumbra’s obligation under this warranty is limited to the repair or replacement of this device and Penumbra shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this device. Penumbra neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. Penumbra assumes no liability with respect to devices reused, reprocessed, or resterilized and makes no warranties, expressed or implied, including but not limited to merchantability or fitness for intended use, with respect to such devices.

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