Esperance **pHLO** Aspiration Catheter

DIRECTIONS FOR USE

ENGLISH

DIRECTIONS FOR USE

Esperance pHLO Aspiration Catheter

phenox Limited Kamrick Court

Ballybrit Business Park,
Galway, Ireland H91 XY38.

Phone: +353 91 740100

Email: info@phenox.ie

Website: www.phenox.ie

Contents

1 x Esperance pHLO Aspiration Catheter

1 x introducer sheath

The Product Family Includes the Following:

Esperance pHLO Aspiration Catheter 125cm and 135cm

Device Description

The Esperance pHLO Aspiration Catheter (Esperance pHLO AC) is a single lumen catheter. The flexibility decreases from distal to proximal, ranging from highly flexible at the distal part to semi-rigid at the proximal end. The Esperance pHLO AC has a radiopaque marker on the distal end to facilitate fluoroscopic visualization. A luer hub is on the proximal end to facilitate attachment of accessories, the passage of devices into the inner lumen of the catheter, and connection to ancillary devices such as a vacuum source. The distal 100cm of the catheter external surface has a hydrophilic coating applied for increased lubricity during use. The Esperance pHLO AC has an atraumatic straight tip.

The Esperance pHLO AC is supplied sterile in a packing hoop, within a sealed pouch and shelf carton configuration. An introducer sheath is supplied on a backing card inside the pouch. The packaging is designed to facilitate ease of handling and aseptic presentation of the device.


Figure 1: Esperance pHLO AC

Table 1: Dimensions of Esperance pHLO AC

Dimension		Esperance pHLO AC 71 125	Esperance pHLO AC 71 135
A	Effective length	125 cm	135 cm
B	Total length	132 cm	142 cm
C	Hydrophilic coating length	100 cm	100 cm
D	Outer diameter	0.082 inch / 6F / 2.1mm	0.082 inch / 6F / 2.1mm
E	Inner diameter	0.071 inch / 1.8 mm	0.071 inch / 1.8 mm
	Distance of distal tip to radiopaque marker	0.6 mm	0.6 mm

Intended Purpose

Esperance pHLO AC is a catheter intended for aspiration in neurovascular vessels as a standalone device or in conjunction with a compatible retrieval device, combining aspiration and mechanical thrombectomy. It can also be used to facilitate the introduction of other suitable devices as well as diagnostic and therapeutic agents subject to use by their IFUs.

Indications for Use

Esperance pHLO AC is indicated for use in the revascularization of patients with acute ischemic stroke secondary to large intracranial vessel occlusive disease (within the internal carotid, middle cerebral-M1 and -M2 segments, basilar and vertebral arteries). Esperance pHLO AC serves also for the introduction of suitable devices (e.g., microcatheter) as well as the administration of diagnostic and therapeutic agents subject to use by their IFUs.

Contraindications

Esperance pHLO AC is not intended for use in peripheral or coronary vasculature.

Esperance pHLO AC is contraindicated for use in patients with known generalized vessel wall weakness.

Patient Target Group

The Esperance pHLO AC is intended for patients with acute ischemic stroke secondary to large intracranial vessel occlusive disease (within the internal carotid, middle cerebral-M1 and -M2

segments, basilar and vertebral arteries) and/or undergoing diagnostic and/or therapeutic procedures in the neurovasculature.

Warnings

- Esperance pHLO AC may only be used by physicians thoroughly trained in percutaneous, intravascular techniques and procedures. Esperance pHLO AC may only be used in a clinical environment capable of supporting neurovascular procedures.
- The safety and effectiveness of the device has not been established for paediatric and neonatal use.
- The Esperance pHLO AC must only be manipulated under fluoroscopy inside the body. Limit the exposure to X-ray radiation doses to patients and physicians / other users by using sufficient shielding, reducing fluoroscopy times, and modifying X-ray technical factors when possible.
- Do not advance or retract the catheter if resistance is felt. Pause the procedure to investigate. Determine the cause of the resistance and take appropriate action. If the cause cannot be determined, withdraw the Esperance pHLO AC. Excessive force against resistance may result in damage of the catheter, damage or dislodgment of other interventional devices, or vessel injury.
- The safety and effectiveness of the device has not been established or is unknown in vascular regions other than those specifically indicated.
- Esperance pHLO AC is coated with a hydrophilic coating at the distal end. Hydrate the catheter for at least 30 seconds before removing from the hoop. Do not allow the catheter to dry once hydrated.
- Do not use if packaging is open or damaged, as sterility cannot be assumed otherwise.
- Do not use automated high-pressure contrast injection equipment with the Esperance pHLO AC because it may damage the device.
- Do not re-use. Esperance pHLO AC is supplied sterile and intended for single use only. Discard the catheter after one procedure. Structural integrity and/or function may be impaired through re-use or cleaning. Catheters are extremely difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused. Re-sterilization or Re-Use may result in ineffective catheter coating lubrication, which may result in high friction and the inability to access the target neurovasculature.

Precautions

- Exercise care in handling the flexible distal portion of the Esperance pHLO AC to reduce the chance of accidental damage.
- Only use the device if DFU and label can be read clearly and can be fully understood and if the label properly adheres to the packaging.
- Inspect catheter prior to use to verify that its condition is suitable for the specific procedure. Inspect for any damage, surface irregularities, bends or kinks. Do not use kinked, damaged, or irregular devices.
- Only insert the catheter into an appropriately sized vessel.
- Sufficient irrigation of the catheter is necessary to remove enclosed air bubbles.
- Use prior to the "Use By" date.
- Only select compatible guidewires, guide catheters and devices or agents to be used with the Esperance pHLO AC. For additional precautions associated with the use of these devices / agents, consult the corresponding DFU.
- Do not expose the catheter to organic solvents, such as disinfection agents, alcohol or DMSO.
- Do not wipe the device.
- Do not excessively torque the device.
- The introducer sheath is not intended for use in the human body.
- If there is suspicion of vasospasm in the affected vascular region, all measures required for reversion should be taken, e.g., with appropriate drug therapy.
- Appropriate premedication with antiplatelet agents and systemic heparinization is to be considered according to the type and duration of the intervention.
- In case of any unusual catheter behaviour, stop the procedure and replace the Esperance pHLO AC.
- Care must be taken to avoid complete blood flow blockage by the catheter.
- If flow through the catheter is restricted, do not attempt to clear the lumen by infusion. The catheter must be removed and replaced.
- Do not wedge the tip of Esperance pHLO AC in the vessel, as this may result in perforation or dissection of the vessel wall. Excessive aspiration with the distal tip of the Esperance pHLO AC covered by the vessel wall may cause vessel injury. Carefully investigate location of the distal tip under fluoroscopy prior to aspiration.
- Maintain a constant infusion of appropriate flush solution.
- Use caution when manipulating, advancing and/or withdrawing Esperance pHLO AC through needles, metal cannulas, stents, or other devices with sharp edges, or through tortuous or calcified blood vessels. Manipulation, advancement, and/or withdrawal past sharp or beveled edges may result in damage and/or separation of the outer coating and/or flexible distal portion, which may lead to clinical adverse events, resulting in coating material remaining in the vasculature or device damage.
- Ensure aspiration is performed for only the minimum time necessary.
- Limit the usage of Esperance pHLO AC to arteries larger than the catheters outer diameter.
- Return all damaged packaging to the manufacturer or distributor.
- Maximum allowable vacuum is -98 kPa (-14.2 psi / -29 inHg) with reference to the atmospheric pressure at sea level.
- For additional precautions associated with the use of a syringe, an aspiration tubing set and an aspiration pump, consult DFUs of the corresponding devices.

Storage

Keep away from heat. Store in a cool, dry place.

Disposal

The device, retrieved clots, blood and if applicable its packaging has to be disposed as biohazard waste.

Potential Adverse Events

The following and further adverse events may occur during or despite the use of Esperance pHLO AC:

- Air-, catheter-, plaque-, thrombus- emboli
- Aneurysm rupture / perforation
- Allergic or toxic reaction
- Adverse response by patient (other than above)
- Emboli (thromboembolic/device components)
- Extracranial Hemorrhage / Hematoma
- Extracranial / intracranial complications during procedure (other than above)
- Harm of 3rd person via contaminated products / infection
- Hemorrhagic stroke
- Infection
- Inflammation
- Ischemic stroke
- Neurovascular / Neurological complications
- Patient Injury
- Severe harm (permanent disability / Death)
- Transient ischemic attack (TIA) / seizure
- Vasospasm
- Vessel occlusion/thrombosis
- Vessel perforation
- Vessel dissection
- Vessel damage (other than above)

Compatibility

Agents and devices used with the Esperance pHLO AC must be compatible with the catheter's dimensions according to the manufacturer's description.

Esperance pHLO AC is not compatible with glue, glue mixture or non-adhesive liquid embolic agents

Esperance pHLO AC is compatible with the following products (not supplied - selected based on compatibility and physician's experience and preference):

- Guidewires with a maximum outer diameter of 0.035 inch / 0.89 mm
- 8F guide catheters/introducer sheaths (ID min: 0.086")
- 9F balloon guide catheters (ID min: 0.086")
- Non-liquid interventional devices which are declared compatible to a catheter with an inner diameter of 0.071 inch according to the respective DFU
- Thrombectomy devices with their microcatheters, which are declared compatible by their respective DFU for withdrawal back into an aspiring catheter (ID: 0.071 inch) during the recanalization procedure
- Aspiration pumps with maximum vacuum of -98 kPa (-14.2 psi / -29 inHg) with reference to the atmospheric pressure at sea level.
- Tubing sets with a standard Luer-connector at one end, to connect the Esperance pHLO AC with an aspiration pump
- Syringes compatible with standard Luer-connector
- Rotating hemostatic valves (RHV) with standard Luer adapter

Esperance pHLO AC is supplied with an Introducer Sheath which facilitates the introduction of the Esperance pHLO AC device into a guide catheter.

Directions for Use

Preparation of procedure

1. Refer to Warnings, Precautions and Potential Adverse Event's prior to use.
2. Introduce a suitable catheter sheath (introducer sheath) into an appropriate vessel according to standard practice / the respective DFU.
3. Prepare an 8F guide catheter (or other suitable catheter) for use per its DFU. Connect the RHV and continuous pressure infusion to the guide catheter.
4. Insert the guide catheter into an appropriate vessel by the help of a suitable guidewire according to standard practice and their DFUs.
5. Alternatively, in case a long sheath is used, the procedure may be conducted without an additional guide catheter.

Preparation of Esperance pHLO AC

6. Confirm the appropriate length Esperance pHLO AC is selected.
7. Carefully open Esperance pHLO AC packaging. Peel open pouch from the chevron end and

remove hoop from pouch.

8. Flush the dispenser hoop with sterile saline through the Luer fitting connected to the end of the dispenser hoop before removing the aspiration catheter. Hydrate for at least 30 seconds.
9. Carefully withdraw the aspiration catheter from the dispenser hoop.
10. Inspect the aspiration catheter thoroughly to ensure that it is not damaged. If the aspiration catheter is damaged, replace the device.

Introduction of Esperance pHLO AC

11. Connect the RHV to the aspiration catheter.
12. Flush the inner lumen of the aspiration catheter with sterile saline.
Warning: Do not allow the catheter to dry once hydrated.
13. Connect the continuous pressure infusion line to the Esperance pHLO AC.
14. Prepare the guidewire according to the respective DFU. Note: for navigating the aspiration catheter only a larger diameter guidewire should be used (0.035").
15. Carefully insert the guidewire into the Esperance pHLO AC.
16. Introduce the Esperance pHLO AC - guidewire assembly through the guide catheter's RHV using the peel-away introducer sheath.
17. Remove the peel-away introducer sheath from the Esperance pHLO AC -guidewire assembly by retracting the introducer sheath from the RHV and peeling it off.
Warning: Removing the peel-away introducer sheath from the Esperance pHLO AC without a guidewire inserted may cause damage to the Esperance pHLO AC.
18. Tighten the O-ring valve around the Esperance pHLO AC just enough to minimize back flow, but not so tight as to inhibit Esperance pHLO AC movement.
Warning: Excessive tightening of the hemostatic valve around the Esperance pHLO AC may lead to damage of the device.
19. Maintain continuous flush in the Esperance pHLO AC and guide catheter.
20. Advance the guidewire and Esperance pHLO AC to the distal tip of the guide catheter.
21. Carefully advance the guidewire and Esperance pHLO AC to the intended vessel under fluoroscopy. Advance the guidewire and follow with the Esperance pHLO AC in an alternate fashion.
Alternatively, Esperance pHLO AC can be used in conjunction with a suitable microcatheter for vessel navigation. If a microcatheter is used, prepare it per its DFU prior to use. Follow the microcatheter-guidewire assembly with Esperance pHLO AC.
Warning: Do not advance or withdraw the device when excessive resistance is met until the cause of resistance is determined.
22. Remove the guidewire (and the microcatheter if required) from the Esperance pHLO AC after the target vessel has been reached.

Steps 23 to 31: Esperance pHLO AC for direct thrombus aspiration by syringe

Steps 32 to 40: Esperance pHLO AC for direct thrombus aspiration by aspiration pump

Steps 41 to 51: Esperance pHLO AC for thrombus aspiration by syringe or aspiration pump facilitated by thrombectomy device

Steps 52 to 54: Esperance pHLO AC for introduction of suitable devices

Esperance pHLO AC for direct thrombus aspiration by syringe

23. Make sure the distal tip of the Esperance pHLO AC is engaged with the emboli or thrombi under fluoroscopy.
24. At the proximal end of the Esperance pHLO AC, remove the perfusion line (of the continuous pressure infusion) and attach a stopcock (valve) with 30cc or 60cc locking syringe to the sidearm of the RHV.
25. Tighten the RHV.
26. With the stopcock closed, pull and lock the plunger of the syringe.
27. To begin aspiration, open the stopcock.
28. If the syringe does not aspirate any blood, thrombi, or emboli, or slow aspiration is observed, carefully investigate the cause of the restriction and reposition the distal tip if necessary.
Warning: Excessive aspiration with the distal tip of the Esperance pHLO AC covered by the vessel wall may cause vessel injury
If blood, thrombi or emboli is stuck in the Esperance pHLO AC remove the entire device out of the patient body. Maintain aspiration to make sure the thrombi or emboli remains fully engaged with the distal tip of the Esperance pHLO AC. With the thrombi or emboli fully engaged pull back the Esperance pHLO AC out of the patient body and clear the inner lumen of Esperance pHLO AC.
Warning: Do not attempt to clear inner lumen of the Esperance pHLO AC by infusion while keeping the device in the patient body
29. After aspiration is completed, remove the Esperance pHLO AC from the patient body.
30. Perform an angiographic run to confirm vessel recanalization.
31. If re-access to the vasculatures with the same Esperance pHLO AC device is needed, flush and clean the inner lumen of the Esperance pHLO AC by infusion of saline. Inspect the device for any damage.
Warning: Do not use the device if any damage or irregularities are observed.
Follow steps 11 through 31 for the re-use of Esperance pHLO AC with syringe.

Esperance pHLO AC for direct thrombus aspiration by aspiration pump

32. Make sure the distal tip of the Esperance pHLO AC is engaged with the emboli or thrombi under fluoroscopy.
33. Connect the aspiration tubing set Luer connector (distal) to the sidearm of the RHV.
34. Connect the funnel adapter (proximal) to a compatible aspiration pump. This is achieved by connecting the funnel adapter to the patient connecting port of the canister to the corresponding aspiration pump. Maximum allowable vacuum is -98 kPa (-14.2 psi / -29 inHg) with reference to the atmospheric pressure at sea level.
35. Consult the DFU of the vacuum source for operating instructions of said device.
36. Start aspiration by opening the flow switch / stopcock / clamp or any other type of lock on the aspiration tubing set.
37. If the aspiration tubing set does not aspirate any blood, thrombi, or emboli, or slow aspiration is observed, carefully investigate the cause of the restriction and reposition the distal tip if necessary.

Warning: Excessive aspiration with the distal tip of the Esperance pHLO AC covered by the vessel wall may cause vessel injury

If blood, thrombi or emboli is stuck in the Esperance pHLO AC remove the entire device out of the patient body. Maintain aspiration to make sure the thrombi or emboli remains fully engaged with the distal tip of the Esperance pHLO AC. With the thrombi or emboli fully engaged pull back the Esperance pHLO AC out of the patient body and clear the inner lumen of Esperance pHLO AC.

Warning: Do not attempt to clear inner lumen of the Esperance pHLO AC by infusion while keeping the device in the patient body

38. After aspiration is completed, remove the Esperance pHLO AC from the patient body. To stop aspiration, close the flow switch / stopcock / clamp or any other type of lock on the aspiration tubing set.
39. Perform an angiographic run to confirm vessel recanalization.
40. If re-access to the vasculatures with the same Esperance pHLO AC device is needed, flush and clean the inner lumen of the Esperance pHLO AC by infusion of saline. Inspect the device for any damage.

Warning: Do not use the device if any damage or irregularities are observed.

Follow steps 11 through 22 and 32 through 40 for the re-use of Esperance pHLO AC with aspiration pump.

Esperance pHLO AC for thrombus aspiration by syringe or aspiration pump facilitated by thrombectomy device

41. Make sure the distal tip of Esperance pHLO AC is located in the selected vessel under fluoroscopy.
42. Prepare, introduce and navigate (if required) a suitable microcatheter with its guidewire into the target vessel according their DFUs.
43. Navigate microcatheter and guidewire through the thrombus. Make sure the microcatheter tip is located distally to the thrombus. Remove the guidewire.
44. Select and prepare a compatible thrombectomy device and introduce this device into the microcatheter following their DFUs.
45. Deploy the thrombectomy device inside the thrombus by withdrawing the microcatheter according to its DFU.
46. Navigate the Esperance pHLO AC to engage with the proximal end of the thrombus.
Thrombus aspiration by syringe: Follow steps 24 through 27
Thrombus aspiration by aspiration pump: Follow steps 33 through 36
47. Slowly withdraw the thrombectomy device and the microcatheter as a unit (according their DFUs) under continuous aspiration to / into the Esperance pHLO AC. Or, alternatively, to perform empty catheter technique, remove the microcatheter before beginning aspiration to increase the lumen available for aspiration.
Warning: Never push the thrombectomy device out of the Esperance pHLO AC distally into the vessel.
48. After the thrombectomy device - microcatheter unit is completely withdrawn and removed from the Esperance pHLO AC as well as aspiration is completed, remove the Esperance pHLO AC from the patient body.
Alternatively, all devices (Esperance pHLO AC & thrombectomy device & microcatheter) can be withdrawn as one unit and removed from the patient body.
49. If the thrombectomy device with its microcatheter, blood, thrombi or emboli is stuck in the Esperance pHLO AC remove all devices (Esperance pHLO AC & thrombectomy device & microcatheter) out of the patient body. Maintain aspiration to make sure the thrombi or emboli remain fully engaged within the Esperance pHLO or with the distal tip of the Esperance pHLO AC. With the thrombi or emboli fully engaged pull back the Esperance pHLO AC out of the patient body and clear the inner lumen of the Esperance pHLO AC.
Warning: Do not attempt to clear inner lumen of the Esperance pHLO AC by infusion while keeping the device in the patient body.
50. Perform an angiographic run to confirm vessel recanalization.
51. If re-access to the vasculatures with the same devices is needed, flush and clean them. Inspect the devices for any damage according their DFUs.
Warning: Do not use the devices if any damage or irregularities are observed.
Follow steps 11 through 22 and 41 through 52 for the re-use of Esperance pHLO AC with syringe or aspiration pump, as well as thrombectomy device and its microcatheter.

Esperance pHLO AC for introduction of suitable devices

52. Make sure the distal tip of Esperance pHLO AC is located in the selected vessel under fluoroscopy.
53. Introduce only suitable devices.
54. Remove the devices before complete removal of Esperance pHLO AC.

Performance Characteristics

The Esperance pHLO AC is used as a standalone or combinational device to perform aspiration or aspiration with mechanical thrombectomy, in neurovascular vessels. The Esperance pHLO AC used in combination with standard guide catheters and accessories, facilitates the introduction of suitable devices and/or agents to enable access/diagnosis/treatment of conditions in the neurovasculature.

Symbols Glossary as per ISO 15223-1

Reference Number or Symbol	Title of symbol	Description of Symbol
 5.1.1	Manufacturer	Indicates the medical device manufacturer, as defined in the EU directives 90/385/EEC, 98/42/EEC and 98/79/EC and Medical Device Regulation 2017/745.
 5.1.3	Date of manufacture	Indicates the date when the medical device was manufactured.
 5.1.4	Use-by date	Indicates the date after which the medical device is not to be used.
 5.1.5	Batch code	Indicates the manufacturers batch code so that the batch or lot can be identified.
 5.1.6	Catalogue number	Indicates the manufacturers catalogue number so that the medical device can be identified.
 5.2.3	Sterilized using ethylene oxide	Indicates medical device has been sterilized using ethylene oxide.
 5.2.6	Do not re-sterilize	Indicates a medical device that is not to be re-sterilized.
 5.2.8	Do not use if package is damaged and consult instructions for use	Indicates that a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information.
 5.2.14	Single Sterile Barrier System with Protective Packaging Outside	Indicates a single sterile barrier system with protective packaging outside. Used in conjunction with symbol 5.2.3.
 5.3.2	Keep away from sunlight	Indicates a medical device that needs protection from light sources.
 5.3.4	Keep dry	Indicates a medical device that needs to be protected from moisture.
 5.4.2	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
 5.4.3	Consult instructions for use	Indicates the need for the user to consult the instructions for use.
 5.4.4	Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
 5.6.3	Non-pyrogenic	Indicates a medical device that is non-pyrogenic.
 5.7.7	Medical device	Indicates the item is a medical device
 5.7.10	Unique Device Identifier	Indicates a carrier that contains Unique Device Identifier information.
 0344	CE Mark with Notified Body Number/ ID	Indicates that a device has received CE approval and the number of the Notified Body.
 CONT	Content	Indicates the contents included in the packaging
	Inner diameter of Catheter	Indicates the inner diameter of the catheter

Summary of Safety and Clinical Performance

The summary of Safety and Clinical Performance (SSCP) for this device is available to view at <https://phenox.net/international/mdr>.

Once EUDAMED is available, the SSCP will also be available to view on <https://ec.europa.eu/tools/eudamed/#/screen/search-device> with the following UDIs:

- UDI-DI 539153091 035 4
- UDI-DI 539153091 036 1

Specification of the Clinical Benefit to be Expected

This device aids in the revascularization of patients with acute ischemic stroke secondary to large intracranial vessel occlusive disease and/or introduction of suitable devices (e.g., microcatheter) as well as the administration of diagnostic and therapeutic agents.

For Use by Professional Clinicians Trained in the Use of these device types only

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the authority per complaints guidance in your jurisdiction.

Limits of Liability

phenox Limited shall not be liable for damage caused by use other than that intended or reuse of the product.

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