INTENDED PURPOSE

The Element™ Vascular Access System is designed to provide a conduit for introduction of interventional devices to the peripheral anatomy.

DEVICE DESCRIPTION

The Element Vascular Access System includes an Element Sheath and a Dilator. The Element Vascular Access System should be used with a 0.038 in [0.97 mm] guidewire to access the desired anatomy. The Element Vascular Access System may be provided with an introducer.

Element Sheath:

The Element Sheath is a single lumen, variable stiffness sheath with a radiopaque markerband on the distal end and a Valve System on the proximal end. The Element Sheath dimensions are included on the label. The Element Sheath has a hydrophilic coating on the distal segment of the shaft. The hydrophilic coating length is noted in Table 1.

Dilator

The Dilator is a single lumen, radiopaque catheter with a tapered distal end and a Luer hub on the proximal end. The Dilator is compatible with the Element Sheath. The Dilator facilitates the percutaneous entry of the Element Sheath by forming an atraumatic transition from the skin through the subcutaneous tissue to the vessel.

Table 1 – Hydrophilic Coating	
Element Sheath Nominal Effective Length (cm)	Hydrophilic Coating Length (cm)
≤ 33	7
> 33	15

INDICATIONS FOR USE

The Element Vascular Access System is indicated for the introduction of interventional devices into the peripheral vasculature.

INTENDED PATIENT POPULATION

The intended patient population is patients who require peripheral intervention.

CONTRAINDICATIONS

There are no known contraindications.

WARNINGS

- The Element Vascular Access System should only be used by physicians who have received appropriate training in interventional techniques.
- The device is intended for single use only. Do not resterilize or reuse. Resterilization and/or reuse may compromise the structural integrity of the device and/or lead to device failure which in turn, may result in patient injury, illness, or death.
- Do not advance or withdraw the Element Vascular Access System against resistance without careful assessment of the cause using fluoroscopy. If the cause cannot be determined, withdraw the device. Unrestrained moving or torquing the device against resistance may result in damage to the vessel or device.
- Do not use automated high-pressure contrast injection equipment with the Element Vascular Access System because it may damage the device.
- Failure to abide by the warnings in this labeling might result in damage to the device coating, which may necessitate intervention or result in serious adverse events.

PRECAUTIONS

- Do not use kinked or damaged devices. Do not use open or damaged packages. Return all damaged devices and packaging to the manufacturer / distributor.
- Use prior to the “Use By” date.
- Use the Element Vascular Access System in conjunction with fluoroscopic visualization.
- Maintain a constant infusion of an appropriate flush solution.
- If flow through the device becomes restricted, do not attempt to clear the lumen by infusion. Remove and replace the device.
- As in all fluoroscopy procedures, consider all necessary precautions to limit patient radiation exposure by using sufficient shielding, reducing fluoroscopy times and modifying radiation technical factors whenever possible.
- Avoid using alcohol, antiseptic solutions, or other solvents to pre-treat the device because this may cause unpredictable changes in the coating which could affect the device safety and performance.

POTENTIAL ADVERSE EVENTS*

- Possible complications include, but are not limited to, the following:
- Access site complications such as hematoma, hemorrhage, spasm, occlusion, or compartment syndrome
 - Acute vessel occlusion
 - Air embolism
 - Allergic reaction and anaphylaxis from contrast media or device material
 - Arrhythmia
 - Arteriovenous fistula
 - Cardiac Injury
 - Cardio-Respiratory Arrest
 - Death
 - Emboli including foreign body embolization
 - Hemoptysis
 - Infection
 - Inflammation
 - Intracranial hemorrhage
 - Ischemia
 - Nerve injury
 - Neurological deficits including stroke
 - Pneumothorax
 - Pseudoaneurysm
 - Radiation exposure that may lead to cataracts, skin reddening, burns, alopecia, or neoplasia from x-ray exposure
 - Renal impairment or acute renal failure from contrast media
 - Respiratory Failure
 - Valvular Damage
 - Vessel dissection
 - Vessel perforation
 - Vessel spasm
 - Vessel thrombosis

*If a serious incident related to the device occurs, contact your Penumbra, Inc. representative and the competent authority in your respective country/region.

Table 2 – Device Compatibility and Selection Information	
Element Vascular Access System	0.038" (0.97 mm) Guidewire

DEVICE PREPARATION and Use

1. Refer to **Warnings, Precautions, and Potential Adverse Events** prior to use.
2. Refer to **Table 2** for device compatibility information.
3. Select an appropriately sized Element Sheath based on the anatomy and length.
4. Gently remove the Element Sheath and packaging card from the pouch by grasping the Element Sheath Valve System and packaging card and slowly pulling them out of the pouch.
5. Remove the Element Sheath from the packaging card by removing the Valve System from the card tabs before gently removing the Element Sheath shaft.

6. Gently remove the Dilator from the packaging card by removing the hub from the card tabs before gently removing the Dilator shaft.

7. Inspect the Element Sheath for kinks or other damage.

8. Flush the Element Sheath lumen and hydrate the outer surface with heparinized saline.

9. Insert the Dilator through the Valve System of the Element Sheath, snapping it into place at the hub.

10. Flush the lumen of the Dilator with heparinized saline.

11. Introduce the cannula of an angiographic needle into the vessel using aseptic technique. Holding the needle in place, insert the flexible end of a 0.038 in [0.97 mm] guidewire through the needle and into the vessel. Gently advance the guidewire to the desired depth.

12. Holding the guidewire in place, withdraw the needle and apply pressure to the puncture site until the Element Sheath is inserted into the vasculature.

13. Thread the Element Sheath/Dilator assembly over the guidewire, grasping the Element Sheath close to the skin to prevent buckling. Using a rotating motion, advance the assembly through the tissue into the vessel.

14. Under fluoroscopy, advance the Element Sheath until the desired position is achieved.

15. Detach the Dilator from the Element Sheath by releasing the snap-fit ring at the hub. Withdraw the guidewire and Dilator. To avoid damage to the Element Sheath, do not withdraw the Dilator until the Element Sheath is in the vessel.

16. Aspirate from the Valve System to remove any potential air. After aspiration and flushing, consider establishing a heparinized solution or suitable isotonic solution via the Valve System.

17. Insert the 0.038 in [0.97 mm] guidewire into the Element Sheath, and advance the guidewire until the distal tip of the Element Sheath.

18. Extend the 0.038 in [0.97 mm] guidewire 5 cm to 6 cm distal to the tip of the Element Sheath.

19. Advance the Element Sheath and guidewire system, and then select the appropriate vessel with the 0.038 in [0.97 mm] guidewire.

20. Advance the Element Sheath and guidewire to the vascular site and remove the guidewire.

21. After the procedure, remove the Element Sheath when clinically indicated by placing compression on the vessel above the puncture site and slowly withdraw the Element Sheath. Discard the Element Sheath appropriately.

22. Important: Upon removing the Element Sheath, aspirate via the Valve System to collect any fibrin that may have been deposited within or at the tip of the sheath.

CLINICAL BENEFITS AND PERFORMANCE CHARACTERISTICS

The Element Vascular Access System aids in the use of interventional devices to the peripheral anatomy. The detailed listing of all known benefits is included in the Summary of Safety and Clinical Performance (SSCP).

SAFETY AND CLINICAL PERFORMANCE

A copy of the Summary of Safety and Clinical Performance can be viewed by searching the device name on the EUDAMED website: <https://ec.europa.eu/tools/eudamed>
Basic UDI DI: 081454801ACCESSY8

STORAGE AND DISPOSAL

- Store in a cool, dry place.

• Dispose of device in accordance with hospital, administrative, and/or local government policy.

SYMBOLS GLOSSARY

		Attention, see Instructions for Use	
Rx Only		Prescription only – US Federal Law restricts this device to use by or on the order of a physician	
		Nonpyrogenic	
		Sterilized using ethylene oxide. Single sterile barrier system with protective packaging outside.	
		Do not use if package is damaged.	
MD	Medical Device		Date of Manufacture
	Manufacturer	EC REP	EU Representative
LOT	Lot number		Do not reuse
	Use by	REF	Catalogue number
	Do not resterilize		Not made with natural rubber latex

WARRANTY

Penumbra Inc. (Penumbra) warrants that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular use. Handling, storage, cleaning, and sterilization of this device as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond Penumbra's control directly affect the device and the results obtained from its use. Penumbra's obligation under this warranty is limited to the repair or replacement of this device and Penumbra shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this device. Penumbra neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. Penumbra assumes no liability with respect to devices reused, reprocessed, or resterilized and makes no warranties, expressed or implied, including but not limited to merchantability or fitness for intended use, with respect to such devices.

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