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INSTRUCTIONS FOR USE – SMR STEMLESS SYSTEM

WARNING: Before using the Instruction for Use, the surgeon shall verify if a more recent version of the document has been released on the eIFU website (eifu.limacorporate.com). In case a more recent version is available, it shall be downloaded and read since it replaces the present Instruction for Use.

Before using a LimaCorporate product, the surgeon shall carefully study the following recommendations, warnings and instructions, as well as the most updated available product-specific information (e.g.: product literature, surgical technique).

1. PRODUCT INFORMATION

The SMR Stemless System is intended for partial or total, primary or revision, shoulder joint replacement. The components are intended for use in uncemented applications as specified in Table 1 of section 2.1. The main goal of the SMR Stemless Anatomic joint prosthesis is to reproduce the articular anatomy, in part or totally, while in the SMR Stemless Reverse shoulder the geometry of the natural gleno-humeral joint is inverted. The joint prosthesis is intended to reduce pain and give articular mobility to the patient. The degree of pain reduction and mobility depend in part on the pre-operative situation, intra-operative options and the post-op rehabilitation. SMR Stemless System components include stemless core, adaptors with screw and reverse liners and it is compatible with some components of the LimaCorporate SMR Shoulder System (see section 3.1). SMR Stemless shoulder components should not be used with components from other systems or other manufacturers.

1.1. MATERIALS

The materials of the SMR Stemless System are reported in the following table:

Components	Material
SMR Stemless Core	Ti6Al4V and Ti6Al4V 3D printed
SMR Stemless Adaptors with screw	Ti6Al4V
SMR Stemless Reverse Liners	CoCrMo
SMR Stemless Adaptor Screw	Ti6Al4V

Ti6Al4V and Ti6Al4V 3D printed (ISO 5832-3) – Ti6Al4V (ISO 5832-3) – CoCrMo (ISO 5832-12 - ASTM F1537)

All materials are purchased according to International Standards and are widely used in the orthopaedic field. The risks related to the use of these materials have been considered and the related risk/benefit ratio has been evaluated as acceptable by the manufacturer. Some patients may be vulnerable (e.g. have hypersensitive or allergic reactions) to the implant materials/substances; this shall be appropriately considered by the Surgeon.

1.1.1. MEDICINAL SUBSTANCES AND TISSUE

The device does not contain or incorporate any: • medicinal substance, including a human blood or plasma derivative; • tissues or cells, or their derivatives, of human origin; • tissues or cells of animal origin, or their derivatives, as referred to in Regulation (EU) No.722/2012.

1.2. HANDLING, STORAGE AND SAFE DISPOSAL

All devices are provided sterile and shall be stored at ambient temperature (indicative range 0-50°C / 32-122°F) in their protective closed packaging in controlled rooms, protected from exposure to light, heat and sudden changes in temperature.

Once the package is opened check that both the model and size of the implant correspond exactly to the description printed on the labels. Avoid any contact between the implant and objects or substances that can alter the sterile condition or surface integrity. A careful visual examination of each implant is recommended before use in order to verify that the implant is not damaged. **Components removed from the package shall not be used if they are dropped or suffer other accidental impacts. Devices shall not be modified in any way.**

The device's code and lot number shall be recorded in the patient's case history by using the labels included in the components packaging.

Re-use of previously implanted devices must be absolutely avoided. Risks associated with the reuse of single use devices are: • infection; • early or late failure of the device or device fixation; • lack of appropriate coupling between modular junctions (e.g. taper connections); • device wear and wear debris associated complications; • transmission of diseases (e.g. HIV, hepatitis); • immune system response / rejection.

DISPOSAL: the disposal of medical devices is to be performed by the hospitals in conformity with applicable laws.

1.3. DEVICE PERFORMANCE, EXPECTED CLINICAL BENEFITS, EXPECTED LIFETIME

The main goal of a joint prosthesis is to reproduce the articular anatomy. The joint prosthesis is intended to reduce pain and give articular mobility to the patient. The degree of pain reduction and mobility depend, in part, on the pre-operative situation, intra-operative options and the post-op rehabilitation.

LimaCorporate annually updates the Summary of Safety and Clinical Performance (SSCP) to be uploaded on EUDAMED (<https://ec.europa.eu/tools/eudamed>).

The expected lifetime of SMR Stemless System is:

at least 93%* at 7 years.

*** with a variability of $\pm 2\%$**

This value is based on the state-of-the-art reference for similar device (stemless shoulder replacement) identified through arthroplasty registries.

NOTE: The expected lifetime of the medical devices, when placed on the market and – moreover – when put into use (namely when implanted in a patient's body) may depend on a variety of factors:

(a) intrinsic characteristic (feature) of the medical device in itself (for example: design of the medical device; materials employed in the medical device manufacturing; manufacturing process of the medical device in itself; surface texture, osteo-conductive and osteo-inductive properties of the medical device);

(b) external factors related to surgery process (for example: right choice of the medical device, stability of the initial fixation of the medical device; accuracy in positioning of the medical device);

(c) external factors related to patient (for example: aetiology, pathology, physical characteristics, age, previous surgical history, daily activities performed);

(d) risk factors mentioned in the IFU leaflet, when not already included in the factors referred to in previous paragraphs from (a) to (c);

(e) factors related to all the possible complications which can affect the duration of the implant.

The expected lifetime was calculated by means of an analytical and rationalized estimation, taking into account:

i. the clinical data collected on the medical device;

ii. the survivorship data reported in arthroplasty registries referring to similar medical devices (identified state of the art benchmark);

iii. the expected (as predictable) behavior of the medical device at higher follow-up.

Due to the analytical evaluations, the expected lifetime value is approximated to the unit.

The value shall be interpreted duly considering the specified value variability, taking into necessary consideration all relevant factors mentioned above, in paragraphs (a) to (e).

2. INTENDED PURPOSE, INDICATIONS, CONTRAINDICATIONS AND RISK FACTORS

SMR Stemless System is intended for cementless shoulder joint replacement and for long-term implantation inside the human body.

2.1. INDICATIONS

The SMR Stemless System is intended for partial or total, primary or revision, shoulder joint replacement in skeletally mature patients. The patient's joint must be anatomically and structurally suited to receive the selected implants and a functional deltoid muscle is necessary to use the device.

The SMR Stemless Anatomic System is indicated for partial or total, primary shoulder joint replacement, in patients suffering from pain and disability due to: • non-inflammatory degenerative joint disease including osteoarthritis; • inflammatory degenerative joint disease such as rheumatoid arthritis; • avascular necrosis of the humeral head; • cuff tear arthropathy (CTA Heads only, not applicable to Sub-Scapularis Sparing technique).

The SMR Stemless Reverse implant is indicated for primary reverse total shoulder replacement or for revision when converting an anatomic SMR Stemless arthroplasty to a reverse total shoulder arthroplasty in case of cuff tear arthropathy or in a grossly rotator cuff deficiency joint with severe arthropathy.

In particular, the SMR Stemless Reverse implant is indicated for patients suffering from pain and disability due to: • rotator cuff tear arthropathy; • osteoarthritis with rotator cuff tear; • rheumatoid arthritis with rotator cuff tear; • massive irreparable rotator cuff tear.

Revision surgery is intended as conversion surgery from anatomic to reverse, where the Stemless Core is stable, well positioned and osseointegrated. Other revisions of the humeral prosthesis part should be treated with traditional shoulder prostheses.

The Modular SMR Stemless System allows the assembly of components in various humeral and glenoid constructs. The constructs are intended for cemented (only glenoids) or uncemented use as specified in the following table.

In the Anatomic shoulder the humeral construct consists of the stemless core and adaptor with screw from the SMR Stemless System coupled to SMR Shoulder System humeral heads. In the Reverse shoulder the humeral construct consists of the stemless core and the reverse liner.

The glenoid construct, anatomic or reverse uses SMR Shoulder System components. The Anatomic glenoid construct consists of an all polyethylene glenoid or a metal back assembled with a liner while the Reverse glenoid construct consists of the metal back, the connector and the glenosphere. On the glenoid side, the fixation of the all polyethylene glenoid or the metal back determines if the construct is cemented or uncemented.

Table 1

System		Components	Use	
Anatomic	Reverse		Cemented	Not Cemented
•	•	SMR Stemless Core		X
•		SMR Stemless Adaptors with screw		X
	•	SMR Stemless Reverse Liners		X
•		SMR Stemless Adaptor Screw		X

2.2. CONTRAINDICATIONS

Absolute contraindications include: • local or systemic general infection; • septicaemia; • persistent acute or chronic local or systemic osteomyelitis; • neurologically confirmed nerve lesion compromising shoulder joint function; • in case of reverse prosthesis, deltoid muscle insufficiency; • poor meta-epiphyseal bone stock compromising stability of the implant (acute fracture of the humeral head, meta-epiphyseal pseudoarthrosis, osteoporosis, extended bone loss after previous prosthetic or non-prosthetic surgery); • tumor.

Relative contraindications include: • humeral head fracture sequelae; • vascular or nerve diseases affecting the concerned limb; • metabolic disorders which may impair fixation and stability of the implant; • any concomitant disease and dependence that might affect the implanted prosthesis; • metal hypersensitivity to implant materials.

2.3. RISK FACTORS

The following factors may result in poor results with this prosthesis: • overweight¹; • strenuous physical activities (active sports, heavy physical work); • incorrect implant positioning; • muscle deficiencies; • multiple joint disabilities; • refusal to modify postoperative physical activities; • patient history of infections or falls; • systemic diseases and metabolic disorders; • local or disseminated neoplastic diseases; • drug therapies that adversely affect bone quality, healing, or resistance to infection; • drug use or alcoholism; • marked osteoporosis or osteomalacia; • immunocompromised patients (HIV, tumour, infections); • severe deformity leading to impaired anchorage or improper positioning of implants; • use or combination with products, prosthesis or instruments of another manufacturer.

¹ According to the definition of World Health Organization (WHO), Body Mass Index (BMI) greater than or equal to 25 Kg/m²

3. WARNINGS

3.1. PRE-OPERATIVE PLANNING

Lima Corporate products shall be implanted only by surgeons familiar with the joint replacement procedures described in the specific surgical techniques. The components of the system must not be used with components of other manufacturer, unless specified by LimaCorporate. The components forming original LimaCorporate systems must be rigorously coupled following the given information included in this document and in the available product-specific information (e.g.: current version of surgical technique). The use of single components, or components belonging to other systems, are subject to the approval by LimaCorporate. The manufacturer and the dealer are not liable for possible coupling incompatibility. The surgeon is solely responsible for the choice and use of the implant. Pre-operative planning, through radiographic templates in different formats, provides essential information regarding the type and size of components to be used and the correct combination of devices required based on the anatomy and specific conditions of each patient. Inadequate pre-operative planning can lead to improper selection of the implants and/or incorrect implant positioning.

In selecting patients for surgery, the following factors can be critical to the eventual success of the procedure:

- **Partial Shoulder replacement:** in cases of a deficient and unreconstructable rotator cuff, a CTA-head is indicated;
- **Total Shoulder replacement:** the rotator cuff must be intact or reconstructable. In cases of a deficient and unreconstructable rotator cuff, a hemiprosthesis with a CTA head or a Reverse Total Shoulder Arthroplasty is indicated;
- **Reverse Shoulder replacement:** the bone stock of the glenoid and humerus must be able to support the implant. In cases of significant bone loss or in which adequate fixation on the glenoid side cannot be obtained, a hemiarthroplasty with a CTA-head should be performed.

The SMR Shoulder System components which are compatible with SMR Stemless System components in anatomical configuration are: • SMR Humeral Heads • SMR CTA Heads • SMR Cemented Glenoids • SMR Cemented Glenoids 3 pegs • SMR Metal back glenoid • SMR Axioma Metal back glenoid • SMR TT Baseplate • SMR Axioma TT Baseplate • SMR Axioma TT Glenoid Peg • SMR Metal Back Liner • SMR Axioma Metal Back Liner • SMR Bone Screws • SMR TT Hybrid Glenoid.

The SMR Shoulder System components which are compatible with SMR Stemless System components in reverse configuration are: • SMR Reverse HP Genospheres • SMR Connectors with screws • SMR Metal back glenoid • SMR Axioma Metal back glenoid • SMR TT Baseplate • SMR Axioma TT Baseplate • SMR Axioma TT Glenoid Peg • SMR Bone Screws • SMR Cortical Bone Screws • SMR Glenoid Plates • SMR TT Hybrid Glenoid Reverse Baseplate + Screw • SMR TT Augmented 360 Baseplate • SMR Lateralized Connector with Screw.

NOTE on SMR Humeral Heads: • SMR Humeral Heads diameter 38 mm are not allowed to be coupled to SMR Stemless • SMR Humeral Heads Dia 50mm H22, Dia 52mm H22, Dia 52mm H23, Dia 54mm H23, Dia 54mm H24 in combination with both neutral and eccentric adaptors are not allowed to be coupled to SMR Stemless.

• SMR Humeral Heads Dia 40mm H12, Dia 40mm H13, Dia 42mm H12, Dia 42mm H13, Dia 44mm H13 both neutral and eccentric versions are not allowed to be coupled to SMR Stemless.

NOTE on SMR CTA Heads: • If the CTA Head has to be used, make sure that an Eccentric Adaptor is used (the coupling with the Neutral Adaptor is not allowed) and that the eccentricity is in the cranial direction only.

NOTE on SMR Reverse HP Glenospheres: • Glenosphere 38 mm (UHMWPE X-LIMA + Ti6Al4V) is not allowed to be coupled to SMR Stemless

NOTE: For SMR Shoulder System components refer to the proper Instructions for Use leaflet.

NOTE: the following SMR Shoulder System products are NOT available in Canada: • SMR Cemented Glenoids 3 pegs • SMR Axioma Metal back glenoid • SMR Axioma TT Baseplate • SMR Axioma Metal Back Liner • SMR Glenoid Plates • SMR Humeral Heads (CoCrMo) diameter 38 mm (neutral and eccentricity 2, 4mm), diameter 40mm H12 (neutral and eccentricity 2, 4mm), diameter 42mm H12 (neutral and eccentricity 2, 4mm), diameter 44mm H13 (neutral and eccentricity 2, 4, 6 and 7mm) and diameter 46mm H14 and diameter 48mm H15 • SMR Humeral Heads (Ti6Al4V) diameter 38 mm (neutral and eccentricity 2, 4mm), diameter 40mm H12 (neutral and eccentricity 2, 4mm), diameter 42mm H12 (neutral and eccentricity 2, 4mm), diameter 44mm H13 (neutral and eccentricity 2, 4, 6 and 7mm) and diameter 46mm H14 and diameter 48mm H15, diameter 50mm H16, diameter 52 H17 and diameter 54mm H18.

LimaCorporate specialized technical staff is available to provide advice regarding pre-operative planning, the surgical technique, and product and instrumentation assistance both prior to and during surgery.

The patient shall be warned that the prosthesis does not replace normal healthy bone, that the prosthesis can break or become damaged as a result of certain activity or trauma, that it has a finite expected implant life, and may need to be replaced at some time in the future: the longevity and the performance of the implant over time, in fact, can be influenced by the natural/physiological progression of the patient pathology, the insurgence of co-morbidities and post-operative complications, even if they are not directly implant-related (i.e. infection pain, stiffness, reduced range of motion).

The possible impact of the factors mentioned in section 2 and 3.4 should be considered preoperatively and the patient informed as to what steps he/she can take to reduce the possible effects of these factors.

Implants are single use devices; do not re-use implants that were previously implanted in another patient. Do not re-use an implant that has previously come into contact with the body fluid or tissue of another person. Do not use any component from a package that has been previously opened or appears to be damaged. Do not use implants after the expiration date printed on the label.

Do not use this product for other than labelled indications (off-label use); the risk of implant failure is higher with inaccurate component alignment or positioning. Use only instruments and provisionals specifically designed for use with these devices to help ensure accurate surgical implantation, soft tissue balancing, and evaluation of shoulder function.

Surgical instruments are subject to wear with normal usage. After extensive use or excessive loads, instruments are susceptible to fracture. Surgical instruments shall be used only for their specific purpose. Before use, the functionality of surgical instruments shall be checked since the use of damaged instruments may lead to early failure of the implants. Damaged instruments shall be replaced before surgery.

3.2. INTRAOPERATIVE

The use of trial devices is recommended to check the correct site preparation, size and positioning of the implants to be used. It is recommended that additional implants be available during surgery for use in those cases requiring prostheses of different sizes or when the preoperatively selected prostheses cannot be used. The correct selection as well as the correct seating/placement of the implant is extremely important. Improper selection, positioning, alignment and fixation of the implant components may result in unusual stress conditions which may negatively affect system performance and survival rate of the implant.

The components forming original Limacorporate systems must be assembled according to the surgical technique and used only for the labelled indications.

Use only instruments and prosthesis trials specifically designed for use with the implants being used. The use of instruments from other manufacturers or the use of instruments designed for use with other systems can lead to inappropriate preparation of the implant site, incorrect positioning, alignment and fixation of the devices followed by loosening of the system, loss of functionality, reduction of the durability of the implant, and the need for further surgery.

Care must be taken to protect the surfaces involved in the coupling between components (tapers); articular surfaces of the implants should be protected from scratches or any other damage. All component coupling surfaces shall be clean and dry before assembly. The stability of component couplings should be verified as described in the surgical technique.

Care must be taken when handling Trabecular Titanium devices; they should not get in contact with fabrics that can release particles inside the Trabecular structure.

3.3. POST-OPERATIVE CARE

Adequate post-operative care shall be provided by the surgeon or other suitably qualified medical staff. Regular postoperative x-ray follow up is recommended to detect any changes in position of the implant or surrounding tissues.

Complications and/or failure of prosthetic implants are more likely to occur in patients with unrealistic functional expectations, heavy patients, physically active patients, and/or with patients who fail to follow through with the required rehabilitation program.

Excessive physical activity or trauma to the replaced joint can lead to premature failure of the arthroplasty

through loosening, fracture, or abnormal wear of the prosthetic implants. The patient shall be cautioned by the surgeon to govern activities accordingly and advised that the implants may fail due to excessive joint wear.

The surgeon shall make the patient aware of the limitation of limb function after shoulder arthroplasty and that caution is needed, especially in the post-operative period. In particular the following precautions should be presented to the patient by the surgeon:

- avoid repeated high weight lifting;
- keep body weight under control;
- avoid sudden peak loads (consequences of activities such as contact sports, playing tennis) or movements which can lead to sudden stops or twisting;
- avoid positions that can increase the risk of dislocation, such as for example excessive degrees of abduction and/or internal rotation.

The patient shall be informed also about the possibility that the implant or its components may wear out, fail or need to be replaced. The implant may not last the rest of the patient's life, or any particular length of time. Because prosthetic implants are not as strong, reliable, or durable as natural, healthy tissue/bones, all such devices may need to be replaced at some point.

Appropriate post-operative rehabilitation instructions and care must be provided to avoid a negative influence on the outcome of the surgical procedure.

3.4. POSSIBLE SIDE EFFECTS

Adverse effects that most commonly and frequently occur in orthopaedic prosthesis procedure include:

- loosening of the prosthetic components, for example due to osteolysis;
- prosthesis dislocation and instability;
- damage to the prosthetic implant;
- instability of the system because of inadequate soft tissue balancing;
- dissociation or malalignment due to incorrect coupling of the devices;
- infection;
- local hypersensitivity;
- local pain;
- periprosthetic fractures, including intraoperative fractures;
- temporary or permanent nerve damage;
- fractures of the devices;
- fretting of modular junctions;
- excessive wear of UHMWPE components due to damaged articular surfaces or the presence of particles;
- heterotopic ossification;
- additional surgery.

Some adverse effects can lead to death.

General complications include venous thrombosis with/without pulmonary embolism, cardiovascular or pulmonary disturbances, haematomas, blood loss, post-operative anaemia, systemic allergic reactions, systemic pain, temporary or permanent paralysis.

4. STERILITY

4.1. Implants

All implantable components of the System are provided sterile with a Sterility Assurance Level (SAL) of 10⁻⁶. Metal components are sterilized by radiations or EtO and UHMWPE components by EtO. Do not use any component from a package that has been previously opened or appears to be damaged. **Do not use implants after the expiration date printed on the label.**

4.2. Instruments

Instruments are supplied non-sterile and must be cleaned, disinfected and sterilized before use according

to appropriate validated methods (refer to the "Instrument Care, Cleaning, Disinfection and Sterilization" brochure for validated sterilization parameters; this brochure is available on request or downloadable from www.limacorporate.com in the Products section). Users should validate their specific cleaning, disinfection and sterilization processes and equipment.

5. INFORMATION USEFUL TO INFORM THE PATIENT

Refer to previous sections for any information that allows the patient to be informed of any warnings, precautions, contra-indications, measures to be taken and limitations of use regarding the device.

The surgeon shall notify the patient about any necessary follow-up specific for the patient's condition.

5.1. MAGNETIC RESONANCE IMAGING (MRI)

LimaCorporate evaluated the interaction between the MRI and components of LimaCorporate shoulder portfolio, by performing non-clinical testing. Stemless Systems is included in this evaluation and have been classified as MR conditional.

A patient with the tested devices can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 Tesla and 3 Tesla, with
 - Maximum spatial field gradient of 8,500 G/cm (85 T/m)
 - Maximum force product of 154,000,000 G²/cm (154 T²/m)
 - Theoretically estimated maximum whole body averaged (WBA) specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode)

Under the scan conditions defined above, the SMR Shoulder System is expected to produce a maximum temperature rise of less than:

- 6.0°C (2 W/kg, 1.5 Tesla) RF-related temperature increase with a background temperature increase of: ≈0.4°C (2 W/kg, 1.5 Tesla, SMR anatomical shoulder system) after 5 minutes of continuous scanning; ≈0.2°C (2 W/kg, 1.5 Tesla, SMR reverse shoulder system) after 4.3 minutes of continuous scanning;
- 6.0°C (2 W/kg, 3 Tesla) RF-related temperature increase with a background temperature increase of: ≈0.6°C (2 W/kg, 3 Tesla, SMR anatomical shoulder system) after 12 minutes of continuous scanning; ≈0.9°C (2 W/kg, 3 Tesla, SMR reverse shoulder system) after 12.5 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device ranges between approximately from 72.9 mm to 98.9 mm for the Shoulder tested components when imaged with a gradient echo pulse sequence and a 3 Tesla MR system.

6. INCIDENT REPORTING

According to the definition of incident/serious incident given by Regulation EU 2017/745², users and/or patients are required to: • **report to the manufacturer and to the relevant Competent Authority any serious incident that has occurred in relation to the device;** • **report to the manufacturer any incident and/or event of malfunction of the device or changes in its performance that may affect safety.**

In non-UE market, refer to incidents' definitions according to applicable laws.

² an incident is defined to be any malfunction or deterioration in the characteristics and/or performance of the device, including use-error due to ergonomic features, as well as any inadequacy in the information supplied by the manufacturer and any undesirable side-effect. An incident is considered to be serious when it directly or indirectly led, might have led or might lead to any of the following: (a) the death of a patient, user or other person, (b) the temporary or permanent serious deterioration of a patient's, user's or other person's state of health, (c) a serious public health threat.

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Limacorporate S.p.A.