PINCE D'EXTRACTION 2 GRIFFES / EXTRACTION FORCEPS 2 PINCERS / EXTRAKTIONSZANGE 2 ZANGEN / PINZAS DE EXTRACCIÓN 2 PINZAS / PINZE DI ESTRAZIONE 2 PINZE / EXTRACTIETANG 2 TANGEN / EKSTRAKTIONSTANG 2 TÆNGER / EKSTRAKSJONSTANG, 2 HAKER / POISTOPIHDIT, 2 NIPISTINTÄ / PINÇA DE EXTRAÇÃO DE DUAS TENAZES / ملقط استخراج ذو كلابين / EXTRAČNÍ KLEŠTĚ S 2 ÚCHOPOVÝMI ČELISTMI / HVATAJKE ZA VAĐENJE S 2 KLJIJEŠTA / FORSEP EKSTRAKSI 2 PENJEPIT / FORCEPS DE EXTRAGERE CU 2 GHEARE / щипцы для экстракции, 2 БРАНШИ

Classe I stérile / Class I sterile / Klasse I steril / Clase I estéril / Klasse I steriel / Luokan I steriili / Classe I estéril / التعقيم من الفئة I / Sterilní třída I / Klasa I sterilno / Steril Kelas I / Clase I, steril / Класс I, стерильное изделие

RÉFÉRENCE / REFERENCE / REFERENZ / REFERENCIA / RIFERIMENTO / REFERENTIE / REFERENCE / REFERANSE / VIITE / REFERÈNCIA / الرقم المرجعي / REFERENCE / REFERENSI / REFERINȚĂ / ССЫЛКА	DIAMÈTRE DE LA GAINE / SHEATH DIAMETER / MANTELDURCHMESSER / DIÁMETRO DE LA VAINA / DIAMETRO DELLA GUAINA / DIAMETER VAN DE MANTEL / SKEDDE DIAMETER / SKILDE DIAMETER / YAIPAN HALKAISIJA / DIÁMETRO DA BAINHA / قطر الغند / PRŮMĚR PLÁŠTĚ / PROMJER PLAŠTA / DIAMETER SELUBUNG / DIAMETRUL TECII / ДИАМЕТР ОБОЛОЧКИ (MM)	LONGUEUR / LENGTH / LÄNGE / LONGITUD / LUNGHEZZA / LENGTE / LÆNGDE / LENGDE / PITUUS / COMPRIMENTO / الطول / DÉLKA / DULJINA / PANJANG / LUNGIME / ДЛИНА (CM)
AS2290718	2,3	180
		

SIGNIFICATION DES SYMBOLES UTILISÉS POUR LES ÉTIQUETTES / SIGNIFICATION OF SYMBOLS USED ON LABELS / BEDEUTUNG DER AUF DEN ETIKETTEN VERWENDETEN SYMBOLE / SIGNIFICADO DE LOS SÍMBOLOS UTILIZADOS PARA LAS ETIQUETAS / SIGNIFICATO DEI SIMBOLI USATI PER LE ETICHETTE / BETEKENIS VAN DE VOOR DE ETIKETTEN GEBRUIKTE SYMBOLEN / BETYDNINGEN AF DE SYMBOLER, DER ANVENDES PÅ ETIKETTERNE / BETYDNING AV SYMBOLENE BRUKT FOR ETIKETTENE / ETIKETEISSÄ KÄYTETTYJEN SYMBOLEIDEN MERKITYS / SIGNIFICADO DOS SÍMBOLOS UTILIZADOS NAS ETIQUETAS / معنى الرموز المستخدمة على الملصقات / VÝZNAM SYMBOLŮ POUŽÍVANÝCH NA ŠTÍTCÍCH / ZNAČENJE SIMBOLA KOJI SE UPOTREBLJAVAJU NA ETIKETAMA / MAKNA SIMBOL YANG DIGUNAKAN DI LABEL / SEMNIFICAȚIA SIMBOLURILOR UTILIZATE PE ETICHETE / ЗНАЧЕНИЕ СИМВОЛОВ, ИСПОЛЬЗУЕМЫХ НА ЭТИКЕТКАХ

	Date de péremption / Expiry date / Verwendbar bis / Fecha de caducidad / Data di scadenza / Vervaldatum / Udløbsdato / Utløpsdato / Viimeinen voimassaolopäivä / Data de validade / تاريخ انتهاء الصلاحية / Datum použitelnosti / Rok valjanosti / Tanggal kedaluwarsa / Data de expirare / Срок годности		Date de fabrication / Manufacturing date / Herstellungsdatum / Fecha de fabricación / Data di fabbricazione / Fabrikatiedatum / Fremstillingsdato / Produktionsdato / Valmistuspäivämäärä / Data de fabric / تاريخ التصنيع / Datum výroby / Datum proizvodnje / Tanggal produksi / Data fabricației / Дата производства
	Numéro de lot / Lot number / Lotnummer / Número de lote / Numero di lotto / Lotnummer / Partiets nummer / Batchnummer / Erän numero / Número de lote / رقم الدفعة / Číslo šarže / Broj serije / Nomor lot / Număr de lot / Номер партии		Fabricant / Manufacturer / Hersteller / Fabricante / Produttore / Fabrikant / Producent / Maker / Valmistaja / Fabricante / الشركة المصنعة / Výrobce / Proizvođač / Prodsusen / Producător / Изготовитель
	Référence du produit / Product reference / Artikelnummer / Referencia del producto / Riferimento del prodotto / Productreferentie / Produktreferenz / Produktreferanse / Tuoteviite / Referència do produto / مرجع المنتج / Odkaz na produkt / Referenca proizvoda / Referensi produk / Referință produs / Код изделия по каталогу		Dispositif médical / Medical device / Medizinprodukt / Dispositivo médico / Dispositivo medico / Medisch hulpmiddel / Medicinsk udstyr / Medisinsk enhet / Lääkinnällinen laite / Aparelho medico / Zdravotnické zařízení / Medicinski uređaj / Alat medis / Aparat medical / Медицинское устройство
	Stérilisation à l'oxyde d'éthylène / Sterilized with Ethylene oxide / Sterilisiert mit Ethylenoxid / Esterilización por óxido de etileno / Sterilizzazione con ossido di etilene / Sterilisatie met ethyleenoxide / Sterilisering med ethylenoxid / Sterilisering med etylenoxid / Etyleenioksidisterilointi / Esterilização por óxido de etileno / التعقيم بأكسيد الإيثيلين / Sterilizace etylenoxidem / Sterilizacija etilen-oksidadom / Sterilisasi etilen oksida / Sterilizare cu oxid de etilenă / Стерилизация оксидом этилена		Identifiant unique de dispositif / Unique device identifier / Eenduitige Identifizierung eines Medizinprodukts / Identificador único del dispositivo / Identificazione unica del dispositivo / Unieke hulpmiddel identificatie / Unik enhedsidentifikation / Unik enhetsidentifikator / Yksilöllinen laitetunniste / Identificador exclusivo do dispositivo / معرف الجهاز الوحيد / Jedinečný identifikátor zařízení / Jedinstveni identifikator uređaja / Pengidentifikasi perangkat unik / Identificator unic de dispozitiv / Уникальный идентификатор устройства
	Ne pas utiliser si l'emballage est endommagé / Do not use if packaging is damaged / Bei beschädigter Verpackung nicht verwenden / No utilizar si el embalaje está dañado / Non utilizzare se l'imballaggio è danneggiato / Niet gebruiken indien de verpakking beschadigd is / Må ikke anvendes, hvis emballagen er beskadiget / Må ikke brukes hvis emballasjen er skadet / Älä käyttää, jos pakkaus on vahingoittunut / Não utilizar se a embalagem estiver danificada / لا تستخدم المنتج إذا كانت العبوة تالفة / Nepoužívejte, pokud je obal poškozený / Nemojte upotrebljavati ako je pakiranje oštećen / Jangan digunakan jika kemasan rusak / Nu utilizați dacă ambalajul este deteriorat / Не использовать, если повреждена упаковка		Système de barrière stérile unique avec emballage de protection intérieure / Single sterile barrier system with inner protective packaging / Einzelnes Sterilbarriersystem mit Schutzverpackung im Inneren / Sistema de barrera estéril único con embalaje protector interior / Sistema singolo di barriera sterile con imballaggio protettivo interno / Enkelvoudig steriel barrièresysteem met binnenvpakking / Enkelt sterilt barrièresystem med indre beskyttelsesemballage / Systeme de barrière stérile unique avec emballage de protection intérieure / Yksittäinen steriili sulkujärjestelmä, jossa on sisäinen suojarakkaus / Sistema de barreira estéril exclusiva com embalagem protetora interna / نظام حاجز معقم وحيد مع عبوة واقية داخلية / Unikátní sterilní bariérový systém s vnitřním ochranným obalem / Jedinstveni sustav sterilne barijere s unutarnjim zaštitnim pakiranjem / Sistem penghalang steril yang unik dengan kemasan pelindung bagian dalam / Sistem unic de barieră sterilă cu ambalaj interior de protecție / Уникальная стерильная барьерная система с внутренней защитной упаковкой

	Ne pas réutiliser / Do not reuse / Nicht wiederverwenden / No reutilizar / Non riutilizzare / Niet hergebruiken / Må ikke genbruges / Ikke gjenbruk / Älä käyttää uudelleen / Não reutilizar / ممنوع إعادة الاستخدام / Nepoužívejte opakovaně / Nemojte ponovno upotrebljavati / Jangan digunakan ulang / Nu reutilizați / Не использовать повторно		Crainť l'humidit� / Avoid humidity / Trocken aufbewahren / Evitar la humedad / Teme l'umidit� / Vermijd vocht / Frygter fugt / Frykt for fuktighet / Pelk� kosteutta / Evitar a humidade / تجنب الرطوبه / Vyhnete se vlhkosti / Izbjegavajte vlagu / Jangan disimpan di tempat lembap / Evitați umiditatea / Избегайте высокой влажности
	Ne pas rest�riliser / Do not resterilize / Nicht erneut sterilisieren / No reesterilizar / Non risterilizzare / Niet opnieuw steriliseren / M� ikke resteriliseres / Ikke re-steriliser / �l� steriloj uudelleen / N�o reesterilizar / ممنوع إعادة التعقيم / Opakovan� nesterilizujte / Nemojte ponovno sterilizirati / Jangan disterilisasi ulang / Nu resterilizați / Не стерилизовать повторно		Conservier a l'abri de la lumi�re du soleil / Store away from sunlight / Vor Sonnenlicht gesch�tzt lagern / Almacenar lejos de la luz solar / Conservare al riparo dalla luce solare / Uit de buurt van zonlicht bewaren / Opbevares v�ek fra sollys / Oppbevares unna sollys / S�ilytett�v� poissa auringonvalolta / Armazena longe da luz solar / يخزن بعيداً عن ضوء الشمس / Skladujte mimo dosah slune�nho z�ařeni / �uvati dalje od sun�eve svjetlosti / Simpan jauh dari sinar matahari / A se p�stra departe de lumina soarelui / Хранить вдали от солнечных лучей
	Consulter les instructions d'utilisation / Read the instructions for use / Gebrauchsanweisung beachten / Consulte las instrucciones de uso / Consultare le istruzioni per l'uso / Raadpleeg de gebruiksaanwijzing / Se brugsanvisningen / Se bruksanvisningen / Tutustu k�ytt�ohjeisiin / Ler as instru�es de utiliza�o / اقرأ تعليمات الاستخدام / Pře�t�te si n�vod k pouz�it� / Pro�tajte upute za uporabu / Baca petunjuk penggunaan / Cititi instrucțiunile de utilizare / Ознакомьтс с инструкцией по использованию		Distributeur / Dealer / H�ndler / Distribuidor / Distributore / Distributeur / Distribut�r / Fordeler / Jakelija / Distribuidor / الموزع / Distributor / Distributer / Distributor / Distributor / Дистрибутор
	Mandataire suisse / Swiss authorised representative / Schweizer Bevollm�chtigter / Mandatario suizo / Mandatario svizzero / Zwitserser agent / Schweizisk agent / Sveitsisk representant / Sveitsin agentti / Agente Su�o / مندوب سويسري / Šv�carsk� agent / ťivcarski predstavnik / perwakilan Swiss / Agent elvețian / Швейцарский агент		Unit� d'emballage / Packing unit / Verpackungseinheit / Unidad de embalaje / Unit� di imballaggio / Verpakkingseenheid / Emballeringsenhet / Embalasje enhet / Pakkausyksikk� / Unidade de embalagem / وحدة التعبئة والتغليف / Jednotka baleni / Jedinica pakiranja / Satuan kemasan / Unitate de ambalare / Упаковочная единица
	Attention / Warning / Warnung / Atenci�n / Avvertimento / Waarschuwing / Advarsel / Advarsel / Varoitus / Aviso / تحذير / Varov�ni / Upozorenje / Peringatan / Avvertizare / предупреждение		

La pince d'extraction est  galement pr sente dans le kit d'extraction de ballon intra-gastrique AS2710718 / The extraction forceps are also included into the intragastric extraction kit AS2710718 / Die Extraktionszange ist auch im Set zur intragastrischen Extraktion AS2710718 enthalten / Las pinzas de extracci n tambi n se incluyen en el kit de extracci n intrag strica AS2710718 / La pinza di estrazione   inclusa nel kit di estrazione intragastrica AS2710718 / De extractietang is ook opgenomen in de intragastrische extractiekit AS2710718 / Ekstraktionstangen er ogs  inkluderet i det intragastriske ekstraktionss t AS2710718 / Ekstraksjonstangen er ogs  inkludert i det intragastriske ekstraksjonsettet AS2710718 / Ulosvetopihdit sis ltyv t my s mahansis iseen ulosvetopakkausseen AS2710718 / As pinças de extrac o tambi m est o inclu das no kit de extrac o intra-g strica AS2710718 / AS2710718 ملىط الاستخراج بالون المعدة / Sou ast  sady pro intragastrickou extrakci AS2710718 jsou tak  extrak n  klešťe / Pinceta za ekstrakciju takoder je uklju ena u komplet za intragastricnu ekstrakciju AS2710718 / Tang ekstraksi juga disertakan ke dalam kit ekstraksi intragastrik AS2710718 / Forcepsul de extracție este, de asemenea, inclus  n kitul de extracție intragastric  AS2710718 / Щипцы для извлечения также входят в набор для внутрижелудочного извлечения AS2710718 :

R�F�RENCE / REFERENCE / REFERENZ / REFERENCIA / RIFERIMENTO / REFERENTIE / REFERENCE / REFERANSE / VIITE / REFER�NCIA / الرقم المرجعي / REFERENCE / REFERENSI / REFERIN�A / ССЫЛКА	D�nominati�n / Denominati�n / Benennung / Denominaci�n / Denominazione / Denominatie / Ben�vnelse / Val�r / Nimellisarvo / Denomina�o / هنة / Denominace / Vjeroispovijest / Denominasi / Denumire / Деноминация
AS2290718	Pince d'extraction 2 griffes / Extraction forceps 2 pincers / Extraktionszange 2 Zangen / Pinzas de extracci�n 2 pinzas / Pinze di estrazione 2 pinze / Extractietang 2 tangen / Ekstraktionstang 2 t�nger / Ekstraksjonstang, 2 haker / Poistopihdit, 2 nipistint� / Pinça de extrac�o de duas tenazes / ملىط استخراج ذو كلابين / Extrak�n� klešťe s 2 �chopov�ymi �elistmi / Hvatalijke za vađenje s 2 klijesta / Forsep ekstraksi 2 penjepit / Forceps de extragere cu 2 gheare / Щипцы для экстракции, 2 бранши
AS2640718	Aiguille de vidange pour ballon gastrique / Emptying needle for gastric balloon / Entleerungsnadel f�r einen magenballon / Aguja de vaciado para bal�n g�strico / Ago di svuotamento per pallone gastrico / Ledigingsnaald voor maagballon / Dr�nkanyle til gastrisk ballon / Dreneringsn�l for gastrisk ballong / Intragastrisen pallon tyhjennysneula / Agulha de drenagem para bal�o g�strico / ابرة التصريف الخاصة بالون المعدة / Dren�zn� jehla pro intragastrick�y balonek / igla za praņjenje gasti�nog balona / Jarum drainase untuk balon lambung / Ac de drenaj pentru balon gastric / дренажная игла для желудочного баллона

Fabricant / Manufacturer / Hersteller / Fabricante / Produttore / Fabrikant / Producent / Maker / Valmistaja / Fabricante / الشركة المصنعة / V yrobce / Proizvođa  / Produsen / Produc tor / Изготовитель:

CE0459

PRINCE MEDICAL, Zac De la Sente du Moulin, 60530 Ercuis - France  
T l. : + 33 (0)3 44 26 54 76 – quality.pm@omerin.com

PRINCE MEDICAL  
DEVICES

Ann e de marquage CE / Year of CE mark / Jahr des CE-Zeichens / A o de marcado CE / Anno di marchio CE / Jaar van CE-markering /  r for CE-m rke /  r for CE-merking / CE-merkin vuosi / Ano da marca o CE / CE  sne  lamate / Rok zisk n  ozna eni CE / Godina oznake CE / Tahun tanda CE / Anul marcatului CE / Год нанесения маркировки CE : 2008

CH REP

NMS Bio Medical S.A., Impasse des Saules 37, CH-1724 Le Mouret. Tel. 026 413 01 50

**1 – DESCRIPTION**

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**1 – DESCRIPTION**

The extraction forceps include:

- A white handle in styrene butadiene copolymer (SBC).
- A blue mobile slider in SBC.
- A polyvinyl chloride (PVC) antikink.
- A stainless steel spring sheath.
- A fluorinated ethylene propylene (FEP) sheath.
- A stainless steel cable.
- 2 pincers in stainless steel.

**2 – INDICATIONS**

The extraction forceps are used for the extraction of intragastric balloon after emptying (realized with the drainage needle AS2640718 of PRINCE MEDICAL range).

**\*DO NOT USE IT FOR ANY OTHER PURPOSES\***

The extraction forceps are delivered STERILE and are for SINGLE USE only.

**3 – CONTRAINDICATIONS**

The extraction forceps must not be used in case of:

- Allergy to one of the components.
- Gastric and / or esophageal lesion.
- Inability to receive endoscopic treatment.

**4 – POSSIBLE UNDESIRABLE SIDE EFFECTS**

Like any invasive medical device insertion related to a body orifice, the use of the extraction forceps may induce the following undesirable side effects (non-exhaustive list):

- Symptoms related to allergy to any of the components.
- Symptoms related to the endoscopic technique.

**5 – RECOMMENDATIONS**

- Read the instructions in their entirety before use.
- This medical device must only be used by trained and qualified medical staff or authorised person, familiar with endoscopy techniques and gastroenterological for treatment against overweight and obesity techniques and acquainted with the pathologies in question as well as their possible complications (under the responsibility of a practitioner).
- **The extraction forceps must be used with a flexible endoscope.**
- The diameter of the operating channel of the endoscope must be at least 0.3 mm bigger than the sheath diameter of the used forceps.
- The progression of the extraction forceps should be carefully done by a series of gentle pushes of low amplitude to prevent the risk of pincers twisting and/or perforation of the endoscope.

**NEVER FORCE.** Check for any obstruction in the operating channel (to prevent tearing of the sheath...)

- It is **NECESSARY** to insert the extraction forceps in the endoscope with the pincers in retracted position in their protective sheath to avoid damaging the endoscope or perforating the patient's mucous membranes.
- The position of the mobile slider allows to distinguish the position of the pincers. When the slider is at the top, the pincers are closed; the slider is at the bottom, the pincers are deployed.
- Deployment of the pincers is done with a constant force for easy sliding.

**DO NOT USE** the device if the packaging has been opened or damaged or if it shows any defect caused by transport or incorrect storage conditions or incorrect handling which might adversely affect its use.

**6 – STORAGE CONDITIONS**

In order to prevent any damage, products should be stored in their original packaging, in a cold and dry place.

**7 – PREPARATION OF MATERIALS – USE**

- Make sure that the packaging is intact and the shelf-life has not expired as these aspects guarantee the sterility and the quality of the product.

- To respect hygiene and safety rules, the wearing of gloves is recommended.
- Open the sachet by peeling off the paper.
- Aseptically remove the medical device from the double packaging.
- Make sure that the device is intact before use.
- Check that the pincers are fully withdrawn inside the sheath before inserting the extraction forceps into the endoscope.
- When the endoscope is in place, after intragastric balloon emptying, insert the extraction forceps into the operating channel of the endoscope.
- When the forceps are in the endoscope's field of vision, open the extraction forceps' pincers pushing the handle slider.
- **Sudden opening or opening without full visualisation of the pincers may cause mucosal trauma.**
- After bringing the balloon to the middle of the pincers, close them by pulling on the slider of the handle until the balloon is fully held.
- Bring the pincers against the operating channel of the endoscope and remove the "endoscope + extraction forceps holding the balloon" assembly from the patient.
- **Never remove the recovered balloon, whatever its size, through the operating channel of the endoscope.**

- Release the balloon from the pincers, retract them into the sheath and take the medical device out of the operating channel of the endoscope.

**\*HANDLE THE EXTRACTION FORCEPS WITH CARE\***

## 8 – DISPOSAL OF THE PRODUCT AFTER USE

After use, dispose of the product and its packaging according to your institution's standard operating procedure on bio-hazardous material and waste handling.

*We remind you that these products are strictly disposable devices which cannot be reused or resterilized.*

*In case of reuse or resterilization, their performances could be affected with an important risk of contamination.*

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