Edwards

KONECT RESILIA Aortic Valved Conduit, Model 11060A Accessory Care and Sterilization Instructions for Tray Model TRAY1190

Directory	
English (en).....	1
Français (fr).....	5
Deutsch (de).....	8
Español (es).....	11
Italiano (it).....	14
Nederlands (nl).....	17
Dansk (da).....	20
Svenska (sv).....	23
Ελληνικά (el).....	26
Português (pt).....	29
Česky (cs).....	32
Magyar (hu).....	35
Polski (pl).....	38
Slovensky (sk).....	41
Norsk (no).....	44
Suomi (fi).....	47
Български (bg).....	50
Română (ro).....	53
Eesti (et).....	56
Lietuvių (lt).....	59
Latviešu (lv).....	62
Türkçe (tr).....	65
Hrvatski (hr).....	68
Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Legenda simbola.....	71

English

KONECT RESILIA Aortic Valved Conduit, Model 11060A Accessory Care and Sterilization Instructions for Tray Model TRAY1190

Rx only

1.0 Tray Model TRAY1190

The Instructions For Use supplied with KONECT RESILIA aortic valved conduit (AVC), Model 11060A provide information on the use of this tray.

2.0 Device Description

TRAY1190 is an accessory to KONECT RESILIA aortic valved conduit (AVC), Model 11060A. It does not have direct patient contact. The product is comprised of a tray base, lid, inserts and latch and latch pin. The material for tray base, lid and inserts

Edwards, Edwards Lifesciences, the stylized E logo, KONECT, and RESILIA are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

is made out of Radel polyphenylsulfone. The lid is secured onto the tray by means of two stainless steel latches and latch pins. This device is reusable, supplied non-sterile, and must be cleaned and sterilized prior to its use. Tray holds the sizers during sterilization, storage, and transfer to the surgical site.

3.0 Intended Use and Indications for Use

The Tray Model TRAY1190 is intended to hold the six Sizers Model 1190 during sterilization and storage and allows for multiple sterilization cycles (re-use) by steam (moist heat).

TRAY1190 is designed to contain, store and hold Sizer Model 1190 during product sterilization and transfer.

4.0 Contraindications

There are no known contraindications with the use of TRAY1190.

5.0 How Supplied

The tray is provided non-sterile and must be cleaned and sterilized before each use. It cannot be sterilized in its original packaging.

6.0 Intended User

Primary users are staff responsible for preparation of KONECT RESILIA aortic valved conduit devices prior to implant; operating room (OR) nurses and the cardiac surgeons who perform valve sizing and implantation during valve and/or aortic root repair/replacement procedures.

7.0 Cleaning Instructions

All trays and lids must be cleaned and sterilized prior to each use in accordance with appropriate hospital procedures and the following recommended instructions and parameters.

Accessories should be retired from use when signs of cracking, crazing, or other deterioration is evident. Contact your Edwards Lifesciences sales representative to obtain appropriate replacements.

CAUTION: Examine tray for signs of wear, such as dullness, cracking or crazing. Replace tray if any deterioration is observed.

7.1 Instructions for Manual Cleaning:

1. Rinse the tray under cold, running tap water to remove gross soil.
2. Prepare a neutral pH cleaning solution (such as Enzol or equivalent) per the manufacturer's recommendation.
3. Fully immerse the tray in the cleaning solution and soak for a minimum of 5 minutes.
4. After the 5-minute soak time, use a soft-bristled brush and scrub the tray for a minimum of 2 minutes to remove any visible soil.
5. Remove the tray from the cleaning solution and rinse with cold, running tap water, then in reverse osmosis/deionized water for a minimum of 1 minute.
6. If soil is visible, repeat cleaning process.

7.2 Instructions for Automated Cleaning:

1. Rinse the tray under cold, running tap water to remove gross soil.
2. Prepare a neutral pH cleaning solution (such as Enzol or equivalent) per the manufacturer's recommendation.
3. Fully immerse the tray in the cleaning solution and soak for a minimum of 1 minute. Avoid contact between the sizers.
4. Load the tray with instruments as indicated by the tray markings.
5. Place the loaded tray into the automated washer with the lid detached.
6. Set motor speed to HIGH, if applicable, and select the cycle and cycle parameters per the following recommendations:

PHASE	RECIRCULATION TIME (MINUTES)	TEMPERATURE	DETERGENT TYPE AND CONCENTRATION
Pre-wash	2:00	Cold tap water	N/A
*Enzyme wash	2:00	Hot tap water	Neutral pH enzymatic detergent (Enzol or equivalent)
Wash	2:00	50-60 °C (122-140 °F)	Alkaline detergent (Neodisher Mediclean Forte or equivalent)

PHASE	RECIRCULATION TIME (MINUTES)	TEMPERATURE	DETERGENT TYPE AND CONCENTRATION
Rinse	1:00	60-70 °C (140-158 °F)	N/A
Drying	7:00	90-115 °C (194-239 °F)	N/A

*The enzyme wash phase is optional and may vary based on the washer model.

8.0 Disinfection

8.1 Instructions for Thermal Disinfection:

1. Load tray with sizers into washer/disinfector with the lid detached.
2. Set the motor speed to HIGH, if applicable, and select the following recommended parameters:

PHASE	RECIRCULATION TIME (MINUTES)	WATER TEMPERATURE	DETERGENT TYPE AND CONCENTRATION
Thermal rinse	5:00	90 °C (194 °F)	N/A

8.2 Instructions for High Level Disinfection:

1. Equilibrate a bath of high-level disinfectant (Cidex OPA or equivalent) to a minimum of 20 °C (68 °F), or follow manufacturer's recommendations.
2. Fully immerse the tray in the disinfectant solution, and ensure all air bubbles are removed from the surface of the tray by wiping with a sterile lint-free cloth.
3. Allow the tray to soak for a minimum of 12 minutes.
4. Thoroughly rinse the tray by fully immersing in sterile purified water or equivalent, agitating, and soaking for a minimum of 1 minute.
5. Dry the tray using a sterile lint-free cloth.

9.0 Sterilization Instructions

CAUTION: Do not sterilize any of the accessories in their original packaging. Each institution should use procedures that include biological indicators to determine the effectiveness of the sterilization procedure.

CAUTION: Do not stack trays during sterilization.

The user is responsible for the qualification of any deviations from the recommended method of cleaning or sterilization.

Sterilizer Type	Configuration	Preconditioning	Minimum Temperature	Exposure Time	Minimum Dry Time
Gravity Displacement	Double Wrapped	N/A	132-137 °C (270-279 °F)	15-18 minutes	40 minutes
Pre-vacuum	Double Wrapped	4 pulses	132-137 °C (270-279 °F)	3-18 minutes	30 minutes

Examine all instruments prior to use. If particles are visible, repeat the cleaning and sterilization process.

10.0 Device Disposal

Devices that have signs of wear, cracking or crazing may be handled and disposed of in the same manner as hospital waste and biohazardous materials. There are no special risks related to the disposal of these devices.

11.0 Basic Unique Device Identification-Device Identifier (UDI-DI)

The Basic UDI-DI is the access key for device-related information entered in the Eudamed.

The following table contains the Basic UDI-DI:

Product	TRAY
Models	TRAY1190
Basic UDI-DI	0690103D002OBT000UR

For a patient/user/third party in the European Economic area; if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and your national competent authority, which can be found at http://ec.europa.eu/growth/sectors/medical-devices/contacts_en.

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

	ISO Reg. No. ¹	Italiano	Nederlands	Dansk	Svenska
#	N/A	Quantità	Hoeveelheid	Mængde	Antal
LOT	2492	Numero di lotto	Lotnummer	Partinummer	Lotnummer
Rx only	N/A	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	1641	Consultare le istruzioni per l'uso sul sito Web	Raadpleeg de gebruiksaanwijzing op de website	Se brugsanvisningen på webstedet	Se bruksanvisningen på webbplatsen
	0434A	Attenzione	Let op	Forsigtig	Var försiktig
UDI	N/A	Identificatore univoco del dispositivo	Unieke hulpmiddelidentificatie	Unik udstyrsidentifikation	Unik produktidentifisering
	2606	Non utilizzare se la confezione è danneggiata	Niet gebruiken als de verpakking is beschadigd	Må ikke anvendes, hvis emballagen er beskadiget	Använd inte om förpackningen är skadad
EC REP	N/A	Rappresentante autorizzato nella Comunità Europea	Geautoriseerde vertegenwoordiger in de Europese Gemeenschap	Autoriseret repræsentant i Det Europæiske Fællesskab	Auktoriserad representant inom Europeiska gemenskapen
	3082	Produttore	Fabrikant	Producent	Tillverkare
	2497	Data di produzione	Fabricagedatum	Fremstillingsdato	Tillverkningsdatum
	2609	Non sterile	Niet-steriel	Ikke-steril	Icke-steril
MD	N/A	Dispositivo medico	Medisch hulpmiddel	Medicinsk udstyr	Medicinteknisk produkt
	N/A	Importatore	Importeur	Importør	Importör

	ISO Reg. No. ¹	Ελληνικά	Português	Česky	Magyar
REF	2493	Αριθμός καταλόγου	Número de catálogo	Katalogové číslo	Katalógusszám
#	N/A	Ποσότητα	Quantidade	Množství	Mennyiség

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




	ISO Reg. No. ¹	Polski	Slovensky	Norsk	Suomi
Rx only	N/A	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	1641	Należy zapoznać się z instrukcją użycia na stronie internetowej	Prečítajte si návod na použitie na webovej stránke	Se bruksanvisningen på nettstedet	Katso käyttöohjeet verkkosivustolta
	0434A	Przestroga	Upozornenie	Forsiktig	Tärkeä huomautus
	N/A	Unikalny identyfikator wyrobu medycznego	Jedinečný identifikátor zariadenia	Unik utstyrsidentifikator	Laitteen yksilöivä tunnus
	2606	Nie stosować, jeśli opakowanie jest uszkodzone	Nepoužívajte, ak je obal poškodený	Skal ikke brukes hvis pakningen er skadet	Älä käytä, jos pakkaus on vahingoittunut
	N/A	Autoryzowany przedstawiciel we Wspólnocie Europejskiej	Autorizovaný zástupca v krajinách Európskeho spoločenstva	Autorisert representant i Det europeiske fellesskap	Valtuutettu edustaja Euroopan yhteisössä
	3082	Producent	Výrobca	Produsent	Valmistaja
	2497	Data produkcji	Dátum výroby	Produksjonsdato	Valmistuspäivämäärä
	2609	Niejałowy	Nesterilné	Usteril	Epästeriili
	N/A	Wyrób medyczny	Zdravotnícka pomôcka	Medisinsk utstyr	Lääkinnällinen laite
	N/A	Importer	Dovozca	Importør	Maahantuoja


	ISO Reg. No. ¹	Български	Română	Eesti	Lietuvių
	2493	Каталожен номер	Număr de catalog	Kataloogi number	Katalogo numeris
	N/A	Количество	Cantitate	Kogus	Kiekis
	2492	Партиден номер	Număr de lot	Partii number	Partijos numeris
Rx only	N/A	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli
 ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos
 ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien
 selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu
 skaidrojums ■ Sembol Açıklamaları ■ Legenda simbola

	ISO Reg. No. ¹	Български	Română	Eesti	Lietuvių
	1641	Направете справка с инструкциите за употреба в уебсайта	Consultați instrucțiunile de utilizare de pe site-ul web	Palun lugege kasutusjuhiseid veebisaidilt	Žr. naudojimo instrukcijas interneto svetainėje
	0434A	Внимание	Atenție	Ettevaatust!	Perspėjimas
	N/A	Уникален идентификатор на изделието	Identificator unic al dispozitivului	Seadme kordumatu identifitseerimistunnus	Unikalasis priemonės identifikatorius
	2606	Да не се използва, ако опаковката е повредена	A nu se utiliza dacă ambalajul este deteriorat	Ärge kasutage, kui pakend on kahjustatud	Nenaudoti, jei pakuotė pažeista
	N/A	Упълномощен представител в Европейската общност	Reprezentant autorizat în Comunitatea Europeană	Volitatud esindaja Euroopa Ühenduses	Igaliotasis atstovas Europos Bendrijoje
	3082	Производител	Producător	Tootja	Gamintojas
	2497	Дата на производство	Data fabricației	Tootmiskuupäev	Pagaminimo data
	2609	Нестерилно	Nesteril	Mittesteriilne	Nesterilus
	N/A	Медицинско изделие	Dispozitiv medical	Meditsiiniseade	Medicinos priemonė
	N/A	Вносител	Importator	Importija	Importuotojas

	ISO Reg. No. ¹	Latviešu	Türkçe	Hrvatski
	2493	Kataloga numurs	Katalog Numarası	Kataloški broj
	N/A	Daudzums	Miktar	Količina
	2492	Partijas numurs	Lot Numarası	Broj serije
Rx only	N/A	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	1641	Skatīt lietošanas instrukciju tīmekļa vietnē	Web sitesindeki kullanım talimatlarına başvurun	Pogledajte upute za upotrebu na internetskoj stranici

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Legenda simbola

	ISO Reg. No. ¹	Latviešu	Türkçe	Hrvatski
	0434A	Uzmanību!	Dikkat	Oprez
	N/A	Unikāls ierīces identifikators	Benzersiz Cihaz Tanımlayıcısı	Jedinstveni identifikator proizvoda
	2606	Nelietot, ja iepakojums ir bojāts	Ambalaj hasarlıysa kullanmayın	Nemojte upotrebljavati ako je pakiranje oštećeno
	N/A	Pilnvarotais pārstāvis Eiropas Kopienā	Avrupa Topluluğu'nda Yetkili Temsilci	Ovlašteni predstavnik u Europskoj zajednici
	3082	Ražotājs	Üretici	Proizvođač
	2497	Izgatavošanas datums	Üretim tarihi	Datum proizvodnje
	2609	Nesterils	Steril değildir	Nesterilno
	N/A	Medicīniska ierīce	Tıbbi cihaz	Medicinski uređaj
	N/A	Importētājs	İthalatçı	Uvoznik

¹ ISO 7000 Graphical symbols for use on medical equipment - Registered symbols ■ ¹ ISO 7000 Symboles graphiques à utiliser sur l'équipement médical — Symboles enregistrés ■ ¹ ISO 7000 Grafische Symbole zur Verwendung auf Medizinprodukten – Registrierte Symbole ■ ¹ ISO 7000 Símbolos gráficos de uso en equipos médicos. Símbolos registrados ■ ¹ ISO 7000 Segni grafici utilizzabili sulle apparecchiature - Simboli registrati ■ ¹ ISO 7000 Grafische symbolen voor gebruik op medische apparatuur – Geregistreerde symbolen ■ ¹ ISO 7000 Grafiske symboler til brug på medicinsk udstyr – Registrerede symboler ■ ¹ ISO 7000 Grafiska symboler för användning på medicinsk utrustning – registrerade symboler ■ ¹ ISO 7000 Γραφικά σύμβολα για χρήση σε ιατροτεχνολογικό εξοπλισμό - Καταχωρισμένα σύμβολα ■ ¹ ISO 7000 Símbolos gráficos para utilização em equipamento médico - Símbolos registados ■ ¹ ISO 7000 Grafické symboly pro použití na zdravotnickém vybavení – registrované symboly ■ ¹ ISO 7000 Grafikus szimbólumok orvostechikai eszközökön való használatra – Regisztrált szimbólumok ■ ¹ ISO 7000 Symbole graficzne umieszczone na urządzeniach medycznych — Zarejestrowane symbole ■ ¹ ISO 7000 Grafické symboly používané na zdravotnických pomôckach – registrované symboly ■ ¹ ISO 7000 Grafiske symboler til bruk på medisinsk utstyr – registrerte symboler ■ ¹ ISO 7000 -standardin mukaiset lääkinnällisissä laitteissa käytettävät graafiset symbolit – rekisteröidyt symbolit ■ ¹ ISO 7000 Графични символи при използване на медицинска апаратура – Регистрирани символи ■ ¹ ISO 7000 Simboluri grafice pentru utilizare pe echipamente medicale – Simboluri înregistrate ■ ¹ ISO 7000 Meditsiiniseadmetel kasutatavad graafilised sümbolid – registreeritud sümbolid ■ ¹ ISO 7000 grafiniai simboliai, skirti naudoti ant medicinos įrangos – registruotieji simboliai ■ ¹ ISO 7000 grafiskie simboli izmantošanai uz medicīniskā aprīkojuma — reģistrētie simboli ■ ¹ ISO 7000 Tıbbi ekipmanlarda kullanılan grafik semboller — Tescilli semboller ■ ¹ ISO 7000 Grafički simboli koji se upotrebljavaju na medicinskoj opremi – Registrirani simboli

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2025-06
10046356002 A
© Copyright 2024, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA
Made in USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU

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Edwards Lifesciences GmbH