INTENDED PURPOSE

The INDIGO Aspiration System is designed to remove thrombus from the vasculature using mechanical aspiration.

DEVICE DESCRIPTION

The INDIGO® Aspiration System is comprised of several devices:

- INDIGO Aspiration Catheter
- Penumbra Aspiration Pump
- INDIGO Aspiration Pump Canister
- INDIGO Aspiration Tubing
- INDIGO Separator™

The INDIGO Aspiration Catheter targets aspiration from the Penumbra Aspiration Pump directly to the thrombus. The INDIGO Separator may be used to clear the lumen of the INDIGO Aspiration Catheter should it become blocked with thrombus. The INDIGO Aspiration Catheter is introduced through a guide catheter or vascular sheath into the peripheral vasculature and guided over a guidewire to the site of the primary occlusion. The INDIGO Aspiration Catheter is used with the Penumbra Aspiration Pump to aspirate thrombus from an occluded vessel. As needed, an INDIGO Separator may be deployed from the INDIGO Aspiration Catheter to assist with thrombus removal. The INDIGO Separator is advanced and retracted through the INDIGO Aspiration Catheter at the proximal margin of the primary occlusion to facilitate clearing of the thrombus from the INDIGO Aspiration Catheter tip. The devices are visible under fluoroscopy. For the aspiration source, the INDIGO Aspiration Catheter is used in conjunction with the Penumbra Aspiration Pump, which is connected using the INDIGO Aspiration Tubing and the INDIGO Aspiration Pump Canister. The INDIGO Aspiration Catheter may be provided with a steam shaping mandrel and rotating hemostasis valve, and introducer. The INDIGO Separator may be provided with an introducer and torque device.

INDIGO System Aspiration Catheters have a hydrophilic coating on the distal segment of the catheter shaft (see **Table 1**)

Table 1 – Hydrophilic Coating Length	
INDIGO Aspiration Catheter	Length (cm)
CAT3	75
CAT5	30
CAT6	11
CAT7/ CAT7D	14
CAT8/ CATD	7
CAT12	22

INDICATION FOR USE

INDIGO Aspiration Catheters and Separators

As part of the INDIGO Aspiration System, the INDIGO Aspiration Catheters and Separators are indicated for the removal of fresh, soft emboli and thrombi from vessels of the peripheral arterial and venous systems, and for the treatment of pulmonary embolism.

INDIGO Aspiration Tubing

As part of the INDIGO Aspiration System, the INDIGO Sterile Aspiration Tubing is indicated to connect the INDIGO Aspiration Catheters to the Penumbra Aspiration Pump.

Penumbra Aspiration Pump

The Penumbra Aspiration Pump is indicated as a vacuum source for Penumbra Aspiration Systems.

INTENDED PATIENT POPULATION

The intended patient population for the INDIGO Aspiration System is patients with soft emboli or thrombi in peripheral arteries and veins, and/or for patients with pulmonary embolism.

CONTRAINDICATIONS

Not for use in the coronaries or the neurovasculature.

WARNINGS

- The INDIGO Aspiration System should only be used by physicians who have received appropriate training in interventional techniques.
- Do not advance, retract, or use any component of the INDIGO System against resistance without careful assessment of the cause using fluoroscopy. If the cause cannot be determined, withdraw the device or system as a unit. Unrestrained torquing or forced insertion of the catheter or separator against resistance may result in damage to the device or vessel.
- Do not use the INDIGO Aspiration System with a pump other than the Penumbra Aspiration Pump.
- Placing guidewire too distal in the pulmonary vasculature or excessive manipulation of aspiration/guiding catheter in the smaller, peripheral, and segmental pulmonary artery branches can result in vessel perforation.

PRECAUTIONS

- The device is intended for single use only. Do not resterilize or reuse. Resterilization or reuse of the device may compromise the structural integrity of the device and/or lead to device failure which in turn, may result in patient injury, illness or death.
- Do not use kinked or damaged devices. Do not use open or damaged packages. Return all damaged devices and packaging to the manufacturer / distributor.
- Use prior to the “Use By” date.
- Use the INDIGO Aspiration System in conjunction with fluoroscopic visualization.
- Maintain a constant infusion of appropriate flush solution.
- When performing aspiration, ensure that the INDIGO Aspiration Tubing is open for only the minimum time needed to remove thrombus. Excessive aspiration or failure to close the INDIGO Aspiration Tubing when aspiration is complete is not recommended.
- Hemoglobin and hematocrit levels should be monitored in patients with >700 mL blood loss from the clot aspiration procedure.
- The INDIGO Separator is not intended for use as a guidewire. If repositioning of the INDIGO Aspiration Catheter is necessary during the revascularization procedure, such repositioning should be performed over an appropriate guidewire using standard catheter and guidewire techniques.
- Do not use automated high-pressure contrast injection equipment with the INDIGO Aspiration Catheter because it may damage the device.
- The INDIGO Separator (SEP3) contains CMR substances (CAS Number: 7440-48-4 and 1344-37-2). The INDIGO Separators (SEP5, and SEP6) contain a CMR substance (CAS Number: 7440-48-4).

INDIGO POTENTIAL ADVERSE EVENTS

Possible complications include, but are not limited to, the following:

- acute vessel occlusion
- air embolism
- allergic reaction and anaphylaxis from contrast media or device material
- anemia
- arrhythmia
- arteriovenous fistula
- cardiac injury, cardiac perforation, cardiac tamponade
- cardio-respiratory arrest
- compartment syndrome
- death
- emboli
- emergent surgery
- foreign body embolization
- hematoma or hemorrhage at access site
- hemoptysis
- hemorrhage
- hypotension/hypertension
- infarction leading to organ damage
- infection
- ischemia
- myocardial infarction
- neurological deficits including stroke
- pneumothorax
- pseudoaneurysm
- renal impairment or acute renal failure from contrast media
- residual thrombus due to inability to completely remove thrombus or control blood flow
- respiratory failure
- valvular damage
- vessel spasm, thrombosis, dissection (intimal disruption), or perforation

*If a serious incident related to the device occurs, contact your Penumbra, Inc. representative and the competent authority in your respective country/region.

PROCEDURE

1. Refer to **Warnings, Precautions, and Potential Adverse Events** prior to use.
2. Refer to **Table 2** for device compatibility information.
3. As each device of the INDIGO Aspiration System is used, remove the device from the packaging and inspect for damage or for kinks.
4. Prepare INDIGO Aspiration System devices for use by flushing the packaging and device with heparinized saline (includes INDIGO Aspiration Catheter and INDIGO Separator only).
5. Prepare a guide catheter or vascular sheath according to the manufacturer's *Instructions for Use*.
6. Using conventional catheterization techniques under fluoroscopic guidance, place the guide catheter or vascular sheath into the appropriate artery or vein that is proximal to the thrombus occlusion site over an appropriate guidewire.

Table 2 – Device Compatibility and Selection Information			
INDIGO Aspiration Catheter	Vascular Sheath / Guide Catheter	Guidewire	INDIGO Separator
CAT3	5F / 7F	0.014"	SEP3
CAT5	6F / 8F	0.014"	SEP5
CAT6	6F / 8F	0.038"	SEP6
CAT7	7F / 9F	0.038"	SEP7
CAT7D	7F / 9F	0.038"	SEP7D
CAT8	8F / 10F	0.038"	SEP8
CATD	8F / 10F	0.038"	SEPD
CAT12	12F / 14F	0.038"	SEP12

The INDIGO Separator is intended to be used only as part of the INDIGO Aspiration System

INDIGO ASPIRATION CATHETER PREPARATION AND USE

1. Confirm vessel diameter and select an appropriately sized INDIGO Aspiration Catheter.
2. If provided with a shaping mandrel, the INDIGO Aspiration Catheter tip may be steam shaped. Insert the shaping mandrel into the INDIGO Aspiration Catheter tip, gently form the desired shape, and position in a flow of steam for ten seconds while rotating gently. To avoid damage from overheating, do not place the catheter tip directly at the steam outlet. Allow tip to cool prior to removing mandrel. Inspect INDIGO Aspiration Catheter for any damage that may have resulted from steam shaping; If any damage is found, do not use the INDIGO Aspiration Catheter.
3. Attach the rotating hemostasis valve provided with the INDIGO Aspiration Catheter.
4. If using a guide catheter, insert the INDIGO Aspiration Catheter into a rotating hemostasis valve connected to the proximal hub of the guide catheter. If a vascular sheath is being used, insert the INDIGO Aspiration Catheter through the valve of the vascular sheath using the introducer (if provided). After inserting the INDIGO Aspiration Catheter, retract the introducer.
5. Advance the INDIGO Aspiration Catheter into the target vessel over the guidewire. Position the INDIGO Aspiration Catheter proximal to the thrombus. Remove the guidewire from the INDIGO Aspiration Catheter.


NOTE: Access to the site of occlusion may be aided by the coaxial use of a compatible inner catheter.

6. Attach the INDIGO Aspiration Tubing to the Penumbra Aspiration Pump Canister and turn on the Penumbra Aspiration Pump (Refer to the Penumbra Aspiration Pump *Operation Manual*). Confirm that the Penumbra Aspiration Pump reaches maximum vacuum level. Ensure that the INDIGO Aspiration Tubing is switched to the OFF position.

INDIGO SEPARATOR PREPARATION AND USE


1. Select an appropriate size INDIGO Separator corresponding to the Aspiration Catheter (see Table 2).
2. Advance the INDIGO Aspiration Catheter to gently embed the distal tip into the thrombus.
3. Insert the INDIGO Separator through the rotating hemostasis valve and into the proximal hub of the INDIGO Aspiration Catheter. If provided, an introducer may be used.
4. Advance the INDIGO Separator through the Aspiration Catheter until the Separator cone is located at the distal tip of the Aspiration Catheter.
5. Connect the Aspiration Tubing to the side port of the rotating hemostasis valve as shown in **Figure 1**.

Figure 1 – Assembled Aspiration System with INDIGO Separator

6. To begin aspiration, switch the valve on the INDIGO Aspiration Tubing to the ON position and advance/retract the INDIGO Separator to assist with aspiration and removal of the thrombus as shown in **Figure 2**.

Figure 2 – Thrombus Aspiration Using the INDIGO Separator

7. To stop aspiration, switch the valve on the INDIGO Aspiration Tubing to the OFF position and turn off the Penumbra Aspiration Pump.
8. Remove the INDIGO Separator.
9. If necessary, additional INDIGO Separators and Aspiration Catheters may be used to further remove thrombus at the discretion of the physician.
10. Using a 5cc or larger syringe, aspirate approximately 5cc of blood from the INDIGO Aspiration Catheter to remove any thrombus that may remain in the catheter.
11. Obtain a post-treatment angiogram by injecting contrast media through the guide catheter or vascular sheath.

CLINICAL EXPERIENCE:**Penumbra EXTRACT-PE Trial Clinical Data Summary:**

The Penumbra EXTRACT-PE trial was a prospective, multicenter, single arm trial to determine the safety and efficacy of the Indigo® Aspiration System for mechanical thrombectomy in subjects with acute pulmonary embolism (PE). The trial enrolled 119 subjects at 22 centers in the U.S. An independent imaging core lab and a clinical events committee (CEC) reviewed safety endpoints data.

INCLUSION CRITERIA:

1. Clinical signs and symptoms consistent with acute PE with duration of 14 days or less. Evidence of PE must be from CTA.
2. Systolic BP ≥ 90 mmHg with evidence of dilated RV with an RV/LV ratio > 0.9
3. Patient is 18 years of age or older

EXCLUSION CRITERIA:

1. tPA use within 14 days prior to baseline CTA
2. Systolic BP < 90 mmHg for 15 mins or the requirement of inotropic support to maintain systolic BP ≥ 90 mmHg
3. Pulmonary hypertension with peak PA > 70 mmHg by right heart catheterization

- History of severe or chronic pulmonary hypertension
- FiO2 requirement >40% or >6 LPM to keep oxygen saturation >90%
- Hematocrit <28%
- Platelets <100,000/ μ L
- Serum creatinine >1.8 mg/dL
- INR >3
- aPTT (or PTT) >50 seconds on no anticoagulation
- History of heparin-induced thrombocytopenia (HIT)
- Contraindication to systemic or therapeutic doses of anticoagulants
- Major trauma <14 days
- Presence of intracardiac lead
- Cardiovascular or pulmonary surgery within last 7 days
- Cancer requiring active chemotherapy
- Known serious, uncontrolled sensitivity to radiographic agents
- Life expectancy <90 days
- Female who is pregnant
- Intracardiac Thrombus
- Patients on ECMO
- Current participation in another investigational study

RESULTS:

There were 119 subjects that were enrolled and included in the Intent-To-Treat (ITT) analysis population. The Modified ITT (mITT) population was a subset of the ITT population and was the primary population for all efficacy parameters. The mITT population excluded subjects who received adjunctive treatments or thrombolytics intra-procedure through 48 hours post-procedure for the purpose of reducing clot burden in the pulmonary artery. There were 110 subjects in the mITT population. Safety endpoints are evaluated based on ITT population.

There were 110 patients who completed the study 30 day follow-up, 3 patients died, 1 withdrew from the study, and 5 were lost to follow-up.

DEMOGRAPHICS AND MEDICAL HISTORY:

Of the 119 subjects, the median age was 59.8 years old (range 21, 92), 44.5% (53/119) were female, a previous deep vein thrombosis (DVT) was reported in 60.5% (72/119), any cancer in 22.7% (27/119), and a previous pulmonary embolism (PE) was reported in 17.6% (21/119) of subjects.

PROCEDURAL CHARACTERISTICS:

The study procedure was initiated in 119 subjects and completed in 118 subjects; one subject had the procedure aborted prior to aspiration. Conscious sedation was used in 97.5% of subjects and 78.2% had a right femoral access site. The median procedure time (venous puncture to Indigo device removal) was 66 min; IQR [46.0, 94.0]. The median time from the first Indigo device insertion to the last Indigo device removal was 37 minutes; IQR [23.5, 60.0]. The median ICU stay was 1.0 day [1.0, 2.0]. The treatment included the following pulmonary embolism locations:

Treatment Location	All Subjects (N=119)
Main PA with Unilateral	0.8% (1/119)
Main PA with Bilateral	23.5% (28/119)
Unilateral Only	5.9% (7/119)
Bilateral Only	69.7% (83/119)

DEVICE USAGE:

There were 129 Penumbra aspiration catheters and 99 Penumbra separators used during the study procedure.

Penumbra Devices Used – By Subject	All Subjects (N=119)
Penumbra Catheters Used per Subject	
1	91.6% (109/119)
2	8.4% (10/119)
Catheters	
Aspiration Catheter 8, TORQ, 85 cm	13.4% (16/119)
Aspiration Catheter 8, XTORQ, 115 cm	83.2% (99/119)
Aspiration Catheter 8, XTORQ, 115 cm Kit	6.7% (8/119)
Penumbra Separators Used per Subject	
0	19.3% (23/119)
1	78.2% (93/119)
2	2.5% (3/119)
Separators (Separator 8, 150 cm)	80.7% (96/119)

PRIMARY EFFICACY ENDPOINT:

The primary efficacy endpoint was the reduction in RV/LV ratio from baseline to 48 hours assessed by CTA and evaluated by independent core laboratory. The lower limit of the 95% confidence interval (CI) of the change in RV/LV ratio at 48 hours was >0.20. The mean baseline RV/LV ratio was 1.46 (0.29) and it was reduced to 1.04 (0.17) at 48 hours. The endpoint was met with absolute reduction in RV/LV ratio of 0.42 $\pm$ 0.25 (95%CI 0.37, 0.46), representing a 26.9% reduction (p<0.0001).

Primary Effectiveness Endpoint	mITT (N=110)	95% CI ^[3]
Absolute Reduction in RV/LV Ratio from Baseline to 48 Hours ^[1]		
N ^[2]	106	
Mean (SD)	0.42 (0.25)	(0.37, 0.46)
Median [IQR]	0.41 [0.25, 0.54]	
Range (Min, Max)	-0.19, 1.23	

^[1] p <0.0001 based on a one-sample, two-sided asymptotic test with a null proportion of 0.20 and the standard error computed from the sample proportion.

^[2] Data were analyzed based on mITT with evaluable 48 hours CTA

^[3] Confidence intervals have not been adjusted for multiplicity

SAFETY

The primary safety endpoint was major adverse events, a composite of device-related death, major bleeding, and device-related SAEs (defined as clinical deterioration, pulmonary vascular injury, and cardiac Injury) within 48 hours. The safety hypothesis stated that the 48-hour rate of major adverse events would not equal 40%. The primary safety endpoint was met and the rate was 1.7% (95% CI 0.0%, 4.0%), p<0.0001). Two subjects (1.7%) experienced 3 events:

1 subject experienced Sustained Ventricular Tachycardia (Device-related Death/Major Bleeding/Clinical Deterioration/Pulmonary Vascular Injury) and Hemoptysis (Major Bleeding/Pulmonary Vascular Injury); 1 subject experienced a groin access site bleeding (Major Bleeding).

PRIMARY SAFETY ENDPOINT (As-TREATED POPULATION):

Primary Safety Analysis	All Subjects (N=119)	95% CI ^[2]
Composite of major adverse events within 48 Hours ^[1]	1.7% (2/119)	(0.0%, 4.0%)

^[1] p <0.0001 based on a one-sample, two-sided asymptotic test with a null proportion of 0.40 and the standard error computed from the sample proportion.

^[2] Confidence intervals have not been adjusted for multiplicity

SECONDARY ENDPOINTS (As-TREATED POPULATION):

Secondary safety endpoints included device-related death within 48 hours, major bleeding within 48 hours, clinical deterioration within 48 hours, pulmonary vascular injury within 48 hours, cardiac injury within 48 hours, any-cause mortality within 30 days, device-related SAEs within 30 days, and symptomatic

PE recurrence within 30 days. Serious Adverse Events (SAEs) and Adverse events (AEs) were collected through the 30-day follow-up. The key events are included below.

A total of three deaths occurred during the trial. One death was adjudicated as procedure and device related. Two additional non-procedure/device related deaths occurred due to progression of pre-existing diseases (cancer progression and a recent brain resection due to a previous ischemic stroke). A total of 7 subjects experienced 8 secondary endpoint events.

SECONDARY SAFETY ENDPOINTS:

Secondary Safety Analysis*	All Subjects (N=119)
Device-related Death within 48 Hours	0.8% (1/119)
Major Bleeding within 48 Hours	1.7% (2/119)
Clinical Deterioration within 48 Hours	1.7% (2/119)
Pulmonary Vascular Injury within 48 Hours	1.7% (2/119)
Cardiac Injury within 48 Hours	0.0%
Any-cause Mortality within 30 Days	2.5% (3/119)
Device-related SAEs within 30 Days	1.7% (2/119)
Symptomatic PE Recurrence within 30 Days	0.0%
Device and procedure related Adverse Events within 30 days from procedure*	
Procedure-related SAEs	2.5% (3/119)
Device-related SAEs	1.7% (2/119)
Procedure-related Non-Serious AEs	9.2% (11/119)
Device-related Non-Serious AEs	4.2% (5/119)

***One subject can meet several endpoint events and can appear in multiple categories**

CLINICAL BENEFITS AND PERFORMANCE CHARACTERISTICS

Clinical benefits of the INDIGO Aspiration System are the removal or reduction of thrombi and emboli from the arterial or venous peripheral vasculature or for patients with pulmonary embolism. The detailed listing of all known benefits is included in the Summary of Safety and Clinical Performance (SSCP).

SAFETY AND CLINICAL PERFORMANCE

A copy of the Summary of Safety and Clinical Performance can be viewed by searching the device name on the EUDAMED website:


<https://ec.europa.eu/tools/eudamed>

Basic UDI DI: 081454801INDIGOSYSTEMYF

STORAGE AND DISPOSAL

- Store in a cool, dry place.
- Dispose of device in accordance with hospital, administrative, and/or local government policy.

SYMBOLS GLOSSARY

		Attention, see instructions for use	
		Prescription only – US Federal Law restricts this device to use by or on the order of a physician	
		Nonpyrogenic	
		Sterilized using ethylene oxide . Single sterile barrier system with protective packaging outside.	
		Medical Device	
	Manufacturer		EU Representative
	Lot number		Do not reuse
	Use by		Catalogue number
	Do not resterilize		Not made with natural rubber latex
	Do not use if package is damaged.		Contains Hazardous Substances CAS Number: 7440-48-4 / 1344-37-2

WARRANTY

Penumbra, Inc. (Penumbra) warrants that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular use. Handling, storage, cleaning, and sterilization of this device as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond Penumbra's control directly affect the device and the results obtained from its use. Penumbra's obligation under this warranty is limited to the repair or replacement of this device and Penumbra shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this device. Penumbra neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. Penumbra assumes no liability with respect to devices reused, reprocessed, or resterilized and makes no warranties, expressed or implied, including but not limited to merchantability or fitness for intended use, with respect to such device.

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