INTENDED PURPOSE

Lightning™ Aspiration Tubing is a sterile aspiration tubing component of the INDIGO® Aspiration System and is designed to serve as a conduit to assist in thrombus removal and restoration of blood flow in the peripheral vasculature and for the treatment of pulmonary embolism.

DEVICE DESCRIPTION

Lightning Aspiration Tubing facilitates the transfer of vacuum between INDIGO Aspiration Catheter and Penumbra Aspiration Pump while providing intermittent or continuous aspiration. Intended users for this device are physicians who have received appropriate training in surgical procedures and/or interventional techniques.

INDICATION FOR USE

As part of the INDIGO Aspiration System, the INDIGO Sterile Aspiration Tubing is indicated to connect the INDIGO Aspiration Catheters to the Penumbra Aspiration Pump.

INTENDED PATIENT POPULATION

The intended patient population for the INDIGO Aspiration System is patients with soft emboli or thrombi in peripheral arteries and veins, and/or for patients with pulmonary embolism.

CONTRAINDICATIONS

There are no known contraindications.

WARNINGS

- Do not use the INDIGO Aspiration System with a pump other than a Penumbra Aspiration Pump.
- Use of Lightning Aspiration Tubing adjacent to other equipment should be avoided because it could result in improper operation. If such use is necessary, Lightning Aspiration Tubing and the other equipment should be observed to verify that they are functioning properly.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 12 inches (30 cm) to any part of Lightning Aspiration Tubing. Otherwise, this could result in degradation of the performance of this equipment.

PRECAUTIONS

- The device is intended for single use only. Do not resterilize or reuse. Resterilization or reuse of the device may compromise the structural integrity of the device and/or lead to device failure which in turn, may result in patient injury, illness or death.
- Do not use kinked or damaged devices. Do not use open or damaged packages. Return all damaged devices and packaging to the manufacturer/distributor.
- Use prior to the "Use By" date.
- When performing aspiration, ensure that the INDIGO Aspiration Tubing is open for only the minimum time needed to remove the thrombus. Excessive aspiration or failure to close the INDIGO Aspiration Tubing when aspiration is complete is not recommended.
- Do not use in the presence of a flammable anesthetic mixture with air or nitrous oxide.
- Do not use in oxygen rich environment.

POTENTIAL ADVERSE EVENTS

Possible complications include, but are not limited to, the following:

- acute vessel occlusion
- air embolism
- allergic reaction and anaphylaxis from contrast media or device material
- anemia
- arrhythmia
- arteriovenous fistula
- cardiac injury, cardiac perforation, cardiac tamponade
- cardio-respiratory arrest
- compartment syndrome
- death
- emboli
- emergent surgery
- foreign body embolization
- hematoma or hemorrhage at access site
- hemoptysis
- hemorrhage
- hypotension/hypertension
- infarction leading to organ damage
- infection
- ischemia
- myocardial infarction
- neurological deficits including stroke
- pneumothorax
- pseudoaneurysm
- renal impairment or acute renal failure from contrast media
- residual thrombus due to inability to completely remove thrombus or control blood flow
- respiratory failure
- valvular damage
- vessel spasm, thrombosis, dissection (intimal disruption), or perforation

*If a serious incident related to the device occurs, contact your Penumbra, Inc. representative and the competent authority in your respective country/region.


PROCEDURE

- Refer to Warnings, Precautions, and Potential Adverse Events prior to use.
- Refer to the product labeling for the size and maximum length of PVC tubing.
- Remove the device from the packaging and inspect for damage or for kinks.

INDIGO LIGHTNING ASPIRATION TUBING PREPARATION AND USE


- Attach Lightning Aspiration Tubing via the suction connector to a compatible canister on the Penumbra Aspiration Pump, plug the Lightning Aspiration Tubing power adapter into Pump Power Port, and turn on the Penumbra Aspiration Pump (refer to the Penumbra Aspiration Pump *Operation Manual*). Confirm that the Aspiration of the Penumbra Aspiration Pump reaches maximum vacuum level. Ensure the switch on Lightning Aspiration Tubing is switched to the OFF position.
- Connect Lightning Aspiration Tubing to the side port of the rotating hemostasis valve as shown in **Figure 1** (refer to the INDIGO Aspiration Catheter *Instructions for Use*).

Figure 1: Assembled Lightning Aspiration Tubing with INDIGO Aspiration System

- To begin aspiration, move the switch of Lightning Aspiration Tubing to the ON position, and use the INDIGO Aspiration Catheter according to the associated instructions for use.
- The lights on Lightning Aspiration Tubing will illuminate to indicate the following conditions:
 - Solid Green = Lightning Aspiration Tubing powered in OFF position, ready for use
 - Flashing Green = Lightning Aspiration Tubing powered in ON position
 - Flashing Yellow = Lightning Aspiration Tubing powered in ON position, and 20+ Seconds of Low/No Flow Detected
 - Solid Red = System Error
 - Flashing Red = No Vacuum Detected
- If a solid red light is observed, stop using the device and replace as needed.
- To stop aspiration, move the switch on Lightning Aspiration Tubing to the OFF position, and turn off the Penumbra Aspiration Pump.

Figure 2: Lightning Aspiration Tubing

INSTRUCTIONS FOR USE

LIGHTNING™
ASPIRATION TUBING

Penumbra 

INDIGO®
Aspiration System

OPERATING CONDITIONS

- Temperature: 65°F - 75°F (18°C - 24°C)
- Humidity: <75% RH
- Pressure: Sea Level – 6000 ft (0 m – 1828 m)

TRANSPORTATION CONDITIONS

- Temperature: 14°F - 131°F (-10°C - 55°C)
- Humidity: <85% RH

STORAGE

- Store in a cool, dry place.

MAINTENANCE PROCEDURE

Lightning Aspiration Tubing is a single use device and has no serviceable parts. For service or service questions on your Lightning Aspiration Tubing device, please call Penumbra Customer Service.

DISPOSAL PROCEDURES

Used Lightning Aspiration Tubing devices should be disposed of using standard hospital biohazard procedures.

APPROVALS

Conforms to

- AAMI Std. ES60601-1
- IEC Std. 60601-1-6
- ISO Std. 10079-1
- Certified to CSA Std. C22.2# 60601-1

Lightning Aspiration Tubing complies with all current international safety and EMC requirements for medical electrical equipment as follows:

- IEC 60601-1-2 including:
 - o CISPR 11 Class A
 - o EN/IEC 61000-3-2
 - o EN/IEC 61000-3-3
 - o EN/IEC 61000-4-2
 - o EN/IEC 61000-4-3
 - o EN/IEC 61000-4-4
 - o EN/IEC 61000-4-5
 - o EN/IEC 61000-4-6
 - o EN/IEC 61000-4-8
 - o EN/IEC 61000-4-11

LIGHTNING ASPIRATION TUBING ELECTROMAGNETIC COMPATIBILITY

Guidance and manufacturer's declaration – electromagnetic immunity			
Lightning Aspiration Tubing is intended for use in the electromagnetic environment specified below. The customer or the user of Lightning Aspiration Tubing should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Voltage Dips 30% reduction, 25/30 periods at 0°	Voltage Dips 30% reduction, 25/30 periods at 0°	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Lightning Aspiration Tubing requires continued operation during power mains interruptions, it is recommended that Lightning Aspiration Tubing be powered from an uninterruptible power supply or a battery.
	Voltage Dips > 95% reduction, 0.5 period at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	Voltage Dips > 95% reduction, 0.5 period at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	
	Voltage Dips > 95% reduction, 1 period at 0°	Voltage Dips > 95% reduction, 1 period at 0°	
	Voltage Interruptions > 95% reduction, 250/300 periods	Voltage Interruptions > 95% reduction, 250/300 periods	
(50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Guidance and manufacturer's declaration – electromagnetic immunity			
Lightning Aspiration Tubing is intended for use in the electromagnetic environment specified below. The customer or the user of Lightning Aspiration Tubing should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz (6 Vrms in ISM radio Bands within 150 kHz to 80 MHz)	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of Lightning Aspiration Tubing, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.7 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	3 V/m	
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which Lightning Aspiration Tubing is used exceeds the applicable RF compliance level above, Lightning Aspiration Tubing should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating Lightning Aspiration Tubing. ^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Recommended separation distances between portable and mobile RF communications equipment and Lightning Aspiration Tubing			
Lightning Aspiration Tubing is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of Lightning Aspiration Tubing can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and Lightning Aspiration Tubing as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.7 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Immunity to RF Wireless Communications Equipment						
Test Frequency (MHz)	Band ^{a)} (MHz)	Service ^{a)}	Modulation ^{b)}	Maximum Power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
385	380 – 390	TETRA 400	Pulse modulation ^{b)} 18 Hz	1.8	0.3	27
450	430 – 470	GMRS 460, FRS 460	FM ^{c)} ±5 kHz deviation 1 kHz sine	2	0.3	28
710	704 – 787	LTE Band 13, 17	Pulse modulation ^{b)} 217 Hz	0.2	0.3	9
745						
780						
810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation ^{b)} 18 Hz	2	0.3	28
870						
930						
1720	1700 – 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation ^{b)} 217 Hz	2	0.3	28
1845						
1970						
2450	2400 – 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^{b)} 217 Hz	2	0.3	28
5240	5100 – 5800	WLAN 802.11 a/n	Pulse modulation ^{b)} 217 Hz	0.2	0.3	9
5500						
5785						
a) For some services, only the uplink frequencies are included. b) The carrier shall be modulated using a 50% duty cycle square wave signal. c) As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.						

Equipment Classifications

Protection from electric shockClass II and externally powered
Degree of protection against electric shock.....Type BF Applied Part
Degree of protection against ingress of liquids or particlesIP21
Mode of Operation.....Continuous Operation
Flow Rate0-18 LPM

CLINICAL BENEFITS AND PERFORMANCE CHARACTERISTICS

Clinical benefits of the INDIGO Aspiration System are the removal or reduction of thrombi and emboli from the arterial or venous peripheral vasculature or for patients with pulmonary embolism. The detailed listing of all known benefits is included in the Summary of Safety and Clinical Performance (SSCP).

SAFETY AND CLINICAL PERFORMANCE

A copy of the Summary of Safety and Clinical Performance can be viewed by searching the device name on the EUDAMED website:
<https://ec.europa.eu/tools/eudamed>
Basic UDI DI: 081454801INDIGOSYSTEMYF

SYMBOL DEFINITIONS

	Manufacturer		European Union Authorized Representative
	Sterilized using ethylene oxide. Single sterile barrier system with protective packaging outside.		Catalogue Number
	Lot number		Medical Device
	Type BF Applied Part		Use by
	Do not reuse		Refer to Instruction Manual / Booklet
	Do not resterilize		Prescription Only – US Federal Law restricts this device to use by or on the order of a physician
	Not made with natural rubber latex		Nonpyrogenic
	High Vacuum/ Low Flow (Intermittent suction)		US and Canada Certification
	Do not use if package is damaged.		

WARRANTY

Penumbra, Inc. (Penumbra) warrants that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular use. Handling, storage, cleaning, and sterilization of this device as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond Penumbra’s control directly affect the device and the results obtained from its use. Penumbra’s obligation under this warranty is limited to the repair or replacement of this device and Penumbra shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this device. Penumbra neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. Penumbra assumes no liability with respect to devices reused, reprocessed, or re-sterilized and makes no warranties, expressed or implied, including but not limited to merchantability or fitness for intended use, with respect to such devices.

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