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INSTRUCTIONS FOR USE – FEMORAL HEADS

Before using a LimaCorporate product, the surgeon should carefully study the following recommendations, warnings and instructions, as well as the available product-specific information (e.g.: product literature, surgical technique).

1. PRODUCT INFORMATION

LimaCorporate femoral heads are indicated for both partial and total hip replacement. Femoral modular heads are used in total hip arthroplasty in combination with cemented or uncemented femoral and acetabular components. Cephalic and Lock Bipolar heads are indicated for partial hip replacement in combination with cemented or uncemented femoral stems.

BioloX@Delta revision heads are designed to allow surgeons to utilize a ceramic head in revision arthroplasty. Its use is recommended when a revision for a fractured ceramic head is required.

The main goal of a joint prosthesis is to reproduce the articular anatomy, in part or totally. The joint prosthesis is intended to reduce pain and give articular mobility to the patient. The degree of pain reduction and mobility depend, in part, on the pre-operative situation, intra-operative options and the post-op rehabilitation.

1.1. MATERIALS

Device	Material
Ceramic Femoral Modular Heads	BioloX@Delta
BioloX@Delta Revision Heads	BioloX@Delta + Ti6Al4V
Metallic Femoral Modular Heads	CoCrMo, AISI 316L, FeCrNiMnMoNbN
Cephalic Heads	AISI 316L
Lock Bipolar Heads	CoCrMo + UHMWPE, AISI 316L + UHMWPE

Applicable Standards

BioloX@Delta (ISO 6474-2) – **CoCrMo** (ISO 5832-4/ASTM F75, ISO 5832-12/ASTM F1537) – **AISI 316L** (ISO 5832-1) – **FeCrNiMnMoNbN** (ISO 5832-9) – **Ti6Al4V** (ISO 5832-3, ASTM F136) – **UHMWPE** (ISO 5834-2, ASTM F648)

Some patients may be hypersensitive or allergic to the implant materials; this should be appropriately considered by the surgeon.

1.2. HANDLING AND STORAGE

Damaged Femoral Modular Heads must not be used. Particular care must be used for the handling of ceramic femoral heads. A ceramic head that has fallen on the floor must not be implanted.

All devices are provided sterile and should be stored at ambient temperature (indicative range 0-50°C, 32-122°F) in their protective closed packaging controlled rooms, protected from exposure to light, heat and sudden changes in temperature.

Once the package is opened, ensure that both the model and size of the implant correspond exactly to the description printed on the labels. Avoid any contact between the implant and objects or substances that can alter the sterile condition or the surface integrity. Careful visual examination of each implant is recommended before use to verify the absence of damages. **Components removed from the package must not be used if they are dropped or suffer other accidental trauma. Devices must not be modified in any way.**

The device's code and lot number should be recorded in the patient's case history by using the labels included with the component packaging.

The disposal of medical devices is to be performed by the hospitals in conformity with applicable laws. **Re-use of previously implanted devices must be absolutely avoided. Risks associated with reuse of single use devices are:** • infection; • early or late failure of the device or device fixation; • lack of appropriate coupling between modular junctions (e.g.: taper connections); • device wear and wear debris associated complications; • transmission of diseases (e.g.: HIV, hepatitis); • immune system response / rejection.

2. GENERAL INFORMATION ON INDICATIONS, CONTRAINDICATIONS AND RISK FACTORS

2.1. INDICATIONS

Femoral Modular heads are intended for use in total hip replacement with compatible cemented or uncemented femoral and acetabular components. Cephalic and Bipolar heads are indicated for partial hip replacement with compatible cemented or uncemented femoral components. LimaCorporate femoral heads can be used in both primary and revision surgery.

Hip arthroplasty is intended for reduction or relief of pain and/or improved hip function in skeletally mature patients with the following conditions: • non-inflammatory degenerative joint disease including osteoarthritis, avascular necrosis and dysplasia; • advanced articular destruction generated by primary degenerative or rheumatoid arthritis, or post-traumatic arthrosis or post-traumatic arthritis; • correction of functional deformity; • treatment of non-union, femoral neck fracture and trochanteric fractures of the proximal femur with head involvement, unmanageable by other techniques; • revision procedures where other treatments or devices have failed.

2.2. CONTRAINDICATIONS

All LimaCorporate Femoral Heads have the same contraindications.

Absolute contraindications include: • local or systemic infection; • acute or chronic infections; • septicæmia; • persistent acute or chronic osteomyelitis; • confirmed nerve or muscle lesion compromising hip joint function. Relative contraindications include: • serious muscular, neurological or vascular diseases affecting the concerned limb; • poor bone stock (for example due to osteoporosis) compromising the stability of the implant; • metabolic disorders which may impair fixation and stability of the implant; • any concomitant disease and dependence that might adversely affect the implanted prosthesis; • metal hypersensitivity to implant materials; • patients with renal impairment.

2.3. RISK FACTORS

The following risk factors may result in poor results with this prosthesis: • overweight; • heavy physical activities, associated with frequent shocks, in which the prosthesis can be excessively overloaded (i.e. physical activity, heavy work, sport competitions, etc.); • incorrect implant positioning (e.g. varus positioning); • wrong size of components; • medical disabilities which can lead to an unnatural gait and loading of the hip joint; • muscle deficiencies; • multiple joint disabilities; • refusal to modify postoperative physical activities; • patient's history of infections or falls; • systemic diseases and metabolic disorders; • local or disseminated neoplastic diseases; • drug therapies that adversely affect bone quality, healing, or resistance to infection; • drug addiction and/or drug, alcohol or medicine abuse; • marked osteoporosis or osteomalacia; • patient's resistance to disease generally weakened (HIV, tumour, infections); • severe deformity leading to impaired anchorage or improper positioning of implants; • local bone tumor; • errors of operative technique.

3. WARNINGS

3.1. PRE-OPERATIVE PLANNING

LimaCorporate products should be implanted only by surgeons familiar with the joint replacement procedures described in the specific surgical techniques.

ALLOWED COMBINATIONS

LimaCorporate femoral heads must not be used with components from other systems or other manufacturers with exception of the following possible combination: Femoral Heads made of • BioloX Delta dia. 28mm taper 12/14 (S, M & L); • CoCrMo dia. 28mm taper 12/14 (S, M & L); • AISI 316/L dia. 28mm taper 12/14 (S, M & L); • FeCrNiMnMoNbN dia. 22mm (-2, 0 & +4); • FeCrNiMnMoNbN dia. 28mm taper 12/14 (S, M & L); • Revision BioloX Delta dia. 28mm taper 12/14 (S, M & L), when used in combination with the following LimaCorporate femoral stems: Modulus, Modulus R, Revision, Friendly short, H-MAX C, H-MAX S, MINIMA S, MASTER SL, can be also coupled with Ø44mm to 68mm Symbol Cup (DM CEM, DM HA, DMR HA) in combination with dual mobility PE Liner (for Ø22.2mm head: Ø44mm to 68mm; for Ø28mm head: Ø46mm to 68mm) manufactured by Dedienne Santé. The cup/liner assemblies come in corresponding sizes; the liner/head assemblies come in corresponding sizes, considering the internal diameter of the liner.

The surgeon must consult also SYMBOL CUP DM instruction for use before using the mixed allowed coupling between LimaCorporate and Dedienne Santé components.

NOT ALLOWED COMBINATIONS

USE OF CERAMIC HEADS:

Ceramic heads, with the exception of BioloX® Delta revision heads, may only be used for brand new, unused stem tapers and never for used ones.

USE OF SKIRTED FEMORAL HEADS:

In primary surgery, due to reduced flexion - extension range of motion, skirted femoral heads (28 - 32 mm #XL, #XXL, #XXXL and 36 mm #XXL, #XXXL) should not be used with the protruded cups or liners.

In revision surgery where the hip femoral stem is left in place and the acetabulum is reconstructed, surgeons may consider the use of a skirted femoral head with protruded cups or liners to achieve appropriate joint stability. However, this may lead to impingement between the head and the cup liner that could potentially cause damage to the implants requiring further surgery. If the surgeon believes that the use of a skirted head with a protruded cup or liner is necessary to achieve joint stability, the patient should be warned of the chance of impingement of the components as well as the possible consequences, including the risks and possible complications. In such circumstances, the patient should also be advised to limit flexion - extension joint movements to help minimize the potential for impingement and associated complications.

Pre-operative planning, through radiographic templates in different formats, provides essential information

regarding the type and size of components to be used and the correct combination of devices required based on the anatomy and specific conditions of each patient. Inadequate pre-operative planning can lead to improper selection of the implants and/or incorrect implant positioning.

Complications or failures of the total hip replacement are more likely to occur in heavy and highly active patients and high offset combinations. The surgeon should perform a careful evaluation of the patient's clinical condition and level of physical activity before performing hip replacement.

LimaCorporate specialized technical staff is available to provide advice regarding pre-operative planning, the surgical technique, and product and instrumentation assistance both prior to and during surgery.

The patient should be warned that the prosthesis does not replace normal healthy bone, that the prosthesis can break or become damaged as a result of certain activity or trauma, that it has a finite expected implant life, and may need to be replaced at some time in the future.

The possible impact of the factors mentioned in sections 2 and 3.4 should be considered preoperatively and the patient informed as to what steps he/she can take to reduce the possible effects of these factors.

Implants are single use devices; do not re-use implants that were previously implanted in another patient. Do not re-use an implant that has previously come into contact with the body fluid or tissue of another person.

Surgical instruments are subject to wear with normal usage. After extensive use or excessive loads, instruments are susceptible to fracture. Surgical instruments should be used only for their specific purpose. Before use, the functionality of surgical instruments should be checked since the use of damaged instruments may lead to early failure of the implants. Damaged instruments should be replaced before surgery.

3.2. INTRAOPERATIVE

In case of revision surgery it is essential to examine the taper of the femoral component, which remains in-situ, for any type of damages. The operating surgeon has to make sure that this damage is acceptable before using BioloX® Delta revision femoral heads.

Inspection of the stem taper and decision criteria:

- Acceptable condition: used stem tapers displaying fine marks from head-stem disassembly.
- Unacceptable stem taper deformations: see Fig 1 The BioloX® Delta revision heads must not be used with stem tapers under these conditions.

Consult the surgical technique for more details on the femoral heads assembling procedures.

Unacceptable taper damage
(DO NOT USE with BIOLOX® delta ceramic femoral heads for revision)

Slanted taper

Taper with broad truncation

Crushed taper

Fig.1

The use of trial devices is recommended to check the correct site preparation, size and positioning of the implants to be used. It is recommended that additional implants be available during surgery for use in those cases requiring prostheses of different sizes or when the preoperatively selected prostheses cannot be used. The correct selection as well as the correct seating/placement of the implant is extremely important. Improper selection, positioning, alignment and fixation of the implant components may result in unusual stress conditions which may negatively affect system performance and survival rate of the implant.

In the extremely rare occurrence of the fracture of a ceramic head, all ceramic particles must be removed during revision surgery and the use of a revision ceramic head is suggested. If a polyethylene cup or liner is being used it must be removed, even if it is firmly in place.

The components forming original LimaCorporate systems must be assembled according to the surgical technique and used only for the labelled indications.

Use only instruments and prosthesis trials specifically designed for use with the implants being used. The use of instruments from other manufacturers or the use of instruments designed for use with other systems can lead to inappropriate preparation of implant site, incorrect positioning, alignment and fixation of the devices followed by loosening of system, loss of functionality and reduction of the durability of the implant, and the need for further surgery.

Care must be taken to protect the surfaces involved in the coupling between components (tapers): do not remove the plastic protective cap, which protect the stem taper from damage, until immediately before the trial head is put on. Articular surfaces of the implants should be protected from scratches or any other damage. All component coupling surfaces should be clean and dry before assembly. The stability of component couplings should be verified as described in the surgical technique.

When placing a ceramic head on the stem: • Place the head on the stem taper by twisting lightly and using axial manual pressure until it sits firmly; • Place the plastic impactor on the pole of the head and with a light tap of the hammer in an axial direction firmly and definitively fix it on the stem taper. The surface structure of the metal taper becomes distorted plastically by the tapping of the impactor, causing an optimal distribution of pressure and a torsion-resistant fixation.

Caution: never use a metal hammer on the ceramic heads. Use only the plastic head impactor provided for this purpose in the instrument set.

A ceramic femoral head that has once been tapped onto a stem taper and then removed may not be placed on again.

3.3. POST-OPERATIVE CARE

Adequate post-operative care should be provided by the surgeon or other suitably qualified medical staff. Regular postoperative x-ray follow up is recommended to detect any changes in position of the implant or surrounding tissue.

The surgeon make the patient aware of the limitations of the limb function after hip arthroplasty and that the reconstructed joint must be protected from full load bearing for a period of time. Excessive physical activity

or trauma to the replaced hip joint can lead to premature failure of the hip arthroplasty through loosening, fracture, or abnormal wear of the prosthetic implants. The patient should be cautioned by the surgeon to govern activities accordingly and advised that the implants may fail due to excessive joint wear.

In particular the following precautions should be presented to the patient by the surgeon: • avoid repeated high weight lifting; • avoid excessive or repeated stair climbing; • keep body weight under control, overweight conditions may adversely affect the outcomes of joint replacement; • avoid sudden peak loads (consequences of activities such as running and skiing) or movements which can lead to sudden stops or twisting; • avoid positions that can increase the risk of dislocation.

Lack of appropriate post-operative rehabilitation instructions and care can negatively affect the outcome of the surgical procedure.

3.4. POSSIBLE ADVERSE EFFECTS AND COMPLICATIONS

Adverse effects that most commonly and frequently occur in hip arthroplasty include: • Excessive wear of UHMWPE components; • loosening of the prosthetic components; • prosthesis dislocation and instability; • damage to the prosthetic implant; • instability of the system because of inadequate soft tissues balancing; • dissociation due to a incorrect coupling of the devices; • infection; • subsidence; • local hypersensitivity; • local pain; • periprosthetic fractures; • temporary or permanent nerve damage; • fractures of the devices; • limb length discrepancies; • additional surgery; • squeaking and clicking (for ceramic heads); • pseudotumor. Some adverse effects can lead to death.

General complications include venous thrombosis with/without pulmonary embolism, cardiovascular or pulmonary disturbances, haematomas, systemic allergic reactions, and systemic pain.

In very isolate cases, fracture of ceramic heads can occur. The cause of this can be an overload of the prosthesis, for example through incorrect placement of the ceramic head on the stem taper or a wrong or missing fit between the head and the taper. In case a ceramic component breaks, a pairing of metal head with a polyethylene liner and of metal-on-metal components is contraindicated in revisions.

4. STERILITY

a. Implants

LimaCorporate Femoral Heads are provided sterile with a Sterility Assurance Level (SAL) of 10^{-6} . Metal and ceramic components are sterilized by radiation or EtO, and UHMWPE components by EtO.

Do not use any component from a package that has been previously opened or appears to be damaged.

Do not use implants after the expiration date printed on the label.

b. Instruments

Instruments are supplied non sterile and must be cleaned, disinfected and sterilized before use according to appropriate validated methods (refer to the "Instrument Care, Cleaning, Disinfections and Sterilization" brochure for validated sterilization parameters; this brochure is available on request or downloadable from

www.limacorporate.com in the Products section). Users should validate their specific cleaning, disinfection and sterilization processes and equipment.

5. MAGNETIC RESONANCE IMAGING (MRI)

There are inherent risks associated with the use of metallic implants in the MR environment; including component migration, heat induction and signal interference or distortion near the component(s).

Heat induction of metallic implants is a risk related to the component geometry and material, as well as the MR power, duration and pulse sequence. Since MR equipment is not standardized, the severity and likelihood of occurrence are unknown for implants.

The femoral heads have not been evaluated for safety and compatibility in the MR environment. The femoral heads have not been tested for heating or migration in the MR environment. Since these devices have not been tested, LimaCorporate cannot make a recommendation for the use of MRIs with these implants, neither for safety considerations nor imaging accuracy. These systems include components which are passive metallic devices, and as with all passive devices, there is potential for reciprocal interference with certain imaging modalities; including image distortion for MR and X-ray scanner in CT.

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Limacorporate S.p.A.