MONOLITHIC CEMENTED STEMS

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INSTRUCTIONS FOR USE – MONOLITHIC CEMENTED STEMS

Before using a LimaCorporate product, the surgeon should carefully study the following recommendations, warnings and instructions, as well as the available product-specific information (e.g. product literature, surgical technique).

1. PRODUCT INFORMATION

Monolithic cemented stems are intended for partial or total primary hip joint replacement. Monolithic cemented stems are intended for use in cemented applications.

The main goal of a joint prosthesis is to reproduce the articular anatomy, in part or totally. The joint prosthesis is intended to reduce pain and give articular mobility to the patient. The degree of pain reduction and mobility depend, in part, on the pre-operative situation, intra-operative options and the post-op rehabilitation. When used in total hip arthroplasty, monolithic cemented stems are intended for use with femoral heads and compatible acetabular cups. When used in partial hip arthroplasty, they are intended for use with femoral heads intended for partial hip replacement or bipolar heads.

1.1. MATERIALS

Components	Material	Accessories	Accessories Material
SL Self Locking femoral stem	Ti6Al4V	Cement Plugs	UHMWPE/UHMWPE + AISI 316L
LOGICA MIRROR femoral stem with Centralizer	FeCrNiMnMoNbN and UHMWPE	Cement Plugs	UHMWPE/UHMWPE + AISI 316L
CL TRAUMA femoral stem	AISI 316L	Cement Plugs	UHMWPE/UHMWPE + AISI 316L

FRIENDLY femoral stem	FeCrNiMnMoNbN	FRIENDLY accessories (plugs, hoods and proximal centralizer)	PMMA and UHMWPE
FRIENDLY SHORT femoral stem	FeCrNiMnMoNbN	FRIENDLY SHORT accessories (plugs and hood)	PMMA and UHMWPE
H-MAX C femoral stems (Standard-Lateralizing)	FeCrNiMnMoNbN	Cement Plugs	UHMWPE/UHMWPE + AISI 316L
Material Standards			
Ti6Al4V (ISO 5832-3 - ASTM F1472) – AISI 316L (ISO 5832-1) – FeCrNiMnMoNbN (ISO 5832-9 – ASTM F1586) – UHMWPE (ISO 5834-2 – ASTM F648) – PMMA (Biocompatibility USP Class VI)			

Note: for Accessories instruction for use, see the specific IFU related (Instruction for use of Cement plugs and Centralizers)

Some patients may be hypersensitive or allergic to the implant materials; this should be appropriately considered by the Surgeon. Follow the instructions of the manufacturer of bone cement (PMMA) for preparing the bone cement applying bone cement for prosthesis fixation.

1.2. HANDLING AND STORAGE

All devices are provided sterile and should be stored at ambient temperature (indicative range 0-50° C/ 32-122°F) in their protective closed packaging in controlled rooms, protected from exposure to light, heat and sudden changes in temperature.

Once the package is opened, ensure that both the model and size of the implant correspond to the description printed on the labels. Avoid any contact between the implant and objects or substances that can alter the sterile condition or the surface integrity. Careful visual examination of each implant is recommended before use in order to verify that the implant is not damaged.

Components removed from the package shall not be used if they are dropped or suffer other accidental trauma. Devices shall not be modified in any way.

The device's code and lot number should be recorded in the patient's case history by using the labels included in the component packaging. The disposal of medical devices is to be performed by the hospitals in conformity with applicable laws.

Re-use of previously implanted devices must be absolutely avoided. Risks associated with reuse of single use devices are: • infection; early or late failure of the device or device fixation; • lack of appropriate coupling between modular junctions (e.g. taper connections); • device wear and wear debris associated complications; • transmission of diseases (e.g. HIV, hepatitis); • immune system response / rejection.

2. GENERAL INFORMATION ON INDICATIONS, CONTRAINDICATIONS AND RISK FACTORS

2.1. INDICATIONS

Monolithic cemented stems are indicated for use in partial or total hip arthroplasty and they are intended for cemented use. When used in total hip arthroplasty, monolithic cemented stems are intended for use with modular heads and compatible acetabular cups. When used in partial hip arthroplasty, they are intended for use with femoral heads intended for partial hip arthroplasty or bipolar heads.

Hip arthroplasty is intended for reduction or relief of pain and/or improved hip function in skeletally mature patients with the following conditions: • non-inflammatory degenerative joint disease including osteoarthritis, avascular necrosis and dysplasia; • rheumatoid arthritis; • treatment of femoral head and neck fractures.

2.2. CONTRAINDICATIONS

Absolute contraindications include: • local or systemic infection; • septicaemia; • persistent acute or chronic osteomyelitis; • confirmed nerve or muscle lesion compromising hip joint function.

Relative contraindications include: • vascular or nerve diseases affecting the concerned limb; • poor bone stock (for example due to osteoporosis or extended previous revision surgery) compromising the stability of the implant; • metabolic disorders which may impair fixation and stability of the implant; • any concomitant disease and dependence that might affect the implanted prosthesis; • metal hypersensitivity to implant materials.

Absolute contraindications include also: • patients with mass higher than 71Kg for: - H-MAX C sizes #9 (Standard and Lateralizing) and #10 (Standard and Lateralizing); - LOGICA MIRROR size #1; - FRIENDLY sizes #1 (132° version) and #2 (125° and 132° versions). • patients with mass higher than 66Kg for: - LOGICA MIRROR size #0; - FRIENDLY sizes #1 (125° version). • patients with mass higher than 63Kg for SL Self Locking size #02.

2.3. RISK FACTOR

The following risk factors may result in poor results with this prosthesis: • overweight; • strenuous physical activities (active sports, heavy physical work); • incorrect implant positioning (e.g varus positioning); • medical disabilities which can lead to an unnatural gait and loading of the hip joint; • muscle deficiencies; • multiple joint disabilities; • refusal to modify postoperative physical activities; • patient history of infections or falls; • systemic diseases and metabolic disorders; • local or disseminated neoplastic diseases; • drug therapies that adversely affect bone quality, healing, or resistance to infection; • drug use or alcoholism; • marked osteoporosis or osteomalacia; • patient's resistance to disease generally weakened (HIV, tumour, infections); • severe deformity leading to impaired anchorage or improper positioning of implants.

3. WARNINGS

3.1. PRE-OPERATIVE PLANNING

LimaCorporate products should be implanted only by surgeons familiar with the joint replacement procedures described in the specific surgical techniques.

Pre-operative planning, through radiographic templates in different formats, provides essential information regarding the type and size of components to be used and the correct combination of devices required based on the anatomy and specific conditions of each patient. Inadequate pre-operative planning can lead to improper selection of the implants and/or incorrect implant positioning.

The surgeon should carefully plan the surgery considering the following:

1. Small sized stems: the Small sized stems are designed for patients with a small intramedullary canal and/or metaphyseal region of the femur. The reduced size (M/L width) of these stems results in a corresponding reduction in the fatigue strength of the implant;

2. High Offset combinations: (use of lateralized stems): the lateralized stems are designed to restore the functional offset of the hip joint to be comparable to that of the contralateral hip, however greater neck lengths are accompanied by a higher risk of failure (e.g. breakage due to fatigue).

Complications or failures of the total hip replacement are more likely to occur in heavy and highly active patients and high offset combinations. The surgeon should perform a careful evaluation of the patient's clinical condition and level of physical activity before performing hip replacement.

NOTE: monolithic cemented stems have been tested according to ISO standards for fatigue strength. Combinations with high offset femoral heads #XXL or #XXXL has not been assessed for all cemented stems with exception Friendly Short stem in standard version. If the surgeon decides to use a #XXL or #XXXL femoral head to achieve appropriate joint stability he should carefully consider the risks introduced by high offset combinations.

COMBINATION ALLOWED/NOT ALLOWED

- Femoral stems made from AISI 316L stainless steel (ISO 5832-1) are not suitable for coupling with CoCrMo (ISO 5832-12, ASTM F1537) femoral heads.
- FRIENDLY SHORT stems in lateralized version is limited in use with heads up to L size.
- FRIENDLY SHORT stem #0 in standard version is limited in use with heads up to XXL size.

LimaCorporate monolithic cemented stems must not be used with components from other systems or other manufacturers with exception of the following possible combination:

For **Friendly short** and **H-MAX C** stems, femoral heads manufactured by LimaCorporate made of • Biolox Delta dia. 28mm taper 12/14 (S, M & L); • CoCrMo dia. 28mm taper 12/14 (S, M & L); • AISI 316/L dia. 28mm taper 12/14 (S, M & L); • FeCrNiMnMoNbN dia. 22mm (-2, 0 & +4); • FeCrNiMnMoNbN dia. 28mm taper 12/14 (S, M & L); • Revision Biolox Delta dia. 28mm taper 12/14 (S, M & L), can be also coupled with Ø44mm to 68mm Symbol Cup (DM CEM, DM HA, DMR HA) in combination with dual mobility PE Liner (for Ø22.2mm head: Ø44mm to 68mm; for Ø28mm head: Ø46mm to 68mm) manufactured by Dedienne Santé.

The cup/liner assemblies come in corresponding sizes; the liner/head assemblies come in corresponding sizes, considering the internal diameter of the liner.

The surgeon must consult also SYMBOL CUP DM instruction for use before using the mixed allowed coupling between LimaCorporate and Dedienne Santé components.

LimaCorporate specialized technical staff is available to provide advice regarding pre-operative planning, the surgical technique, and product and instrumentation assistance both prior to and during surgery.

The patient should be warned that the prosthesis does not replace normal healthy bone, that the prosthesis can break or become damaged as a result of certain activity or trauma, that it has a finite expected implant life, and may need to be replaced at some time in the future.

The possible impact of the factors mentioned in sections 2 and 3.4 should be considered preoperatively and the patient informed as to what steps he/she can take to reduce the possible effects of these factors.

Implants are single use devices; do not re-use implants that were previously implanted in another patient. Do not re-use an implant that has previously come into contact with the body fluid or tissue of another person.

Surgical instruments are subject to wear with normal usage. After extensive use or excessive loads, instruments are susceptible to fracture. Surgical instruments should be used only for their specific purpose. Before use, the functionality of surgical instruments should be checked since the use of damaged instruments may lead to early failure of the implants. Damaged instruments should be replaced before surgery.

3.2. INTRAOPERATIVE

The use of trial devices is recommended to check the correct site preparation, size and positioning of the

implants to be used. It is recommended that additional implants be available during surgery for use in those cases requiring prostheses of different sizes or when the preoperatively selected prostheses cannot be used. The correct selection as well as the correct seating/placement of the implant is extremely important. Improper selection, positioning, alignment and fixation of the implant components may result in unusual stress conditions which may negatively affect system performance and the survival rate of the implant.

The components forming original LimaCorporate systems must be assembled according to the surgical technique and used only for the labelled indications. Use only instruments and prosthesis trials specifically designed for use with the implants being used. The use of instruments from other manufacturers or the use of instruments designed for use with other systems can lead to inappropriate preparation of the implant site, incorrect positioning, alignment and fixation of the devices followed by loosening of the system, loss of functionality, reduction of the durability of the implant, and the need for further surgery. Care must be taken to protect the surfaces involved in the coupling between components (tapers); articular surfaces of the implants should be protected from scratches or any other damage. All component coupling surfaces should be clean and dry before assembly. The stability of component couplings should be verified as described in the surgical technique.

3.3. POST-OPERATIVE CARE

Adequate post-operative care should be provided by the surgeon or other suitably qualified medical staff. Regular postoperative x-ray follow up is recommended to detect any changes in position of the implant or surrounding tissues.

The surgeon should make the patient aware of the limitations of limb function after hip arthroplasty and that the reconstructed joint must be protected from full load bearing for a period of time. Excessive physical activity or trauma to the replaced hip joint can lead to premature failure of the hip arthroplasty through loosening, fracture, or abnormal wear of the prosthetic implants. The patient should be cautioned by the surgeon to govern activities accordingly and advised that the implants may fail due to excessive joint wear.

In particular the following precautions should be presented to the patient by the surgeon:

- avoid repeated high weight lifting;
- avoid excessive or repeated stair climbing;
- keep body weight under control (overweight conditions or increase of weight over the limitation in mass may adversely affect the outcomes of joint replacement);
- avoid sudden peak loads (consequences of activities such as running and skiing) or movements which can lead to sudden stops or twisting;
- avoid positions that can increase the risk of dislocation.

Appropriate post-operative rehabilitation instructions and care must be provided to avoid a negative influence on the outcome of the surgical procedure.

3.4. POSSIBLE ADVERSE EFFECTS

Adverse effects that most commonly and frequently occur in hip arthroplasty include:

- excessive wear of UHMWPE components;
- loosening of the prosthetic components;
- prosthesis dislocation and instability;
- damage to the prosthetic implant;
- instability of the system because of inadequate soft tissue balancing;
- dissociation due to incorrect coupling of the devices;
- infection;
- subsidence;
- local hypersensitivity;
- local

pain; • periprosthetic fractures; • temporary or permanent nerve damage; • fractures of the devices; • limb length discrepancies; • blood loss; • additional surgery.

Some adverse effects can lead to death. General complications include venous thrombosis with/without pulmonary embolism, cardiovascular disturbances, haematomas, systemic allergic reactions and systemic pain.

4. STERILITY

a. Implants

All implantable components of monolithic cemented stems are provided sterile with a Sterility Assurance Level (SAL) of 10^{-6} . Metal components are sterilized by radiation or EtO. UHMWPE components by EtO.

Do not use any component from a package that has been previously opened or appears to be damaged.

Do not use implants after the expiration date printed on the label.

b. Instruments

Instruments are supplied non sterile and must be cleaned, disinfected and sterilized before use according to appropriate validated methods. (Refer to the "Instrument Care, Cleaning, Disinfections and Sterilization" brochure for validated sterilization parameters; this brochure is available on request or downloadable from www.LimaCorporate.com in the Products section). Users should validate their specific cleaning, disinfection and sterilization processes and equipment.

5. MAGNETIC RESONANCE IMAGING (MRI)

There are inherent risks associated with the use of metallic implants in the MR environment; including component migration, heat induction and signal interference or distortion near the component(s). Heat induction of metallic implants is a risk related to the component geometry and material, as well as the MR power, duration and pulse sequence. Since MR equipment is not standardized, the severity and likelihood of occurrence are unknown for these implants.

The Monolithic cemented stems have not been evaluated for safety and compatibility in the MR environment. The Monolithic cemented stems have not been tested for heating or migration in the MR environment. Since these devices have not been tested, LimaCorporate cannot make a recommendation for the use of MRIs with these implants, neither for safety considerations nor imaging accuracy.

These systems include components which are passive metallic devices and as with all passive devices, there is potential for reciprocal interference with certain imaging modalities, including image distortion in MR and X-ray scanner in CT.

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Limacorporate S.p.A.