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## INSTRUCTIONS FOR USE – PRIMA GLENOID SYSTEM

Before using a LimaCorporate product, the surgeon shall carefully study the following recommendations, warnings, and instructions, as well as the most updated available product-specific information (e.g.: product literature, surgical technique).

### 1. PRODUCT INFORMATION

The PRIMA Glenoid System is intended for total primary or revision shoulder joint replacement in reverse shoulder configuration. The components are intended for use in uncemented applications. The joint prosthesis is intended to reduce pain and give articular mobility to the patient. The degree of pain reduction and mobility depend in part on the pre-operative situation, intra-operative options and the post-op rehabilitation. PRIMA Glenoid System components include modular and monoblock glenoid baseplates, glenoid peg, glenosphere, glenosphere connector and screw, peripheral screw, central compressive screw and screw locking cap. PRIMA Glenoid System components can be used in combination with components of the LimaCorporate SMR Shoulder System, SMR Stemless System and PRIMA Humeral System; PRIMA Glenoid System components shall not be used alone as single components, or with components from other systems or other manufacturers.

#### 1.1. MATERIALS

The materials of the PRIMA Glenoid System are reported in the following table.

Components	Material	Raw Material Reference <sup>a</sup>
PRIMA TT Glenoid Monoblock Reverse TT Baseplate	Ti6Al4V 3D printed	M3
PRIMA TT Glenoid Modular Reverse TT Baseplate	Ti6Al4V 3D printed	M3
PRIMA TT Glenoid TT Peg	Ti6Al4V 3D printed	M3
PRIMA TT Glenoid Glenosphere	CoCrMo	M4
	UHMWPE X-Lima, Ti6Al4V	M8, M1
	LimaVit™ (Vitamin E highly crosslinked UHMWPE), Ti6Al4V	M9, M1

PRIMA TT Glenoid Glenosphere Connector Screw	Ti6Al4V, UHMWPE	M8, M1
PRIMA TT Glenoid Glenosphere Connector with Screw	Ti6Al4V, UHMWPE	M8, M1
PRIMA TT Glenoid Central Compressive Cortical Screw	Ti6Al4V	M1
PRIMA TT Glenoid Central Compressive Cancellous Screw	Ti6Al4V	M1
PRIMA TT Glenoid Peripheral Cortical Screw	Ti6Al4V	M1
PRIMA TT Glenoid Locking Cap	Ti6Al4V	M1
<b>Material Standards</b>		
<b>Ti6Al4V</b> (ISO 5832-3 – ASTM F1472) – <b>Ti6Al4V 3D printed</b> (ISO 5832-3*) – <b>CoCrMo</b> (ISO 5832-12 – ASTM F1537) – <b>UHMWPE</b> (ISO 5834-2 – ASTM F648) – <b>LimaVit™</b> (Vitamin E highly crosslinked UHMWPE) (ISO 5834-2 - ASTM F648 - ASTM F2695 – ASTM F2565)		

# For the chemical composition of the raw materials please refer to Section “**Raw Material Compositional Limits (%)**” considering the raw material reference indicated in the table.

\* Mechanical properties and chemical composition are in compliance with ISO 5832-3.

All materials are purchased according to International Standards and are widely used in the orthopaedic field. CoCrMo, used to manufacture some components of PRIMA Glenoid System, includes cobalt (which is a CMR class 1B) in percentage above 0.1% (w/w). Nevertheless, chemical tests show that the release of cobalt in the human body is below 0.1% (w/w). Additionally, biocompatibility tests and long-term clinical experience of the use of this material in orthopaedic applications have shown that an acceptable level of biological response can be expected. Therefore, the risks related to the use of CoCrMo and other materials have been considered and the related risk/benefit ratio has been evaluated as acceptable by the manufacturer. Some patients may be vulnerable (e.g. have hypersensitive or allergic reactions) to the implant materials/substances; this shall be appropriately considered by the Surgeon.

#### **1.1.1. MEDICINAL SUBSTANCES AND TISSUE**

The device does not contain or incorporate any: • medicinal substance, including a human blood or plasma derivative; • tissues or cells, or their derivatives, of human origin; • tissues or cells of animal origin, or their derivatives, as referred to in Regulation (EU) No.722/2012.

## 1.2. HANDLING, STORAGE AND SAFE DISPOSAL

All devices are provided sterile and shall be stored at ambient temperature (indicative range 0-50°C / 32-122°F) in their protective closed packaging in controlled rooms, protected from exposure to light, heat and sudden changes in temperature.

Once the package is opened check that both the model and size of the implant correspond exactly to the description printed on the labels. Avoid any contact between the implant and objects or substances that can alter the sterile condition or surface integrity. A careful visual examination of each implant is recommended before use in order to verify that the implant is not damaged. **Components removed from the package shall not be used if they are dropped or suffer other accidental impacts. Devices shall not be modified in any way.**

The device's code and lot number shall be recorded in the patient's case history by using the labels included in the components packaging.

**Re-use of previously implanted devices must be absolutely avoided. Risks associated with the reuse of single use devices are:** • infection; • early or late failure of the device or device fixation; • lack of appropriate coupling between modular junctions (e.g. taper connections); • device wear and wear debris associated complications; • transmission of diseases (e.g. HIV, hepatitis); • immune system response / rejection.

**DISPOSAL:** the disposal of medical devices is to be performed by the hospitals in conformity with applicable laws.

## 1.3. DEVICE PERFORMANCE, EXPECTED CLINICAL BENEFITS, EXPECTED LIFETIME

The main goal of a joint prosthesis is to reproduce the articular anatomy. The joint prosthesis is intended to reduce pain and give articular mobility to the patient. The degree of pain reduction and mobility depend, in part, on the pre-operative situation, intra-operative options and the post-op rehabilitation. LimaCorporate annually updates the Summary of Safety and Clinical Performance (SSCP) to be uploaded on EUDAMED (<https://ec.europa.eu/tools/eudamed>).

**The expected lifetime of PRIMA Glenoid System is at least 95%\* at 7 years.**

**\* with a variability of  $\pm 2\%$**

**These values are based on the state-of-the-art reference for similar device (reverse shoulder replacement, stemmed or stemless) identified through arthroplasty registries.**

NOTE: The expected lifetime of the medical devices, when placed on the market and – moreover – when put into use (namely when implanted in a patient's body) may depend on a variety of factors:

- (a) intrinsic characteristic (feature) of the medical device in itself (for example: design of the medical device; materials employed in the medical device manufacturing; manufacturing process of the medical device in itself; surface texture, osteo-conductive and osteo-inductive properties of the medical device);
- (b) external factors related to surgery process (for example: right choice of the medical device, stability of the initial fixation of the medical device; accuracy in positioning of the medical device);
- (c) external factors related to patient (for example: aetiology, pathology, physical characteristics, age, previous surgical history, daily activities performed);
- (d) risk factors mentioned in the IFU leaflet, when not already included in the factors referred to in previous paragraphs from (a) to (c);
- (e) factors related to all the possible complications which can affect the duration of the implant.

The expected lifetime was calculated by means of an analytical and rationalized estimation, taking into account:

- i. the clinical data collected on the medical device;
- ii. the survivorship data reported in arthroplasty registries referring to similar medical devices (identified state of the art benchmark);
- iii. the expected (as predictable) behavior of the medical device at higher follow-up.

Due to the analytical evaluations, the expected lifetime value is approximated to the unit.

The value shall be interpreted duly considering the specified value variability, taking into necessary consideration all relevant factors mentioned above, in paragraphs (a) to (e).

## **2. INTENDED PURPOSE, INDICATIONS, CONTRAINDICATIONS AND RISK FACTORS**

Prosthetic implants are designed for long-term implantation inside the human body. Their intended purpose is the shoulder joint replacement.

### **2.1. INDICATIONS**

The PRIMA Glenoid System is indicated for primary, fracture or revision total shoulder replacement in a grossly rotator cuff deficient joint with severe arthropathy (disabled shoulder) in skeletally mature patients. The patient's joint must be anatomically and structurally suited to receive the selected implants and a functional deltoid muscle is necessary to use the device. The PRIMA Glenoid System components are intended for uncemented use with the addition of screw fixation.

## **2.2. CONTRAINDICATIONS**

Absolute contraindications include: • local or systemic infection; • septicaemia; • persistent acute or chronic osteomyelitis; • confirmed nerve lesion compromising shoulder joint function; • deltoid muscle insufficiency.

Relative contraindications include: • vascular or nerve diseases affecting the concerned limb; • poor bone stock (for example due to osteoporosis or extended previous revision surgery) compromising the stability of the implant; • metabolic disorders which may impair fixation and stability of the implant; • any concomitant disease and dependence that might affect the implanted prosthesis; • metal hypersensitivity to implant materials.

In cases of bone tumors, use an appropriate system designed to treat cases requiring large bone resections. The use of primary or revision implants not designed and intended for use in cases of bone resection may result in a poor outcome and / or failure of the implant or implant fixation.

## **2.3. RISK FACTORS**

The following risk factors may result in poor results with this prosthesis: • overweight<sup>1</sup>; • strenuous physical activities (active sports, heavy physical work); • incorrect implant positioning; • wrong size of components; • muscle deficiencies; • multiple joint disabilities; • refusal to modify postoperative physical activities; • patient's history of infections or falls; • systemic diseases and metabolic disorders; • local or disseminated neoplastic diseases; • drug therapies that adversely affect bone quality, healing, or resistance to infection; • drug use or alcoholism; • marked osteoporosis or osteomalacia; • patient's resistance to disease generally weakened (HIV, tumour, infections); • severe deformity leading to an impaired anchorage or improper positioning of implants; • use in combination with products, prosthesis or instruments of another manufacturer; • errors of operative technique.

## **3. WARNINGS**

### **3.1. PRE-OPERATIVE PLANNING**

LimaCorporate products shall be implanted only by surgeons familiar with the joint replacement procedures described in the specific surgical techniques.

The components of the system must not be used with components of other manufacturer, unless specified by LimaCorporate. The components forming original LimaCorporate systems must be rigorously coupled following the given information included in this document and in the available product-specific information (e.g.: current version of surgical technique). The use of single components, or components belonging to other systems, are subject to the approval by LimaCorporate. The manufacturer and the dealer are not liable for possible coupling incompatibility. The surgeon is solely responsible for the choice and use of the implant.

Pre-operative planning, through radiographic templates in different formats, provides essential information regarding the type and size of components to be used and the correct combination of required devices

---

<sup>1</sup> According to the definition of World Health Organization (WHO), Body Mass Index (BMI) greater than or equal to 25 Kg/m<sup>2</sup>

based on the anatomy and specific conditions of each patient. Inadequate pre-operative planning can lead to improper selection of the implants and/or incorrect implant positioning.

In selecting patients for surgery, the following factors can be critical to the eventual success of the procedure: the bone stock of the glenoid and humerus must be able to support the implant. In cases with significant bone loss and in which adequate fixation on the glenoid side cannot be obtained, a hemiarthroplasty with a CTA-head should be performed. PRIMA Glenoid System components can be used in combination with the humeral components of the SMR Shoulder System and the PRIMA Humeral System to be only used in reverse shoulder configuration. Instructions for use of SMR Shoulder System, SMR Stemless System and PRIMA Humeral System shall be consulted for warnings that apply to these systems.

## COMBINATIONS ALLOWED/NOT ALLOWED

**NOTE:** Regular sizes of monoblock and modular baseplates must be fixed to the glenoid bone by means of two peripheral screws with minimum length 20mm and one peripheral screw with minimum length 18mm. Alternatively, the same baseplates can be fixed to the glenoid bone by means of two peripheral screws with minimum length 20mm and one central compressive screw. Bone screws must be positioned according to the surgical technique.

**NOTE:** Lateralized and Full Wedge sizes of monoblock and modular baseplates must be fixed to the glenoid bone by means of two peripheral screws with minimum length 26mm and one peripheral screw with minimum length 18mm. Alternatively, the same baseplates can be fixed to the glenoid bone by means of two peripheral screws with minimum length 26mm and one central compressive screw. Bone screws must be positioned according to the surgical technique.

**NOTE:** When using the modular baseplate, the central compressive screw can be used only with modular TT Peg size Short and Medium; the 25mm long central compressive screws can be used only with the modular TT Peg size Short.

**NOTE:** The following table identifies the allowed (✓) / not allowed (X) combinations between the baseplates (both modular and monoblock, any diameter) and the glenosphere connectors:

		GLENOSPHERE CONNECTOR SIZES		
		Low Lateralization	Medium Lateralization	High Lateralization
BASEPLATE SIZES	Regular	✓	✓	✓
	Full Wedge 10°	✓	✓	X
	Full Wedge 15° / 15°X	✓	X	X
	Full Wedge 20° / 20°X	✓	X	X
	Lateralized +2mm	✓	✓	X
	Lateralized +4mm	✓	X	X

LimaCorporate specialized technical staff is available to provide advice regarding pre-operative planning, the surgical technique, and product and instrumentation assistance both prior to and during surgery.

The patient shall be warned that the prosthesis does not replace normal healthy bone, that the prosthesis can break or become damaged as a result of certain activity or trauma, that it has a finite expected implant life, and may need to be replaced at some time in the future: the longevity and the performance of the implant over time,

in fact, can be influenced by the natural/physiological progression of the patient pathology, the insurgence of co-morbidities and post-operative complications, even if they are not directly implant-related (i.e. infection pain, stiffness, reduced range of motion).

The possible impact of the factors mentioned in sections "2. INTENDED PURPOSE, INDICATIONS, CONTRAINDICATIONS AND RISK FACTORS" and "3.4. POSSIBLE SIDE EFFECTS" shall be considered preoperatively and the patient informed as to what steps he/she can take to reduce the possible effects of these factors.

**Implants are single use devices.**

**Do not re-use implants that were previously implanted in another patient.**

**Do not re-use an implant that has previously come into contact with the body fluid or tissue of another person.**

**Do not use any component from a package that has been previously opened or appears to be damaged.**

**Do not use implants after the expiration date printed on the label.**

Surgical instruments used to perform the surgery are subject to wear with normal usage. After extensive use or excessive loads, instruments are susceptible to fracture. Surgical instruments shall be used only for their specific purpose. Before use, the functionality of surgical instruments shall be checked since the use of damaged instruments may lead to early failure of the implants. Damaged instruments shall be replaced before surgery.

### **3.2. INTRAOPERATIVE**

The use of trial devices is recommended to check the correct site preparation, size and positioning of the implants to be used. It is recommended that additional implants be available during surgery for use in those cases requiring prostheses of different sizes or when the preoperatively selected prostheses cannot be used. The correct selection as well as the correct seating/placement of the implant is extremely important. Improper selection, positioning, alignment and fixation of the implant components may result in unusual stress conditions which may negatively affect system performance and the survival rate of the implant.

The components forming original LimaCorporate systems must be assembled according to the surgical technique and used only for the labelled indications.

Use only instruments and trial components specifically designed for use with the implants being used. The use of instruments from other manufacturers or the use of instruments designed for use with other systems can lead to inappropriate preparation of the implant site, incorrect positioning, alignment and fixation of the devices followed by loosening of the system, loss of functionality, reduction of the durability of the implant, and the need for further surgery.

Care must be taken to protect the surfaces involved in the coupling between components (tapers); the articular surfaces of the implants shall be protected from scratches or any other damage. All component coupling surfaces shall be clean and dry before assembly. The stability of component couplings shall be verified as described in the surgical technique. Care must be taken when handling Trabecular Titanium devices; they must not get in contact with fabrics that can release particles inside the trabecular structure.

### **3.3. POST-OPERATIVE CARE**

Adequate post-operative care shall be provided by the surgeon or other suitably qualified medical staff. Regular postoperative x-ray follow-up is recommended to detect any changes in position or condition of the implant or surrounding tissues. Complications and/or failure of prosthetic implants are more likely to occur in patients with unrealistic functional expectations, heavy patients, physically active patients, and/or with patients who fail to follow through with the required rehabilitation program. Excessive physical activity or trauma to the replaced joint can lead to premature failure of the arthroplasty through loosening, fracture, or abnormal wear of the prosthetic implants. The patient shall be cautioned by the surgeon to govern activities accordingly and advised that the implants may fail due to excessive joint wear.

The surgeon shall make the patient aware of the limitation of limb function after shoulder arthroplasty and that caution is needed, especially in the post-operative period, since the reconstructed joint must be protected from full load bearing. In particular the following precautions shall be presented to the patient by the surgeon:

- avoid repeated high weightlifting;
- keep body weight under control;
- avoid sudden peak loads (consequences of activities such contact sports, playing tennis) or movements which can lead to sudden stops or twisting;
- avoid positions that can increase the risk of dislocation, such as for example degrees of abduction an/or internal rotation.

The patient shall be informed also about the possibility that the implant or its components may wear out, fail or need to be replaced. The implant may not last the rest of the patient's life, or any particular length of time. Because prosthetic implants are not as strong, reliable, or durable as natural, healthy tissue/bones, all such devices may need to be replaced at some point.

Appropriate post-operative rehabilitation instructions and care must be provided to avoid a negative influence on the outcome of the surgical procedure.

### **3.4. POSSIBLE SIDE EFFECTS**

Adverse effects that most commonly and frequently occur in orthopaedic prosthesis procedure include:

- loosening of the prosthetic components, for example due to osteolysis;
- prosthesis dislocation and instability;
- breakage or damage of the device;
- instability of the system because of inadequate soft tissue balancing;
- dissociation due to incorrect coupling of the devices;
- infection;
- local hypersensitivity;
- local pain;
- periprosthetic fractures, including intraoperative fractures;
- bone resorption;
- temporary or permanent nerve or muscle damage;
- fretting of modular junctions;
- excessive wear of UHMWPE components due to damaged articular surfaces or the presence of particles;
- metallosis;
- heterotopic ossification;
- scapular notching;
- loss of range of motion;
- blood loss;
- additional surgery.

Some adverse effects can lead to death. General complications include venous thrombosis with/without pulmonary embolism, cardiovascular or pulmonary disturbances, haematomas, systemic allergic reactions, systemic pain, blood loss, donor site morbidity, temporary or permanent paralysis.

## 4. STERILITY

### a. Implants

All implantable components of the system are provided sterile with a Sterility Assurance Level (SAL) of  $10^{-6}$ . Metal components are sterilized by radiation (E-Beam or Gamma rays) or EtO while UHMWPE components by EtO. Do not resterilize implantable components. Do not use any component from a package that has been previously opened or appears to be damaged. **Do not use implants after the expiration date printed on the label.**

### b. Instruments

Instruments are supplied non-sterile and must be cleaned, disinfected and sterilized before use according to appropriate validated methods (refer to the "Instrument Care, Cleaning, Disinfection and Sterilization" brochure for validated sterilization parameters; this brochure is available on request or downloadable from [www.limacorporate.com](http://www.limacorporate.com) in the Products section). Users shall validate their specific cleaning, disinfection and sterilization processes and equipment.

## 5. INFORMATION USEFUL TO INFORM THE PATIENT

Refer to previous sections for any information that allows the patient to be informed of any warnings, precautions, contra-indications, measures to be taken and limitations of use regarding the device.

**The surgeon shall notify the patient about any necessary follow-up specific for the patient's condition.**

### 5.1. MAGNETIC RESONANCE IMAGING (MRI)

There are inherent risks associated with the use of metallic implants in the MR environment, including component migration, heat induction and signal interference or distortion near the component(s). Heat induction of metallic implants is a risk related to the component geometry and material, as well as the MR power, duration and pulse sequence. Since MR equipment is not standardized, the severity and likelihood of occurrence are unknown for these implants.

The PRIMA Glenoid System has not been evaluated for safety and compatibility in the MR environment. The devices have not been tested for heating or migration in the MR environment. Since these devices have not been tested, LimaCorporate cannot make a recommendation for the use of MRIs with these implants, neither for safety considerations nor imaging accuracy. The risks associated with a passive implant in an MR environment have been evaluated and are known to include heating, migration and image artifacts at or near the implant site.

## 6. INCIDENT REPORTING

According to the definition of incident/serious incident given by Regulation EU 2017/745<sup>2</sup>, users and/or patients are required to: • **report to the manufacturer and to the relevant Competent Authority any serious incident that has occurred in relation to the device;** • **report to the manufacturer any incident and/or event of malfunction of the device or changes in its performance that may affect safety.**

In non-UE market, refer to incidents' definitions according to applicable laws.

---

<sup>2</sup> an incident is defined to be any malfunction or deterioration in the characteristics and/or performance of the device, including use-error due to ergonomic features, as well as any inadequacy in the information supplied by the manufacturer and any undesirable side-effect. An incident is considered to be serious when it directly or indirectly led, might have led or might lead to any of the following: (a) the death of a patient, user or other person, (b) the temporary or permanent serious deterioration of a patient's, user's or other person's state of health, (c) a serious public health threat.

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**Limacorporate S.p.A.**