# SMR SHOULDER SYSTEM

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## **INSTRUCTIONS FOR USE – SMR SHOULDER SYSTEM**

Before using a LimaCorporate product, the surgeon shall carefully study the following recommendations, warnings, and instructions, as well as the most updated available product-specific information (e.g.: product literature, surgical technique).

### **1. PRODUCT INFORMATION**

The main goal of the SMR Anatomic joint prosthesis is to reproduce the articular anatomy, in part or totally, while in the SMR Reverse Shoulder the geometry of the natural gleno-humeral joint is inverted. The joint prosthesis is intended to reduce pain and give articular mobility to the patient. The degree of pain reduction and mobility depends in part on the pre-operative situation, intra-operative options and on the post-op rehabilitation. The SMR Shoulder System is intended for partial or total primary or revision shoulder joint replacement. SMR Shoulder System components include humeral stems, humeral bodies, adaptor tapers, humeral heads and glenoid components. The components are intended to be used in cemented or uncemented applications as specified in Section 2.1.

SMR System components shall not be used alone as single components, or with components from other manufacturer.

NOTE: some components of SMR System are allowed for the use with SMR Stemless System and PRIMA Humeral System, manufactured by Limacorporate S.p.A.

#### **1.1. MATERIALS**

The materials of the SMR System are reported in the table in Section 2.1.

All materials are purchased according to International Standards and are widely used in the orthopaedic field. CoCrMo, used to manufacture some components of SMR System, includes cobalt (which is a CMR class 1B) in percentage above 0.1% (w/w). Nevertheless, chemical tests show that the release of cobalt in the human body is below 0.1% (w/w). Additionally, biocompatibility tests and long-term clinical experience of the use of this material in orthopaedic applications has shown that an acceptable level of biological response can be expected. Therefore, the risks related to the use of CoCrMo and other materials have been considered and the related risk/benefit ratio has been evaluated as acceptable by the manufacturer. Some patients may be vulnerable (e.g. have hypersensitive or allergic reactions) to the implant materials/substances; this shall be appropriately considered by the Surgeon. Follow the instructions of the manufacturer of bone cement (PMMA) for preparing the bone cement applying bone cement for prosthesis fixation.

### 1.1.1. MEDICINAL SUBSTANCES AND TISSUE

The device does not contain or incorporates any: • medicinal substance, including a human blood or plasma derivative; • tissues or cells, or their derivatives, of human origin; • tissues or cells of animal origin, or their derivatives, as referred to in Regulation (EU) No.722/2012.

### 1.2. HANDLING, STORAGE AND SAFE DISPOSAL

All devices are provided sterile and shall be stored at ambient temperature (indicative range 0-50°C / 32-122°F) in their protective closed packaging in controlled rooms, protected from exposure to light, heat and sudden changes in temperature.

Once the package is opened check that both the model and size of the implant correspond exactly to the description printed on the labels. Avoid any contact between the implant and objects or substances that can alter the sterile condition or surface integrity. A careful visual examination of each implant is recommended before use in order to verify that the implant is not damaged. **Components removed from the package shall not be used if they are dropped or suffer other accidental impacts. Devices shall not be modified in any way.**

The device's code and lot number shall be recorded in the patient's case history by using the labels included in the components packaging.

**Re-use of previously implanted devices must be absolutely avoided. Risks associated with the reuse of single use devices are:** • infection; • early or late failure of the device or device fixation; • lack of appropriate coupling between modular junctions (e.g. taper connections); • device wear and wear debris associated complications; • transmission of diseases (e.g. HIV, hepatitis); • immune system response / rejection.

**DISPOSAL:** the disposal of medical devices is to be performed by the hospitals in conformity with applicable laws.

### 1.3. DEVICE PERFORMANCE, EXPECTED CLINICAL BENEFITS, EXPECTED LIFETIME

The main goal of a joint prosthesis is to reproduce the articular anatomy. The joint prosthesis is intended to reduce pain and give articular mobility to the patient. The degree of pain reduction and mobility depend, in part, on the pre-operative situation, intra-operative options and the post-op rehabilitation.

LimaCorporate annually updates the Summary of Safety and Clinical Performance (SSCP) to be uploaded on EUDAMED (<https://ec.europa.eu/tools/eudamed>).

**The expected lifetime of SMR Reverse System is at least 96%\* at 7 years.**  
**The expected lifetime of SMR Total Anatomic System is at least 92%\* at 7 years.**  
**The expected lifetime of SMR Hemi Anatomic System is at least 91%\* at 7 years.**  
**\* with a variability of  $\pm 2\%$**

**These values are based on the state-of-the-art reference for similar device identified through arthroplasty registries.**

NOTE: The expected lifetime of the medical devices, when placed on the market and – moreover – when put into use (namely when implanted in a patient's body) may depend on a variety of factors:

- (a) intrinsic characteristic (feature) of the medical device in itself (for example: design of the medical device; materials employed in the medical device manufacturing; manufacturing process of the medical device in itself; surface texture, osteo-conductive and osteo-inductive properties of the medical device);
- (b) external factors related to surgery process (for example: right choice of the medical device, stability of the initial fixation of the medical device; accuracy in positioning of the medical device);
- (c) external factors related to patient (for example: aetiology, pathology, physical characteristics, age, previous surgical history, daily activities performed);
- (d) risk factors mentioned in the IFU leaflet, when not already included in the factors referred to in previous paragraphs from (a) to (c);
- (e) factors related to all the possible complications which can affect the duration of the implant.

The expected lifetime was calculated by means of an analytical and rationalized estimation, taking into account:

- i. the clinical data collected on the medical device;
- ii. the survivorship data reported in arthroplasty registries referring to similar medical devices (identified state of the art benchmark);
- iii. the expected (as predictable) behavior of the medical device at higher follow-up.

Due to the analytical evaluations, the expected lifetime value is approximated to the unit.

The value shall be interpreted duly considering the specified value variability, taking into necessary consideration all relevant factors mentioned above, in paragraphs (a) to (e).

## **2. INTENDED PURPOSE, INDICATIONS, CONTRAINDICATIONS AND RISK FACTORS**

Prosthetic implants are designed for long-term implantation inside the human body. Their intended purpose is the shoulder joint replacement.

### **2.1. INDICATIONS**

The SMR System is intended for partial or total, primary or revision shoulder joint replacement in skeletally mature patients. The SMR Anatomic System is indicated for partial or total, primary or revision shoulder joint replacement in patients suffering from disability due to: • Non-inflammatory degenerative joint disease including

osteoarthritis and avascular necrosis; • Inflammatory degenerative joint diseases such as rheumatoid arthritis; • Treatment of acute fractures of the humeral head that cannot be treated with other fracture fixation methods; • Revision of a failed primary implant (not applicable to Sub-Scapularis Sparing technique; in case of SMR Cementless Finned Short Stems only if sufficient bone stock remains); • Cuff tear arthropathy (CTA heads only, not applicable to Sub-Scapularis Sparing technique); • Glenoid arthrosis without excessive glenoid bone loss: A1, A2 and B1 according to Walch classification (SMR TT Hybrid Glenoid only);

The Large Resection Stems are indicated for oncology applications.

The intended use of the Glenoid Bone Graft Instruments is to enable the preparation of the bone graft from the patient's humeral head or allograft from the bone bank, to restore the glenoid anatomy in case of glenoid deficiency (e.g. glenoid type B2 or C according to Walch's classification).

The SMR Reverse System is indicated for primary, fracture or revision total shoulder replacement in grossly rotator cuff deficient joint with severe arthropathy (disabled shoulder).

The patient's joint must be anatomically and structurally suited to receive the selected implants and functional deltoid muscle is necessary to use the device. The modular SMR System allows the assembly of components in various humeral and glenoid constructs.

In the Anatomic shoulder the humeral construct consists of the humeral stem, the humeral body, the adaptor taper and the humeral head. In the Reverse shoulder the humeral construct consists of the humeral stem, the reverse humeral body and the reverse liner. On the humeral side the fixation of the humeral stem determines if the construct is cemented or uncemented. The Anatomic glenoid construct consists of an all polyethylene glenoid or a metal back assembled with a liner while the Reverse glenoid construct consists of the metal back, the connector and the glenosphere. On the glenoid side, the fixation of the all-polyethylene glenoid or metal back determines if the construct is cemented or uncemented. The peripheral pegs of the SMR TT Hybrid Glenoid are intended for fixation by means of bone cement, while the central peg is intended for cementless fixation. The SMR TT Hybrid Glenoid Reverse Baseplate is indicated for conversion to reverse shoulder replacement of a previously implanted SMR TT Hybrid Glenoid; this device cannot be used for primary surgeries.

More information on the cemented or uncemented use of the components of SMR System is available in the following table.

| System |   | Components                              | Material | Raw Material Reference (**) | Use |    | Available in Canada |
|--------|---|-----------------------------------------|----------|-----------------------------|-----|----|---------------------|
| A      | R |                                         |          |                             | C   | NC |                     |
| •      | • | SMR Stems (Cemented, Cemented Revision) | Ti6Al4V  | M1                          | X   |    | •                   |

|   |   |                                                    |              |         |   |   |   |
|---|---|----------------------------------------------------|--------------|---------|---|---|---|
| • | • | SMR Stems (Cementless Finned, Cementless Revision) | Ti6Al4V      | M1      |   | X | • |
| • | • | SMR Short Stems (Cementless Finned)                | Ti6Al4V      | M1      |   | X |   |
| • | • | SMR Large Resection Stems                          | Ti6Al4V      | M1      | X |   | • |
| • | • | SMR Modular Spacers                                | Ti6Al4V      | M1      | X |   | • |
| • |   | SMR Humeral Bodies (Trauma, Finned)                | Ti6Al4V      | M1      | X | X | • |
| • | • | SMR Reverse Humeral Body                           | Ti6Al4V      | M1      | X | X | • |
|   |   |                                                    | Ti6Al4V + HA | M1, M13 | X | X | • |
| • | • | SMR 140° Reverse Humeral Body                      | Ti6Al4V      | M1      | X | X |   |
|   |   |                                                    | Ti6Al4V + HA | M1, M13 | X | X |   |
|   | • | SMR Humeral Extension                              | Ti6Al4V      | M1      | X | X | • |
| • |   | SMR Humeral Heads (Standard, CTA)                  | CoCrMo       | M4      | X | X | • |
|   |   |                                                    | Ti6Al4V      | M1      | X | X | • |
| • |   | SMR Adaptor Tapers (Neutral, Eccentric)            | Ti6Al4V      | M1      | X | X | • |

|   |   |                                               |                                                |        |   |   |   |
|---|---|-----------------------------------------------|------------------------------------------------|--------|---|---|---|
| • |   | SMR CTA Head Adaptor for Reverse Humeral Body | Ti6Al4V                                        | M1     | X | X | • |
|   | • | SMR Glenospheres                              | CoCrMo                                         | M4     |   | X | • |
|   |   |                                               | Ti6Al4V                                        | M1     |   | X | • |
|   | • | SMR Reverse HP Glenosphere                    | UHMWPE X-Lima + Ti6Al4V                        | M1, M8 |   | X | • |
|   | • | SMR Connectors with Screw                     | Ti6Al4V + UHMWPE                               | M1, M8 |   | X | • |
|   | • | SMR Reverse Liners                            | UHMWPE                                         | M8     | X | X | • |
|   |   |                                               | UHMWPE X-Lima                                  | M8     | X | X |   |
|   |   |                                               | LimaVit™ (Vitamin E highly crosslinked UHMWPE) | M9     | X | X |   |
|   | • | SMR Reverse HP Liners                         | CoCrMo                                         | M4     | X | X | • |
|   |   |                                               | Alumina                                        | M10    | X | X | • |
| • |   | SMR Cemented Glenoids                         | UHMWPE                                         | M8     | X |   | • |

|   |   |                                                    |                                   |                          |                 |   |   |   |
|---|---|----------------------------------------------------|-----------------------------------|--------------------------|-----------------|---|---|---|
| • |   | SMR Cemented Glenoids<br>3 Pegs                    |                                   | UHMWPE *                 | M8              | X |   | • |
|   |   |                                                    |                                   | UHMWPE<br>X-Lima         | M8              | X |   |   |
| • | • | SMR Metal Back Glenoid                             |                                   | Ti6Al4V +<br>PoroTi + HA | M1, M12,<br>M13 |   | X | • |
|   |   |                                                    |                                   | Ti6Al4V +<br>PoroTi FDA  | M1, M12         |   | X | • |
| • | • | SMR TT Baseplate                                   |                                   | Ti6Al4V                  | M1              |   | X | • |
| • | • | SMR Axioma TT Baseplate                            |                                   | Ti6Al4V                  | M1              |   | X |   |
| • | • | SMR Glenoid Peg TT                                 |                                   | Ti6Al4V 3D<br>printed    | M3              |   | X | • |
| • |   | SMR Metal Back Liner                               |                                   | UHMWPE                   | M8              |   | X | • |
| • |   | SMR Axioma Metal Back<br>Liner                     |                                   | UHMWPE                   | M8              |   | X |   |
|   | • | SMR Glenoid Plate                                  |                                   | Ti                       | M7              |   | X |   |
| • |   | SMR TT<br>Hybrid<br>Glenoid                        | Peripheral pegs<br>(subcomponent) | UHMWPE + Ta              | M8, M11         | X |   | • |
|   |   |                                                    | Central Peg<br>(subcomponent)     | Ti6Al4V 3D<br>printed    | M3              |   | X |   |
|   | • | SMR TT Hybrid Glenoid<br>Reverse Baseplate + Screw |                                   | Ti6Al4V +<br>UHMWPE      | M1, M8          |   | X | • |

|                                                                                                                                                                                                                                                                                                                                                                                                                                                        |   |                                |         |    |  |   |   |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---|--------------------------------|---------|----|--|---|---|
|                                                                                                                                                                                                                                                                                                                                                                                                                                                        | • | SMR TT Augmented 360 Baseplate | Ti6Al4V | M1 |  | X | • |
| Material Standards                                                                                                                                                                                                                                                                                                                                                                                                                                     |   |                                |         |    |  |   |   |
| <b>Ti6Al4V</b> (ISO 5832-3 – ASTM F1472) – <b>Ti6Al4V 3D printed</b> (ISO 5832-3) – <b>CoCrMo</b> (ISO 5832-12 – ASTM F1537) – <b>Ti</b> (ASTM F67) – <b>UHMWPE</b> (ISO 5834-2 – ASTM F648) – <b>LimaVit™</b> (Vitamin E highly crosslinked UHMWPE) (ISO 5834-2 - ASTM F648 - ASTM F2695 – ASTM F2565) – <b>Alumina</b> (ISO 6474) – <b>Ta</b> (ISO13782 - ASTM F560) – <b>PoroTi Titanium Coating</b> (ASTM F1580) – <b>HA coating</b> (ISO 13779-6) |   |                                |         |    |  |   |   |

A = Anatomic; R = Reverse; C = Cemented; NC = Not Cemented

**\* Not available in EU.**

**\*\* For the chemical composition of the raw materials please refer to Section “Raw Material Compositional Limits (%)” considering the raw material reference indicated in the table.**

## 2.2. CONTRAINDICATIONS

Absolute contraindications include: • local or systemic infection; • septicemia; • persistent acute or chronic osteomyelitis; • confirmed nerve lesion compromising shoulder joint function; • deltoid muscle insufficiency.

Relative contraindications include: • vascular or nerve diseases affecting the concerned limb; • poor bone stock (for example due to osteoporosis or extended previous revision surgery) compromising the stability of the implant; • metabolic disorders which may impair fixation and stability of the implant; • any concomitant disease and dependence that might affect the implanted prosthesis; • metal hypersensitivity to implant materials.

Relative contraindications related to the Sub-Scapularis Sparing Technique include: • revision surgery; • obesity; • medial glenoid erosion; • severe rotator cuff tendinopathy; • internal rotation contracture.

In cases of bone tumors, use an appropriate system designed to treat cases requiring large bone resections (SMR Large Resections stems). The use of primary or revision implants not designed and intended for use in cases of bone resection may result in a poor outcome and / or failure of the implant or implant fixation.

## 2.3. RISK FACTORS

The following risk factors may result in poor results with this prosthesis: • overweight<sup>1</sup>; • strenuous physical activities (active sports, heavy physical work); • incorrect implant positioning; • wrong size of components; • muscle deficiencies; • multiple joint disabilities; • refusal to modify postoperative physical activities; • patient's history of infections or falls; • systemic diseases and metabolic disorders; • local or disseminated neoplastic diseases; • drug therapies that adversely affect bone quality, healing, or resistance to infection; • drug use or alcoholism; • marked osteoporosis or osteomalacia; • patient's resistance to disease generally weakened (HIV, tumour, infections); • severe deformity leading to an impaired anchorage or improper positioning of implants; • use in combination with products, prosthesis or instruments of another manufacturer; • errors of operative technique.

## 3. WARNINGS

### 3.1. PRE-OPERATIVE PLANNING

LimaCorporate products shall be implanted only by surgeons familiar with the joint replacement procedures described in the specific surgical techniques.

### COMBINATIONS ALLOWED/NOT ALLOWED

Note: With **CTA heads** the use of **Trauma Humeral Bodies** is recommended to avoid possible impingement between the head and the body when using the Finned Humeral Body and the eccentric tapers in the lower position.

NOTE: the **size large Metal Back Glenoid** is not suitable for coupling with the 36 mm Glenospheres and the 40 mm Glenospheres.

NOTE: The Dia. 50, 52 and 54 mm **Humeral Heads** with + 3mm increased height cannot be coupled to the Long Adaptor Tapers (both concentric and eccentric).

The Dia. 52 and 54 mm Humeral Heads with + 2mm increased height cannot be coupled to the Long Adaptor Tapers (both concentric and eccentric).

NOTE: The SMR Metal Backs are compatible with bone screws dia.6,5mm and dia. 5 mm with minimum length 20 mm.

NOTE: The SMR Glenoid Plates are compatible with SMR Reverse HP Glenospheres only.

---

<sup>1</sup> According to the definition of World Health Organization (WHO), Body Mass Index (BMI) greater than or equal to 25 Kg/m<sup>2</sup>

## NOTES ON GLENOID BONE GRAFT

The Glenoid Bone Graft Instruments are intended to be used only with Axioma TT Metal Back and L1 TT Metal Back.

For glenoid bone grafting techniques, it is important to use a peg size Medium, Long or X-Long in order to enable a minimum section of the peg into the native bone, providing component stability. The following table identifies the allowed (✓) / not allowed (X) combinations between the bone graft thickness and the peg dimensions:

|                      |                  | TT METAL-BACK - PEG SIZE |        |      |        |             |        |      |        |
|----------------------|------------------|--------------------------|--------|------|--------|-------------|--------|------|--------|
|                      |                  | SMALL-R                  |        |      |        | SMALL / STD |        |      |        |
|                      |                  | Short                    | Medium | Long | X-Long | Short       | Medium | Long | X-Long |
| Bone Graft Thickness | 5 mm GRAFT       | X                        | ✓      | ✓    | ✓      | X           | ✓      | ✓    | ✓      |
|                      | 10 mm GRAFT      | X                        | X      | ✓    | ✓      | X           | ✓      | ✓    | ✓      |
|                      | 15 mm GRAFT      | X                        | X      | ✓    | ✓      | X           | X      | ✓    | ✓      |
|                      | 15° SLOPED GRAFT | X                        | ✓      | ✓    | ✓      | X           | ✓      | ✓    | ✓      |
|                      | 20° SLOPED GRAFT | X                        | ✓      | ✓    | ✓      | X           | ✓      | ✓    | ✓      |

## NOTES ON SUB-SCAPULARIS SPARING TECHNIQUE (not approved in Canada)

The Sub-Scapularis Sparing Technique can be used to implant only the SMR Glenoid Peg TT, SMR TT Metal Back glenoid, SMR Axioma TT Metal Back Glenoid, the SMR Metal Back Glenoid, SMR Axioma Metal Back Glenoid, Bone screws and SMR Metal Back Liner in primary SMR Total Anatomic configuration replacements. In case it is found intraoperatively that the reverse configuration is required, the Sub-Scapularis Sparing Technique can be used to implant only SMR Glenoid Peg TT, SMR TT Metal Back glenoid, SMR Axioma TT Metal Back Glenoid, the SMR Metal Back Glenoid, SMR Axioma Metal Back Glenoid, Bone screws. The standard surgical technique for the SMR Reverse Shoulder is to be used to implant the other components.

The Sub-Scapularis Sparing Technique is not to be used for the implantation of SMR CTA heads or for cases in which bone grafting of glenoid defects is necessary.

**NOTES ON SMR TT HYBRID GLENOID AND SMR TT HYBRID GLENOID REVERSE BASEPLATE**

The following table identifies the allowed (✓) / not allowed (X) combinations between the SMR TT Hybrid Glenoid and the SMR Humeral Heads:

|                     |              | SMR HUMERAL/RESURFACING HEADS DIAMETER (MM) |    |    |    |    |    |    |    |
|---------------------|--------------|---------------------------------------------|----|----|----|----|----|----|----|
|                     |              | 40                                          | 42 | 44 | 46 | 48 | 50 | 52 | 54 |
| Hybrid Glenoid Size | SMALL LOW    | ✓                                           | ✓  | ✓  | ✓  | ✓  | X  | X  | X  |
|                     | SMALL        | ✓                                           | ✓  | ✓  | ✓  | ✓  | ✓  | ✓  | ✓  |
|                     | STANDARD LOW | ✓                                           | ✓  | ✓  | ✓  | ✓  | ✓  | ✓  | X  |
|                     | STANDARD     | ✓                                           | ✓  | ✓  | ✓  | ✓  | ✓  | ✓  | ✓  |
|                     | LARGE LOW    | ✓                                           | ✓  | ✓  | ✓  | ✓  | ✓  | ✓  | ✓  |
|                     | LARGE        | ✓                                           | ✓  | ✓  | ✓  | ✓  | ✓  | ✓  | ✓  |

If a SMR TT Hybrid Glenoid is in place and revision to a reverse prosthesis is required, the patient can be revised by removing the polyethylene baseplate, leaving the metal peg in place and by connecting it to the SMR TT Hybrid Glenoid Reverse Baseplate.

The SMR TT Hybrid Glenoid Reverse Baseplate is not suitable if one of the following conditions arises:

- Peg not properly seated into the cavity, that comes out from the glenoid surface;
- Peg not stable in glenoid bone;
- Bone loss;
- Bone graft is required.

The SMR TT Hybrid Glenoid Reverse Baseplate is intended to be coupled with the SMR TT Hybrid Glenoid metal peg only.

The SMR TT Hybrid Glenoid Reverse Baseplate is intended to be used in combination with the following SMR Glenospheres:

| Component                                    | Material                | Product code |
|----------------------------------------------|-------------------------|--------------|
| Glenosphere dia. 36 mm                       | CoCrMo                  | 1374.09.111  |
|                                              | Ti6Al4V                 | 1374.15.111  |
| Eccentric Glenosphere dia. 36 mm             | CoCrMo                  | 1376.09.031  |
|                                              | Ti6Al4V                 | 1376.15.031  |
| Reverse HP Glenosphere dia. 40 mm            | UHMWPE X-LIMA + Ti6Al4V | 1374.50.400  |
| Reverse HP Glenosphere dia. 44 mm            |                         | 1374.50.440  |
| Reverse HP Corrective Glenosphere dia. 44 mm |                         | 1374.50.444  |

The SMR TT Hybrid Glenoid Reverse Baseplate must not be used in combination with the SMR Glenospheres with connector and screw, belonging to the SMR Shoulder System.  
The SMR TT Hybrid Glenoid Reverse Baseplate cannot be used with SMR Glenoid Plates.

#### **SMR TT AUGMENTED 360 BASEPLATE**

The SMR TT Augmented 360 Baseplate is not suitable when bone grafting technique is needed.

The SMR TT Augmented 360 Baseplate is compatible only with glenosphere connectors size Small-R and TT Glenoid Peg sizes Small-R Medium, S-R Long and S-R X-Long.

The SMR TT Augmented 360 Baseplate is compatible with cortical bone screws dia. 4,5mm (product code 8431.15.XXX) with minimum length 26mm and must be fixed to the glenoid bone with a minimum of three bone screws.

The SMR TT Augmented 360 Baseplate is compatible with bone screws dia. 6,5mm with minimum length 25mm and must be fixed to the glenoid bone with a minimum of two bone screws.

The SMR TT Augmented 360 Baseplate is compatible with cortical bone screws dia. 5mm (product code 8432.15.XXX) and must be fixed to the glenoid bone with a minimum of three bone screws (two bone screws with minimum length 26mm and one with minimum length 18mm). Bone screws must be positioned according to the surgical technique.

The SMR TT Augmented 360 Baseplate cannot be used with SMR Glenoid Plates and SMR Lateralized Connector with Screw.

### SMR LATERALIZED CONNECTORS WITH SCREWS

SMR Lateralized Connectors are not indicated for use with glenoid bone grafting techniques and cannot be used with SMR TT Augmented 360 Baseplate.

The following table identifies the allowed (✓) / not allowed (X) combinations between the SMR Lateralized Connectors and the SMR Glenospheres:

|                                       | LATERALIZED CONNECTOR |         |         |         |         |
|---------------------------------------|-----------------------|---------|---------|---------|---------|
|                                       | LAT+1mm               | LAT+2mm | LAT+3mm | LAT+4mm | LAT+5mm |
| Glenosphere Ø36                       | ✓                     | ✓       | ✓       | ✓       | ✓       |
| Glenosphere Ø36 Eccentric             | X                     | X       | X       | X       | X       |
| Reverse HP Glenosphere Ø40            | ✓                     | ✓       | ✓       | ✓       | X       |
| Reverse HP Glenosphere Ø44            | ✓                     | ✓       | X       | X       | X       |
| Reverse HP Glenosphere Ø44 Corrective | X                     | X       | X       | X       | X       |

### SMR 140° REVERSE HUMERAL BODY

SMR 140° Reverse humeral bodies are indicated for use only with 36mm glenospheres and the related 36mm reverse humeral liner.

**REVERSE LINER LimaVit™ (Vitamin E highly crosslinked UHMWPE)** cannot be coupled with Ti6Al4V glenospheres.

The components of the system must not be used with components of other manufacturer, unless specified by LimaCorporate. The components forming original LimaCorporate systems must be rigorously coupled following the given information included in this document and in the available product-specific information (e.g.: current version of surgical technique). The use of single components, or components belonging to other systems, are subject to the approval by LimaCorporate. The manufacturer and the dealer are not liable for possible coupling incompatibility. The surgeon is solely responsible for the choice and use of the implant.

Pre-operative planning, through radiographic templates in different formats, provides essential information regarding the type and size of components to be used and the correct combination of required devices based on the anatomy and specific conditions of each patient. Inadequate pre-operative planning can lead to improper selection of the implants and/or incorrect implant positioning. In selecting patients for surgery, the following factors can be critical to the eventual success of the procedure:

- **Partial Shoulder replacement:** in cases with a deficient and unreconstructable rotator cuff, a CTA-head is indicated;
- **Total Shoulder replacement:** the rotator cuff must be intact or reconstructable. In cases with a deficient and unreconstructable rotator cuff, a hemiprosthesis with a CTA head or a Reverse Total Shoulder Arthroplasty is indicated;
- **Reverse Shoulder replacement:** the bone stock of the glenoid and humerus must be able to support the implant. In cases with significant bone loss and in which adequate fixation on the glenoid side cannot be obtained, a hemiarthroplasty with a CTA-head should be performed.
- **Bone Grafting:** once removed the graft should be inspected to ensure the bone quality is adequate for use with the glenoid bone grafting technique. The glenoid bone grafting technique should never be used with poor quality bone, as it may compromise bone healing.

LimaCorporate specialized technical staff is available to provide advice regarding pre-operative planning, the surgical technique, and product and instrumentation assistance both prior to and during surgery.

The patient shall be warned that the prosthesis does not replace normal healthy bone, that the prosthesis can break or become damaged as a result of certain activity or trauma, that it has a finite expected implant life, and may need to be replaced at some time in the future: the longevity and the performance of the implant over time, in fact, can be influenced by the natural/physiological progression of the patient pathology, the insurgence of co-morbidities and post-operative complications, even if they are not directly implant-related (i.e. infection pain, stiffness, reduced range of motion).

The possible impact of the factors mentioned in sections "GENERAL INFORMATION ON INDICATIONS, CONTRAINDICATIONS AND RISK FACTORS" and "POSSIBLE ADVERSE EFFECTS" shall be considered preoperatively and the patient informed as to what steps he/she can take to reduce the possible effects of these factors.

**Implants are single use devices.**

**Do not re-use implants that were previously implanted in another patient.**

**Do not re-use an implant that has previously come into contact with the body fluid or tissue of another person.**

**Do not use any component from a package that has been previously opened or appears to be damaged.**

**Do not use implants after the expiration date printed on the label.**

Surgical instruments used to perform the surgery are subject to wear with normal usage. After extensive use or excessive loads, instruments are susceptible to fracture. Surgical instruments shall be used only for their specific purpose. Before use, the functionality of surgical instruments shall be checked since the use of damaged instruments may lead to early failure of the implants. Damaged instruments shall be replaced before surgery.

### **3.2. INTRAOPERATIVE**

The use of trial devices is recommended to check the correct site preparation, size and positioning of the implants to be used. It is recommended that additional implants be available during surgery for use in those cases requiring prostheses of different sizes or when the preoperatively selected prostheses cannot be used. The correct selection as well as the correct seating/placement of the implant is extremely important. Improper selection, positioning, alignment and fixation of the implant components may result in unusual stress conditions which may negatively affect system performance and the survival rate of the implant.

The components forming original LimaCorporate systems must be assembled according to the surgical technique and used only for the labelled indications.

Use only instruments and trial components specifically designed for use with the implants being used. The use of instruments from other manufacturers or the use of instruments designed for use with other systems can lead to inappropriate preparation of the implant site, incorrect positioning, alignment and fixation of the devices followed by loosening of the system, loss of functionality, reduction of the durability of the implant, and the need for further surgery.

Care must be taken to protect the surfaces involved in the coupling between components (tapers); the articular surfaces of the implants shall be protected from scratches or any other damage. All component coupling surfaces shall be clean and dry before assembly. The stability of component couplings shall be verified as described in the surgical technique. Care must be taken when handling Trabecular Titanium devices; they must not get in contact with fabrics that can release particles inside the trabecular structure.

### **3.3. POST-OPERATIVE CARE**

Adequate post-operative care shall be provided by the surgeon or other suitably qualified medical staff. Regular postoperative x-ray follow-up is recommended to detect any changes in position or condition of the implant or surrounding tissues. Complications and/or failure of prosthetic implants are more likely to occur in patients with unrealistic functional expectations, heavy patients, physically active patients, and/or with patients who fail to follow through with the required rehabilitation program. Excessive physical activity or trauma to the replaced joint can lead to premature failure of the arthroplasty through loosening, fracture, or abnormal wear of the prosthetic implants. The patient shall be cautioned by the surgeon to govern activities accordingly and advised that the implants may fail due to excessive joint wear.

The surgeon shall make the patient aware of the limitation of limb function after shoulder arthroplasty and that caution is needed, especially in the post-operative period, since the reconstructed joint must be protected from full load bearing. In particular the following precautions shall be presented to the patient by the surgeon:

- avoid repeated high weightlifting;
- keep body weight under control;
- avoid sudden peak loads (consequences of activities such contact sports, playing tennis) or movements which can lead to sudden stops or twisting;
- avoid positions that can increase the risk of dislocation, such as for example degrees of abduction and/or internal rotation.

The patient shall be informed also about the possibility that the implant or its components may wear out, fail or need to be replaced. The implant may not last the rest of the patient's life, or any particular length of time.

Because prosthetic implants are not as strong, reliable, or durable as natural, healthy tissue/bones, all such devices may need to be replaced at some point.

Appropriate post-operative rehabilitation instructions and care must be provided to avoid a negative influence on the outcome of the surgical procedure.

### **3.4. POSSIBLE SIDE EFFECTS**

Adverse effects that most commonly and frequently occur in orthopaedic prosthesis procedure include:

- loosening of the prosthetic components, for example due to osteolysis;
- prosthesis dislocation and instability;
- breakage or damage of the device;
- stem subsidence;
- instability of the system because of inadequate soft tissue balancing;
- dissociation due to incorrect coupling of the devices;
- infection;
- local hypersensitivity;
- local pain;
- periprosthetic fractures, including intraoperative fractures;
- bone resorption;
- temporary or permanent nerve or muscle damage;
- fretting of modular junctions;
- excessive wear of UHMWPE components due to damaged articular surfaces or the presence of particles;
- metallosis;
- heterotopic ossification;
- scapular notching;
- rotator cuff failure;
- glenoid erosion (only in hemi procedures);
- loss of range of motion;
- blood loss;
- additional surgery.

Some adverse effects can lead to death. General complications include venous thrombosis with/without pulmonary embolism, cardiovascular or pulmonary disturbances, haematomas, systemic allergic reactions, systemic pain, blood loss, donor site morbidity, temporary or permanent paralysis.

## **4. STERILITY**

### **a. Implants**

All implantable components of the system are provided sterile with a Sterility Assurance Level (SAL) of  $10^{-6}$ . Metal components are sterilized by radiation (E-Beam or Gamma rays) or EtO while UHMWPE components by EtO. HA coated implants are sterilized only by radiation (E-Beam or Gamma rays). Do not use any component from a package that has been previously opened or appears to be damaged. **Do not use implants after the expiration date printed on the label.**

### **b. Instruments**

Instruments are supplied non-sterile and must be cleaned, disinfected and sterilized before use according to appropriate validated methods (refer to the "Instrument Care, Cleaning, Disinfection and Sterilization" brochure for validated sterilization parameters; this brochure is available on request or downloadable from [www.limacorporate.com](http://www.limacorporate.com) in the Products section). Users shall validate their specific cleaning, disinfection and sterilization processes and equipment.

## 5. INFORMATION USEFUL TO INFORM THE PATIENT

Refer to previous sections for any information that allows the patient to be informed of any warnings, precautions, contra-indications, measures to be taken and limitations of use regarding the device. The surgeon shall notify the patient about any necessary follow-up specific for the patient's condition.

### 5.1. MAGNETIC RESONANCE IMAGING (MRI)

Non-clinical testing has demonstrated that LimaCorporate **SMR System** is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 Tesla and 3 Tesla, with
- Maximum spatial field gradient of 8,500 G/cm (85 T/m)
- Maximum force product of 154,000,000 G<sup>2</sup>/cm (154 T<sup>2</sup>/m)
- Theoretically estimated maximum whole body averaged (WBA) specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode)

#### RF Heating

Under the scan conditions defined above, the SMR Shoulder System is expected to produce a maximum temperature rise of less than:

- 6.0°C (2 W/kg, 1.5 Tesla) RF-related temperature increases with a background temperature increase of: ≈0.4°C (2 W/kg, 1.5 Tesla, SMR anatomical shoulder system) after 5 minutes of continuous scanning; ≈0.2°C (2 W/kg, 1.5 Tesla, SMR reverse shoulder system) after 4.3 minutes of continuous scanning;
- 6.0°C (2 W/kg, 3 Tesla) RF-related temperature increases with a background temperature increase of: ≈0.6°C (2 W/kg, 3 Tesla, SMR anatomical shoulder system) after 12 minutes of continuous scanning; ≈0.9°C (2 W/kg, 3 Tesla, SMR reverse shoulder system) after 12.5 minutes of continuous scanning.

#### MR Artifact at 3T

In non-clinical testing, the image artifact caused by the device extends approximately 98.9 mm (SMR anatomical shoulder system) / 72.9 mm (SMR reverse shoulder system) from the SMR Shoulder System when imaged with a gradient echo pulse sequence and a 3 Tesla MR system.

**NOTE:** The SMR System has been tested for MR environment in combination with bone screws dia 6.5mm lengths 20+35 mm (codes 8420.15.010, 8420.15.020, 8420.15.030, 8420.15.040).

The use of SMR TT Augmented 360 Baseplate, SMR Cementless Finned Short Stem, SMR Lateralized Connectors with Screw, SMR Reverse LimaVit™ (Vitamin E highly crosslinked UHMWPE) Liners and SMR 140° Reverse Humeral Body and the use of the SMR System in combination with other bone screws have not been evaluated for safety and compatibility in the MR environment.

These devices have not been tested for heating or migration in the MR environment. The risks associated with a passive implant in an MR environment have been evaluated and are known to include heating, migration and image artifacts at or near implant side.

## 6. INCIDENT REPORTING

According to the definition of incident/serious incident given by Regulation EU 2017/745<sup>2</sup>, users and/or patients are required to: **• report to the manufacturer and to the relevant Competent Authority any serious incident that has occurred in relation to the device; • report to the manufacturer any incident and/or event of malfunction of the device or changes in its performance that may affect safety.**

In non-UE market, refer to incidents' definitions according to applicable laws.

---

<sup>2</sup> an incident is defined to be any malfunction or deterioration in the characteristics and/or performance of the device, including use-error due to ergonomic features, as well as any inadequacy in the information supplied by the manufacturer and any undesirable side-effect. An incident is considered to be serious when it directly or indirectly led, might have led or might lead to any of the following: (a) the death of a patient, user or other person, (b) the temporary or permanent serious deterioration of a patient's, user's or other person's state of health, (c) a serious public health threat.

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