INSTRUCTIONS FOR USE – REVISION STEM

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INSTRUCTION FOR USE – REVISION STEM INTENDED FOR UNCEMENTED FIXATION

Before using a LimaCorporate product, the surgeon should carefully study the following recommendations, warnings and instructions, as well as the available product-specific information (e.g. product literature, surgical technique).

1. PRODUCT INFORMATION

The LimaCorporate Revision hip stem is indicated for cementless revision hip arthroplasty of either an uncemented or cemented previous femoral hip implant. The main goal of a joint prosthesis is to reproduce the articular anatomy, in part or totally. The joint prosthesis is intended to reduce pain and give articular mobility to the patient. The degree of pain reduction and mobility depends in part on the pre-operative situation, intra-operative options and post-op rehabilitation. Revision is a modular system which includes femoral stems and modular necks.

1.1. MATERIALS

Device	Material
Femoral stems	Ti6Al4V
Modular necks	Ti6Al4V
	UHMWPE
Applicable Standards	
Ti6Al4V (ISO 5832-3 - ASTM F1472) – UHMWPE (ISO 5834-2- ASTM F648)	

Some patients may be hypersensitive or allergic to the implant materials; this should be appropriately considered by the surgeon.

1.2. HANDLING AND STORAGE

All devices are provided sterile and should be stored at ambient temperature (indicative range 0-50°C / 32-122°F) in their protective closed packaging in controlled rooms, protected from exposure to light, heat and sudden changes in temperature. Once the package is opened, ensure that both the model and size of the

implant correspond exactly to the description printed on the labels. Avoid any contact between the implant and objects or substances that can alter the sterile condition or the surface integrity. Careful visual examination of each implant is recommended before use in order to verify that the implant is not damaged. **Components removed from the package should not be used if they are dropped or suffer other accidental trauma. Devices should not be modified in any way.** The device's code and lot number should be recorded in the patient's case history by using the labels included in the component packaging.

The disposal of medical devices and their packaging is to be performed by the hospitals in conformity with applicable laws.

Re-Use of previously implanted devices must be absolutely avoided. Risks associated with reuse of single use devices are: • infection; • early or late failure of the device or device fixation; • lack of appropriate coupling between modular junctions (e.g. taper connections); • device wear and wear debris associated complications; • transmission of diseases (e.g. HIV, hepatitis); • immune system response/rejection.

2. GENERAL INFORMATION ON INDICATIONS, CONTRAINDICATIONS AND RISK FACTORS

2.1. INDICATIONS

The LimaCorporate Revision hip stem is indicated for cementless revision hip arthroplasty of either an uncemented or cemented previous femoral hip implant. In particular, it is indicated in case of: • Revision procedures due to trauma or peri-prosthetic fracture; • Revision procedures due to aseptic or septic loosening; • Revision procedures in which other treatments or devices have failed.

2.2. CONTRAINDICATIONS

Absolute contraindications include: • local or systemic infection; • septicaemia; • persistent acute or chronic osteomyelitis; • confirmed nerve or muscle lesion compromising hip joint function.

Relative contraindications include: • vascular or nerve diseases affecting the concerned limb; • poor bone stock compromising the fixation and the stability of the implant; • metabolic disorders which may impair fixation and stability of the implant and bone formation; • any concomitant disease and dependence that might affect the implanted prosthesis; • metal hypersensitivity to implant materials.

2.3. RISK FACTORS

The following risk factors may result in poor results with this prosthesis: • overweight; • heavy physical activities in which the prosthesis can be excessively overloaded (i.e. physical labor, heavy work, active sports, etc); • fretting of modular junctions; • incorrect implant positioning (e.g. varus positioning); • wrong size of components; • muscle deficiencies; • multiple joint disabilities; • medical disabilities which can lead to an unnatural gait and loading of the hip joint; • refusal to modify postoperative physical activities; • patient's history of infections or falls; • systemic diseases and metabolic disorders; • local or disseminated neoplastic diseases; • drug therapies that adversely affect bone quality, healing, or resistance to infection; • drug addiction

and/or drug, alcohol or medicine abuse; • marked osteoporosis or osteomalacia; • patient's resistance to disease generally weakened (HIV, tumour, infections); • severe deformity leading to impaired anchorage or improper positioning of implants; • errors of operative technique.

3. WARNINGS

3.1. PRE-OPERATIVE PLANNING

LimaCorporate products should be implanted only by surgeons familiar with the joint replacement procedures described in the specific surgical techniques.

COMBINATIONS ALLOWED/NOT ALLOWED

Allowed combination between LimaCorporate Femoral heads and Revision femoral necks: only the sizes S, M and L can be coupled with Revision necks.

Use of femoral heads with greater neck lengths may result in failure of the hip stem (e.g. breakage due to fatigue).

The combinations of diameters 14 mm and 15 mm Revision stems with Lateralized Revision necks are not allowed.

Revision stems must not be used with components from other systems or other manufacturers with exception of the following possible combination:

For Revision stems, femoral heads manufactured by LimaCorporate made of • Biolox Delta dia. 28mm taper 12/14 (S, M & L); • CoCrMo dia. 28mm taper 12/14 (S, M & L); • AISI 316/L dia. 28mm taper 12/14 (S, M & L); • FeCrNiMnMoNbN dia. 22mm (-2, 0 & +4); • FeCrNiMnMoNbN dia. 28mm taper 12/14 (S, M & L); • Revision Biolox Delta dia. 28mm taper 12/14 (S, M & L), can be also coupled with Ø44mm to 68mm Symbol Cup (DM CEM, DM HA, DMR HA) in combination with dual mobility PE Liner (for Ø22.2mm head: Ø44mm to 68mm; for Ø28mm head: Ø46mm to 68mm) manufactured by Dedienne Santé.

The cup/liner assemblies come in corresponding sizes; the liner/head assemblies come in corresponding sizes, considering the internal diameter of the liner.

The surgeon must consult also SYMBOL CUP DM instruction for use before using the mixed allowed coupling between LimaCorporate and Dedienne Santé components.

Pre-operative planning, through radiographic templates in different formats (the templates and x-rays must have the same magnification), provides useful, but indicative information regarding the type and size of components to be used and the correct combination of devices required based on the anatomy and specific conditions of each patient. Inadequate pre-operative planning can lead to improper selection of the implants and/or incorrect implant positioning.

The surgeon should carefully plan surgery considering the following:

1. **Smaller sized femoral implants** (e.g. 14, 15 and 16 mm stems): the smaller sized femoral implants are designed for patients with a small intramedullary canal and/or diaphyseal region of the femur. The reduced size (diameter) of these stems results in a corresponding reduction in the fatigue strength of the implant;
2. **High Offset combinations** (use of lateralized modular necks): the lateralized necks are designed to restore the functional offset of the hip joint, however greater neck lengths are accompanied by a higher risk of failure (e.g. breakage due to fatigue).

Complications or failures of the total hip replacement are more likely to occur in heavy and highly active patients and high offset combinations.

The surgeon should perform a careful evaluation of the patient's clinical condition and level of physical activity before performing hip replacement. **Patients who are overweight (BMI >25 kg/m²) or have high activity levels may not be candidates for a modular hip replacement.** In case of bone loss or insufficient femoral bone stock, bone grafting or other adjunctive reinforcement procedures is advisable to provide proximal support to the stem. This is necessary because, without proximal support, the stem is vulnerable to fracture. If proximal support cannot be achieved, an alternative surgical option should be considered. If proximal support is weak, the patient should be warned of an increased risk for potential fatigue fracture. LimaCorporate specialized technical staff is available to provide advice regarding pre-operative planning, the surgical technique, and product and instrumentation assistance both prior to and during surgery. The patient should be warned that the prosthesis does not replace normal healthy bone, that the prosthesis can break or become damaged as a result of certain activity or trauma, that it has a finite expected implant life, and may need to be replaced at some time in the future. The possible impact of the factors mentioned in sections 2 and 3.4 should be considered preoperatively and the patient informed as to what steps he/she can take to reduce the possible effects of these factors. **Implants are single use devices; do not re-use implants that were previously implanted in another patient. Do not re-use an implant that has previously come into contact with the body fluid or tissue of another person.** Surgical instruments are subject to wear with normal usage. After extensive use or excessive loads, instruments are susceptible to fracture. Surgical instruments should be used only for their specific purpose. Before use, the functionality of surgical instruments should be checked since the use of damaged instruments may lead to early failure of the implants. Damaged instruments should be replaced before surgery.

3.2. INTRAOPERATIVE

The use of trial devices is recommended to check the correct site preparation, size and positioning of the implants to be used. It is recommended that additional implants are available during surgery for use in those cases requiring prostheses of different sizes or when the preoperatively selected prostheses cannot be used. The correct selection as well as the correct seating/placement of the implant is extremely important. Improper selection, positioning, alignment and fixation of the implant components may result in unusual stress conditions which may negatively affect system performance and survival rate of the implant. The components forming original LimaCorporate systems must be assembled according to the surgical technique and used only

for the labelled indications. Use only instruments and prosthesis trials specifically designed for use with the implants being used. The use of instruments from other manufacturers or the use of instruments designed for use with other systems can lead to inappropriate preparation of implant site, incorrect positioning, alignment and fixation of the devices followed by loosening of the system, loss of functionality, reduction of the durability of the implant, and the need for further surgery. Care must be taken to protect the surfaces involved in the coupling between components (tapers); articular surfaces of the implants should be protected from scratches or any other damage. All component coupling surfaces should be clean and dry before assembly. The stability of component couplings should be verified as described in the surgical technique.

3.3. POST-OPERATIVE CARE

Adequate post-operative care should be provided by the surgeon or other suitably qualified medical staff. Regular postoperative x-ray follow up is recommended to detect any changes in position of the implant or surrounding tissues. The surgeon should make the patient aware of the limitations of the limb function after hip arthroplasty and that the reconstructed joint must be protected from full load bearing for a period of time. Revision surgeries normally need a long rehabilitation period to restore joint functions. Excessive physical activity or trauma to the replaced hip joint can lead to premature failure of hip arthroplasty through loosening, fracture, or abnormal wear of the prosthetic implants. The patient should be cautioned by the surgeon to govern activities accordingly and advised that the implants may fail due to excessive joint wear.

In particular the following precautions should be presented to the patient by the surgeon: • avoid repeated high weight lifting; • avoid excessive or repeated stair climbing; • keep body weight under control (overweight conditions may adversely affect the outcomes of joint replacement); • avoid sudden peak loads (consequences of activities such as running and skiing) or movements which can lead to sudden stops or twistings; • avoid positions that can increase the risk of dislocation.

Lack of appropriate post-operative rehabilitation instructions and care can negatively influence the outcome of the surgical procedure.

3.4. POSSIBLE ADVERSE EFFECTS

Adverse effects can occur in hip arthroplasty include: • loosening of the prosthetic components; • prosthesis dislocation and instability; • fretting of modular junctions; • excessive wear of UHMWPE components; • breakage or damage of the prosthetic implant; • instability of the system because of inadequate soft tissues balancing; • dissociation due to incorrect coupling of the devices; • infection; • subsidence; • local hypersensitivity; • local pain; • periprosthetic fractures; • temporary or permanent nerve damage; • metallosis; • limb length discrepancies; • blood loss; • additional surgery.

Some adverse effects can lead to death. General complications include venous thrombosis with/without pulmonary embolism, cardiovascular or pulmonary disturbances, hematomas, systemic allergic reactions, systemic pain, temporary or permanent paralysis.

4. STERILITY

a. Implants

All implantable components of the Revision system are provided sterile with a Sterility Assurance Level (SAL) of 10^{-6} . Metal components are sterilized by irradiation or EtO and UHMWPE components by EtO. HA coated implants are sterilized only by irradiation. Do not use any component from a package that has been previously opened or appears to be damaged. **Do not use implants after the expiration date printed on the label.**

b. Instruments

Instruments are supplied non sterile and must be cleaned, disinfected and sterilized before use according to appropriate validated methods. (refer to the "Instrument Care, Cleaning, Disinfections and Sterilization" brochure for validated sterilization parameters; this brochure is available on request or downloadable from www.limacorporate.com in the Products section). Do not expose polyethylene to steam sterilization. Users should validate their specific cleaning, disinfection and sterilization processes and equipment.

5. MAGNETIC RESONANCE IMAGING (MRI)

There are inherent risks associated with the use of metallic implants in the MR environment; including component migration, heat induction and signal interference or distortion near the component(s). Heat induction of metallic implants is a risk related to the component geometry and material, as well as the MR power, duration and pulse sequence. Since MR equipment is not standardized, the severity and likelihood of occurrence are unknown for these implants. LimaCorporate performed tests to evaluate displacement force, torque, heating and artifacts produced by the interaction between the MRI and some configurations of LimaCorporate systems. The results indicate that the application of a MR up to 3T is not expected to raise concerns of migration and heating of the devices. Prosthetic systems include components which are passive metallic devices and as with all passive devices, there is potential for reciprocal interference with certain imaging modalities, including image distortion in MR and CT scans.

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Limacorporate S.p.A.