DEVICE DESCRIPTION

The Neuron™ MAX System is a three-component system comprised of the Neuron™ MAX Delivery Catheter, Neuron™ 6F Select Catheter, and a Dilator. The Neuron MAX Delivery Catheter can be used individually with a 0.038in [0.97 mm] guidewire or together with the Neuron 6F Select Catheter to access the desired anatomy.

Neuron MAX Delivery Catheter:

The Neuron MAX Delivery Catheter is a single lumen, braid-reinforced, variable stiffness catheter with a radiopaque zone on the distal end and a Luer hub on the proximal end. The Neuron MAX Delivery Catheter dimensions are included on the individual device label. The Neuron MAX Delivery Catheter is compatible with introducer sheaths appropriately sized for the outer diameter of Neuron MAX.

Neuron 6F Select Catheter:

The Neuron 6F Select Catheter is a single lumen, braid-reinforced, variable stiffness catheter with a radiopaque distal end and a Luer hub on the proximal end. The Neuron 6F Select Catheter is available in four tip shapes (SIMMONS, H1, BER, or SIM-V). The Neuron 6F Select Catheter tip shape and dimensions are included on the individual device label. The Neuron 6F Select Catheter is compatible with the Neuron MAX Delivery Catheter.

Dilator:

The Dilator is a single lumen, radiopaque catheter with a tapered distal end and a Luer hub on the proximal end. The Dilator is compatible with the Neuron MAX Delivery Catheter. The Dilator facilitates the percutaneous entry of the Neuron MAX Delivery Catheter by forming an atraumatic transition from the skin through the subcutaneous tissue to the vessel.

INDICATION FOR USE

The Neuron MAX System is indicated for the introduction of interventional devices into the peripheral, coronary, and neuro vasculature.

CONTRAINDICATIONS

There are no known contraindications.

WARNINGS

- The Neuron MAX System should only be used by physicians who have received appropriate training in interventional techniques.

PRECAUTIONS

- The device is intended for single use only. Do not resterilize or reuse. Resterilization and/or reuse may result in ineffective catheter coating lubrication, which may result in high friction and the inability to access the target vasculature location; and/or may compromise the structural integrity of the device.
- Do not use kinked or damaged devices. Do not use open or damaged packages. Return all damaged devices and packaging to the manufacturer/distributor.
- Use prior to the "Use By" date.
- Use the Neuron MAX System in conjunction with fluoroscopic visualization.
- Do not advance or withdraw the Neuron MAX System against resistance without careful assessment of the cause using fluoroscopy. If the cause cannot be determined, withdraw the device. Unrestrained moving or torquing the device against resistance may result in damage to the vessel or device.
- Maintain a constant infusion of an appropriate flush solution.
- If flow through the device becomes restricted, do not attempt to clear the lumen by infusion. Remove and replace the device.

POTENTIAL ADVERSE EVENTS

Possible complications include, but are not limited to, the following:

- | | | |
|-----------------------|---|--|
| • acute occlusion | • false aneurysm formation | • ischemia |
| • air embolism | • hematoma or hemorrhage at puncture site | • neurological deficits including stroke |
| • death | • infection | • vessel spasm, thrombosis, dissection, or perforation |
| • distal embolization | • intracranial hemorrhage | |
| • emboli | | |

DEVICE PREPARATION AND USE (SHEATH CONFIGURATION)

1. Select an appropriately sized Neuron MAX Delivery Catheter based on the anatomy and length.
2. Gently remove the Neuron MAX Delivery Catheter and packaging card from the pouch by grasping the Neuron MAX Delivery Catheter hub and packaging card and slowly pulling them out of the pouch.
3. Remove the Neuron MAX Delivery Catheter from the packaging card by removing the hub from the card tabs before gently removing the Neuron MAX Delivery Catheter shaft.
4. Gently remove the Dilator from the packaging card by removing the hub from the card tabs before gently removing the Dilator shaft.
5. Inspect the Neuron MAX Delivery Catheter for kinks or other damage.
6. Flush the lumen of the Neuron MAX Delivery Catheter with heparinized saline.
7. Insert the Dilator through the hemostasis valve of the Neuron MAX Delivery Catheter, snapping it into place at the hub.
8. Flush the lumen of the Dilator with heparinized saline.
9. Introduce the cannula of an angiographic needle into the vessel using aseptic technique. Holding the needle in place insert the flexible end of a 0.038in [0.97 mm] guidewire through the needle and into the vessel. Gently advance the guidewire to the desired depth.
10. Holding the guidewire in place, withdraw the needle and apply pressure to the puncture site until the Neuron MAX Delivery Catheter is inserted into the vasculature.
11. Thread the Neuron MAX Delivery Catheter/Dilator assembly over the guidewire, grasping the Delivery Catheter close to the skin to prevent buckling. Using a rotating motion, advance the assembly through the tissue into the vessel.
12. Under fluoroscopy, advance the Neuron MAX Delivery Catheter until the desired position is achieved.
13. Detach the Dilator from the Delivery Catheter by releasing the snap-fit ring at the hub. Withdraw the guidewire and Dilator. To avoid damage to the Delivery Catheter, do not withdraw the Dilator until the Delivery Catheter is in the vessel.
14. Aspirate from the Hemostasis Valve Adapter (HVA) to remove any potential air. After aspiration and flushing, consider establishing a heparinized solution or suitable isotonic solution via the HVA.
15. If using, select appropriate Neuron 6F Select Catheter shape based on the target vessel and its surrounding anatomy.
16. If using, insert the Neuron 6F Select Catheter into the Neuron MAX Delivery Catheter, and advance the Neuron 6F Select Catheter until the distal tip of the Neuron 6F Select Catheter is at the distal tip of the Neuron MAX Delivery Catheter.
17. If a 0.038in [0.97 mm] guidewire is used instead of the Neuron 6F Select Catheter, insert the 0.038in [0.97 mm] guidewire into the Neuron MAX Delivery Catheter, and advance the guidewire until the distal tip of the wire is at the distal tip of the Neuron MAX Delivery Catheter.
18. Advance the Neuron MAX Delivery Catheter and Neuron 6F Select Catheter or 0.038in [0.97 mm] guidewire simultaneously into the introducer sheath, then extend the Neuron 6F Select Catheter or 0.038in [0.97 mm] guidewire 5cm to 6cm distal to the tip of the Neuron MAX Delivery Catheter.
19. If desired, a 0.038in [0.97 mm] guidewire may be introduced through the Neuron 6F Select Catheter.
20. Advance the Neuron MAX Delivery Catheter and Neuron 6F Select Catheter or guidewire system, and then select the appropriate vessel with the Neuron 6F Select Catheter or 0.038in [0.97 mm] guidewire.
21. If using, remove the Neuron 6F Select Catheter, and replace with a 0.038in [0.97 mm] or smaller guidewire.
22. Advance the Neuron MAX Delivery Catheter and guidewire to the vascular site and remove the guidewire.
23. After the procedure, remove the Neuron MAX Delivery Catheter when clinically indicated by placing compression on the vessel above the puncture site and slowly withdraw the Neuron MAX Delivery Catheter. Discard the Neuron MAX Delivery Catheter appropriately.
24. Important: Upon removing the Neuron MAX Delivery Catheter, aspirate via the HVA to collect any fibrin that may have been deposited within or at the tip of the catheter.

DEVICE PREPARATION AND USE (GUIDE CATHETER CONFIGURATION)

1. Select an appropriately sized Neuron MAX Delivery Catheter based on the anatomy and length.
2. If using, select appropriate Neuron 6F Select Catheter shape based on the target vessel and its surrounding anatomy.
3. Gently remove the Neuron MAX Delivery Catheter and packaging card from the pouch by grasping the Neuron MAX Delivery Catheter hub and packaging card and slowly pulling them out of the pouch.
4. Remove the Neuron MAX Delivery Catheter from the packaging card by removing the hub from the card tabs before gently removing the Neuron MAX Delivery Catheter shaft.
5. Remove the Rotating Hemostasis Valve (RHV) from the card tabs.
6. Inspect the Neuron MAX Delivery Catheter for kinks or other damage.
7. Connect the RHV to the hub of the Neuron MAX Delivery Catheter. Flush the lumen with heparinized saline.
8. If using, insert the Neuron 6F Select Catheter into the Neuron MAX Delivery Catheter, and advance the Neuron 6F Select Catheter until the distal tip of the Neuron 6F Select Catheter is at the distal tip of the Neuron MAX Delivery Catheter.
9. If a 0.038in [0.97 mm] guidewire is used instead of the Neuron 6F Select Catheter, insert the 0.038in [0.97 mm] guidewire into the Neuron MAX Delivery Catheter, and advance the guidewire until the distal tip of the wire is at the distal tip of the Neuron MAX Delivery Catheter.
10. Place an introducer sheath appropriately sized for the Neuron MAX's outer diameter.
11. Advance the Neuron MAX Delivery Catheter and Neuron 6F Select Catheter or 0.038in [0.97 mm] guidewire simultaneously into the introducer sheath, and extend the Neuron 6F Select Catheter or 0.038in [0.97 mm] guidewire 5 cm to 6 cm distal to the tip of the Neuron MAX Delivery Catheter.
12. If desired, a 0.038in [0.97 mm] guidewire may be introduced through the Neuron 6F Select Catheter.
13. Advance the Neuron MAX Delivery Catheter and Neuron 6F Select Catheter or guidewire system, and then select the appropriate vessel with the Neuron 6F Select Catheter or 0.038in [0.97 mm] guidewire.
14. After gaining access to the desired vessel with the Neuron 6F Select Catheter or 0.038in [0.97 mm] guidewire, advance the Neuron MAX Delivery Catheter over the Neuron 6F Select Catheter or 0.038in [0.97 mm] guidewire into the desired vessel.
15. Under fluoroscopy, advance the Neuron MAX Delivery Catheter until the desired position is achieved.
16. If used, remove the Neuron 6F Select Catheter, and replace with a 0.038in [0.97 mm] or smaller guidewire.
17. Advance the Neuron MAX Delivery Catheter and guidewire to the vascular site and remove the guidewire.

PACKAGING

The Neuron MAX System components are affixed to a paper packaging card, whereas the Delivery Catheter and the Dilator are inside a protective HDPE tube, which is then placed in a heat-sealed PET/PE/ Tyvek® pouch to maintain sterility post EtO sterilization.

The Neuron MAX System will remain sterile unless the pouch is opened, damaged, or the "Use By" date has passed.




Materials

The Neuron MAX System and the packaging do not contain latex.

STORAGE

Store in a cool, dry place.

SYMBOLS GLOSSARY

	Attention, see Instructions For Use
	Prescription only – US Federal Law restricts this device to use by or on the order of a physician
	Nonpyrogenic
	Sterilized using ethylene oxide
	Medical Device
 Manufacturer	 EU Representative
 Lot number	 Do not reuse
 Use by	 Catalogue number
 Do not resterilize	 Not made with natural rubber latex

WARRANTY

Penumbra Inc. (Penumbra) warrants that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular use. Handling, storage, cleaning, and sterilization of this device as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond Penumbra's control directly affect the device and the results obtained from its use. Penumbra's obligation under this warranty is limited to the repair or replacement of this device and Penumbra shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this device. Penumbra neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. Penumbra assumes no liability with respect to devices reused, reprocessed, or resterilized and makes no warranties, expressed or implied, including but not limited to merchantability or fitness for intended use, with respect to such devices.

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