DEVICE DESCRIPTION

Neuron™ Intracranial Access System is a two-catheter system comprised of the Neuron™ Delivery Catheter and the Neuron™ Select™ Catheter. The Neuron Delivery Catheter can be used individually with a 0.038in [0.97 mm] guidewire or together with the Neuron Select Catheter to access the desired anatomy.

Neuron Delivery Catheter:

The Neuron Delivery Catheter is a single lumen, coil-reinforced, variable stiffness catheter with a radiopaque markerband on the distal end and a Luer hub on the proximal end. The Neuron Delivery Catheter dimensions are included on the individual device label. The Neuron Delivery Catheter is compatible with introducer sheaths having an inner diameter of 6F or greater.

Neuron Select Catheter:

The Neuron Select Catheter is a single lumen, hypo-tube and braid-reinforced, variable stiffness catheter with a radiopaque distal end and a Luer hub on the proximal end. The Neuron Select Catheter is available in four tip shapes (straight, Simmons, H1, or BER). The Neuron Select Catheter tip shape and dimensions are included on the individual device label. The Neuron Select Catheter is compatible with the Neuron Delivery Catheter.

INDICATION FOR USE

The Neuron Intracranial Access System is indicated for the introduction of interventional devices into the peripheral, coronary, and neuro vasculature.

CONTRAINDICATIONS

There are no known contraindications.

WARNINGS

- The Neuron Intracranial Access System should only be used by physicians who have received appropriate training in interventional techniques.

PRECAUTIONS

- The device is intended for single use only. Do not resterilize or reuse. Resterilization and/or Reuse may result in ineffective catheter coating lubrication, which may result in high friction and the inability to access the target vasculature location; and/or may compromise the structural integrity of the device.
- Do not use kinked or damaged devices. Do not use open or damaged packages. Return all damaged devices and packaging to the manufacturer/distributor.
- Use prior to the "Use By" date.
- Use the Neuron Intracranial Access System in conjunction with fluoroscopic visualization.
- Do not advance or withdraw the Neuron Intracranial Access System against resistance without careful assessment of the cause using fluoroscopy. If the cause cannot be determined, withdraw the device. Unrestrained moving or torquing the device against resistance may result in damage to the vessel or device.
- Maintain a constant infusion of an appropriate flush solution.
- If flow through the device becomes restricted, do not attempt to clear the lumen by infusion. Remove and replace the device.

POTENTIAL ADVERSE EVENTS

Possible complications include, but are not limited to, the following:

- acute occlusion
- air embolism
- death
- distal embolization
- emboli
- false aneurysm formation
- hematoma or hemorrhage at puncture site
- infection
- intracranial hemorrhage
- ischemia
- neurological deficits including stroke
- vessel spasm, thrombosis, dissection, or perforation

DEVICE PREPARATION AND USE

1. Select an appropriate sized Neuron Delivery Catheter based on the anatomy and length.
2. If used, select appropriate Neuron Select Catheter shape based on the target vessel and its surrounding anatomy.
3. Gently remove the Neuron Delivery Catheter and packaging card from the pouch by grasping the Neuron Delivery Catheter hub and packaging card and slowly pulling them out of the pouch.
4. Remove the Neuron Delivery Catheter from the packaging card by removing the hub from the card tabs before gently removing the Neuron Delivery Catheter shaft.
5. Remove the Rotating Hemostasis Valve (RHV) from the card tabs.
6. Inspect the Neuron Delivery Catheter for kinks or other damage.
7. Connect a rotating hemostasis valve to the hub of the Neuron Delivery Catheter. Flush the lumen with heparinized saline.
8. If used, insert the Neuron Select Catheter into the Neuron Delivery Catheter, and advance the Neuron Select Catheter until the distal tip of the Neuron Select Catheter is at the distal tip of the Neuron Delivery Catheter.
9. If a 0.038in [0.97 mm] guidewire is used instead of the Neuron Select Catheter, insert the 0.038in [0.97 mm] wire into the Neuron Delivery Catheter, and advance the guidewire until the distal tip of the wire is at the distal tip of the Neuron Delivery Catheter.
10. Place an introducer sheath with minimum 6F inner diameter in the primary access artery.
11. Advance the Neuron Delivery Catheter and Neuron Select Catheter or 0.038in [0.97 mm] guidewire simultaneously into the introducer sheath, then extend the Neuron Select Catheter or 0.038in [0.97 mm] wire 5cm to 6cm distal to the tip of the Neuron Delivery Catheter.
12. If desired, a 0.018in [0.46 mm] guidewire may be introduced through the Neuron Select Catheter.
13. Advance the Neuron Intracranial Access System, and then select the appropriate vessel with the Neuron Select Catheter or 0.038in [0.97 mm] guidewire.
14. After gaining access to the desired vessel with the Neuron Select Catheter or 0.038in [0.97 mm] guidewire, advance the Neuron Delivery Catheter over the Neuron Select Catheter or 0.038in [0.97 mm] wire into the desired vessel.
15. If used, remove the Neuron Select Catheter, and replace with an 0.038in [0.97 mm] or smaller guidewire.
16. Advance the Delivery Catheter and guidewire to the vascular site and remove the guidewire.

PACKAGING

The Neuron Intracranial Access System components are affixed to a paper packaging card, whereas the Delivery Catheter is inside a protective HDPE tube, which is then placed in a heat-sealed PET/PE/Tyvek® pouch to maintain sterility post EtO sterilization. Additionally, Neuron Delivery Catheters with shaped tips have a stainless steel support mandrel inserted into the distal tip to aid in securing the catheter to the support card.

The Neuron Intracranial Access System will remain sterile unless the pouch is opened, damaged, or the "Use By" date has passed.


MATERIALS

The Neuron Intracranial Access System and the packaging do not contain latex.

STORAGE

Store in a cool, dry place.

SYMBOLS GLOSSARY

	Attention, see Instructions For Use
	Prescription only – US Federal Law restricts this device to use by or on the order of a physician
	Nonpyrogenic
	Sterilized using ethylene oxide
	Medical Device
 Manufacturer	 EU Representative
 Lot number	 Do not reuse
 Use by	 Catalogue number
 Do not resterilize	 Not made with natural rubber latex

WARRANTY

Penumbra Inc. (Penumbra) warrants that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular use. Handling, storage, cleaning, and sterilization of this device as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond Penumbra's control directly affect the device and the results obtained from its use. Penumbra's obligation under this warranty is limited to the repair or replacement of this device and Penumbra shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this device. Penumbra neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. Penumbra assumes no liability with respect to devices reused, reprocessed, or resterilized and makes no warranties, expressed or implied, including but not limited to merchantability or fitness for intended use, with respect to such devices.

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Penumbra, Inc.