DEVICE DESCRIPTION

Benchmark™ Intracranial Access System is a two-catheter system comprised of the Benchmark™ Delivery Catheter and the Penumbra Select™ Catheter. The Benchmark Delivery Catheter can be used individually with a 0.038in [0.97 mm] guidewire or together with the Penumbra Select Catheter to access the desired anatomy.

Benchmark Delivery Catheter:

The Benchmark Delivery Catheter is a single lumen, coil-reinforced, variable stiffness catheter with a radiopaque distal end and a Luer hub on the proximal end. The Benchmark Delivery Catheter dimensions are included on the individual device label. The Benchmark Delivery Catheter is compatible with introducer sheaths having an inner diameter of 6F or greater.

Penumbra Select Catheter:

The Penumbra Select Catheter is a single lumen, braid-reinforced, variable stiffness catheter with a radiopaque distal end and a Luer hub on the proximal end. The Penumbra Select Catheter is available in tip shapes (Simmons, H1, or BER). The Penumbra Select Catheter tip shape and dimensions are included on the individual device label. The Penumbra Select Catheter is compatible with the Benchmark Delivery Catheter.

INDICATION FOR USE

The Benchmark Intracranial Access System is indicated for the introduction of interventional devices into the peripheral, coronary, and neuro vasculature.

CONTRAINDICATIONS

There are no known contraindications.

WARNINGS

- The Benchmark Intracranial Access System should only be used by physicians who have received appropriate training in interventional techniques.

PRECAUTIONS

- The device is intended for single use only. Do not resterilize or reuse. Resterilization and/or Reuse may result in ineffective catheter coating lubrication, which may result in high friction and the inability to access the target vasculature location; and/or may compromise the structural integrity of the device.
- Do not use kinked or damaged devices. Do not use open or damaged packages. Return all damaged devices and packaging to the manufacturer/distributor.
- Use prior to the “Use By” date.
- Use the Benchmark Intracranial Access System in conjunction with fluoroscopic visualization.
- Do not advance or withdraw the Benchmark Intracranial Access System against resistance without careful assessment of the cause using fluoroscopy. If the cause cannot be determined, withdraw the device. Unrestrained moving or torquing the device against resistance may result in damage to the vessel or device.
- Maintain a constant infusion of an appropriate flush solution.
- If flow through the device becomes restricted, do not attempt to clear the lumen by infusion. Remove and replace the device.

POTENTIAL ADVERSE EVENTS

Possible complications include, but are not limited to, the following:

- | | | |
|-----------------------|----------------------------|--|
| • acute occlusion | • false aneurysm formation | • ischemia |
| • air embolism | • hematoma or hemorrhage | • neurological deficits including stroke |
| • death | • at puncture site | • vessel spasm, thrombosis, dissection, or perforation |
| • distal embolization | • infection | |
| • emboli | • intracranial hemorrhage | |

DEVICE PREPARATION AND USE

1. Select an appropriate sized Benchmark Delivery Catheter based on the anatomy and length.
2. If used, select appropriate Penumbra Select Catheter shape based on the target vessel and its surrounding anatomy.
3. Gently remove the Benchmark Delivery Catheter and packaging card from the pouch by grasping the Benchmark Delivery Catheter hub and packaging card and slowly pulling them out of the pouch.
4. Remove the Benchmark Delivery Catheter from the packaging card by removing the hub from the card tabs before gently removing the Benchmark Delivery Catheter shaft.
5. Inspect the Benchmark Delivery Catheter for kinks or other damage. If any damage is observed, discard the Benchmark Delivery Catheter.
6. Connect a rotating hemostasis valve to the hub of the Benchmark Delivery Catheter. Flush the lumen with heparinized saline.
7. If used, insert the Penumbra Select Catheter into the Benchmark Delivery Catheter, and advance the Penumbra Select Catheter until the distal tip of the Penumbra Select Catheter is at the distal tip of the Benchmark Delivery Catheter.
8. If a 0.038in [0.97 mm] guidewire is used instead of the Penumbra Select Catheter, insert the 0.038in [0.97 mm] wire into the Benchmark Delivery Catheter, and advance the guidewire until the distal tip of the wire is at the distal tip of the Benchmark Delivery Catheter.
9. Place an introducer sheath with minimum 6F inner diameter in the primary access artery.
10. Advance the Benchmark Delivery Catheter and Penumbra Select Catheter or 0.038in [0.97 mm] guidewire simultaneously into the introducer sheath, then extend the Penumbra Select Catheter or 0.038in [0.97 mm] wire distal to the tip of the Benchmark Delivery Catheter.
11. Advance the Benchmark Intracranial Access System, and then select the appropriate vessel with the Penumbra Select Catheter or 0.038in [0.97 mm] guidewire.
12. After gaining access to the desired vessel with the Penumbra Select Catheter or 0.038in [0.97 mm] guidewire, advance the Benchmark Delivery Catheter over the Penumbra Select Catheter or 0.038in [0.97 mm] wire into the desired vessel.
13. If used, remove the Penumbra Select Catheter, and replace with an 0.038in [0.97 mm] or smaller guidewire.
14. Advance the Delivery Catheter and guidewire to the vascular site and remove the guidewire.

PACKAGING

The Benchmark Intracranial Access System components are affixed to a paper packaging card, whereas the Delivery Catheter is inside a protective HDPE tube, which is then placed in a heat-sealed PET/PE/Tyvek® pouch to maintain sterility post EtO sterilization. Additionally, Benchmark Delivery Catheters with shaped tips have a stainless steel support mandrel inserted into the distal tip to aid in securing the catheter to the support card.

The Benchmark Intracranial Access System will remain sterile unless the pouch is opened, damaged, or the “Use By” date has passed.




MATERIALS

The Benchmark Intracranial Access System and the packaging do not contain latex.

STORAGE

Store in a cool, dry place.

SYMBOLS GLOSSARY

	Attention, see Instructions For Use		
	Prescription only – US Federal Law restricts this device to use by or on the order of a physician		
	Nonpyrogenic		
	Sterilized using ethylene oxide		
	Medical Device		
	Manufacturer		EU Representative
	Lot number		Do not reuse
	Use by		Catalogue number
	Do not resterilize		Not made with natural rubber latex

WARRANTY

Penumbra Inc. (Penumbra) warrants that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular use. Handling, storage, cleaning, and sterilization of this device as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond Penumbra's control directly affect the device and the results obtained from its use. Penumbra's obligation under this warranty is limited to the repair or replacement of this device and Penumbra shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this device. Penumbra neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. Penumbra assumes no liability with respect to devices reused, reprocessed, or resterilized and makes no warranties, expressed or implied, including but not limited to merchantability or fitness for intended use, with respect to such devices.

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Penumbra, Inc.