PHYSICA SYSTEM

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INSTRUCTIONS FOR USE – PHYSICA SYSTEM

Before using a LimaCorporate product, the surgeon shall carefully study the following recommendations, warnings and instructions, as well as the most updated available product-specific information (e.g.: product literature, surgical technique).

1. PRODUCT INFORMATION

The main goal of a knee joint prosthesis is to reproduce the articular anatomy. The joint prosthesis is intended to reduce pain and give articular mobility to the patient. The degree of pain reduction and mobility depends in part on the pre-operative situation, intra-operative options and on the post-op rehabilitation. Physica system is a modular knee replacement system. The CR (cruciate retaining) and KR (kinematic retaining) designs are intended to be used in patients with a preserved and well-functioning posterior cruciate ligament. The PS (posterior stabilized) is intended to be used in situations where the posterior cruciate ligament is absent or cannot be preserved. All three versions (CR, KR and PS) require preserved and well-functioning collateral ligaments. Physica system includes femoral components, tibial implants with tibial stems, tibial liners and patellar components. These devices are intended for cemented use only. The components shall not be used alone as single components, or with components from other manufacturer.

1.1. MATERIALS

Components	Material	Material standard
Femoral component	CoCrMo	ISO 5832-4 – ASTM F75
Pegs for PS femoral component	CoCrMo	ISO 5832-12 – ASTM F1537
Tibial plate	Ti6Al4V	ISO 5832-3 – ASTM F1472
Tibial liner	UHMWPE	ISO 5834-2 – ASTM F648
Tibial stem	Ti6Al4V	ISO 5832-3 – ASTM F1472
Patella	UHMWPE	ISO 5834-2 – ASTM F648
Tibial plate plug	UHMWPE	ISO 5834-2 – ASTM F648

All materials are purchased according to International Standards and are widely used in the orthopaedic field. CoCrMo, used to manufacture some components of Physica System, includes cobalt (which is a CMR class 1B) in percentage above 0.1% (w/w). Nevertheless, chemical tests show that the release of cobalt in the human body is below 0.1% (w/w). Additionally, biocompatibility tests and long-term clinical experience of the use of this material in orthopaedic applications has shown that an acceptable level of biological response can be expected. Therefore, the risks related to the use of CoCrMo and other materials have been considered and the related risk/benefit ratio has been evaluated as acceptable by the manufacturer. Some patients may be vulnerable (e.g. have hypersensitive or allergic reactions) to the implant materials/substances, this shall be appropriately considered by the Surgeon. Follow the instructions of the manufacturer of bone cement (PMMA) for preparing the bone cement applying bone cement for prosthesis fixation.

1.1.1. MEDICINAL SUBSTANCES AND TISSUE

The device does not contain or incorporates any: • medicinal substance, including a human blood or plasma derivative; • tissues or cells, or their derivatives, of human origin; • tissues or cells of animal origin, or their derivatives, as referred to in Regulation (EU) No.722/2012.

1.2. HANDLING, STORAGE AND SAFE DISPOSAL

All devices are provided sterile and shall be stored at ambient temperature (indicative range 0-50°C / 32-122°F) in their protective closed packaging in controlled rooms, protected from exposure to light, heat and sudden changes in temperature.

Once the package is opened check that both the model and size of the implant correspond exactly to the description printed on the labels. Avoid any contact between the implant and objects or substances that can alter the sterile condition or surface integrity. A careful visual examination of each implant is recommended before use in order to verify that the implant is not damaged. **Components removed from the package shall not be used if they are dropped or suffer other accidental impacts. Devices shall not be modified in any way.**

The device's code and lot number shall be recorded in the patient's case history by using the labels included in the components packaging.

Re-use of previously implanted devices must be absolutely avoided. Risks associated with the reuse of single use devices are: • infection; • early or late failure of the device or device fixation; • lack of appropriate coupling between modular junctions (e.g. taper connections); • device wear and wear debris associated complications; • transmission of diseases (e.g. HIV, hepatitis); • immune system response / rejection.

DISPOSAL: the disposal of medical devices is to be performed by the hospitals in conformity with applicable laws.

1.3. DEVICE PERFORMANCE, EXPECTED CLINICAL BENEFITS, EXPECTED LIFETIME

The main goal of a joint prosthesis is to reproduce the articular anatomy. The joint prosthesis is intended to reduce pain and give articular mobility to the patient. The degree of pain reduction and mobility depend, in part,

on the pre-operative situation, intra-operative options and the post-op rehabilitation.
LimaCorporate annually updates the Summary of Safety and Clinical Performance (SSCP) to be uploaded on EUDAMED (<https://ec.europa.eu/tools/eudamed>).

The expected lifetime of Physica CR and Physica KR is at least 95%* at 15 years.

The expected lifetime of Physica PS is at least 94%* at 15 years.

*** with a variability of $\pm 2\%$**

NOTE: The expected lifetime of the medical devices, when placed on the market and – moreover – when put into use (namely when implanted in a patient's body) may depend on a variety of factors:

- (a) intrinsic characteristic (feature) of the medical device in itself (for example: design of the medical device; materials employed in the medical device manufacturing; manufacturing process of the medical device in itself; surface texture, osteo-conductive and osteo-inductive properties of the medical device);
- (b) external factors related to surgery process (for example: right choice of the medical device, stability of the initial fixation of the medical device; accuracy in positioning of the medical device);
- (c) external factors related to patient (for example: aetiology, pathology, physical characteristics, age, previous surgical history, daily activities performed);
- (d) risk factors mentioned in the IFU leaflet, when not already included in the factors referred to in previous paragraphs from (a) to (c);
- (e) factors related to all the possible complications which can affect the duration of the implant.

The expected lifetime was calculated by means of an analytical and rationalized estimation, taking into account:

- i. the clinical data collected on the medical device;
- ii. the survivorship data reported in arthroplasty registries referring to similar medical devices (identified state of the art benchmark);
- iii. the expected (as predictable) behavior of the medical device at higher follow-up.

Due to the analytical evaluations, the expected lifetime value is approximated to the unit.

The value shall be interpreted duly considering the specified value variability, taking into necessary consideration all relevant factors mentioned above, in paragraphs (a) to (e).

2. INTENDED PURPOSE, INDICATIONS, CONTRAINDICATIONS AND RISK FACTORS

Physica System is intended for the cemented knee joint replacement and for long-term implantation inside the human body.

2.1. INDICATIONS

Physica system is indicated for use in knee arthroplasty in skeletally mature patients with the following conditions: • Non-inflammatory degenerative joint disease: including osteoarthritis, traumatic arthritis or avascular necrosis; • Inflammatory degenerative joint disease including rheumatoid arthritis.

In presence of one of the previous conditions affecting the joint cartilages (e.g., osteoarthritis, rheumatoid arthritis, etc.), additional indications for Physica PS components, in patients with preserved and well-functioning collateral ligaments, are: • Absent or not-functioning posterior cruciate ligament; • Severe antero-posterior instability of the knee joints.

2.2. CONTRAINDICATIONS

Common contraindications for Physica CR, KR and PS versions include: • Severe instability of the knee joint secondary to the absence of collateral ligament integrity and/or function; • Local or systemic infection; • Important bone loss on femoral or tibial joint side; • Progressive tumour diseases; • Known incompatibility or allergy to the product materials; • Septicaemia; • Persistent acute or chronic osteomyelitis; • Open epiphyses (immature patient with active bone growth).

Specific contraindications for Physica CR and KR versions include: • Important joint instability; • Deficiency of posterior cruciate ligament.

The relative contraindications are: • Vascular or nerve diseases affecting the concerned limb; • Bone stock compromised by disease, infection or prior implantation that cannot provide adequate support and/or fixation to the prosthesis; • Metabolic disorders which may impair fixation and stability of the implant; • Any concomitant disease and dependence that might affect the implanted prosthesis; • Metal hypersensitivity to implant materials; • Important osteoporosis, haemophilic disease; • Internistic problems with high risk for surgery; • Skeletal immaturity.

2.3. RISK FACTORS

The following risk factors may result in poor results with this prosthesis: • overweight¹; • strenuous physical activities (i.e. active sports, heavy physical work); • fretting of modular junctions; • incorrect implant positioning; • insufficient bone to support the femoral and/or tibial components; • medical disabilities which can lead to an unnatural gait and loading of the knee joint; • muscle deficiencies; • multiple joint disabilities; • refusal to modify postoperative physical activities; • patient's history of infections or falls; • systemic diseases and metabolic disorders; • local or disseminated neoplastic diseases; • drug therapies that adversely affect bone quality, healing, or resistance to infection; • drug use or alcoholism; • marked osteoporosis or osteomalacia; • patient's resistance to disease generally weakened (HIV, tumour, infections); • severe deformity leading to an impaired anchorage or improper positioning of implants.

3. WARNINGS

3.1. PRE-OPERATIVE PLANNING

LimaCorporate products shall be implanted only by surgeons familiar with the joint replacement procedures described in the specific surgical techniques.

¹ According to the definition of World Health Organization (WHO), Body Mass Index (BMI) greater than or equal to 25 Kg/m²

COMBINATIONS ALLOWED/NOT ALLOWED

- The permitted combinations of femur and tibial liner sizes for Physica System are shown in this table:

		Tibial Liner									
		#1	#2	#3	#4	#5	#6	#7	#8	#9	#10
Femoral Component	#1										
	#2										
	#3										
	#4										
	#5										
	#6										
	#7										
	#8										
	#9										
	#10										

- Physica tibial liners match exactly the tibial plate sizes: they are mated one-to-one (e.g. size 1 tibial liner is only used with size 1 tibial plate);
- Physica CR femoral component can be coupled with CR and LMC (UHMWPE and LimaVit) tibial liner;
- Physica KR femoral component can be coupled with KR (UHMWPE and LimaVit) tibial liner;
- Physica PS femoral component can be coupled with PS tibial liner.

For more information on the allowed combinations of Physica System with Physica TT tibial plate, LMC (UHMWPE and LimaVit) and KR (LimaVit) tibial liner, Patellar prostheses (LimaVit) and AMF TT cones please refer to the specific IFU.

The components of the system must not be used with components of other manufacturer, unless specified by LimaCorporate. The components forming original LimaCorporate systems must be rigorously coupled following the given information included in this document and in the available product-specific information (e.g.: current version of surgical technique). The use of single components, or components belonging to other systems, are subject to the approval by LimaCorporate. The manufacturer and the dealer are not liable for possible coupling incompatibility. The surgeon is solely responsible for the choice and use of the implant.

Pre-operative planning, through radiographic templates in different formats, provides essential information regarding the type and size of components to be used and the correct combination of required devices based on the anatomy and specific conditions of each patient. Inadequate pre-operative planning can lead to improper selection of the implants and/or incorrect implant positioning. **Complications or failures of the total knee replacement are more likely to occur in heavy and highly active patients.** The surgeon shall perform a careful evaluation of the patient's clinical condition and level of physical activity before performing knee replacement. An inadequate clinical condition and level of physical activities may adversely affect the outcome of the joint replacement.

LimaCorporate specialized technical staff is available to provide advice regarding pre-operative planning, the surgical technique, and product and instrumentation assistance both prior to and during surgery.

The patient shall be warned that the prosthesis does not replace normal healthy bone, that the prosthesis can break or become damaged as a result of certain activity or trauma, that it has a finite expected implant life, and may need to be replaced at some time in the future: the longevity and the performance of the implant over time, in fact, can be influenced by the natural/physiological progression of the patient pathology, the insurgence of co-morbidities and post-operative complications, even if they are not directly implant-related (i.e. infection pain, stiffness, reduced range of motion).

The possible impact of the factors mentioned in sections "GENERAL INFORMATION ON INDICATIONS, CONTRAINDICATIONS AND RISK FACTORS" and "POSSIBLE ADVERSE EFFECTS" shall be considered preoperatively, and the patient informed as to what steps he/she can take to reduce the possible effects of these factors.

Implants are single use devices.

Do not re-use implants that were previously implanted in another patient.

Do not re-use an implant that has previously come into contact with the body fluid or tissue of another person.

Do not use any component from a package that has been previously opened or appears to be damaged.

Do not use implants after the expiration date printed on the label.

Do not use this product for other than labelled indications (off-label use); the risk of implant failure is higher with inaccurate component alignment or positioning. Use only instruments and provisionals specifically designed for use with these devices to help ensure accurate surgical implantation, soft tissue balancing, and evaluation of knee function.

Surgical instruments used to perform the surgery are subject to wear with normal usage. After extensive use or excessive loads, instruments are susceptible to fracture. Surgical instruments shall be used only for their specific purpose. Before use, the functionality of surgical instruments shall be checked since the use of damaged instruments may lead to early failure of the implants. Damaged instruments shall be replaced before surgery.

3.2. INTRAOPERATIVE

The use of trial devices is recommended to check the correct site preparation, size and positioning of the implants to be used. It is recommended that additional implants be available during surgery for use in those cases requiring prostheses of different sizes or when the preoperatively selected prostheses cannot be used. The correct selection as well as the correct seating/placement of the implant is extremely important. Improper selection, positioning, alignment and fixation of the implant components may result in unusual stress conditions which may negatively affect system performance and the survival rate of the implant.

The components forming original LimaCorporate systems must be assembled according to the surgical technique and used only for the labelled indications.

Use only instruments and trial components specifically designed for use with the implants being used. The use of instruments from other manufacturers or the use of instruments designed for use with other systems can lead to inappropriate preparation of the implant site, incorrect positioning, alignment and fixation of the devices followed by loosening of the system, loss of functionality, reduction of the durability of the implant, and the need for further surgery.

Care must be taken to protect the surfaces involved in the coupling between components (taper coupling between the tibial plate and the tibial stem and coupling between the tibial plate and the tibial liner); the articular surfaces of the implants shall be protected from scratches or any other damage. All component coupling surfaces shall be clean and dry before assembly. The stability of component couplings shall be verified as described in the surgical technique.

3.3. POST-OPERATIVE CARE

Adequate post-operative care shall be provided by the surgeon or other suitably qualified medical staff. Regular postoperative x-ray follow-up is recommended to detect any changes in position or condition of the implant or surrounding tissues. Complications and/or failure of prosthetic implants are more likely to occur in patients with unrealistic functional expectations, heavy patients, physically active patients, and/or with patients who fail to follow through with the required rehabilitation program. Excessive physical activity or trauma to the replaced joint can lead to premature failure of the arthroplasty through loosening, fracture, or abnormal wear of the prosthetic implants. The patient shall be cautioned by the surgeon to govern activities accordingly and advised that the implants may fail due to excessive joint wear.

The surgeon shall make the patient aware of the limitation of limb function after knee arthroplasty and that caution is needed, especially in the post-operative period. In particular the following precautions shall be presented to the patient by the surgeon: • avoid repeated high weight lifting; • keep body weight under control (overweight conditions may adversely affect the outcomes of joint replacement); • avoid sudden peak loads (consequences of activities such as running and skiing) or movements which can lead to sudden stops or twisting; • avoid positions that can increase the risk of dislocation.

The patient shall be informed also about the possibility that the implant or its components may wear out, fail or need to be replaced. The implant may not last the rest of the patient's life, or any particular length of time. Because prosthetic implants are not as strong, reliable, or durable as natural, healthy tissue/bones, all such

devices may need to be replaced at some point.

Appropriate post-operative rehabilitation instructions and care must be provided to avoid a negative influence on the outcome of the surgical procedure.

3.4. POSSIBLE SIDE EFFECTS

Adverse effects that most commonly and frequently occur in orthopaedic prosthesis procedure include:

- loosening of the prosthetic components;
- prosthesis dislocation and instability;
- damage to the prosthetic implant;
- instability of the system because of inadequate soft tissue balancing;
- ligaments instability;
- dissociation or malalignment due to incorrect coupling of the devices;
- dislocation, luxation or subluxation;
- infection;
- tissue reaction to the implant material or abrasion;
- local hypersensitivity;
- bone fracture or lysis;
- patella erosion;
- local pain;
- joint stiffness or arthrofibrosis;
- periprosthetic fractures;
- temporary or permanent nerve damage;
- flexion contracture;
- limited range of motion or incorrect rotation;
- lengthening or shortening of the leg;
- fracture of the devices;
- excessive wear of UHMWPE components due to damaged articular surfaces or presence of particles;
- additional surgery.

Some adverse effects can lead to death. General complications include venous thrombosis with/without pulmonary embolism, cardiovascular or pulmonary disturbances, haematomas, blood loss, post-operative anaemia, systemic allergic reactions, systemic pain, temporary or permanent paralysis.

4. STERILITY

a. Implants

All implantable components of the system are provided sterile with a Sterility Assurance Level (SAL) of 10^{-6} . Metal components are sterilized by radiation or EtO and UHMWPE components by EtO. Do not use any component from a package that has been previously opened or appears to be damaged. **Do not use implants after the expiration date printed on the label.**

b. Instruments

Instruments are supplied non-sterile and must be cleaned, disinfected and sterilized before use according to appropriate validated methods (refer to the "Instrument Care, Cleaning, Disinfection and Sterilization" brochure for validated sterilization parameters; this brochure is available on request or downloadable from www.limacorporate.com in the Products section). Users shall validate their specific cleaning, disinfection and sterilization processes and equipment.

5. INFORMATION USEFUL TO INFORM THE PATIENT

Refer to previous sections for any information that allows the patient to be informed of any warnings, precautions, contra-indications, measures to be taken and limitations of use regarding the device.

The surgeon shall notify the patient about any necessary follow-up specific for the patient's condition.

5.1. MAGNETIC RESONANCE IMAGING (MRI)

Non-clinical testing has demonstrated that LimaCorporate Physica System is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 Tesla (1.5 T) and 3 Tesla (3 T)
- Maximum spatial field gradient of 3,880 G/cm (38.80 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of:
 - 2.0 W/kg (Normal Level Controlled Operating Mode) at 1.5 T
 - 1.6 W/kg at 3T

RF Heating

Under the scan conditions defined above, LimaCorporate Physica System is expected to produce a maximum temperature rise of less than 6.0°C after 15 minutes of continuous scanning

Caution: The RF heating behaviour does not scale with static field strength. Devices that do not exhibit detectable heating at one field strength may exhibit high values of localized heating at another field strength.

MR Artifact

In non-clinical testing, the image artifact caused by the device extends approximately:

- 5.5 cm from the Physica System when imaged with a gradient echo pulse sequence in a 1.5 T MR system;
- 7.9 cm from the Physica System when imaged with a spin echo pulse sequence in a 3 T MR system.

6. INCIDENT REPORTING

According to the definition of incident/serious incident given by Regulation EU 2017/745², users and/or patients are required to: **• report to the manufacturer and to the relevant Competent Authority any serious incident that has occurred in relation to the device; • report to the manufacturer any incident and/or event of malfunction of the device or changes in its performance that may affect safety.**

In non-UE market, refer to incidents' definitions according to applicable laws.

² an incident is defined to be any malfunction or deterioration in the characteristics and/or performance of the device, including use-error due to ergonomic features, as well as any inadequacy in the information supplied by the manufacturer and any undesirable side-effect. An incident is considered to be serious when it directly or indirectly led, might have led or might lead to any of the following: (a) the death of a patient, user or other person, (b) the temporary or permanent serious deterioration of a patient's, user's or other person's state of health, (c) a serious public health threat.

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Limacorporate S.p.A.