
# neopump<sup>®</sup>-R

ELASTOMERIC  
INFUSION SYSTEM

## Instructions for use Mode d'emploi

**Instructions for use****Description**

Neopump® - R is a disposable sterile device, for the continuous infusion therapy of pharmaceutical products; no electricity/battery required.

Neopump® - R can be used only by health care workers for hospital or domestic treatments.

**Indications**

Neopump® - R is used for the continuous infusion via intravenous, intra-arterial, epidural or subcutaneous routes according to the specific clinical protocol.

It may be used for the following treatments: chemotherapy, acute and chronic pain therapy, chelation therapy, antibiotic and antiviral drug therapy.

**Precautions**

- Carefully read this descriptive folder before using.
- Determine the integrity of the product packaging. If it has been opened or damaged, do not use.
- Single use only: do not re-sterilize or reuse. After use the device is not sterile and the correct performance is no longer guaranteed.
- When preparing and administering the drug, follow the instructions supplied by the producer and the prescribing physician.
- Verify the volume and the nominal flow which should be the same as indicated in the protocol. Take care in cases of the administration of critical drugs of which interruption or infusion unsteadiness could cause serious patient damage.
- When preparing and administering antineoplastic chemotherapy drugs it is necessary to carefully observe the established prevention systems including individual personal protection devices in accordance with normal protocols.
- Do not use during transportation due to possible elastomeric damage.
- Elastomeric devices do not guarantee the infusion as precisely as electronic systems. Infusion can last more or less than the indicated duration.
- An accurate flow rate is guaranteed when the housing is kept at the same height as the access point.
- Nominal flow rate is obtained at specific temperatures: flow restrictor = 32 °C / 89 °F - pump body = 24°C / 76° F
- Nominal flow refers to a saline solution (NaCl 0.9%). The flow rate may increase or decrease with decreasing or increasing viscosity of the solution to be injected.
- Health care workers must inform patients about the correct instructions for use and potential dangers.

**DFO models**

Dedicated setting for Deferoxamine (DFO) infusion. Nominal flow refers to a 3.0 grams +/- 10% DFO solution. The flow rate may increase or decrease with decreasing or increasing viscosity of the solution to be injected.

**Potential dangers**

In case of accidental elastomeric damage take the right precautions, keep the housing in a correct position to avoid drug dispersion, and substitute the device.

In case of incorrect flow rate, verifiable through the graduated scale, seek health care assistance.

**Instruction for use****Preparation**

- Before opening the package, check that the flow rate shown on the package is the same as that required for use
- Connect the syringe to the needleless filling port and inject the medication into the membrane
- Wait until the tube is completely filled and any air has been expelled.

**Patient connection**

- Remove the distal cap
- Verify that the first drop is visible
- Connect the distal connector to the patient

**Waste**

The device must be disposed of in accordance with clinical waste law.

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