ΔVΔNOS*

MIC*
FEEDING PORT
WITH ENFit® CONNECTOR

For use with
14 Fr Tube

Instructions for Use

<https://www.eifu.online/AVA>

AVANOS®

MIC® Feeding Port with ENFit® Connector • For use with 14 Fr

Instructions for Use

Rx Only: Federal Law (USA) restricts this device to sale by or on the order of a physician.

Non sterile contents • One time use only

Contents:

- Instructions
- ENFit® Connector with tethered cap (25 eaches)

⚠ Warning: Do not reuse or reprocess this medical device. Reuse or reprocessing may 1) adversely affect the known biocompatibility characteristics of the device, 2) compromise the structural integrity of the device, 3) lead to the device not performing as intended, or 4) create a risk of contamination and cause the transmission of infectious diseases resulting in patient injury, illness, or death.

Intended Use


The MIC® Feeding Port with ENFit® Connector is a single use medical device intended to be connected to 14 Fr MIC® Tubes. This device is intended to be used for any patient with a 14 Fr MIC® Tube and to be used by trained nurses and physicians in a hospital or alternative care setting.

Read all instructions before beginning procedure.

1. Once the tube has been placed, squarely cut off the proximal end of the tube close to the adapter **①**.
2. Discard the male adapter feeding port head per facility protocol.
3. Press fit the barbed ENFit® Connector with tethered cap into the feeding tube until fully seated **②**.
4. If the tube is torn or ruptured, the ENFit® Feeding Connector breaks, or the connector and/or tube is obstructed, repeat steps 1–3 with a new ENFit® Feeding Connector.

Note: The MIC® Feeding Port with ENFit® Connector is designed for only that portion of the tube that protrudes sufficiently past the skin line. Placement is not possible close to the skin line.

⚠ Warning: For enteral nutrition and/or medications only.

Tube Patency Guidelines

Proper tube flushing is the best way to avoid clogging and maintain tube patency. The following are guidelines to avoid clogging and maintain tube patency.

- Flush the feeding tube with water every 4–6 hours during continuous feeding, anytime the feeding is interrupted, or at least every 8 hours if the tube is not being used.
- Flush the feeding tube after checking gastric residuals.
- Flush the feeding tube before and after medication administration and between medications. This will prevent the medication from interacting with formula and potentially causing the tube to clog.
- Use liquid medication when possible and consult the pharmacist to determine if it is safe to crush solid medication and to mix with water. If safe, pulverize the solid medication into a fine powder form

and dissolve the powder in warm water before administering through the feeding tube. Never crush enteric-coated medication or mix medication with formula.

- Avoid using acidic fluids such as cranberry juice and cola beverages to flush feeding tubes as the acidic quality when combined with formula proteins may actually contribute to tube clogging.

General Flushing Guidelines

- Use a 30 to 60 ml ENFit® syringe. Do not use smaller size syringes as this can increase pressure on the tube and potentially rupture smaller tubes.
- Use room temperature water for tube flushing. Sterile water may be appropriate where the quality of municipal water supplies is of concern. The amount of water will depend on the patient's needs, clinical condition, and type of tube, but the average volume ranges from 10 to 50 ml for adults, and 3 to 10 ml for infants. Hydration status also influences the volume used for flushing feeding tubes. In many cases, increasing the flushing volume can avoid the need for supplemental intravenous fluid. However, individuals with renal failure and other fluid restrictions should receive the minimum flushing volume necessary to maintain patency.
- Do not use excessive force to flush the tube. Excessive force can perforate the tube and can cause injury to the gastrointestinal tract.
- Document the time and amount of water used in the patient's record. This will enable all caregivers to monitor the patient's needs more accurately.

Tube Occlusion

Tube occlusion is generally caused by:

- Poor flushing techniques
- Failure to flush after measurement of gastric residuals
- Inappropriate administration of medication
- Pill fragments
- Thick formulas, such as concentrated or enriched formulas that are generally thicker
- Formula contamination that leads to coagulation
- Reflux of gastric or intestinal contents up the tube

To Unclog a Tube

1. Make sure that the feeding tube is not kinked or clamped off.
2. If the clog is visible above the skin surface, gently massage or milk the tube between fingers to break up the clog.
3. Connect a 30 to 60 ml ENFit® syringe filled with warm water into the appropriate access port of the tube and gently pull back on then depress the plunger to dislodge the clog. Do not use smaller size syringes as this can increase pressure on the tube and potentially rupture smaller tubes.
4. If the clog remains, repeat step #3. Gentle suction alternating with syringe pressure will relieve most obstructions.
5. If this fails, consult with the physician. Do not use cranberry juice, cola drinks, meat tenderizer or chymotrypsin, as they can actually cause clogs or create adverse reactions in some patients.

⚠ Caution: Do not insert foreign objects through the tube.

 Diameter	 MD Medical Device			Rx Only				
 MR MR Safe				 UDI Unique device identifier				

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*Registered Trademark or Trademark of Avanos Medical, Inc., or its affiliates.
© 2018 AVNS. All rights reserved. 2020-10-26 15-M1-980-01 / 50000321
Distributed in the USA by Avanos Medical Sales, LLC,
Avanos Medical, Inc.,
Australia Pty Ltd, 475 Victoria Avenue,