
AVANOS*

GASTROINTESTINAL ANCHOR SET

WITH SAF-T-PEXY* T-FASTENERS

Instructions for Use


Fig 1


Fig 2


Fig 3


Fig 4


Fig 5


Fig 6


Fig 7


Fig 8

AVANOS* Gastrointestinal Anchor Set with SAF-T-PEXY* T-Fasteners • Instructions for Use

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 Single Use Only	 Sterilized Using Ethylene Oxide	 Do not use if package is damaged	 Do not resterilize
Not made with natural rubber latex	Rx Only	 Caution	 Consult instructions for use

Instructions for Use

Rx Only: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Kit Contents:

- 1 Gastrointestinal Anchor Set with SAF-T-PEXY® T-Fasteners
- 2 Hemostat

Intended Use:

The AVANOS® Gastrointestinal Anchor Set with SAF-T-PEXY® T-Fasteners is intended to affix the stomach to the anterior abdominal wall facilitating primary placement of the AVANOS® MIC® and MIC-KEY® brand Enteral Feeding Tubes. It is recommended that these T-Fasteners be used only with the AVANOS® MIC® and MIC-KEY® brand Enteral Feeding Tubes.

Contraindications:

Contraindications include, but are not limited to ascites, colonic interposition, portal hypertension, gastric varices, peritonitis, aspiration pneumonia and morbid obesity (stoma lengths longer than 10 cm).

⚠ Warning: Verify package integrity prior to opening. Do not use if package is damaged or sterile barrier compromised.

⚠ Warning: Do not reuse, reprocess, or resterilize this medical device. Reuse, reprocessing, or resterilization may 1) adversely affect the known biocompatibility characteristics of the device, 2) compromise the structural integrity of the device, 3) lead to the device not performing as intended, or 4) create a risk of contamination and cause the transmission of infectious diseases resulting in patient injury, illness, or death.

Suggested Placement Procedure:

⚠ Caution: Consult Glucagon Instructions for Use for rate of IV injection and recommendations for use with insulin dependent patients.

Note: PO/NG contrast may be administered the night prior or an enema administered just prior to placement to opacify the transverse colon.

1. Place the patient in the supine position.
2. Prep and sedate the patient according to clinical protocol.
3. Ensure that the left lobe of the liver is not over the fundus or the body of the stomach.
4. Identification of the medial edge of the liver by CT scan or ultrasound may be helpful.
5. Glucagon 0.5 to 1.0 mg IV may be administered to diminish gastric peristalsis.
6. If performing endoscopic or radiologic procedures, insufflate the stomach with air using a nasogastric catheter, usually 500 to 1,000 ml or until adequate distention is achieved. It is often necessary to continue air insufflation during the procedure, especially at the time of needle puncture and tract dilation, to keep the stomach distended in order to appose the gastric wall against the anterior abdominal wall. **(Fig 1)**
7. Choose a catheter insertion site in the left sub-costal region, preferably over the lateral aspect or lateral to the rectus abdominis muscle **(Note: the superior epigastric artery courses along the medial aspect of the rectus)** and directly over the body of the stomach toward the greater curvature. Using fluoroscopy, endoscopy, or laparoscopy choose a location that allows as direct a needle path as possible. Obtain a cross table lateral view prior to placement of gastrostomy when interposed colon or small bowel anterior to the stomach is suspected.
8. If performing a laparoscopic procedure, place a skin mark at the tube insertion site and define the gastropexy pattern (before peritoneal insufflation) by placing three skin marks equidistant (approximately 2 cm apart) from the tube insertion site and in a triangle configuration.

(Fig 2) Allow adequate distance between the insertion site and SAF-T-PEXY® placement so as to prevent interference of the anchor set and balloon once inflated. Insufflate the stomach sufficiently to improve gastric access and allow the SAF-T-PEXY® needle insertion and T-fastener deployment. (Skin marks for endoscopic or radiologic procedures are performed during step 10.)

9. Prep and drape according to facility protocol.

Placing the SAF-T-PEXY®:

Note: The SAF-T-PEXY® device contains 3/0 Biosyn® synthetic absorbable suture that in non-clinical studies retained tensile strength to approximately 75% of U.S.P. and E.P. minimum knot strength at 14 days and approximately 40% at 21 days post implantation. Absorption of the suture is essentially complete within 90 to 110 days.

⚠ Warnings:

- The kinetics of gastric wall adhesion to the anterior abdominal wall relative to suture absorption must be considered prior to using the SAF-T-PEXY® device when a compromised healing response is anticipated, especially when adhesion of the gastric wall to the anterior abdominal wall is not expected within 14 days.
- T-Fasteners may migrate and be retained in the gastric mucosa, abdominal musculature, or subcutaneous tissues and in rare circumstances have exited through the skin adjacent to the stoma site.
- If performing an early replacement gastrostomy within the first few weeks after initial SAF-T-PEXY® placement, verify that the adhesion of the gastric wall to the interior abdominal wall is established and maintained. Consider use of an additional T-Fastener to ensure proper early replacement of gastrostomy.

⚠ Warning: The suture locks may pose a choking hazard.

⚠ Caution: The SAF-T-PEXY® needle point is sharp.

⚠ Caution: Maintaining proper stomach insufflation is critical to preventing T-Fastener placement through the posterior gastric wall.

Note: It is recommended to perform a three-point gastropexy that approximates an equilateral triangle to help ensure secure and uniform attachment of the gastric wall to the anterior abdominal wall. An alternate pattern may need to be identified if placing a low volume balloon gastrostomy tube.

1. If performing endoscopic or radiologic procedures, now place a skin mark at the tube insertion site and define the gastropexy pattern by placing three skin marks equidistant (approximately 2 cm apart) from the tube insertion site and in a triangle configuration. Allow adequate distance between the insertion site and SAF-T-PEXY® placement in order to prevent interference of the anchor set and balloon once inflated. **(Fig 2)**
2. If performing a laparoscopic procedure, reconfirm gastropexy placement sites.
3. Administer local anesthetic to the planned puncture sites as per hospital protocol.
4. Carefully remove the preloaded SAF-T-PEXY® device from the protective sheath and maintain slight tension on the trailing suture, noting that the suture is held to the needle by a retaining snap on the side of the needle hub.
5. Attach a Luer slip syringe containing 1–2 ml of sterile water or saline to the needle hub. **(Fig 3)**
6. Under fluoroscopic, endoscopic or laparoscopic guidance, insert the preloaded SAF-T-PEXY® slotted needle with a single sharp thrust through one of the marked corners of the triangle until it is within the gastric lumen. **(Fig 4)** The simultaneous return of air into the syringe and visualization (during fluoroscopic visualization contrast

may be injected upon air return to visualize gastric folds and confirm intraluminal position) confirms correct intragastric position. After confirmation of correct position, remove the syringe from the device.

7. Release the suture strand. Bend the locking tab on the needle hub. **(Fig 5)** Firmly push the inner hub into the outer hub until the locking mechanism clicks into place. **(Fig 6)** This will dislodge the T-Bar from the end of the needle and lock the inner stylet into position. **(Fig 7)**
8. Withdraw the needle while continuing to gently pull the T-Bar until it is flush against the gastric mucosa, avoid having the T-Bar exert excessive tension onto the gastric mucosa. Discard the needle according to facility protocol.
9. Gently slide the suture lock down to the skin surface of the abdominal wall. A small hemostat may be clamped above the suture lock to temporarily hold it in place.
⚠ Caution: Do NOT place hemostat below the suture lock, or between the suture lock and the abdominal wall, as damage to the suture may occur.
10. Repeat the procedure until all three anchor sets have been inserted in the corners of the triangle. After the three SAF-T-PEXY* devices are properly positioned, gently pull on the sutures to appose the stomach to the anterior abdominal wall avoiding exerting excessive tension on the sutures. Close the suture lock with the supplied hemostat until an audible “click” is heard securing the suture. Any excess suture length may be cut and removed. **(Fig 8)**

Note: For additional suture security, a knot may be tied in the suture strand at the surface of the suture lock.

Post Procedure:

1. Inspect the stoma and gastropexy sites daily and assess for signs of infection, including: redness, irritation, edema, swelling, tenderness, warmth, rashes, purulent or gastrointestinal drainage. Assess for any signs of pain, pressure or discomfort.
2. After the assessment, routine care should include cleansing the skin around the stoma site and gastropexy sites with warm water and mild soap, using a circular motion, moving from the tube and external bolsters outward, followed by a thorough rinsing and drying well. The sutures may be left to be absorbed or they may be cut when deemed appropriate by the placing physician. Letting the sutures be absorbed or cutting them will allow the T-bars to pass through the gastrointestinal system. After the sutures dissolve (or are cut) the suture locks may be removed and discarded.

Note: It is recommended not to cut the sutures within two weeks post procedure.

MRI labeling based on the test results : MRI Information

MR Conditional

The T-Bar Assembly (T-Fastener) was determined to be MR-conditional according to the terminology specified in the American Society for Testing and Materials (ASTM) International, Designation: F2503-08, Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment. ASTM International, 100 Barr Harbor Drive, PO Box C700, West Conshohocken, Pennsylvania.

Non-clinical testing demonstrated that the T-Bar Assembly (T-Fastener) is MR Conditional. A patient with this device can be scanned safely under the following conditions:

Static Magnetic Field

-Static magnetic field of 3-Tesla or less

-Maximum spatial gradient magnetic field of 720-Gauss/cm or less

MRI-Related Heating

In non-clinical testing, the T-Bar Assembly (T-Fastener) produced the following temperature rise during MRI performed for 15-min of scanning (i.e., per pulse sequence) in the 3-Tesla (3- Tesla/128-MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR system: Highest temperature change +1.6 °C

Therefore, the MRI-related heating experiments for the T-Bar Assembly (T-Fastener) at 3-Tesla using a transmit/receive RF body coil at an MR system reported whole body averaged SAR of 2.9 -W/kg (i.e., associated with a calorimetry measured whole body averaged value of 2.7-W/kg) indicated that the greatest amount of heating that occurred in association with these specific conditions was equal to or less than +1.6 °C.

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the T-Bar Assembly (T-Fastener). Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	1,064-mm ²	638-mm ²	2,361-mm ²	1,540-mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

Biosyn® is a registered trademark of US Surgical Corporation.

For more information, please call 1-844-4AVANOS (1-844-428-2667) in the United States, or visit our web site at www.avanos.com.

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*Registered Trademark or Trademark of Avanos Medical, Inc.,
or its affiliates. © 2018 AVNS. All rights reserved. 2019-01-28
Distributed in the USA by Avanos Medical Sales, LLC,
Avanos Medical, Inc., 5405 Windward Parkway, Alpharetta, GA 30004 USA.
Sponsored in Australia by Avanos Medical Australia Pty Ltd, 475 Victoria Avenue, Chatswood,