INTENDED PURPOSE

The Penumbra System is designed to remove thrombus and restore blood flow in the neurovasculature using aspiration.

DEVICE DESCRIPTION

The Penumbra System® is comprised of several devices:

- Penumbra Reperfusion Catheter
- Penumbra Aspiration Pump
- Penumbra Aspiration Pump Canister
- Penumbra Aspiration Tubing
- Penumbra Separator

The Reperfusion Catheter targets aspiration from the pump directly to the thrombus. The Separator may be used to clear the lumen of the Reperfusion Catheter should it become blocked with thrombus. The use of the Separator may not be necessary when using a Reperfusion Catheter with an I.D. of .054in [1.37 mm] or larger. The Reperfusion Catheter is introduced through a guide catheter or long sheath and into the intracranial vasculature and guided over a neurovascular guidewire to the site of the primary occlusion. SENDit™ Technology may be used with its prepackaged Reperfusion Catheter to facilitate access to the primary occlusion. The Penumbra Reperfusion Catheter is used with the Aspiration Pump to aspirate thrombus from an occluded vessel. As needed, a Penumbra Separator may be deployed from the Reperfusion Catheter to assist with thrombus removal. The Penumbra Separator is advanced and retracted through the Penumbra Reperfusion Catheter at the proximal margin of the primary occlusion to facilitate clearing of the thrombus from the Reperfusion Catheter tip. For the aspiration source, the Penumbra Reperfusion Catheter is used in conjunction with the Penumbra Aspiration Pump, which is connected using the Penumbra Aspiration Pump Canister and Penumbra Aspiration Tubing. The Penumbra Reperfusion Catheter is provided with a steam shaping mandrel, rotating hemostasis valve, peelable sheath and, optionally, SENDit Technology. The SENDit Technology compatibility is limited to its prepackaged Reperfusion Catheter. The Penumbra Separator is provided with an introducer and torque device. The devices are visible under fluoroscopy.

Penumbra System Reperfusion Catheters have a hydrophilic coating on the distal segment of the catheter shaft (see **Table 1**).

Table 1 – Hydrophilic Coating Length	
Penumbra Reperfusion Catheter	Length (cm)
3MAX	95
041 / 4MAX / RED 43	60
5MAX / ACE / Penumbra JET / RED 62 / RED 68 / RED 72 / RED 78	30

INDICATION FOR USEPenumbra Reperfusion Catheters and Separators

As part of the Penumbra System, the Reperfusion Catheters and Separators are indicated for use in the revascularization of patients with acute ischemic stroke secondary to intracranial large vessel occlusive disease (within the internal carotid, middle cerebral – M1 and M2 segments, basilar, and vertebral arteries) within 8 hours of symptom onset. Patients who are ineligible for intravenous tissue plasminogen activator (IV t-PA) or who fail IV t-PA therapy are candidates for treatment.

Penumbra Aspiration Tubing

As part of the Penumbra System, the Penumbra Sterile Aspiration Tubing is indicated to connect the Penumbra Reperfusion Catheters to the Penumbra Aspiration Pump.

Penumbra Aspiration Pump

The Penumbra Aspiration Pump is indicated as a vacuum source for Penumbra Aspiration Systems.

INTENDED PATIENT POPULATION

The intended patient population are those patients with acute ischemic stroke secondary to large vessel occlusion.

CONTRAINDICATIONS

There are no known contraindications.

WARNINGS

- The device is intended for single use only. Do not resterilize or reuse. Resterilization or reuse of the device may compromise the structural integrity of the device and/or lead to device failure which in turn may result in patient injury, illness, or death.
- Do not use kinked or damaged devices. Do not use open or damaged packages. Return all damaged devices and packaging to the manufacturer/distributor.

- Do not use automated high-pressure contrast injection equipment with the Penumbra Reperfusion Catheter or SENDit Technology because it may damage the device.
- Confirm vessel diameter and select appropriately sized Penumbra Reperfusion Catheter. Do not use in arteries with diameters smaller or equal to the distal outer diameter of the Penumbra Reperfusion Catheters. Refer to the Reperfusion Catheter labeling for dimensional information.
- Do not advance, retract or use any component of the Penumbra System against resistance without careful assessment of the cause using fluoroscopy. If the cause cannot be determined, withdraw the device or system as a unit. Unrestrained torquing or forced insertion of the catheter or separator against resistance may result in damage to the device or vessel.
- Do not use the Penumbra System with a pump other than the Penumbra Aspiration Pump.
- Failure to abide by the warnings in this labeling might result in damage to the device or device coating, which may necessitate intervention or result in serious adverse events.

PRECAUTIONS

- The Penumbra System should only be used by physicians who have received appropriate training in interventional neuro-endovascular techniques and treatment of acute ischemic stroke.
- Use prior to the “Use By” date.
- Use the Penumbra System in conjunction with fluoroscopic visualization.
- As in all fluoroscopy procedures, consider all necessary precautions to limit patient radiation exposure by using sufficient shielding, reducing fluoroscopy times, and modifying radiation technical factors whenever possible.
- Maintain a constant infusion of appropriate flush solution.
- When performing aspiration, ensure that the Penumbra Aspiration Tubing is open for only the minimum time needed to remove thrombus. Excessive aspiration or failure to close the Penumbra Aspiration Tubing when aspiration is complete is not recommended.
- The Penumbra Separator is not intended for use as a neurovascular guidewire. If repositioning of the Penumbra Reperfusion Catheter is necessary during the revascularization procedure, such repositioning should be performed over an appropriate neurovascular guidewire using standard microcatheter and guidewire techniques.
- As in all surgical interventions, monitoring of intra-procedural blood loss is recommended so that appropriate management may be instituted.
- Avoid using alcohol, antiseptic solutions, or other solvents to pre-treat the device because it may cause unpredictable changes in the coating which could affect the device safety and performance.
- The Penumbra Separators (041, 3MAX, 4MAX and 5MAX) contain a CMR substance (CAS Number: 7440-48-4).

POTENTIAL ADVERSE EVENTS

Possible complications include, but are not limited to, the following:

- allergic reaction and anaphylaxis from contrast media or device material
- acute vessel occlusion
- air embolism
- arrhythmia
- arteriovenous fistula
- death
- foreign body embolization
- emboli
- pseudoaneurysm
- hematoma or hemorrhage at access site
- residual thrombus due to inability to completely remove thrombus
- infection
- inflammation
- intracranial hemorrhage
- ischemia
- renal impairment or acute renal failure from contrast media
- neurological deficits including stroke
- vessel spasm, thrombosis, dissection, or perforation.
- radiation exposure that may lead to cataracts, skin reddening, burns, alopecia, or neoplasia from x-ray exposure

*If a serious incident related to the device occurs, contact your Penumbra Inc., representative and the competent authority in your representative country/region.

PROCEDURE

1. Refer to **Warnings**, **Precautions**, and **Potential Adverse Events** prior to use.
2. Refer to **Table 2** for device compatibility information.
3. As each device of the Penumbra System is used, remove the device from the packaging, and inspect for damage or for kinks.
4. Prepare Penumbra System devices for use by flushing the packaging hoop and device with heparinized saline (includes Penumbra Reperfusion Catheter, SENDit Technology, and Penumbra Separator only). Flushing with heparinized saline activates the hydrophilic coating located on the distal segment of the Penumbra Reperfusion Catheter and SENDit Technology.
5. Prepare a guide catheter or long sheath according to the manufacturer's *Instructions for Use*.
6. Place the guide catheter or long sheath into the appropriate artery that is proximal to the thrombus occlusion site.

Penumbra Reperfusion Catheter	Long Sheath / Guide Catheter	Guidewire	Penumbra Separator
3MAX	6F / 8F	0.014"	3MAX
041 / 4MAX / RED 43			041 / 4MAX
5MAX / ACE / Penumbra JET / RED 62 / RED 68 / RED 72 / RED 78			5MAX

The Penumbra Separator is intended to be used only as part of the Penumbra System.

PENUMBRA REPERFUSION CATHETER PREPARATION AND USE

1. Confirm vessel diameter and select an appropriately sized Penumbra Reperfusion Catheter. Ensure the catheter's labeled outer diameter is smaller than the treatment vessel diameter. Refer to the Reperfusion Catheter labeling for dimensional information.
2. The Penumbra Reperfusion Catheter tip may be shaped using the steam shaping mandrel provided. Insert the shaping mandrel into the Reperfusion Catheter tip, gently form the desired shape, and position in a flow of steam for ten seconds while rotating gently. To avoid damage from overheating, do not place the catheter directly at the steam outlet. Allow tip to cool prior to removing mandrel. Inspect Reperfusion Catheter for any damage that may have resulted from steam shaping; if any damage is found, do not use the catheter. **NOTE: Attempting to alter the shape of the Catheter by bending, twisting, or similar methods beyond the instructed methods may compromise the coating integrity. Such damage to the coating may not always be noticeable to the naked eye.**
3. Attach the Penumbra Aspiration Tubing to the Aspiration Pump Canister and turn on the Aspiration Pump (Refer to the Aspiration Pump *Operation Manual*). Confirm that the Penumbra Aspiration Pump reaches maximum vacuum level. Ensure the Penumbra Aspiration Tubing is switched to the OFF position.
4. Insert the Penumbra Reperfusion Catheter into a rotating hemostasis valve connected to the proximal hub of a guide catheter. If a guide catheter is not used, insert the Penumbra Reperfusion Catheter through the valve of the long sheath using the peelable sheath. After inserting the Penumbra Reperfusion Catheter, remove the peelable sheath from the vascular sheath, and peel from the Penumbra Reperfusion Catheter shaft.
5. Using conventional catheterization techniques under fluoroscopic guidance, advance the Penumbra Reperfusion Catheter into the target vessel over an appropriate neurovascular guidewire. Alternatively, if applicable, advance the Penumbra Reperfusion Catheter with SENDit Technology over the appropriate neurovascular guidewire to aide access to the site of occlusion as shown in **Figure 1**. Position the Penumbra Reperfusion Catheter proximal to the thrombus. Remove the guidewire and SENDit Technology (if present) from the Penumbra Reperfusion Catheter. **NOTE: Access to the site of occlusion may be aided by the triaxial use of a compatible inner microcatheter.**
6. Attach the rotating hemostasis valve provided to the Penumbra Reperfusion Catheter.
7. Connect the Aspiration Tubing to the side port of the rotating hemostasis valve as shown in **Figure 2**.
8. To begin aspiration, switch the Penumbra Aspiration Tubing to the ON position. To stop aspiration, switch the Penumbra Aspiration Tubing to the OFF position and turn off the Aspiration Pump.
9. Using a 5cc or 10cc syringe, aspirate approximately 5cc of blood from the Penumbra Reperfusion Catheter to remove any thrombus that may remain in the catheter.
10. Obtain a post-treatment angiogram by injecting contrast media through the guide catheters. Remove the Reperfusion Catheter from the guide catheter if necessary for adequate visualization.

Figure 1 – Reperfusion Catheter with SENDit Technology




Figure 2 – Assembled Aspiration System with Penumbra Separator

PENUMBRA SEPARATOR PREPARATION AND USE

1. If needed, select an appropriate size Penumbra Separator corresponding to the Reperfusion Catheter (see **Table 2**).
2. Using the introducer sheath provided, insert the Penumbra Separator through the rotating hemostasis valve and into the proximal hub of the Penumbra Reperfusion Catheter.
3. To begin aspiration, switch the Penumbra Aspiration Tubing to the ON position. Advance/retract the Penumbra Separator to assist with aspiration and removal of the thrombus from the Reperfusion Catheter as shown in **Figure 3**.

Figure 3 – Thrombus Aspiration Using the Penumbra Separator

4. To stop aspiration, switch the Penumbra Aspiration Tubing to the OFF position, and turn off the Aspiration Pump.
5. Remove the Penumbra Separator, if used.
6. If necessary, additional Penumbra Separators and Reperfusion Catheters may be used to further remove thrombus at the discretion of the physician.

CLINICAL DATA SUMMARY

The Penumbra Pivotal trial was a prospective, single arm trial to assess the safety and effectiveness of the Penumbra System. A total of 125 patients were enrolled at 24 international centers. Main entry criteria were NIH Stroke Scale (NIHSS) score ≥ 8 , presentation within 8 hours of symptom onset, and an occlusion (TIMI 0 or 1) of a treatable intracranial vessel. Patients who presented within 3 hours from symptom onset must have been ineligible for or refractory to tPA therapy. Primary endpoints were revascularization of the target vessel (TIMI 2 or 3) and incidence of procedural serious adverse events (SAE). All patients were followed and assessed up to 90 days from the procedure. All enrolled patients treated by the Penumbra System used the Reperfusion Catheter in the following sizes (026, 032 and 041). The Separators used were in the following sizes (026, 032, 041).

STATISTICAL ANALYSIS

A sample size of 125 was selected for this study in order to generate at least 115 evaluable patients. A patient was considered evaluable if the patient met eligibility requirements for primary endpoint assessment and was treated with the Penumbra System. Since this was an uncontrolled study, a formal hypothesis test was not planned for the overall response rate for establishing a TIMI grade of 2 or 3. Standard descriptive statistics for categorical endpoints were the number and percent of patients with each level of the endpoint. For numeric endpoints, the standard descriptive statistics included the number of non-missing observations (n), the mean, the median, the standard deviation, the minimum value, and the maximum value.

STUDY PROCEDURES

Arterial access to the target vessel was achieved using standard percutaneous technique and a four-vessel digital subtraction angiography was performed to define the angio-architecture of the occluded vascular segment. Each patient was required to have their cerebral large vessel occlusion confirmed by angiography and a site of primary occlusion designated a priori by the investigator. If the occlusion is in a treatable vessel, an appropriate guide catheter was positioned and the Penumbra Reperfusion Catheter advanced to the proximal end of the occluded target vessel. A patient was considered enrolled when the Reperfusion Catheter was deployed from the guide catheter.

KEY INCLUSION CRITERIA

- Clinical signs consistent with acute ischemic stroke, including:
 - Sudden numbness or weakness of face or extremities
 - Sudden confusion, trouble speaking or understanding speech
 - Sudden trouble walking, dizziness, or loss of balance
- From 18 to 79 years of age
- Neurological deficit resulting in a NIH Stroke Scale (NIHSS) score ≥ 8
- Complete occlusion (TIMI Grade 0) or contrast penetration with minimal perfusion but no distal branch filling (TIMI Grade I) in intracranial vessels accessible to the Penumbra Stroke System
- Signed informed consent
- Present within 8 hours of stroke symptom onset
- Ineligible or refractory to intravenous tPA therapy if presenting within 3 hours of symptom onset

KEY EXCLUSION CRITERIA

- Evidence of rapidly improving neurological signs of stroke at time of enrollment
- NIHSS >30 or coma
- Females who are pregnant
- Vessel tortuosity too difficult to allow endovascular access
- Known hemorrhagic diathesis, coagulation deficiency, or on oral anticoagulant therapy with an INR >3.0
- PTT greater than 2 times the lab normal
- Admission platelets $<30,000$
- Pre-existing neurological or psychiatric disease that could confound the study results
- Known severe allergy to contrast media
- Uncontrolled hypertension (defined as systolic blood pressure >185 mmHg or diastolic blood pressure >110 mmHg)
- CT evidence of the following conditions on presentation:
 - Significant mass effect with midline shift
 - Large hypodensity region $>1/3$ of the middle cerebral artery territory
 - Evidence of intracranial hemorrhage
- Angiographic evidence of an arterial stenosis proximal to the occlusion that could prevent thrombus removal
- Angiographic evidence of preexisting arterial injury
- Life expectancy less than 90 days
- Participation in another clinical investigation that could confound the evaluation of the study device

SUMMARY OF REASONS FOR EXCLUSION OF ENROLLED SUBJECTS, ESPECIALLY THOSE BASED ON ANGIOGRAPHIC CRITERIA

No enrolled patient was excluded for angiographic reasons.

PRIMARY EFFECTIVENESS ENDPOINT

The Primary Effectiveness Endpoint was revascularization of the occluded target vessel. Revascularization was defined by a TIMI score of II or III following use of the Penumbra System only and prior to the use of any additional or adjunctive devices or therapies. Angiographic results were adjudicated by an independent Core Laboratory.

Effectiveness	% (n/N) [95% CI]†
Primary Effectiveness Endpoint (TIMI score of II or III) <i>TIMI scores before adjunctive therapies are used</i>	81.6% (102/125) [73.7%, 88.0%]
Post-Device Revascularization Success (TIMI ≥ 2) <i>Any Adjunctive Therapy Use Considered Treatment Failure</i>	81.6% (102/125) [73.7%, 88.0%]
Post-Device Revascularization Success (TIMI ≥ 2) <i>Any IA Lytic Use Considered Treatment Failure</i>	81.6% (102/125) [73.7%, 88.0%]

†Exact Clopper-Pearson confidence intervals

PRIMARY SAFETY ENDPOINT

The Primary Safety Endpoint was incidence of device and/or procedure-related serious adverse events that occurred during the procedure as reported by the investigator. There were 18 procedural events reported in 16 patients (12.8%). Among these, 4 events in 3 patients were rated as serious by the CEC (2.4%).

Composite Adverse Events during the procedure	% (n/N) [95% CI]†
Composite SAE during the procedure	2.4% (3/125)
Primary Safety Endpoint: Composite Procedure/Device related SAE during the procedure	2.4% (3/125)

†Exact Clopper-Pearson confidence intervals

EACH SPECIFIC SAE THROUGH 90 DAYS OCCURRING IN MORE THAN 2% OF PATIENTS

SAE	% (n/N)
Composite SAEs	52.0% (65/125)
Cerebral Edema	11.2% (14/125)
ICH	8.0% (10/125)
Respiratory Dysfunction/Failure	8.0% (10/125)
Worsening Symptoms of Qualifying Stroke	7.2% (9/125)
Infection	5.6% (7/125)
Other Pulmonary Complications	5.6% (7/125)
Pneumonia	4.0% (5/125)
Cardiac Arrhythmia	3.2% (4/125)
Access Site Complication	2.4% (3/125)
Myocardial Infarction	2.4% (3/125)
Peripheral Vascular Complication	2.4% (3/125)
Renal Dysfunction/Failure	2.4% (3/125)
Stroke in New Territory	2.4% (3/125)

CATEGORIZATION OF AEs THROUGH 90 DAYS UNDER ANATOMICAL SYSTEM

Anatomical System	% (n/N)
Nervous System Disorders	62.4% (78/125)
Infections and Infestations	44.0% (55/125)
Cardiac Disorders	28.8% (36/125)
Metabolism and Nutrition Disorders	25.6% (32/125)
General Disorders and Administration Site Conditions	24.8% (31/125)
Respiratory, Thoracic and Mediastinal Disorders	24.8% (31/125)
Injury, Poisoning and Procedural Complications	18.4% (23/125)
Gastrointestinal Disorders	16.0% (20/125)
Blood and Lymphatic System Disorders	14.4% (18/125)
Vascular Disorder	11.2% (14/125)
Psychiatric Disorders	10.4% (13/125)
Renal and Urinary Disorders	8.0% (10/125)
Hepatobiliary Disorders	4.0% (5/125)
Investigations	3.2% (4/125)
Skin and Subcutaneous Tissue Disorders	3.2% (4/125)
Neoplasms Benign, Malignant and Unspecified (Incl Cysts And Polyps)	2.4% (3/125)
Musculoskeletal and Connective Tissue Disorders	0.8% (1/125)

CLINICAL BENEFITS AND PERFORMANCE CHARACTERISTICS

The main clinical benefit of the Penumbra System is the revascularization of patients with acute ischemic stroke secondary to intracranial large vessel occlusive disease leading to improved survival and reduced disability. The detailed listings of all known benefits are included in the Summary of Safety and Clinical Performance (SSCP).

SAFETY AND CLINICAL PERFORMANCE


A copy of the Summary of Safety and Clinical Performance can be viewed by searching the device name on the EUDAMED website: <https://ec.europa.eu/tools/eudamed>

Basic UDI-DI: 081454801PSSYSTEMDS

STORAGE CONDITIONS AND DISPOSAL

- Store in a cool and dry place.
- Dispose of device in accordance with hospital, administrative and/or local/government policy.

SYMBOLS GLOSSARY

	Attention, see Instructions For Use		
	Prescription only – US Federal Law restricts this device to use by or on the order of a physician		
	Nonpyrogenic		
	Sterilized using ethylene oxide. Single sterile barrier system with protective packaging outside.		
	Medical Device		Date of Manufacture
	Manufacturer		EU Representative
	Lot number		Do not reuse
	Use by		Catalogue number
	Do not resterilize		Not made with natural rubber latex
	Do not use if package is damaged		Contains Hazardous Substance CAS Number: 7440-48-4

WARRANTY

Penumbra, Inc. (Penumbra) warrants that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular use. Handling, storage, cleaning, and sterilization of this device as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond Penumbra's control directly affect the device and the results obtained from its use. Penumbra's obligation under this warranty is limited to the repair or replacement of this device and Penumbra shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this device. Penumbra neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. Penumbra assumes no liability with respect to devices reused, reprocessed, or resterilized and makes no warranties, expressed or implied, including but not limited to merchantability or fitness for intended use, with respect to such device.

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