

## EXAMPLE OF LABELS OF MONOLITHIC CEMENTED STEMS

Label of the first code of each kind of device for Monolithic Cemented Stems is following reported. For each label, all the information which changes in the other labels of the same device is showed.

The labels are reported as follows:

Section 1	Example of labels for CL stems
Section 2	Example of labels for SL Self Locking stems
Section 3	Example of labels for LOGICA MIRROR stems
Section 4	Example of labels for FRIENDLY stems
Section 5	Example of labels for H-MAX C stems

## **STERILIZATION – DATA ON THE LABELS**

The information about the sterilization method for the device is reported on the labels. The external cardboard box (first label in each of the following pages) indicates only the possible sterilization methods:



- **STERILE R:** sterilization by **irradiation** (Beta rays or Gamma rays according to the applicable sterilization method).



- **STERILE EO:** sterilization by **ethylene oxide**.

Always refer to device labels for information about the sterilization method.

For the information on applicable sterilization method for each device refer to Design Dossier.

**SECTION 1 – EXAMPLE OF LABELS FOR CL STEMS**

**Data Legend**

(the information which changes in the other labels of the same device is squared and showed in the table below).

Code	Size	Model	Material	Applicable Sterilization
2606.05.010	Size S	CL Trauma	AISI 316L	Beta rays / Gamma rays
2606.05.015	Size M	CL Trauma	AISI 316L	Beta rays / Gamma rays
2606.05.020	Size L	CL Trauma	AISI 316L	Beta rays / Gamma rays

Example of labels for the first code

REF 2606.05.010

Taper 12/14  
CL TRAUMA-AISI316L



Size  
S

CL TRAUMA Hip/Anca-AISI316L  
SMALL Taper/Cono=12/14  
Stem/Stelo  
Cemented Fem. Stem **EN**  
Stelo Fem. Cementato **IT**  
Vástago Fem. Cementado **ES**  
Zementierte Fem. **DE**  
Tige Fém. Cimentée **FR**  
Haste Fem. cimentada **PT**

2022-06-30 Qty 1 STERILE R Limacorporate S.p.a.

REF 2606.05.010

03

**EN** Cemented Fem. stem  
**IT** Stelo Fem. Cementato  
**ES** Vástago Fem. Cementado  
**DE** Zementierte Fem.  
**FR** Tige Fém. Cimentée  
**PT** Haste Fem. Cimentada

Taper 12/14 2022-08-31 Limacorporate S.p.A.

REF 2606.05.010

CL TRAUMA Hip/Anca-AISI316L  
Stem/Stelo  
SMALL Taper/Cono=12/14



03

2022-08-31 LOT 2017 08211 STERILE R 1700254



REF 2606.05.010

CL TRAUMA Hip/Anca-AISI316L  
Stem/Stelo  
SMALL Taper/Cono=12/14















03

2022-08-31 LOT 2017 08211 STERILE R 1700254



Limacorporate S.p.A. - Via Nazionale, 52  
33038 Villanova di San Daniele - Udine - Italy 0123

-  Single use
-  Do not use if package is damaged
-  Keep away from sunlight
-  Temperature limit
-  Use By
-  Reference Number
-  Keep dry
-  Lot Number
-  Caution
-  Sterilized using Ethylene Oxide
-  IFU Consult instructions for use
-  Manufacturer
-  Sterilized using irradiation

Limacorporate S.p.A. - Via Nazionale, 52  
33038 Villanova di San Daniele - Udine - Italy

CE 0123




**SECTION 2 – EXAMPLE OF LABELS FOR SL SELF LOCKING STEMS**
**Data Legend**

(the information which changes in the other labels of the same device is squared and showed in the table below).

**Code**  
REF 3515.15.003

Taper 12/14  
SL - Ti6Al4V



**Size**  
Size  
02

SL SELF LOCKING Stem/stelo  
#02 - (4.5mm)  
Taper/Cono 12/14 - Ti6Al4V

Cemented Fem. Stem **EN**  
Stelo Fem. Cementato **IT**  
Vástago Fem. Cementado **ES**  
Zementierte Fem. **DE**  
Tige Fém. Cimentée **FR**  
Haste Fem. Cimentada **PT**

2023-07-31 Qty 1 STERILE R Limacorporate S.p.a.

REF 3515.15.003 01

**EN** Cemented Fem. stem  
**IT** Stelo Fem. Cementato  
**ES** Vástago Fem. Cementato  
**DE** Zementierte Fem.  
**FR** Tige Fém. Cimentée  
**PT** Haste Fem. Cimentada


Taper 12/14 2022-03-31 Limacorporate S.p.A.

REF 3515.15.003 01

SL SELF LOCKING Stem/stelo  
Taper/Cono 12/14 - Ti6Al4V  
#02 - (4.5mm)

2022-03-31 LOT 2017 01814 STERILE R 1700095




(01)08033390003421(17)220331

REF 3515.15.003 01


SL SELF LOCKING Stem/stelo  
Taper/Cono 12/14 - Ti6Al4V  
#02 - (4.5mm)

2022-03-31 LOT 2017 01814 STERILE R 1700095




(20)01(240)93515150030(10)1701814(31)1700095(92)8139  
(01)08033390003421(17)220331

Limacorporate S.p.A. - Via Nazionale, 52  
33038 Villanova di San Daniele - Udine - Italy 0123



Limacorporate S.p.A. - Via Nazionale, 52  
33038 Villanova di San Daniele - Udine - Italy

**CE** 0123 

Code	Size
3515.15.003	#02 – (4.5mm)
3515.15.005	#0 – (6.5mm)
3515.15.010	#1 – (7.5mm)
3515.15.015	#1 PLUS – (8.75mm)
3515.15.020	#2 – (10mm)
3515.15.030	#3 – (11-5mm)
3515.15.040	#4 – (12.5mm)
3515.15.050	#5 – (13.5mm)
3515.15.060	#6 – (15mm)

Code	Size
3515.15.070	#7 – (17.5mm)
3515.15.080	#8 – (20mm)

**Example of labels for the first code**

**REF** 3515.15.003  
Taper 12/14  
SL - Ti6Al4V



**Size**  
02

SL SELF LOCKING Stem/stelo  
#02 - (4.5mm)  
Taper/Cono 12/14 - Ti6Al4V

Cemented Fem. Stem **EN**  
Stelo Fem. Cementato **IT**  
Vástago Fem. Cementado **ES**  
Zementierte Fem. **DE**  
Tige Fém. Cimentée **FR**  
Haste Fem. Cimentada **PT**

2023-07-31 Qty 1 **STERILE R** Limacorporate S.p.a.

**REF** 3515.15.003 01

**EN** Cemented Fem. stem  
**IT** Stelo Fem. Cementato  
**ES** Vástago Fem. Cementato  
**DE** Zementierte Fem.  
**FR** Tige Fém. Cimentée  
**PT** Haste Fem. Cimentada

Taper 12/14 2022-03-31 Limacorporate S.p.A.

**REF** 3515.15.003 01

SL SELF LOCKING Stem/stelo  
Taper/Cono 12/14 - Ti6Al4V  
#02 - (4.5mm)

2022-03-31 **LOT** 2017 01814 **STERILE R** 1700095



(20)01(240)93515150030(10)1701814(91)1700095(92)8139



(01)08033390003421(17)220331


**REF** 3515.15.003 01

SL SELF LOCKING Stem/stelo  
Taper/Cono 12/14 - Ti6Al4V  
#02 - (4.5mm)

2022-03-31 **LOT** 2017 01814 **STERILE R** 1700095



(20)01(240)93515150030(10)1701814(91)1700095(92)8139



(01)08033390003421(17)220331

Limacorporate S.p.A. - Via Nazionale, 52  
33038 Villanova di San Daniele - Udine - Italy 0123

-  Single use
-  Do not use if package is damaged
-  Keep away from sunlight
-  Temperature limit
-  Use By
-  Reference Number
-  Keep dry
-  Lot Number
-  Caution
-  Sterilized using Ethylene Oxide
-  IFU Consult Instructions for use
-  Manufacturer
-  Sterilized using irradiation

**SECTION 3 – EXAMPLE OF LABELS FOR LOGICA MIRROR STEMS**
**Data Legend**


(the information which changes in the other labels of the same device is squared and showed in the table below).

**Code**  
REF **3715.07.005**

LOGICA Mirror Hip - FeCrNiMnMoNbN  
Bright Stem/Stelo Lucido  
#0 H=130mm Taper/Cono=12/14

**Size**  
#0  
H130

LOGICA Mirror



2024-04-30 Qty 1 STERILE EO Limacorporate S.p.a.

REF **3715.07.005** 03

EN Cemented Fem. Stem  
IT Stelo Fem. Cementato  
ES Vástago Fem. Cementado  
DE Zementierte Fem.  
FR Tige Fem. Cimentée  
PT Haste Fem. Cimentata


2022-05-31 Limacorporate S.p.A.

REF **3715.07.005** 03

LOGICA Mirror Hip - FeCrNiMnMoNbN  
Bright Stem/Stelo Lucido  
#0 H=130mm Taper/Cono=12/14

2022-05-31 LOT 2017 04630 STERILE EO 1700177




(20)03(240)93715070050(10)1704630(91)1700177(92)6801  
(01)08033390003711(17)220531

REF **3715.07.005** 03

LOGICA Mirror Hip - FeCrNiMnMoNbN  
Bright Stem/Stelo Lucido  
#0 H=130mm Taper/Cono=12/14


2022-05-31 LOT 2017 04630 STERILE EO 1700177





(20)03(240)93715070050(10)1704630(91)1700177(92)6801  
(01)08033390003711(17)220531


Limacorporate S.p.A. - Via Nazionale, 52  
33038 Villanova di San Daniele - Udine - Italy 0123


 Single use


 Do not use if package is damaged


 Keep away from sunlight

 Temperature limit


  


 Use By


 REF Reference Number

 Keep dry


  


 LOT Lot Number


 Caution

 STERILE EO Sterilized using Ethylene Oxide


  


 IFU Consult instructions for use


 Manufacturer

 STERILE R Sterilized using irradiation

 Limacorporate S.p.A. - Via Nazionale, 52  
33038 Villanova di San Daniele - Udine - Italy

 0123



Code	Size
3715.07.005	#0 H130
3715.07.010	#1 H138
3715.07.020	#2 H144
3715.07.030	#3 H150
3715.07.040	#4 H156
3715.07.050	#5 H162


**Example of labels for the first code**



REF 3715.07.005

12/14

LOGICA Mirror



#0

H130

LOGICA Mirror Hip - FeCrNiMnMoNbN  
Bright Stem/Stelo Lucido  
#0 H=130mm Taper/Cono=12/14

Cemented Fem. Stem **EN**

Stelo Fem. Cementato **IT**

Vástago Fem. Cementado **ES**

Zementierte Fem. **DE**

Tige Fém. Cimentée **FR**

Haste Fem. cimentada **PT**

REF 3715.07.005

03

**EN** Cemented Fem. Stem

**IT** Stelo Fem. Cementato

**ES** Vástago Fem. Cementado

**DE** Zementierte Fem.

**FR** Tige Fem. Cimentée

**PT** Haste Fem. Cimentada

2024-04-30 Qty 1 STERILE EO Limacorporate S.p.a.

12/14 2022-05-31 Limacorporate S.p.A.

REF 3715.07.005

LOGICA Mirror Hip - FeCrNiMnMoNbN

Bright Stem/Stelo Lucido

#0 H=130mm Taper/Cono=12/14

2022-05-31 LOT 2017 04630 STERILE EO 1700177



(20)03(240)93715070050(10)1704630(91)1700177(92)6801

(01)08033390003711(17)220531



03

REF 3715.07.005

LOGICA Mirror Hip - FeCrNiMnMoNbN

Bright Stem/Stelo Lucido

#0 H=130mm Taper/Cono=12/14

2022-05-31 LOT 2017 04630 STERILE EO 1700177



(20)03(240)93715070050(10)1704630(91)1700177(92)6801

(01)08033390003711(17)220531



03

Limacorporate S.p.A. - Via Nazionale, 52  
33038 Villanova di San Daniele - Udine - Italy 0123



Single use



Do not use if  
package  
is damaged



Keep away  
from sunlight



Temperature limit



Use By



Reference  
Number



Keep dry



Lot Number



Caution



Sterilized  
using  
Ethylene  
Oxide



IFU Consult  
instructions  
for use



Manufacturer



Sterilized  
using  
irradiation



Limacorporate S.p.A. - Via Nazionale, 52  
33038 Villanova di San Daniele - Udine - Italy

CE  
0123




**SECTION 4 – EXAMPLE OF LABELS FOR FRIENDLY STEMS**

**Data Legend**

(the information which changes in the other labels of the same device is squared and showed in the table below).

**Code** REF **4110.07.220**

12/14 FRIENDLY



**Angle**

**Size** #1 H133

FRIENDLY Hip - FeCrNiMnMoNbN  
 Stem/Stelo - Taper/Cono=12/14  
 #1 125degree/gradi H=133mm

Cemented Fem. Stem **EN**  
 Stelo Fem. Cementato **IT**  
 Vástago Fem. Cementado **ES**  
 Zementierte Fem. **DE**  
 Tige Fém. Cimentée **FR**  
 Haste Fem. cimentada **PT**

2022-08-31 Qty 1 STERILE R Limacorporate S.p.a.

REF **4110.07.220** 06

**EN** Cemented Fem. Stem  
**IT** Stelo Fem. Cementato  
**ES** Vástago Fem. Cementado  
**DE** Zementierte Fem.  
**FR** Tige Fem. Cimentée  
**PT** Haste Fem. Cimentata

12/14 2021-06-30 Limacorporate S.p.A.

REF **4110.07.220** 06

FRIENDLY Hip - FeCrNiMnMoNbN  
 Stem/Stelo - Taper/Cono=12/14  
 #1 125degree/gradi H=133mm

2021-06-30 LOT 2015 19164 STERILE R 1600167



(20)06(240)94110072200(10)1519164(91)1600167(92)6966  
 (01)08033390004305(17)210630

REF **4110.07.220** 06














FRIENDLY Hip - FeCrNiMnMoNbN  
 Stem/Stelo - Taper/Cono=12/14  
 #1 125degree/gradi H=133mm

2021-06-30 LOT 2015 19164 STERILE R 1600167




(20)06(240)94110072200(10)1519164(91)1600167(92)6966  
 (01)08033390004305(17)210630

Limacorporate S.p.A. - Via Nazionale, 52  
 33038 Villanova di San Daniele - Udine - Italy 0123

-  Single use
-  Do not use if package is damaged
-  Keep away from sunlight
-  Temperature limit
-  Use By
-  Reference Number
-  Keep dry
-  Lot Number
-  Caution
-  Sterilized using Ethylene Oxide
-  IFU Consult instructions for use
-  Manufacturer
-  Sterilized using irradiation

Limacorporate S.p.A. - Via Nazionale, 52  
 33038 Villanova di San Daniele - Udine - Italy

CE 0123 

Code	Size	Angle
4110.07.220	#1 H133	125°
4110.07.230	#2 H135	125°
4110.07.240	#3 H137	125°
4110.07.320	#1 H135	132°
4110.07.330	#2 H137	132°
4110.07.340	#3 H139	132°

**Example of labels for the first code**



REF 4110.07.220

12/14

FRIENDLY



#1

H133

FRIENDLY Hip - FeCrNiMnMoNbN  
Stem/Stelo - Taper/Cono=12/14  
#1 125degree/gradi H=133mm

Cemented Fem. Stem **EN**  
Stelo Fem. Cementato **IT**  
Vástago Fem. Cementado **ES**  
Zementierte Fem. **DE**  
Tige Fém. Cimentée **FR**  
Haste Fem. cimentada **PT**

2022-08-31 Qty 1 STERILE R Limacorporate S.p.a.

REF 4110.07.220

06

**EN** Cemented Fem. Stem  
**IT** Stelo Fem. Cementato  
**ES** Vástago Fem. Cementado  
**DE** Zementierte Fem.  
**FR** Tige Fem. Cimentée  
**PT** Haste Fem. Cimentata

12/14 2021-06-30 Limacorporate S.p.A.

REF 4110.07.220

FRIENDLY Hip - FeCrNiMnMoNbN

Stem/Stelo - Taper/Cono=12/14

#1 125degree/gradi H=133mm

2021-06-30 LOT 2015 19164 STERILE R 1600167



(20)06(240)94110072200(10)1519164(91)1600167(92)6966



(01)08033390004305(17)210630

06

REF 4110.07.220

FRIENDLY Hip - FeCrNiMnMoNbN

Stem/Stelo - Taper/Cono=12/14

#1 125degree/gradi H=133mm

2021-06-30 LOT 2015 19164 STERILE R 1600167



(20)06(240)94110072200(10)1519164(91)1600167(92)6966



(01)08033390004305(17)210630

Limacorporate S.p.A. - Via Nazionale, 52  
33038 Villanova di San Daniele - Udine - Italy 0123

- Single use
- Do not use if package is damaged
- Keep away from sunlight
- Temperature limit
- Use By
- Reference Number
- Keep dry
- Lot Number
- Caution
- Sterilized using Ethylene Oxide
- IFU Consult instructions for use
- Manufacturer
- Sterilized using irradiation


Limacorporate S.p.A. - Via Nazionale, 52  
33038 Villanova di San Daniele - Udine - Italy

CE  
0123



**SECTION 5 – EXAMPLE OF LABELS FOR H-MAX C STEMS**
**Data Legend**

(the information which changes in the other labels of the same device is squared and showed in the table below).



**Code** REF 4261.07.090

**Model** H-MAX C Lateralized

**Size** #9 H134

**Lat** Lat

H-MAX C Lateral. Hip-FeCrNiMnMoNbN  
 Stem/Stelo - Taper/Cono 12/14  
 #9 H=134mm

Cemented Fem. Stem EN  
 Stelo Fem. Cementato IT  
 Vástago Fem. Cementado ES  
 Zementierte Fem. DE  
 Tige Fem. Cimentée FR  
 Haste Fem. cimentada PT

Cemented Use Only  
 2022-08-31 Qty 1 STERILE R Limacorporate S.p.a.

REF 4261.07.090

H-MAX C Lateral. Hip-FeCrNiMnMoNbN  
 Stem/Stelo - Taper/Cono 12/14  
 #9 H=134mm

2022-08-31 LOT 2017 01741 STERILE R 1700260




(20)01(240)94261070900(10)1701741(91)1700260(92)5358  
 (01)06033390064378(17)220631

REF 4261.07.090


H-MAX C Lateral. Hip-FeCrNiMnMoNbN  
 Stem/Stelo - Taper/Cono 12/14  
 #9 H=134mm

2022-08-31 LOT 2017 01741 STERILE R 1700260




(20)01(240)94261070900(10)1701741(91)1700260(92)5358  
 (01)06033390064378(17)220631

Limacorporate S.p.A. - Via Nazionale, 52  
 33038 Villanova di San Daniele - Udine - Italy 0123



Limacorporate S.p.A. - Via Nazionale, 52  
 33038 Villanova di San Daniele - Udine - Italy

CE 0123




Code	Size	Model
4261.07.090	#9 H134	H-MAX C lateralized
4261.07.100	#10 H139	H-MAX C lateralized
4261.07.110	#11 H144	H-MAX C lateralized
4261.07.120	#12 H149	H-MAX C lateralized
4261.07.130	#13 H154	H-MAX C lateralized
4261.07.140	#14 H159	H-MAX C lateralized
4261.07.150	#15 H164	H-MAX C lateralized
4261.07.160	#16 H169	H-MAX C lateralized
4261.07.170	#17 H174	H-MAX C lateralized
4261.07.180	#18 H179	H-MAX C lateralized
4260.07.090	#9 H134	H-MAX C standard
4260.07.100	#10 H139	H-MAX C standard

Code	Size	Model
4260.07.110	#11 H144	H-MAX C standard
4260.07.120	#12 H149	H-MAX C standard
4260.07.130	#13 H154	H-MAX C standard
4260.07.140	#14 H159	H-MAX C standard
4260.07.150	#15 H164	H-MAX C standard
4260.07.160	#16 H169	H-MAX C standard
4260.07.170	#17 H174	H-MAX C standard
4260.07.180	#18 H179	H-MAX C standard

**Example of labels for the first code**

**REF 4261.07.090**  
12/14  
H-MAX C Lateralized



**#9**  
**H134**  
**Lat**

Cemented Use Only

H-MAX C Lateral. Hip-FeCrNiMnMoNbN  
Stem/Stelo - Taper/Cono 12/14  
#9 H=134mm

Cemented Fem. Stem **EN**  
Stelo Fem. Cementato **IT**  
Vástago Fem. Cementado **ES**  
Zementierte Fem. **DE**  
Tige Fém. Cimentée **FR**  
Haste Fem. cimentada **PT**

2022-08-31 Qty 1 **STERILE R** Limacorporate S.p.a.

**REF 4261.07.090** 01

**EN** Cemented Fem. Stem  
**IT** Stelo Fem. Cementato  
**ES** Vástago Fem. Cementado  
**DE** Zementierte Fem.  
**FR** Tige Fem. Cimentée  
**PT** Haste Fem. Cimentata

12/14 2022-08-31 Limacorporate S.p.a.

**REF 4261.07.090** 01

H-MAX C Lateral. Hip-FeCrNiMnMoNbN  
Stem/Stelo - Taper/Cono 12/14  
#9 H=134mm


2022-08-31 **LOT** 2017 01741 **STERILE R** 1700260



**REF 4261.07.090** 01

H-MAX C Lateral. Hip-FeCrNiMnMoNbN  
Stem/Stelo - Taper/Cono 12/14  
#9 H=134mm

2022-08-31 **LOT** 2017 01741 **STERILE R** 1700260



Limacorporate S.p.A. - Via Nazionale, 52  
33038 Villanova di San Daniele - Udine - Italy 0123

- Single use
- Do not use if package is damaged
- Keep away from sunlight
- Temperature limit
- Use By
- Reference Number
- Keep dry
- Lot Number
- Caution
- Sterilized using Ethylene Oxide
- IFU Consult instructions for use
- Manufacturer
- Sterilized using irradiation

Limacorporate S.p.A. - Via Nazionale, 52  
33038 Villanova di San Daniele - Udine - Italy

**CE** 0123 