

EXAMPLE OF LABELS OF FEMORAL HEADS

Label of the first code of each kind of device for Femoral Heads is following reported. For each label, all the information which changes in the other labels of the same device is showed.

The labels are reported as follows:

Section 1	Example of labels for Heads (BioloX® Delta) - Taper 12/14
Section 2	Example of labels for Heads (CoCrMo) - Taper 12/14
Section 3	Example of labels for Heads (AISI 316/L) - Taper 12/14
Section 4	Example of labels for Heads (FeCrNiMnMoNbN) - Taper 12/14
Section 5	Example of labels for Revision Heads (BioloX® Delta + Ti6Al4V) - Taper 12/14
Section 6	Example of labels for Heads (CoCrMo)
Section 7	Example of labels for Bipolar Heads (AISI 316/L) with Liner (UHMWPE)
Section 8	Example of labels for Bipolar Heads (CoCrMo) with Liner (UHMWPE)
Section 9	Example of labels for Lock Self-Centering Heads (AISI 316/L) with Liner (UHMWPE)
Section 10	Example of labels for Lock Self-Centering Heads (CoCrMo) with Liner (UHMWPE)
Section 11	Example of labels for Femoral Cephalic Heads (AISI 316/L) - Taper 12/14

STERILIZATION – DATA ON THE LABELS

The information about the sterilization method for the device is reported on the labels. The external cardboard box (first label in each of the following pages) indicates only the possible sterilization methods:



- **STERILE R:** sterilization by **irradiation** (Beta rays or Gamma rays according to the applicable sterilization method).



- **STERILE EO:** sterilization by **ethylene oxide**.

Labels provided below are here included as example. For the sterilization method applicable to each product refer to Chapter 1.2.3.1.

SECTION 1 – EXAMPLE OF LABELS FOR HEADS (BioloX® Delta) - TAPER 12/14

Data Legend


(the information which changes in the other labels of the same device is squared and showed in the table below).

Code	Size	Diameter
5010.42.281	S	28
5010.42.282	M	28
5010.42.283	L	28
5010.42.321	S	32
5010.42.322	M	32
5010.42.323	L	32
5010.42.361	S	36
5010.42.362	M	36
5010.42.363	L	36
5010.42.364	XL	36
5010.42.401	S	40
5010.42.402	M	40
5010.42.403	L	40
5010.42.404	XL	40

Example of labels for the first code

REF 5010.42.281

D=28mm
BIOLOX DELTA



Femoral Head/Testa Femorale
BIOLOX DELTA Dia=28mm
SHORT Taper/Cono=12/14

Femoral Head **EN**
Testa Femorale **IT**
Cabeza Femoral **ES**
Prothesenkopf **DE**
Tete Femoral **FR**
Cabeça Femoral **PT**

2022-07-31 **qty:** 1 **STERILE**EO **Limacorporate S.p.a.**

- Single use
- Do not use if package is damaged
- Keep away from sunlight
- Temperature limit
- Use By
- Reference Number
- Keep dry
- Lot Number
- Caution
- Sterilized using Ethylene Oxide
- IFU Consult instructions for use
- Manufacturer
- Sterilized using irradiation

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE
0123



REF 5010.42.281

Femoral Head/Testa Femorale
BIOLOX DELTA Dia=28mm
SHORT Taper/Cono=12/14

2022-07-31 **LOT** 2017 81023 **STERILE**EO 1700239



REF 5010.42.281

Femoral Head/Testa Femorale
BIOLOX DELTA Dia=28mm
SHORT Taper/Cono=12/14

2022-07-31 **LOT** 2017 81023 **STERILE**EO 1700239




Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE
0123



SECTION 2 – EXAMPLE OF LABELS FOR HEADS (CoCrMo) - TAPER 12/14
Data Legend

(the information which changes in the other labels of the same device is squared and showed in the table below).



Code
REF 5010.09.221

Diameter
D=22mm
CoCrMo

Size
S

2022-08-31 QTY: 1 STERILE EO Sterilization Limacorporate S.p.a.

Single use Do not use if package is damaged Keep away from sunlight Temperature limit

Use By REF Reference Number Keep dry

LOT Lot Number Caution Sterilized using Ethylene Oxide Sterilized using irradiation

IFU Consult instructions for use Manufacturer

Limacorporate S.p.A. - Via Nazionale, 52 33038 Villanova di San Daniele - Udine - Italy CE 0123

REF 5010.09.221
Femoral Head/Testa Femorale
CoCrMo Dia=22mm
SHORT Taper/Cono=12/14
Testa Femorale IT
Cabeza Femoral ES
Prothesenkopf DE
Tete Femoral FR
Cabeça Femoral PT

2022-08-31 LOT 2016 10428 STERILE EO 1700245

REF 5010.09.221
Femoral Head/Testa Femorale
CoCrMo Dia=22mm
SHORT Taper/Cono=12/14
2022-08-31 LOT 2016 10428 STERILE EO 1700245

Limacorporate S.p.A. - Via Nazionale, 52 33038 Villanova di San Daniele - Udine - Italy CE 0123


Code	Size	Diameter
5010.09.221	S	22
5010.09.222	M	22
5010.09.223	L	22
5010.09.224	XL	22
5010.09.281	S	28
5010.09.282	M	28
5010.09.283	L	28
5010.09.284	XL	28
5010.09.285	XXL	28
5010.09.286	XXXL	28
5010.09.321	S	32
5010.09.322	M	32
5010.09.323	L	32
5010.09.324	XL	32
5010.09.325	XXL	32
5010.09.326	XXXL	32
5010.09.361	S	36
5010.09.362	M	36
5010.09.363	L	36
5010.09.364	XL	36
5010.09.365	XXL	36

Code	Size	Diameter
5010.09.366	XXXL	36

Example of labels for the first code

REF 5010.09.221

D=22mm
CoCrMo



S

Femoral Head/Testa Femorale
CoCrMo Dia=22mm
SHORT Taper/Cono=12/14

Femoral Head **EN**

Testa Femorale **IT**

Cabeza Femoral **ES**

Prothesenkopf **DE**

Tete Femoral **FR**

Cabeça Femoral **PT**

2022-08-31 **qty:** 1 **STERILE EO** Limacorporate S.p.a.

- Single use
- Do not use if package is damaged
- Keep away from sunlight
- Temperature limit
- Use By
- Reference Number
- Keep dry
- Lot Number
- Caution
- Sterilized using Ethylene Oxide
- Sterilized using irradiation
- IFU Consult instructions for use
- Manufacturer

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE
0123



REF 5010.09.221

Femoral Head/Testa Femorale
CoCrMo Dia=22mm
SHORT Taper/Cono=12/14

2022-08-31 **LOT** 2016 10428 **STERILE EO** 1700245



(20)03(240)95010092210(10)1610428(91)1700245(92)5516

(01)08033390006750(17)220831



REF 5010.09.221

Femoral Head/Testa Femorale
CoCrMo Dia=22mm
SHORT Taper/Cono=12/14

2022-08-31 **LOT** 2016 10428 **STERILE EO** 1700245



(20)03(240)95010092210(10)1610428(91)1700245(92)5516

(01)08033390006750(17)220831




Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE
0123

SECTION 3 – EXAMPLE OF LABELS FOR HEADS (AISI 316/L) - TAPER 12/14

Data Legend

(the information which changes in the other labels of the same device is squared and showed in the table below).



Code	Size	Diameter
5010.05.221	S	22
5010.05.222	M	22
5010.05.223	L	22
5010.05.224	XL	22
5010.05.281	S	28
5010.05.282	M	28
5010.05.283	L	28
5010.05.284	XL	28
5010.05.285	XXL	28
5010.05.286	XXXL	28
5010.05.321	S	32
5010.05.322	M	32
5010.05.323	L	32
5010.05.324	XL	32
5010.05.325	XXL	32
5010.05.326	XXXL	32

Example of labels for the first code

REF 5010.05.221

D=22mm
AISI 316L



Femoral Head/Testa Femorale
AISI 316L Dia=22mm
SHORT Taper/Cono=12/14

Femoral Head EN
Testa Femorale IT
Cabeza Femoral ES
Prothesenkopf DE
Tete Femoral FR
Cabeça Femoral PT

2022-05-31 **qty:** 1 **STERILE** EO Limacorporate S.p.a.

- Single use
- Do not use if package is damaged
- Keep away from sunlight
- Temperature limit
- Use By
- Reference Number
- Keep dry
- Lot Number
- Caution
- Sterilized using Ethylene Oxide
- IFU Consult instructions for use
- Manufacturer
- Sterilized using Irradiation

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy



REF 5010.05.221

Femoral Head/Testa Femorale
AISI 316L Dia=22mm
SHORT Taper/Cono=12/14

2022-05-31 **LOT** 2016 09318 **STERILE** EO 1700162



REF 5010.05.221

Femoral Head/Testa Femorale
AISI 316L Dia=22mm
SHORT Taper/Cono=12/14

2022-05-31 **LOT** 2016 09318 **STERILE** EO 1700162



Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy



SECTION 4 – EXAMPLE OF LABELS FOR HEADS (FeCrNiMnMoNbN) - TAPER

12/14

Data Legend

(the information which changes in the other labels of the same device is squared and showed in the table below).

Code
REF 2416.07.221
D=22mm
ACCENT FeCrNiMnMoNbN

Size
#-2

Sterilization
STERILE EO

ACCENT - FeCrNiMnMoNbN
Femoral Head / Testa Femorale
Dia=22mm Taper/Cono=12/14
Femoral Head **EN**
Testa Femorale **IT**
Cabeza Femoral **ES**
Prothesenkopf **DE**
Tete Femoral **FR**
Cabeça Femoral **PT**

2021-02-28 qty: 1 STERILE EO 1600040
Lmacorporate S.p.a.

Single use Do not use if package is damaged Keep away from sunlight Temperature limit
Use By REF Reference Number Keep dry
LOT Lot Number Caution Sterilized using Ethylene Oxide
IFU Consult instructions for use Manufacturer Sterilized using irradiation

LimaCorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy 0123

Code
REF 2416.07.221
ACCENT - FeCrNiMnMoNbN
Femoral Head / Testa Femorale
Dia=22mm Taper/Cono=12/14
#-2

Sterilization
STERILE EO 1600040

2021-02-28 LOT 2015 11612 STERILE EO 1600040
Lmacorporate S.p.a.

Single use Do not use if package is damaged Keep away from sunlight Temperature limit
Use By REF Reference Number Keep dry
LOT Lot Number Caution Sterilized using Ethylene Oxide
IFU Consult instructions for use Manufacturer Sterilized using irradiation

LimaCorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy 0123


Code	Size
2416.07.221	# -2
2416.07.222	#0
2416.07.223	#+4

Example of labels for the first code

REF 2416.07.221

D=22mm

ACCENT FeCrNiMnMoNbN



#-2

ACCENT - FeCrNiMnMoNbN
Femoral Head / Testa Femorale

#-2 Dia=22mm Taper/Cono=12/14

Femoral Head **EN**

Testa Femorale **IT**


Cabeza Femoral **ES**

Prothesenkopf **DE**

Tete Femoral **FR**

Cabeça Femoral **PT**

2021-02-28 qty:1 **STERILE EO** Limacorporate S.p.a.



Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE
0123



REF 2416.07.221

ACCENT - FeCrNiMnMoNbN

Femoral Head / Testa Femorale

#-2 Dia=22mm Taper/Cono=12/14

2021-02-28 **LOT** 2015 11612 **STERILE EO** 1600040



(20)00(240)92416072210(10)1511612(91)1600040(92)6366



(01)08033390003018(17)210228

REF 2416.07.221

ACCENT - FeCrNiMnMoNbN

Femoral Head / Testa Femorale

#-2 Dia=22mm Taper/Cono=12/14

2021-02-28 **LOT** 2015 11612 **STERILE EO** 1600040



(20)00(240)92416072210(10)1511612(91)1600040(92)6366



(01)08033390003018(17)210228

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy


CE
0123

Data Legend

(the information which changes in the other labels of the same device is squared and showed in the table below).

Code
REF 5010.07.281

D=28mm
FeCrNiMnMoNbN




Size
S

2018-03-08 qty: 1 STERILE EO Limacorporate S.p.a.

Sterilization

- Single use
- Do not use if package is damaged
- Keep away from sunlight
- Temperature limit
- Use By
- REF Reference Number
- Keep dry
- LOT Lot Number
- Caution
- Sterilized using Ethylene Oxide
- IFU Consult instructions for use
- Manufacturer
- Sterilized using irradiation

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE 0123 


REF 5010.07.281

Femoral Head/Testa Femorale
FeCrNiMnMoNbN Dia=28mm

SHORT Taper/Cono=12/14

Femoral Head EN
Testa Femorale IT
Cabeza Femoral ES
Prothesenkopf DE
Tete Femoral FR
Cabeça Femoral PT

2018-03-08 LOT 2013 00361 STERILE EO 1300084




(20)04(240)95010072810(10)1300361(91)1300084(92)5502
(01)08033390006699(17)180308

REF 5010.07.281

Femoral Head/Testa Femorale
FeCrNiMnMoNbN Dia=28mm


SHORT Taper/Cono=12/14

2018-03-08 LOT 2013 00361 STERILE EO 1300084



(20)04(240)95010072810(10)1300361(91)1300084(92)5502
(01)08033390006699(17)180308

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy


CE 0123 

Code	Size
5010.07.281	S
5010.07.282	M
5010.07.283	L
5010.07.284	XL
5010.07.285	XXL
5010.07.286	XXXL

Example of labels for the first code

REF 5010.07.281

D=28mm
FeCrNiMnMoNbN



Femoral Head/Testa Femorale
FeCrNiMnMoNbN Dia=28mm
SHORT Taper/Cono=12/14
Femoral Head EN
Testa Femorale IT
Cabeza Femoral ES
Prothesenkopf DE
Tete Femoral FR
Cabeça Femoral PT

2018-03-08 Qty: 1 STERILE EO Limacorporate S.p.a.

- Single use
- Do not use if package is damaged
- Keep away from sunlight
- Temperature limit
- Use By
- Reference Number
- Keep dry
- Lot Number
- Caution
- Sterilized using Ethylene Oxide
- Sterilized using irradiation
- IFU Consult instructions for use
- Manufacturer

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy



REF 5010.07.281

Femoral Head/Testa Femorale
FeCrNiMnMoNbN Dia=28mm
SHORT Taper/Cono=12/14



2018-03-08 LOT 2013 00361 STERILE EO 1300084



REF 5010.07.281

Femoral Head/Testa Femorale
FeCrNiMnMoNbN Dia=28mm
SHORT Taper/Cono=12/14



2018-03-08 LOT 2013 00361 STERILE EO 1300084




Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy



SECTION 5 – EXAMPLE OF LABELS FOR REVISION HEADS (Biolox ® Delta + Ti6Al4V) - TAPER 12/14
Data Legend

(the information which changes in the other labels of the same device is squared and showed in the table below).



Code	Size	Diameter
5010.42.021	S	28
5010.42.022	M	28
5010.42.023	L	28
5010.42.024	XL	28
5010.42.031	S	32
5010.42.032	M	32
5010.42.033	L	32
5010.42.034	XL	32
5010.42.041	S	36
5010.42.042	M	36
Code	Size	Diameter
5010.42.043	L	36
5010.42.044	XL	36
5010.42.051	S	40
5010.42.052	M	40
5010.42.053	L	40
5010.42.054	XL	40

Example of labels for the first code

REF 5010.42.021

D=28mm
BIOLOX DELTA+Ti6Al4V



S

BIOLOX DELTA+Ti6Al4V
Revision Femoral Head
SHORT - Dia. 28mm - 12/14

Revision Head **EN**
Testa da Revisione **IT**
Cabeza para Revisione **ES**
Revisionkopf **DE**
Tete de Revision **FR**
Cabeça para Revisao **PT**

2021-07-31 qty 1 STERILE EO Limacorporate S.p.a.

- Single use
- Do not use if package is damaged
- Keep away from sunlight
- Temperature limit
- Use By
- Reference Number
- Keep dry
- Lot Number
- Caution
- Sterilized using Ethylene Oxide
- Sterilized using irradiation
- IFU Consult instructions for use
- Manufacturer

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE 0123



REF 5010.42.021

BIOLOX DELTA+Ti6Al4V
Revision Femoral Head
SHORT - Dia. 28mm - 12/14



2021-07-31 LOT 2015 80614 STERILE EO 1600209



REF 5010.42.021

BIOLOX DELTA+Ti6Al4V
Revision Femoral Head
SHORT - Dia. 28mm - 12/14



2021-07-31 LOT 2015 80614 STERILE EO 1600209



Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE 0123

SECTION 6 – EXAMPLE OF LABELS FOR HEADS (CoCrMo)

Data Legend

(the information which changes in the other labels of the same device is squared and showed in the table below).

Code
REF 5012.09.221

Taper 14/16
Head - CoCrMo

Diameter
Dia 22

S Size

2019-09-30 qty: 1 STERILE EO Sterilization Limacorporate S.p.a.

Single use Do not use if package is damaged Keep away from sunlight Temperature limit

Use By REF Reference Number Keep dry

LOT Lot Number Caution Sterilized using Ethylene Oxide

IFU Consult instructions for use Manufacturer Sterilized using Irradiation

Limacorporate S.p.A. - Via Nazionale, 52 33038 Villanova di San Daniele - Udine - Italy CE 0123

Code
REF 5012.09.221

Femoral Head/Testa Femorale
CoCrMo Dia=22mm

SHORT Taper/Cono=14/16

Femoral Head EN
Testa Femorale IT
Cabeza Femoral ES
Prothesenkopf DE
Tete Femoral FR
Cabeça Femoral PT

2019-09-30 LOT 2000 00910 STERILE EO 1400249

(20)02(240)9501209221(10)0000910(91)1400249(92)6167
(01)08033390007573(17)190930

Limacorporate S.p.A. - Via Nazionale, 52 33038 Villanova di San Daniele - Udine - Italy CE 0123

Code	Size	Diameter
5012.09.221	S	22
5012.09.222	M	22
5012.09.223	L	22
5012.09.281	S	28
5012.09.282	M	28
5012.09.283	L	28
5012.09.284	XL	28
5012.09.285	XXL	28
5012.09.286	XXXL	28
5012.09.321	S	32
5012.09.322	M	32
5012.09.323	L	32
5012.09.324	XL	32
5012.09.325	XXL	32
5012.09.326	XXXL	32

Example of labels for the first code

REF 5012.09.221

Taper 14/16
Head - CoCrMo



Dia
22
S

Femoral Head/Testa Femorale
CoCrMo Dia=22mm
SHORT Taper/Cono=14/16
Femoral Head **EN**
Testa Femorale **IT**
Cabeza Femoral **ES**
Prothesenkopf **DE**
Tete Femoral **FR**
Cabeça Femoral **PT**

2019-09-30 **qty:** 1 **STERILE** **EO** Limacorporate S.p.a.

- Single use
- Do not use if package is damaged
- Keep away from sunlight
- Temperature limit
- Use By
- Reference Number
- Keep dry
- Lot Number
- Caution
- Sterilized using Ethylene Oxide
- Sterilized using irradiation
- IFU Consult instructions for use
- Manufacturer

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE
0123



REF 5012.09.221

Femoral Head/Testa Femorale
CoCrMo Dia=22mm
SHORT Taper/Cono=14/16



2019-09-30 **LOT** 2000 00910 **STERILE** **EO** 1400249



(20)02(240)95012092210(10)0000910(91)1400249(92)6167



(01)08033390007573(17)190930

REF 5012.09.221

Femoral Head/Testa Femorale
CoCrMo Dia=22mm
SHORT Taper/Cono=14/16



2019-09-30 **LOT** 2000 00910 **STERILE** **EO** 1400249



(20)02(240)95012092210(10)0000910(91)1400249(92)6167



(01)08033390007573(17)190930

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE
0123

SECTION 7 – EXAMPLE OF LABELS FOR BIPOLAR HEADS (AISI 316/L) WITH LINER (UHMWPE)

Data Legend

(the information which changes in the other labels of the same device is squared and showed in the table below).

Code
REF 5030.05.380

For Heads 22
Head - AISI316L/PE

Diameter
Dia
38

Bipolar Head/Testa Biart.-AISI316L
Liner/Inserto - UHMWPE
Dia=38mm for Heads/per Teste 22mm

Bipolar Head
Testa Biarticolare
Cabeza Biarticular
Biertikularkopf
Tete Bipolaire
Cabeça Biarticular

2021-07-31 qty: 1 STERILE EO Limacorporate S.p.a.

Sterilization

Single use
Do not use if package is damaged
Keep away from sunlight
Temperature limit
Use By
Reference Number
Keep dry
Lot Number
Caution
Sterilized using Ethylene Oxide
Sterilized using irradiation
IFU Consult instructions for use
Manufacturer

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE 0123

Code
REF 5030.05.380

Bipolar Head/Testa Biart.-AISI316L
Liner/Inserto - UHMWPE
Dia=38mm for Heads/per Teste 22mm

2021-07-31 LOT 2016 12220 STERILE EO 1600216

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE 0123

Code	Diameter
5030.05.380	38
5030.05.390	39
5030.05.400	40

Example of labels for the first code

REF 5030.05.380

For Heads 22
Head - AISI316LPE



Dia
38

Bipolar Head/Testa Biart.-AISI316L
Liner/Inserto - UHMWPE
Dia=38mm for Heads/per Teste 22mm

Bipolar Head EN
Testa Biarticolare IT
Cabeza Biarticular ES
Biartikularkopf DE
Tete Bipolaire FR
Cabeça Biarticular PT

2021-07-31 **qty:** 1 **STERILE EO** **Limacorporate S.p.a.**

- | | | | |
|----------------------------------|----------------------------------|---------------------------------|-------------------|
| Single use | Do not use if package is damaged | Keep away from sunlight | Temperature limit |
| Use By | Reference Number | Keep dry | |
| Lot Number | Caution | Sterilized using Ethylene Oxide | |
| IFU Consult instructions for use | Manufacturer | Sterilized using irradiation | |

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE
0123



REF 5030.05.380

Bipolar Head/Testa Biart.-AISI316L
Liner/Inserto - UHMWPE

Dia=38mm for Heads/per Teste 22mm

2021-07-31 **LOT** 2016 12220 **STERILE EO** 1600216



(20)00(240)95030053800(10)1612220(9)11600216(92)6185

(01)08033390007818(17)210731



REF 5030.05.380

Bipolar Head/Testa Biart.-AISI316L
Liner/Inserto - UHMWPE
Dia=38mm for Heads/per Teste 22mm

2021-07-31 **LOT** 2016 12220 **STERILE EO** 1600216



(20)00(240)95030053800(10)1612220(9)11600216(92)6185

(01)08033390007818(17)210731




Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE
0123

SECTION 8 – EXAMPLE OF LABELS FOR BIPOLAR HEADS (CoCrMo) WITH LINER (UHMWPE)

Data Legend

(the information which changes in the other labels of the same device is squared and showed in the table below).



Code 5030.09.380

For Heads 22
Head-CoCrMo/UHMWPE

Diameter Dia 38

Dia=38mm for Heads/per Teste 22mm

Bipolar Head/Testa Biart. - CoCrMo
Liner/Inserito - UHMWPE

Bipolar Head
Testa Biarticolare IT
Cabeza Biarticular ES
Biertikularkopf DE
Tete Bipolaire FR
Cabeça Biarticular PT

2022-07-31 [REF] 5030.09.380 [LOT] 2017 09024 [STERILE EO] 1700240

Limacorporate S.p.a.

Sterilization

Single use
Use By
Lot Number
IFU Consult instructions for use

Do not use if package is damaged
Reference Number
Caution
Manufacturer

Keep away from sunlight
Keep dry
Sterilized using Ethylene Oxide
Sterilized using irradiation

Temperature limit

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE 0123

Code	Diameter
5030.09.380	38
5030.09.390	39
5030.09.400	40

Example of labels for the first code

REF 5030.09.380

For Heads 22
Head-CoCrMo/UHMWPE



Dia
38

Bipolar Head/Testa Biart. - CoCrMo
Liner/Inserto - UHMWPE
Dia=38mm for Heads/per Teste 22mm

Bipolar Head EN
Testa Biarticolare IT
Cabeza Biarticular ES
Biartikularkopf DE
Tete Bipolaire FR
Cabeça Biarticular PT

2022-07-31 qty: 1 **STERILE** EO Limacorporate S.p.a.

REF 5030.09.380

Bipolar Head/Testa Biart. - CoCrMo
Liner/Inserto - UHMWPE
Dia=38mm for Heads/per Teste 22mm

2022-07-31 **LOT** 2017 09024 **STERILE** EO 1700240



(20)00(240)95030093800(10)1709024(91)1700240(92)6188
(01)08033390007849(17)220731



- Single use
- Do not use if package is damaged
- Keep away from sunlight
- Temperature limit
- Use By
- Reference Number
- Keep dry
- Lot Number
- Caution
- Sterilized using Ethylene Oxide
- Sterilized using irradiation
- IFU Consult instructions for use
- Manufacturer

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE
0123



REF 5030.09.380

Bipolar Head/Testa Biart. - CoCrMo
Liner/Inserto - UHMWPE
Dia=38mm for Heads/per Teste 22mm

2022-07-31 **LOT** 2017 09024 **STERILE** EO 1700240



(20)00(240)95030093800(10)1709024(91)1700240(92)6188
(01)08033390007849(17)220731



Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE
0123

SECTION 9 – EXAMPLE OF LABELS FOR LOCK SELF-CENTERING HEADS (AISI 316/L) WITH LINER (UHMWPE)

Data Legend

(the information which changes in the other labels of the same device is squared and showed in the table below).

Code	Diameter
5527.05.410	41
5527.05.420	42
5527.05.430	43
5527.05.440	44
5527.05.450	45
5527.05.460	46
5527.05.470	47
5527.05.480	48
5527.05.490	49
Code	Diameter
5527.05.500	50
5527.05.510	51
5527.05.520	52
5527.05.530	53
5527.05.540	54
5527.05.550	55
5527.05.560	56

5527.05.570

57

Example of labels for the first code

REF 5527.05.410

For Heads 28mm
LOCK-AISI316L/UHMWPE



Dia.
41

LOCK Bipolar Head/Testa B.-AISI316L
Liner/Inserto - UHMWPE
Dia=41mm for Heads/per Teste 28mm

Bipolar Head EN
Testa Biarticolare IT
Cabeza Biarticular ES
Biartikularkopf DE
Tete Bipolaire FR
Cabeça Biarticular PT

2022-08-31 **qty:** 1 **STERILE** EO **Limacorporate S.p.a.**

REF 5527.05.410

LOCK Bipolar Head/Testa B.-AISI316L
Liner/Inserto - UHMWPE
Dia=41mm for Heads/per Teste 28mm

2022-08-31 **LOT** 2017 09351 **STERILE** EO 1700252



(20)02(240)95527054100(10)1709351(91)1700252(92)6214



(01)08033390009515(17)220631



- Single use
- Do not use if package is damaged
- Keep away from sunlight
- Temperature limit
- Use By
- Reference Number
- Keep dry
- Lot Number
- Caution
- Sterilized using Ethylene Oxide
- Sterilized using irradiation
- IFU Consult instructions for use
- Manufacturer

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE
0123



REF 5527.05.410

LOCK Bipolar Head/Testa B.-AISI316L
Liner/Inserto - UHMWPE
Dia=41mm for Heads/per Teste 28mm


2022-08-31 **LOT** 2017 09351 **STERILE** EO 1700252



(20)02(240)95527054100(10)1709351(91)1700252(92)6214



(01)08033390009515(17)220631




Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE
0123

SECTION 10 – EXAMPLE OF LABELS FOR LOCK SELF-CENTERING HEADS (CoCrMo) WITH LINER (UHMWPE)

Data Legend

(the information which changes in the other labels of the same device is squared and showed in the table below).



Code
REF 5527.09.410

For Heads 28mm
LOCK - CoCrMo/UHMWPE

Diameter
Dia. 41

LOCK Bipolar Head/Testa B.-CoCrMo
Liner/Inserto - UHMWPE
Dia=41mm for Heads/per Teste 28mm

Bipolar Head EN
Testa Biarticolare IT
Cabeza Biarticular ES
Biartikularkopf DE
Tete Bipolaire FR
Cabeça Biarticular PT

2022-07-31 qty 1 STERILE EO Limacorporate S.p.a.

STERILIZATION

Single use Do not use if package is damaged Keep away from sunlight Temperature limit
Use By REF Reference Number Keep dry
LOT Lot Number Caution Sterilized using Ethylene Oxide
IFU Consult instructions for use Manufacturer Sterilized using irradiation

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE 0123

REF 5527.09.410
LOCK Bipolar Head/Testa B.-CoCrMo
Liner/Inserto - UHMWPE
Dia=41mm for Heads/per Teste 28mm

2022-07-31 LOT 2017 06601 STERILE EO 1700226

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE 0123

Code	Diameter
5527.09.410	41
5527.09.420	42
5527.09.430	43
5527.09.440	44
5527.09.450	45
5527.09.460	46
5527.09.470	47
5527.09.480	48
5527.09.490	49
5527.09.500	50
5527.09.510	51
5527.09.520	52
5527.09.530	53
5527.09.540	54
5527.09.550	55
5527.09.560	56
5527.09.570	57

Example of labels for the first code

REF 5527.09.410

For Heads 28mm
LOCK - CoCrMo/UHMWPE



Dia.
41

LOCK Bipolar Head/Testa B.-CoCrMo
Liner/Inserto - UHMWPE
Dia=41mm for Heads/per Teste 28mm

Bipolar Head EN
Testa Biarticolare IT
Cabeza Biarticular ES
Biartikularkopf DE
Tete Bipolaire FR
Cabeça Biarticular PT

2022-07-31 **qty:** 1 **STERILE** EO **Limacorporate S.p.a.**

REF 5527.09.410

LOCK Bipolar Head/Testa B.-CoCrMo
Liner/Inserto - UHMWPE
Dia=41mm for Heads/per Teste 28mm

2022-07-31 **LOT** 2017 06601 **STERILE** EO 1700226



(20)02(240)95527094100(10)1706601(91)1700226(92)6231



(01)08033390009690(17)220731



- Single use
- Do not use if package is damaged
- Keep away from sunlight
- Temperature limit
- Use By
- Reference Number
- Keep dry
- Lot Number
- Caution
- Sterilized using Ethylene Oxide
- IFU Consult instructions for use
- Manufacturer
- Sterilized using irradiation

Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE
0123



REF 5527.09.410

LOCK Bipolar Head/Testa B.-CoCrMo
Liner/Inserto - UHMWPE
Dia=41mm for Heads/per Teste 28mm

2022-07-31 **LOT** 2017 06601 **STERILE** EO 1700226



(20)02(240)95527094100(10)1706601(91)1700226(92)6231



(01)08033390009690(17)220731



Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE
0123

SECTION 11 – EXAMPLE OF LABELS FOR FEMORAL CEPHALIC HEADS (AISI 316/L) - TAPER 12/14

Data Legend

(the information which changes in the other labels of the same device is squared and showed in the table below).

Code 5061.05.410

Diameter Dia 41

REF 5061.05.410
Femoral Head/Testa Femorale for Endoprosthesis - AISI316L
Dia=41mm Taper/Cono=12/14

2022-05-31 | 1 | STERILE R | Limacorporate S.p.a.

STERILIZATION

Single use | Do not use if package is damaged | Keep away from sunlight | Temperature limit

Use By | Reference Number | Keep dry

Lot Number | Caution | Sterilized using Ethylene Oxide

IFU Consult instructions for use | Manufacturer | Sterilized using irradiation

Limacorporate S.p.A. - Via Nazionale, 52 33038 Villanova di San Daniele - Udine - Italy

CE 0123

Code	Diameter
5061.05.410	41.0
5061.05.420	42.0
5061.05.430	43.0
5061.05.445	44.5
5061.05.460	46.0
5061.05.475	47.5
5061.05.490	49.0
5061.05.510	51.0
5061.05.525	52.5
5061.05.540	54.0
5061.05.555	55.5
5061.05.570	57.0

Example of labels for the first code

REF 5061.05.410

Taper 12/14
Head - AISI316L



Dia
41

Femoral Head/Testa Femorale
for Endoprosthesis - AISI316L
Dia=41mm Taper/Cono=12/14

Femoral Head EN
Testa Femorale IT
Cabeza Femoral ES
Prothesenkopf DE
Tete Femoral FR
Cabeça Femoral PT

2022-05-31 **qty:** 1 **STERILE R** Limacorporate S.p.a.


REF 5061.05.410

Femoral Head/Testa Femorale
for Endoprosthesis - AISI316L
Dia=41mm Taper/Cono=12/14

2022-05-31 **LOT** 2016 14263 **STERILE R** 1700160



(20)00(240)95061054100(10)1614263(91)1700160(92)6191
(01)08033390008327(17)220531



Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE
0123



REF 5061.05.410

Femoral Head/Testa Femorale
for Endoprosthesis - AISI316L
Dia=41mm Taper/Cono=12/14

2022-05-31 **LOT** 2016 14263 **STERILE R** 1700160



(20)00(240)95061054100(10)1614263(91)1700160(92)6191
(01)08033390008327(17)220531



Limacorporate S.p.A. - Via Nazionale, 52
33038 Villanova di San Daniele - Udine - Italy

CE
0123