


Electrodes réutilisables en acier inoxydable

Référence Produit	Description
	Stainless Steel Reusable Electrode  Please note: These electrodes require connection cable CB112-TC
CSK-TC5	Stainless Steel Reusable Electrode – 5 cm
CSK-TC10	Stainless Steel Reusable Electrode – 10 cm
CSK-TC15	Stainless Steel Reusable Electrode – 15 cm
CSK-TC20	Stainless Steel Reusable Electrode – 20 cm 
CB112-TC	Connection Cable for Reusable Electrodes
	Stainless Steel Reusable Electrode with integrated 3m Cable 
CSK-TC5-3M	Stainless Steel Reusable Electrode – 5cm with integrated 3m cable
CSK-TC10-3M	Stainless Steel Reusable Electrode – 10cm with integrated 3m cable
CSK-TC15-3M	Stainless Steel Reusable Electrode – 15cm with integrated 3m cable
CSK-CASE	Sterilizing Case for CSK Electrode 

Electrodes Jetables

Référence Produit	Description
	Disposable Electrode <i>Package Contents</i> (10) Disposable Electrodes Please note: These electrodes require connection cable CB114-TC 
TCD-5-P	Disposable Electrode – 5 cm
TCD-10-P	Disposable Electrode – 10 cm
TCD-15-P	Disposable Electrode – 15 cm
TCD-20-P	Disposable Electrode – 20 cm
	Disposable Electrode with integrated 3m Cable <i>Package Contents</i> (10) Disposable Electrodes Please note: These electrodes require adapter CB114-ADP or connection cable CB114-TC 
TCD-5-3M-P	Disposable Electrode – 5 cm with integrated 3m cable
TCD-10-3M-P	Disposable Electrode – 10 cm with integrated 3m cable
TCD-15-3M-P	Disposable Electrode – 15 cm with integrated 3m cable
CB114-ADP	Adapter Plug for Disposable Electrodes 
CB114-TC	Connection Cable for Disposable Electrodes 


Electrode et Canule Unifiées

Référence Produit	Déscription
	Unified Disposable Electrode/Cannula <i>Package Contents</i> (10) Unified Disposable Electrode/Cannula Please note: All CU/CUR electrodes require connection cable CB114-TC 
	Unified Disposable Electrode/Cannula Straight Tip 
CU-6522-P	Length: 6 cm – Tip: 5 mm – Gauge: 22ga
CU-61022-P	Length: 6 cm – Tip: 10 mm – Gauge: 22ga
CU-10522-P	Length: 10 cm – Tip: 5 mm – Gauge: 22ga
CU-101018-P	Length: 10 cm – Tip: 10 mm – Gauge: 18ga
CU-101022-P	Length: 10 cm – Tip: 10 mm – Gauge: 22ga
CU-15522-P	Length: 15 cm – Tip: 5 mm – Gauge: 22ga
	Unified Disposable Electrode/Cannula Straight Tip – Echo RF™ 
CU-6520U-P	Length: 6 cm – Tip: 5 mm – Gauge: 20ga
CU-10520U-P	Length: 10 cm – Tip: 5 mm – Gauge: 20ga
CU-101020U-P	Length: 10 cm – Tip: 10 mm – Gauge: 20ga
CU-15520U-P	Length: 15 cm – Tip: 5 mm – Gauge: 20ga
CU-151020U-P	Length: 15 cm – Tip: 10 mm – Gauge: 20ga


Electrode et Canule Unifiées

Référence Produit	Description	
	Unified Disposable Electrode/Cannula Curved Tip	
CUR-5616-P	Length: 5 cm – Tip: 6 mm – Gauge: 16ga	
CUR-51016-P	Length: 5 cm – Tip: 10 mm – Gauge: 16ga	
CUR-6518-P	Length: 6 cm – Tip: 5 mm – Gauge: 18ga	
CUR-101016-P	Length: 10 cm – Tip: 10 mm – Gauge: 16ga	
CUR-101018-P	Length: 10 cm – Tip: 10 mm – Gauge: 18ga	
CUR-101022-P	Length: 10 cm – Tip: 10 mm – Gauge: 22ga	
CUR-151016-P	Length: 15 cm – Tip: 10 mm – Gauge: 16ga	
CUR-151018-P	Length: 15 cm – Tip: 10 mm – Gauge: 18ga	
	Unified Disposable Electrode/Cannula Curved Tip – Echo RF™	
CUR-101020U-P	Length: 10 cm – Tip: 10 mm – Gauge: 20ga	
CUR-151020U-P	Length: 15 cm – Tip: 10 mm – Gauge: 20ga	
CB114-TC	Connection Cable for Disposable Electrodes	

Aiguilles et Electrodes d'injection jetables

Référence Produit	Description
	Voltage Controlled Injection Electrode (Non-temperature Controlled) <i>Package Contents</i> (10) Voltage Controlled Injection Electrodes Please note: These electrodes require connection cable CB116 or CB116-4 
CR-6-P	Length: 6 cm – Tip: 5 mm
CR-6522U-P	Length: 6 cm – Tip: 5 mm – Echo RF
CR-10-P	Length: 10 cm – Tip: 5 mm
CR-10522U-P	Length: 10 cm – Tip: 5 mm – Echo RF
CR-15-P	Length: 15 cm – Tip: 5 mm
	Stimulation Guided Nerve Block Needle <i>Package Contents</i> (10) Stimulation Guided Nerve Block Needles Please note: These needles require connection cable CB116 or CB116-4 
CP-6-P	Length: 6 cm – Tip: 1 mm
CP-10-P	Length: 10 cm – Tip: 1 mm
	Disposable Nerve Block Needle <i>Package Contents</i> (10) Disposable Nerve Block Needles 
CN-6-P	Length: 6 cm
CN-10-P	Length: 10 cm
CN-15-P	Length: 15 cm


Aiguilles et Electrodes d'injection jetables

Référence Produit	Description	
	Disposable RF Injection Electrode Connection Cables	
CB116	Single Injection Electrode Cable	
CB116-4	Four Injection Electrode Cable	
CR-CO	Injection Electrode Organizer Clip	
CB103-1A	Injection Electrode to Alligator Clip Adapter	

Procedure Kits Jetables

PRODUCT NUMBER	DESCRIPTION	
BDK-15	Bipolar Disc Kit <i>Package Contents</i> (2) CC152020 Straight Sharp Cannulae (2) TCD-15 15cm Disposable Electrodes Please note: These electrodes require connection cable CB114-TC	
CB114-TC	Connection Cable for Disposable Electrodes	

Plaque de mise à la terre jetables

Référence Produit	Description
	Disposable Ground Pads 
DGP-PM-10	<i>Package Contents</i> (10) Disposable Ground Pads
DGP-PM-25	<i>Package Contents</i> (25) Disposable Ground Pads
	Disposable Cordless Ground Pads  Please note: These ground pads require connection cable CB105
DGP-PMC-10	<i>Package Contents</i> (10) Disposable Cordless Ground Pads
DGP-PMC-25	<i>Package Contents</i> (25) Disposable Cordless Ground Pads
DGP-PMC-50	<i>Package Contents</i> (50) Disposable Cordless Ground Pads
CB105	Connection Cable for Cordless Ground Pads 
CB103-B	Ground Pad Cable for Non-Disposable Ground Plate 

Kits et bloc de guidage Palisade™ jetables

Référence Produit	Description
	Disposable Palisade Guideblock <i>Package Contents</i> (10) Disposable Palisade Guideblocks 
PG-1018-P	Disposable Palisade Guideblock for 18ga Cannulae
PG-1020-P	Disposable Palisade Guideblock for 20ga Cannulae
	Disposable Palisade Guideblock Kit <i>Package Contents</i> (6) Straight Sharp Cannulae (4) TCD-10 Disposable Electrodes  These electrodes require connection cable CB114-TC Disposable Palisade Guideblock to be ordered separately
PALISADE-18	Included Cannulae: CC101018 - Length: 10 cm – Tip: 10 mm – Gauge: 18ga
PALISADE-20	Included Cannulae: CC101020 - Length: 10 cm – Tip: 10 mm – Gauge: 20ga
CB114-TC	Connection Cable for Disposable Electrodes 



Electrodes Stéréotaxiques Jetables

Référence Produit	Description
	Guide/Reducing Tubes <i>Package Contents</i> (1) Guide/Reducing Tube 
GTRT-11-76	Guide/Reducing Tube for 1.1 mm Electrodes – Length 76 mm
GTRT-16-76	Guide/Reducing Tube for 1.6 mm Electrodes – Length 76 mm
	Disposable Depth Stops 
DS-11-D-P	(10) Disposable Depth Stops for 1.1 mm Electrodes
CB114-TC	Connection Cable for Disposable Electrodes 

Kits Electrodes et Canules

Référence Produit	Description
FLEX-K151320-5P	Flextrode™ Introducer Kit <i>Package Contents</i> (5) Flextrode Disposable Electrodes Tip: 13 mm – Gauge: 20ga (5) Flextrode Cannulae Tip: 15 cm – Gauge: 17ga  Please note: The Flextrode Introducer Kit requires connection cable CB114-TC
	LCED Disposable Cordotomy Kit <i>Package Contents</i> (1) LCED-TC Electrode (1) BD Spinal Needle  Please note: The Cordotomy Kit requires connection cable CB114-TC
LCED	(1) Disposable Cordotomy Electrode Kit
LCED-P	(10) Disposable Cordotomy Electrode Kits
	Reig-Cosman Epidural (RCE) Pulsed RF Electrode <i>Package Contents</i> (10) RCE Disposable Electrodes Length: 40 cm Length: 66 cm  Please note: The Pulsed Electrode Kit requires connection cable CB114-TC
RCE-E401519-P	(10) RCE 40 cm Electrodes
RCE-E661519-P	(10) RCE 66 cm Electrodes
RCE-C916S-P	(10) RCE Introducer Cannulae
CB-114TC	Connection Cable for Disposable Electrodes Please note: The Bipolar Disc Kit requires connection cable CB114-TC

Connexion - Câbles et Adaptateurs

Référence Produit	Description	
	Connection Cables and Adapters for Boston Scientific (Cosman) Generators	
CB112-TC	Connection Cable for Reusable Electrodes	
CB114-TC	Connection Cable for Disposable Electrodes	
CB114-ADP	Adapter Plug for Disposable Electrodes	
CB105	Connection Cable for Cordless Ground Pads	
CB116	Single Injection Electrode Cable	
CB116-4	Four Injection Electrode Cable	
CR-CO	Injection Electrode Organizer Clip	
CB103-1A	Injection Electrode to Alligator Clip Adapter	
CB103-B	Ground Pad Cable for Non-Disposable Ground Plate	
CB103-R	Alligator Clip Connection Cable	
CB113	RFG-1A/1B Bipolar Adapter Cable for Ground Pad Jack (Reusable Electrodes)	
CB130-TC	Neurotherm Simplicity Electrode to Cosman Generator Cable	
CB135-TC	Neurotherm Reusable Electrode to Cosman Generator Cable	
CB-SN	Smith & Nephew Spinecath Electrode to Cosman Generator Cable	

Connexion - Câbles et Adaptateurs

Référence Produit	Description	
	Connection Cables and Adapters for Major Manufacturers Generators <i>Please note: Some adapters may require an additional connection cable for use with some generators</i>	
C125-TC	Adapter for Stryker Generators	
C120-TC	Adapter Cable for Halyard Single-Output Generators (Reusable Electrodes)	
C121-TC	Adapter Cable for Halyard Single-Output Generators (Disposable Electrodes)	
C124-TC	Adapter Cable for Halyard Four-Output Generators (Reusable Electrodes)	
C134-TC	Adapter Cable for Halyard Four-Output Generators (Disposable Electrodes)	
C126-TC	Adapter Cable for Smith & Nephew Generators (Reusable Electrodes)	
C127-TC	Adapter Cable for Smith & Nephew Generators (Disposable Electrodes)	
C115-TC	Adapter Cable for Neurotherm Generators (Reusable Electrodes)	
C123-TC	Adapter Cable for Neurotherm Generators (Disposable Electrodes)	
C110	Adapter Cable for Neurotherm Generators (Injection Electrodes)	
C112-TC	Adapter Cable for Radionics Generators (Reusable Electrodes)	
C114-TC	Adapter Cable for Radionics Generators (Disposable Electrodes)	
C116	Adapter Cable for Radionics Generators (Injection Electrodes)	
C116-4	Adapter Cable for Radionics Generators (Injection Electrodes x 4)	
C105	Ground Pad Adapter Cable for Halyard/Radionics Generators (DGP-PMC)	
C103-RTM	Radionics Thermistor Cable to Alligator Clip	
C104-TM	Radionics Thermistor Cable for Thermistor Electrode	