


Thrombektomie Instrument - Dispositif de thrombectomie - Dispositivo per trombectomia - Instrumento para trombectomia  
- Instrumento para trombectomia - инструмент для тромбэктомии

phenox GmbH  
Lise-Meitner-Allee 31 - 44801 Bochum, Germany  
Tel: +49 (0) 234 / 36919-0 - Fax: +49 (0) 234 / 36919-19  
www.phenox.net

**4 mm** x **20 mm**  
Shaft Ø Working Length

**ID 0.021"** **2 mm**  
Microcatheter Min. Vessel Ø

## LUX



A: Shaft Diameter [mm]	B: Shaft Length [mm]	C: Insertion Wire Length [mm]	D: Marker Distance [mm]	E: Working Length [mm]	Microcatheter Compatibility [Inch]	Min. Vessel Diameter [mm]
4	31	1800	1260	20	ID 0.021"	2

REF PRE-LUX-4-20 LOT 123456 2021-12 CONT 1 x pRESET LUX

pRESET LUX REF PRE-LUX-4-20 LOT 123456



(01)04260123782481(17)211231(10)123456

pRESET LUX phenox GmbH: www.phenox.net



REF PRE-LUX-4-20

LOT 123456 2021-12

(01)04260123782481(17)211231(10)123456

pRESET LUX phenox GmbH: www.phenox.net



REF PRE-LUX-4-20

LOT 123456 2021-12

(01)04260123782481(17)211231(10)123456

pRESET LUX phenox GmbH: www.phenox.net



REF PRE-LUX-4-20

LOT 123456 2021-12

(01)04260123782481(17)211231(10)123456

pRESET LUX phenox GmbH: www.phenox.net



REF PRE-LUX-4-20

LOT 123456 2021-12

(01)04260123782481(17)211231(10)123456

CE 0297



globe icon



Laser-free



MAT

NITI



STERILE EO

Sterilized with Ethylene Oxide Gas



Consult instructions for use



Caution



Do not reuse



Keep dry



Keep away from sunlight



Do not use if the packaging is damaged



Do not sterilize



Non-synaptic

LAB-8849-PRESET-LUX B - 2019/01



REF PRE-LUX-4-20

LOT 123456 2021-12