

p64 MW HPC phenox Flow Modulation Device (pFMD)

P64-MW-HPC-450-21

Flow Modulation Device - Vorrichtung zur Blutstrommodulation - Dispositif de modulation du débit - Dispositivo modulatore di flusso - Dispositivo de modulación de flujo - Dispositivo de modulação de fluxo - Устройство модуляции потока

phenox GmbH
Lise-Meitner-Allee 31 - 44801 Bochum, Germany
Tel: +49 (0) 234 / 36919-0 - Fax: +49 (0) 234 / 36919-19
www.phenox.net

4.5 mm x 21 mm

Max. Vessel Ø


Length in Max. Vessel

ID 0.021" 4 mm

Microcatheter

Min. Vessel

HPC - Hydrophilic Polymer Coating



	Vessel Diameter [mm]	A: Length [mm]		A: Length [mm]	B: Device Diameter [mm]
Min.	4	31.5	Unconstrained	8.9	4.9
Max.	4.5	21	In Microcatheter ID 0.021"	51	0.53

REF **P64-MW-HPC-450-21** LOT 123456 2021-12 CONT 1 x p64 MW HPC

p64 MW HPC

REF **P64-MW-HPC-450-21**

LOT 123456



(01)04260123780128(17)211231(10)123456

p64 MW HPC

phenox GmbH: www.phenox.net



REF **P64-MW-HPC-450-21**

LOT 123456

2021-12

(01)04260123780128(17)211231(10)123456

p64 MW HPC

phenox GmbH: www.phenox.net



REF **P64-MW-HPC-450-21**

LOT 123456

2021-12

(01)04260123780128(17)211231(10)123456

p64 MW HPC

phenox GmbH: www.phenox.net



REF **P64-MW-HPC-450-21**

LOT 123456

2021-12

(01)04260123780128(17)211231(10)123456

p64 MW HPC

phenox GmbH: www.phenox.net



REF **P64-MW-HPC-450-21**

LOT 123456

2021-12

(01)04260123780128(17)211231(10)123456

CE 0297



Phthalate-free



Latex-free



MAT NiTi



Sterilized with Ethylene Oxide Gas



Consult instructions for use



Caution



Do not reuse



Keep dry



Keep away from sunlight



Do not use if the package is damaged



Do not resterilize



Non-pyrogenic

LAB-B865-P64-MW-HPC A - 2019/07

p64 MW HPC phenox Flow Modulation Device (pFMD)

4.5 mm x 21 mm



REF **P64-MW-HPC-450-21**

LOT 123456

2021-12